

Chocolate™

PTA Balloon Catheter

Minimize vessel trauma, dissections and the need for bailout stenting above or below the knee with the Chocolate™ PTA Balloon. The balloon’s unique nitinol-constraining structure creates pillows and grooves that provide a predictable, uniform, and atraumatic dilatation.

- Minimize vessel trauma
- Reduce dissections
- Decrease bailout stenting

Technical specifications

Catheter design	Over-the-wire
Usable catheter lengths	120, 135, and 150 cm
Introducer sheath compatibility	5 and 6 F
Guidewire compatibility	0.014” and 0.018”



Chocolate™

PTA Balloon Catheter

Home

Drug-coated
Balloons

Stent Systems

PTA Catheters

Directional
Atherectomy

Embolic Protection
Devices




Crossing
Catheters

Thrombus
Management

Guidewires

Snares

Ellipsys™



Order information

Product catalogue number	Balloon diameter (mm)	Balloon length (mm)	Catheter length (cm)	Guidewire (in)	Introducer sheath (F)	Nominal pressure (atm)	Rated burst pressure (atm)
CB1415025040 OTW	2.5	40	150	0.014	5	9	14
CB1415025080 OTW	2.5	80	150	0.014	5	9	14
CB1415025120 OTW	2.5	120	150	0.014	5	9	14
CB1415030040 OTW	3.0	40	150	0.014	5	9	14
CB1415030080 OTW	3.0	80	150	0.014	5	9	14
CB1415030120 OTW	3.0	120	150	0.014	5	9	14
CB1413535040 OTW	3.5	40	135	0.014	5	9	14
CB1413535080 OTW	3.5	80	135	0.014	5	9	14
CB1413535120 OTW	3.5	120	135	0.014	5	9	14
CB1413540040 OTW	4.0	40	135	0.014	5	9	14
CB1413540080 OTW	4.0	80	135	0.014	5	9	14
CB1413540120 OTW	4.0	120	135	0.014	5	9	14
CB1812050040 OTW	5.0	40	120	0.018	6	6	12
CB1812050080 OTW	5.0	80	120	0.018	6	6	12
CB1812050120 OTW	5.0	120	120	0.018	6	6	12
CB1812060040 OTW	6.0	40	120	0.018	6	6	12
CB1812060080 OTW	6.0	80	120	0.018	6	6	12
CB1812060120 OTW	6.0	120	120	0.018	6	6	12