


# INTRODUCERS

Comprehensive product line



- Excellent transitions
- Hemostasis valve
- Color-coded hubs
- Variety of sizes

BALTON  
since  
1980

**BALTON**®

CE 0197

# INTRODUCERS

## INTRODUCER

**Set with hemostasis valve to introduce and exchange catheters and endocavitary electrodes**

Set contains:

- sheath with a valve
- angiographic needle
- guide wire
- dilatator

| Catalog Number | Sheath size / length | Colour identification | Dilatator size | Angiographic needle size | Guide wire size / length / type |
|----------------|----------------------|-----------------------|----------------|--------------------------|---------------------------------|
| INT3F          | 3F / 60 mm           | white                 | 3F             | 21G                      | .018" / 35 cm / STRAIGHT        |
| INT4F          | 4F / 70 mm           | red                   | 4F             | 19G                      | .022" / 40 cm / "J"             |
| INT5F          | 5F / 110 mm          | grey                  | 5F             | 18G                      | .035" / 40 cm / "J"             |
| INT5F23        | 5F / 230 mm          | grey                  | 5F             | 18G                      | .035" / 60 cm / "J"             |
| INT6F          | 6F / 110 mm          | green                 | 6F             | 18G                      | .035" / 40 cm / "J"             |
| INT6F23        | 6F / 230 mm          | green                 | 6F             | 18G                      | .035" / 60 cm / "J"             |
| INT7F          | 7F / 110 mm          | orange                | 7F             | 18G                      | .035" / 40 cm / "J"             |
| INT7F23        | 7F / 230 mm          | orange                | 7F             | 18G                      | .035" / 60 cm / "J"             |
| <b>INT8F</b>   | <b>8F</b> / 110 mm   | blue                  | 8F             | 18G                      | .035" / 40 cm / "J"             |
| <b>INT8F23</b> | <b>8F</b> / 230 mm   | blue                  | 8F             | 18G                      | .035" / 60 cm / "J"             |
| INT9F          | 9F / 110 mm          | black                 | 9F             | 18G                      | .035" / 40 cm / "J"             |
| INT9F23        | 9F / 230 mm          | black                 | 9F             | 18G                      | .035" / 60 cm / "J"             |
| INT10F         | 10F / 110 mm         | white                 | 10F            | 18G                      | .035" / 40 cm / "J"             |
| INT10F23       | 10F / 230 mm         | white                 | 10F            | 18G                      | .035" / 60 cm / "J"             |
| INT11F         | 11F / 110 mm         | blue                  | 11F            | 18G                      | .035" / 40 cm / "J"             |
| INT11F23       | 11F / 230 mm         | blue                  | 11F            | 18G                      | .035" / 60 cm / "J"             |

When ordering use Ref. Number, e.g.: INT6F.

Special order can contain introducer with radiopaque marker (available for sizes greater than 4F).

Model reference number: **INT8FM** – introducer with a marker.

### Introducers are intended for introducing catheters and endocavitary electrodes.

The introducer's sheath is delivered with hemostasis valve, which prevents from losing blood during the surgery and it minimizes the embolism manifestation possibility.

The connection of the drain with sheath's head enables the infusion of fluids during the electrode placement procedure.

#### Caution:

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information on herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time. This publication is only a commercial information and cannot be handled the commercial offer in the meaning of the art. 66 §1 of Polish Civil Code.



294 Modlinska Str., 03-152 Warsaw, POLAND  
ph.: (+48) 22 597 44 00, fax: (+48) 22 597 44 44  
e-mail: [balton@balton.pl](mailto:balton@balton.pl), [www.balton.pl](http://www.balton.pl)