

GENTA 1&3



SURGICAL CEMENTS with antibiotic

ORCEM® 1G & ORCEM 3G are radiopaque, sterile and self hardening bone cements.

- They contain an aminoside antibiotic (Gentamicine)
- Both are dedicated to arthroplasty surgery
- They allow immediate and stable bone anchorage of implants

Cements ideally respond to the needs of the surgeon, two viscosities for two applications:

OR**CEM**® 1G: standard viscosity bone cement, manual application.

ORCEM® 3G: low viscosity bone cement, with optimum fluidity to allow an easy application with a syringe or injection gun.

HEADQUARTER & INTERNACIONAL: Avda. Leonardo Da Vinci, 12-14 · Parque Tecnológico · 46980 Paterna

SPAIN SALES DEPARTMENT:C/Ignasi Iglesias, 70 · 08950 ·
Esplugues de Llobregat · Barcelona





GENTA 1&3

New generation bone cements advantages

- Mechanical properties above ISO 5833 standard
- Low exothermicity release in situ
- Less release of monomer
- A stable fixation of implants
- · Barium sulfate, an excellent radiopaque agent
- Easily prepared with a bowl & spatula or in a vacuum mixing system


ORCEM® 1G & 3G bone cements allow local antibiotherapy.


Gentamicin, issued from family of aminosides, has a broad spectrum and is effective against most of the bacteria encountered.


The quantity of gentamicin contained in one dose has been optimised to guarantee a long term release


and then prevent infections.

Mechanical properties higher than ISO 5833 standard:









Indications

Fixation of implants, arthroplasty revision procedures, patients with high risk of bone infection, revision of infected prosthesis.

Storage & Products

Bone cements are delivered sterile, ready to use in the operating theatre. They must be stored in their unopened original packaging in a clean, dark and dry place at a maximum temperature of 25°C. For any further information, please refer to the IFUs.

PACKAGING

Liquid: within glass ampoule Powder: within double pouch

DESIGNATION REFERENCE

| ORCEM® 1G - 40 g of powder and 14.4 g of liquid | F9001011 |
|---|----------|
| ORCEM® 3G - 40 g of powder and 16.4 g of liquid | F9001013 |
| BOWL MIXER for bone cement* | F9000055 |
| SYRINGE MIXER & GUN for bone cement* | F9000921 |







