





Safe and easy insertion.

The introduction of the spray pipe in the endoscope is easy and safe. Thanks to the kink-resistant material the catheter adapts at each turn of the endoscope. To use the staining solution the catheter is equipped with a Y-cannula.

| | Model | Diameter of jaws (mm) | Length (mm) |
|---|-----------------------|-----------------------|-------------|
|  | Straight shape | | |
| | ATE-PSG-Z-10x500 | 1.0 | 500 |
| | ATE-PSG-Z-10x1200 | 1.0 | 1,200 |
| | ATE-PSG-Z-18x1200 | 1.8 | 1,200 |
| | ATE-PSG-Z-18x1800 | 1.8 | 1,800 |
| | ATE-PSG-Z-23x2300 | 2.3 | 2,300 |
|  | Mist | | |
| | ATE-PSG-W-10x500 | 1.0 | 500 |
| | ATE-PSG-W-10x1200 | 1.0 | 1,200 |
| | ATE-PSG-W-18x1200 | 1.8 | 1,200 |
| | ATE-PSG-W-18x1800 | 1.8 | 1,800 |
| | ATE-PSG-W-23x2300 | 2.3 | 2,300 |