





LUBRICANT GEL
WITH LIDOCAINE

OneGel

LUBRICANT GEL WITH LIDOCAINE (STERILE)

USAGE AREAS

- ONEGEL is a sterile, lubricating gel used to lubricate catheters and other medical devices during urethral, colon and rectal applications such as urethral catheterization, endoscopy and cystoscopy.
- It can be also used for rectal and colonic applications as a lubricant gel.
- The slider effect of the gel; Helps prevent trauma between urethral, colon and rectal mucosa and catheter/other medical devices. In this way, the patient is relieved and minimizes spasticity or discomfort caused by iatrogenic injury.



PRODUCT PROPERTIES

- Disposable.
- Soluable in water.
- Easy to use thanks to the syringe and accordion tube.
- Eliminates the risk of damage to tissue.
- Contains antiseptic Chlorhexidine Gluconate to eliminate the risk of infection.
- Lidocaine Hydrochloride is a local anesthetic agent used to prevent pain.
- It provides a painless catheterization with its anesthetic effect.
- Thanks to the lubricating feature provided by hydroxyethylcellulose, the patient relaxes and minimizes the spasticity and discomfort caused by iatrogenic injury.



Contains 2% Lidocaine



Anesthetic



Antiseptic

	ONEGEL 6ML	ONEGEL 11ML	ONEGEL 12,5GR
Product Reference Code	KAF G27-6	KAF G27-11	KAF G27-12,5
Packaging Type	Syringe	Syringe	Accordion Tube
GMDN	37717	37717	37717
Packaging	25 x 6ml / Box	25 x 11ml / Box	25 x 12,5gr / Box





USAGE

- Remove the syringe/accordion tube (6ml/11ml/12.5g) by tearing off the sterile packaging.
- Remove the stopper at the end of the syringe/accordion tube.
- Pour a drop of gel to facilitate application.
- After placing the syringe tip on the area to be lubricated, gently squeeze some ONEGEL Lidocaine Lubricating Gel by pressing the plunger of the syringe/accordion tube.



COMPOSITION





- Pure Water
- Propylene Glycol, Hydroxyethylcellulose (Lubricant)
- Lidocaine Hydrochloride (Local Anesthetic)
- Chlorhexidine Gluconate - 20% Concentration (Antiseptic)
- Methyl Hydroxybenzoate (Preservative)
- Propyl Hydroxybenzoate (Preservative)

It is a Class III Medical Device within
the scope of 93 / 42 / EEC Medical
Device Regulation



STERILE





KAFGRUP
SAĞLIK HİZMETLERİ

KAF GRUP SAĞLIK HİZMETLERİ İNŞ. SAN. VE TİC. LTD. ŞTİ.
www.kafgrup.com | www.wancare.com.tr | info@kafgrup.com | wan.care & kafgrup



REGISTRUL DE STAT AL DISPOZITIVELOR MEDICALE

Nr	Denumire	Den.comerc.	Model	Nr. catalog	Tara	Producatorul	Reprezentant	Ordin	Data
						KAF	sofrag		
DM000368066	GEL LUBRIFIANT UREATRAL CU LIDOCAINĂ	ONEGEL	11 ML	KAF G27-11	Turcia	KAF GRUP SAG. HIZM. INS. SAN. VE TIC. LTD. STI.	SOFRAGRUP S.R.L.	Rg04-000198	24-0
DM000368065	GEL LUBRIFIANT UREATRAL CU LIDOCAINĂ	ONEGEL	12.5 ML	KAF G27-12.5	Turcia	KAF GRUP SAG. HIZM. INS. SAN. VE TIC. LTD. STI.	SOFRAGRUP S.R.L.	Rg04-000198	24-0
DM000368064	GEL LUBRIFIANT UREATRAL CU LIDOCAINĂ	ONEGEL	6 ML	KAF G27-6	Turcia	KAF GRUP SAG. HIZM. INS. SAN. VE TIC. LTD. STI.	SOFRAGRUP S.R.L.	Rg04-000198	24-0

☒ Содержит('Reprezentant'), 'sofrag') И Содержит('Producatorul'), 'KAF')