

# TERUMOBCT

## 40814

Anticoagulant Citrate Dextrose Solution Ph Eur (ACD) Solution A 500 mL


## 40818

Anticoagulant Citrate Dextrose Solution Ph Eur (ACD) Solution A 750 mL

### Ingredient Information




Component	Percent (%w/v)	grams/100 mL	CAS Number
Citric Acid Monohydrate (Acidum Citricum Monohydricum Ph Eur)	0.80%	0.80 g	5949-29-1
Sodium Citrate Dihydrate (Natrii Citras Dihydricum Ph Eur)	2.20%	2.20 g	6132-04-3
Dextrose Monohydrate (Glucosum Monohydricum Ph Eur)	2.45%	2.45 g	14431-43-7
Water for Injection (Aqua ad Injectabilia Ph Eur)	—	—	—

### General Specifications





Catalog Numbers	40814, 40818
Legal Manufacturer	 Terumo BCT, Inc., Lakewood, CO
Manufacturing Location	Larne, UK
Declaration of Conformity	The conformity assessment was done in conformance with the procedures set out in Annex II, Section 3.2, Full Quality Assurance System of the Council Directive 93/42/EEC as amended, under the supervision of BSI.
Referenced Standards	See applicable standards document. <ul style="list-style-type: none"><li>• EN ISO 10993-1:2009; 10993-2:2006; 10993-3:2003/(R)2009; 10993-4:2009; 10993-5:2009; 10993-10:2010; 10993-11:2009; 10993-12:2009; 10993-17:2002; 10993-18:2005</li><li>• EN ISO 13485:2003</li><li>• EN ISO 14971:2012</li><li>• EN ISO 15223-1:2012</li><li>• EN ISO 17665-1:2006</li></ul>
Notified Body Within the European Union	BSI
GMDN Code	46812, Blood Storage Solution, Anticoagulant

### Packaging Specifications

Quantity of Bags Per Case	500 mL bags:	18
	750 mL bags:	12

External Case Dimension	385 mm × 300 mm × 240 mm (15.16 in × 11.81 in × 9.45 in) (L × W × H)	
Case Weight	500 mL bags:	Containing bags: 10.31 kg (22.73 lb)
	750 mL bags:	Containing bags: 10.90 kg (24.03 lb)
Bag Primary Packaging Dimension	500 mL bags:	220 mm × 135 mm (8.66 in × 5.31 in)
	750 mL bags:	290 mm × 135 mm (11.42 in × 5.31 in)
Recycling		Bag overwrap: plastic (0.015 kg, 0.03 lb)
	Box: cardboard (0.60 kg, 1.32 lb)	
	IFU/Insert: paper (variable)	
Lot Number and Expiry Date Location	Refer to the Terumo BCT Support Center (TERUMOBCT.COM) for certificates of analysis.	
		Users can find the lot number adjacent to this symbol in the following locations on the packaging: box label, bag label.
		Users can find the expiration date adjacent to this symbol in the following locations on the packaging: box label, bag label.

## Technical Specifications

Sterilization		Indicates that the product is sterilized using steam or dry heat. This product is steam sterilized.
		Indicates that the product should not be reused.
		Indicates that the product is non-pyrogenic.
Expiry Date	24 months from the date of manufacture. Refer to product labeling for Expiry Date.	
Shipping Conditions	No special shipping conditions required.	
Bag Storage Conditions		Indicates the temperature limitations for product storage when accompanied by temperature(s).
	<ul style="list-style-type: none"> <li>• Store at up to 25 °C.</li> <li>• Protect from freezing.</li> <li>• Avoid excessive heat.</li> </ul>	


## Bag Specifications

Bag Composition	<ul style="list-style-type: none"> <li>• Bag: Polyolefin</li> <li>• Connectivity: <ul style="list-style-type: none"> <li>– Luer (made of polycarbonate) with a frangible connector</li> <li>– Injection port: polycarbonate, polypropylene and elastomer</li> </ul> </li> </ul>	
	Not made with natural rubber latex, PVC, or Di (2-ethylhexyl) phthalate (DEHP).	





Terumo BCT, Inc.  
10811 W. Collins Avenue  
Lakewood, Colorado 80215  
USA  
USA Phone: +1.877.339.4228  
Phone: +1.303.231.4357  
USA Fax: +1.866.715.6768  
Fax: +1.303.542.5215



Terumo BCT Europe N.V.  
Ikaroslaan 41  
1930 Zaventem  
Belgium  
Phone: +32.2.715.05.90  
Fax: +32.2.721.07.70

TERUMOBCT.COM

©2015 Terumo BCT, Inc.

2015-08  
Part No. 777020-370