





WWW.ZARYS.COM



# **since 1989...** ZARYS. We create good solutions



# WWW.ZARYS.COM

ZARYS International Group is a company with 32-year experience in the medical device industry. Since the foundation in 1989 ZARYS has surmounted a variety of business, trading and economic upheavals in becoming the market leader in healthcare distribution in Poland and abroad.

ZARYS offers a wide range of products including dressings, surgical drapes, examination gloves, disposable medical devices for urology, gynecology, oxygen therapy and many other items. The company is still expanding its product list in order to meet new market needs.

What distinguishes ZARYS from the competitors is strong ethics and continuous search for new solutions. In order to meet our clients' expectations we are constantly reviewing our product range. All of our devices are selected very carefully with great attention paid to current needs, quality convenient availability. The company has the pleasure to cooperate with highly qualified health care specialists, who share with us their expertise and understanding of the latest trends in clinical practice.

We offer products under our own brand and distribute the world's top manufacturers' medical devices. Our team believes that maintaining fair and reliable relationships with our international partners is the key to turn a business opportunity into a major success story, therefore we regularly visit our customers and suppliers in their own countries and handle challenges together.

Our values, strategies and comprehensiveness guarantee successful and satisfactory cooperation.

# **EXAMINATION GLOVES**

### nitrile examination gloves

11	mediCARE nitrile lime green
11	mediCARE nitrile pink
12	mediCARE nitrile AMG
13	mediCARE nitrile Hydra+
14	easyCARE nitrile black
14	easyCARE nitrile pink
15	easyCARE nitrile

### latex examination gloves

16	easyCARE latex
16	easyCARE latex PF

#### vinvl examination gloves

17 easyCARE vinyl easyCARE vinyl PF 17

# **BASIC MEDICAL DEVICES** pp.18-71

# **BASIC DISPOSABLE EQUIPMENT**

#### dispoFINE injection needle 21

- 21 dispoGUARD safety injection needle
- 23 duoNEX single-use syringe, 2-part, luer
- 23 dicoNEX SN single-use syringe, 3-part, luer, with needle alongside the syringe
- 24 dicoNEX single-use syringe, 3-part, luer
- dicoNEX single-use syringe, 3-part, luer lock 24
- 25 dicoNEX single-use catheter syringe, 3-part
- 25 dicoNEX single-use catheter syringe, 3-part
- 26 dicoTUBER tuberculin syringe
- dicoSULIN insulin syringe 26
- 27 safeLANCE pressure-activated safety lancet
- 27 dispoSULIN insulin pen needle
- 28 side hole needle **NEW!**
- 29 stomach tube/duodenal tube
- feeding tube graduated, with X-ray contrast 30
- 30 suction catheter
- 31 suction catheter with vacuum control
- suction catheter with frozen surface, phthalate-free 31
- 32 vaginal speculum cusco type
- dicoSPIKE withdrawal cannula 32
- 33 cervical brush standard
- cervical brush special 33
- 33 umbilical cord clamp

### **INFUSION THERAPY**

- SANVIflon I.V. cannula 35
- SANVIFION S I.V. cannula without port 35
- 36 easyFLOW IS infusion set
- easvFLOW IS ECO infusion set 36
- easyFLOW LINE extension tube for infusion pump 37
- 37 easyFLOW TS transfusion set
- 38 3-way stopcock
- 38 3-way stopcock with extension tube
- 39 needle-free valve
- 39 cannula stopper

# **UROLOGICAL CATHETERS AND URINE BAGS**

- 41 foley catheter 2-way, silicone-coated latex
- foley catheter 3-way, silicone-coated latex 42
- 42 foley catheter 2-way, 100% silicone
- 43 foley catheter 3-way, 100% silicone
- foley catheter 2-way with tiemann tip, silicone-coated latex 43
- 44 pezzer catheter, silicone-coated latex
- 44 nelaton catheter
- 45 tiemann catheter 45
- urine bag 46
- urine bag with sample port 46 urine bag with sample port 2W
- 47 samples taking urine bag
- 47 urine bag hanger
- 47 enema bag, sterile
- 48 enema bag, non-sterile
- 48 enema set, non-sterile
- 48 vomit bag. non-sterile

# **OXYGEN THERAPY**

- 50 oxygen mask with tubing
- non-rebreather mask with tubing 50
- 51 nebulizer mask with tubing
- multi-vent mask with tubing 51
- 52 nasal oxygen cannula
- 52 oxygen tubing
- nebulizer with mouthpiece and tubing 52
- 53 nebulizer with T adaptor, mouthpiece,
- corrugated tube and tubing 53 nebulizer with T F/M adaptor, mouthpiece,

## corrugated tube and tubing NEW!

### **ANAESTHESIA PRODUCTS**

- 55 laryngeal mask PVC
- 55 laryngeal mask silicone
- 56 intubation stylet
- 56 tracheal tube introducer curved
- 56 endotracheal tube uncuffed
- 57 endotracheal tube cuffed
- 57 reinforced endotracheal tube cuffed
- 58 endotracheal tube holder
- 58 tracheostomy tube cuffed
- 59 dual-limb breathing circuit
- 59 catheter mount with double swivel elbow connector
- 60 oropharyngeal airway/guedel airway
- 60 air cushion anaesthetic mask

### **GENERAL SURGERY**

- 62 Yankauer surgical suction cannula, standard
- 62 Yankauer surgical suction cannula, bulb tip
- 63 surgical suction tube
- surgical suction set, with Yankauer cannula, standard 64
- 65 surgical suction set, with Yankauer cannula, bulb tip

### **OTHER DISPOSABLE PRODUCTS**

- 67 latex-free tourniquet
- automatic tourniquet 67
- 68 dental brush with suction
- 68 oral swab

70

70

71

71

- 68 tonque depressor
- wash glove MOLTON 69
- wash glove SPUNLACE 69 69 wash glove with film lining ID wristband

electronic thermometer

entrance mat

PE gloves

# DRESSINGS pp.72-107

# TRADITIONAL DRESSING MATERIALS

#### gauze and non-woven dressings

- 75 GAZA lux dressing gauze, zig-zag folded
- 75 GAZA lux cutting gauze
- GAZA lux S cutting gauze 75
- 76 KOMPRI lux gauze swabs BESTSELLER!
- 77 KOMPRI lux S gauze swabs
- SERVI lux gauze lap sponges 78 78
- SERVI lux S gauze lap sponges TUPFER lux gauze balls 79
- 79 TUPFER lux S gauze balls
- NONVI lux non-woven swabs 80
- NONVI lux S non-woven swabs
- 80
- 80 NONVI lux S non-woven swabs with Y-incision
- MULTIabsorb ABD pads, multi-layer 81
- MULTIabsorb S ABD pads, multi-layer 81
- elastoLUMENAL eye pads 82 82 elastoLUMENAL S eye pads

## ADHESIVE TAPES AND PLASTERS WITH DRESSING

#### adhesive tapes

- 84 SENSIplast adhesive fabric tape, white
- 84 PLASTIplast adhesive fabric tape, skin colour
- 85 SOFTplast adhesive non-woven tape
- 85 SILKplast adhesive silk tape
- 86 FILMplast adhesive PE tape
- 86 POREplast adhesive perforated non-woven tape

#### plasters with dressing

- NONVIplast non-woven plaster with dressing 87
- FABRIplast fabric plaster with dressing 87

### SPECIAL PURPOSE DRESSINGS

#### self-adhesive non-woven dressings

- 89 elastopor, elastopor E adhesive non-woven dressing tape
- 90 elastopor STERIL non-woven dressing with absorbent pad BESTSELLER!
- 90 elastopor IV non-woven IV cannula dressing BESTSELLER!
- 91 elastopor STERIL D non-woven dressing with absorbent pad,
- incision and O-hole
- elastopor EYE non-woven eye dressing with absorbent pad 91
- 92 elastoKIDS IV non-woven IV cannula dressing
- elastoKIDS STERIL non-woven dressing with absorbent pad 92

# MEDICAL CLOTHING pp.108-131

### MEDICAL GOWNS

- medical gown with elastic cuffs, non-woven 111
- medical gown with knitted cuffs, non-woven 111

### SURGICAL CLOTHING

- 113 surgical suit, short-sleeved shirt+trousers, non-woven
- 114 surgical shirt. short-sleeved. non-woven
- 114 surgical shirt, long-sleeved, non-woven
- surgical trousers, with ties, non-woven 115

# **CAPS AND PROTECTIVE MASKS**

- medical bouffant cap, non-woven
- DONA medical bouffant cap 117
- LILI medical cap with rubber 117 118 POLA medical cap with rubber
- **SOL** medical cap with ties 118
- ABI medical cap with ties 118
- NINO medical cap 119
- KOKO long medical cap with snood 119
- 119 DORA long medical cap with snood
- BEN medical cap with ties 120
- TOM medical cap with ties 120
- 120 ROS medical cap with ties
- 121 medical mask with ties
- 121 medical mask with earloops

### stabilizing tapes

- 93 elastoNASAL non-woven nasal cannula fixation plaster
- 93 elastoFIXAL adhesive tape for catheter tubing fixation
- 94 elastoFIXAL Set adhesive tape for endotracheal tube
- 94 elastoFIXAL TUBE band for tracheostomy tube fixation NEW!

# STABILIZING DRESSINGS

- 96 elastoBAND BASIC knitted supporting bandage
- 96 elastoBAND FLEX elastic supporting bandage, woven
- 97 elastoNET elastic tubular net bandage

#### immobilizing dressings

- 98 orthoNATURAL orthopedic padding, natural
- orthoSYNTHETIC orthopedic padding, synthetic 98
- 99 elastoSLING triangular sling
- orthoBAND FAST, orthoBAND NORMAL, orthoBAND LONG 100 P.O.P. bandages

### SPECIALIST DRESSINGS

103 | paraffiNET paraffin gauze dressing

#### foil dressings

- elastoDERM foil dressing 104
- 104 elastoDERM F foil dressing with frame
- 105 elastoDERM F-IV IV cannula dressing with frame
- 105 elastoDERM IV IV cannula dressing
- 105 elastoDERM PAD foil dressing with absorbent pad

### ADDITIONAL PRODUCTS

- 107 alkoCLEAN gauze pad moistened with alcohol
- 107 ALKOPAD gauze pad moistened with alcohol NEW!
- ALKOPAD DUO gauze pad moistened with alcohol NEW! 107
- 107 thermCARE rescue blanket, foil, termic

- **PATIENT CLOTHING**
- patient's gown, tied-up, non-woven 123
- 123 patient's gown, Y-neck , non-woven
- 123 patient's gown, V-neck, non-woven
- 124 patient's gown, for delivery, non-woven 124 colonoscopy boxers. SMS non-woven
- 124
- colonoscopy boxers, PP non-woven 125
- gynecological skirt, non-woven examination underpants, single-use 125
- 125
- examination underpants, single-use

127shoe covers127shoe covers, calf-high

129

129

129

130

130

130

131

131 131

**BEDCLOTHES AND** 

SHOE COVERS AND SLIPPERS

**ABSORBENT UNDERPADS** 

disposable medical bed sheet

disposable medical pillow case

disposable medical quilt cover

disposable medical bedding set

disposable absorbent underpad

thermal patient cover soft

thermal patient cover

disposable absorbent underpad for operating table

disposable medical bed sheet, reinforced

# PERSONAL PROTECTIVE EQUIPMENTS pp.132-133

### **BETASAFE**

133 | protective coverall SAFE + against chemical and infective agents type 4b-5b-6b

# SURGICAL SETS, DRAPES AND GOWNS pp.134-187

# SURGICAL DRAPES

### 2-layer surgical drapes

- 141 | sterile surgical drapes with cellulose layer
- 141 sterile surgical drapes
- 141 sterile surgical drapes with adhesive edge
- 142 sterile surgical drapes with adhesive edge
- sterile surgical drapes with central fenestration
   sterile surgical drapes with adhesive central fenestration
- 143 sterile surgical drapes with adhesive central adjustable fenestration
- 143 sterile adhesive epidural drapes
- 143 sterile instrument table covers

#### **3-layer surgical drapes**

- 144 | sterile surgical drapes
- 144 sterile surgical drapes with adhesive edge
- 144 sterile surgical drapes with central fenestration
- 145 sterile surgical drapes with adhesive central
- fenestration 145 sterile surgical drapes with adhesive fenestration, reinforced

#### specialistic surgical drapes

146 146 147 147 147 147 148 148 148	abdominal drape 1-IF abdo-perineal drape 1 laparoscopy drape 3 orthopaedic drape 1 orthopaedic drape 2 orthopaedic drape 3 orthopaedic drape 4 extremity drape 1 extremity drape 2 extremity drape 3
149	extremity drape 5
149	gynaecology drape 1
149	gynaecology drape 2
149	delivery drape 1
150	ophthalmic drape 1-IF
150	ophthalmic drape 2-IF
150	ophthalmic drape 3-IF
150	ophthalmic drape 4-IF

#### surgical covers

152	sterile mayo stand cover, reinforced - red BESTSELLER!
152	sterile set of limb covers
152	sterile limb covers
152	sterile camera cabls covers - perforated,
	with 2 adhesive tapes BESTSELLER!
152	sterile camera cable cover - unperforated,
	without adhesive tape

# **SPORTS LINE** pp.188-193

### yellow**BAND**

190 | cohesive bandage

### yellowTAPE

190 | kinesiotape

### yellowRIGID TAPE

191 | sports tape, line glue coating

- 153 | sterile camera cable cover unperforated,
- with an adhesive tape 153 sterile circular banded covers for medical devices
- 153 sterile rectangular banded covers for medical devices
- 153 sterile c-arm cover set
- 153 sterile armboard cover set
- . . . . . . . . . . .

#### accessories

- 154 | sterile surgical tapes
- 154 sterile adhesive pouches
- 154 sterile adhesive fluid collection pouch
- 154 sterile absorbent drape
- 155 sterile absorbent drape
- 155 sterile incise film **elastoFILM**

# SURGICAL SETS

#### universal sets

unity	i sui sets
157	basic set 1
157	universal set 1
157	universal set 2
157	universal set 4 BESTSELLER!
158	universal set 5 BESTSELLER!
158	universal set 6
158	universal set 7
159	universal set 8
159	paediatric universal set 9
159	universal set 10
160	universal set 5-G
160	universal set 5-GP

### general surgery

- 161 | abdominal set 1-IF
- 161 abdo-perineal set 1
- 162 laparoscopy set 1-IF
- 162 laparoscopy set 2 **BESTSELLER!**
- 163 laparoscopy set 3
- 163 varicose vein set 1
- 163 | varicose vein set 2

#### orthopaedics and traumatology

164   extremity set 1		
164 extremity set 2		
165 arthroscopy set 1		
165	165 arthroscopy set 2	
166	66 arthroscopy set 2-IP	
166	166 arthroscopy set 3	
167	167 arthroscopy set 4	
167	arthroscopy set 5	
168	hip set 1	
168	hip set 2	

- 168 hip set 2 169 hip set 3
- 169 hip set 3169 shoulder set 1
- 170 shoulder set 2

# gynaecology and urology

- 171 gynaecology set 1
- 171 gynaecology set 2 BESTSELLER!172 gynaecology set 3
- 172 cystoscopy set 1
- 172 uro/gynaecology set 2
- 173 uro/gynaecology set 3
- 173 tur set 1
- 173 tur set 1 P1 NEW!
- 174 percutaneous lithotripsy set 1-IP

### obstetrics

174	c-section set 1
175	c-section set 1-IF
175	c-section set 2-IF
176	c-section set 3-IF
176	delivery set 1
176	delivery set 2
177	delivery set 3
177	delivery set 4
177	delivery set 5

#### ophthalmology

178 ophthalmic set 2-IF 178 ophthalmic set 3-IF

#### neurosurgery and otolaryngology

179	otolaryngology set 1
179	otolaryngology set 2
179	craniotomy set 1-IF

#### cardiosurgery and angiography

180	ablation set 1-IF
180	angiography set 2

# SURGICAL GOWNS

182	sterile gown NORMAL
182	surgical gown STANDARD BESTSELLER!
182	surgical gown STANDARD
182	surgical gown STANDARD
183	surgical gown STANDARD
183	surgical gown SOFT
183	surgical gown STANDART PLUS BESTSELLER!
183	surgical gown STANDART PLUS
184	surgical gown STANDART PLUS
184	surgical gown STANDART PLUS
184	surgical gown EXTRA SAFE
184	surgical gown COMFORT
185	surgical gown COMFORT PLUS
	tavila aurainal anuma

### non-sterile surgical gowns

185	L	surgical gown STANDARD
185		surgical gown PREMIUM

# yellow**SPORTS TAPE**

191 | sports tape

### yellow WRAP

191 | foam underwrap

### yellow**FOAM**

192 | foam tape

### scissors


192 | for sports tapes

### yellow**BOX**

192 | dispenser for kinesiotapes

# **Comprehensive range** of medical equipment

Operating theater is the heart of every hospital, where houndreds or even thousands of surgical procedures take place, during which the patient's skin layers are opened and cut. We present a comprehensive offer that allows you to protect the patient in the operating room while ensuring the safety of medical personnel.





# Custom surgical sets

**O** omegapack

# Custom procedure kits

**A** deltaset

())



# Basic medical devices

- devices for injection
- devices for infusion and transfusion
- urological devices
- devices for oxygen therapy and many others

**BETA Safe** 

Personal protective equipment

# Medical clothing

- surgical clothing
- shoe covers
- caps and protective masks

stila

care

patient clothing



# Anaesthesia products

# **Prevention of hospital infections**

- patient hygiene
- prevention of surgical site infections
- patient isolation









# **Medical gloves**

**EXAMINATION GLOVES** are intended for temporary use during medical examinations, diagnostic and therapeutic procedures, for handling medical contaminated material in order to reduce the risk of: transferring dangerous biological material through hands to a patient, acquisition of infections from patients by healthcare workers. Gloves are subject to the restrictive requirements of standards harmonized with the Regulation (EU) 2017/745 concerning medical devices and the Regulation (EU) 2016/425 of the European Parliament and of the Council on personal protective equipment.

**NITRILE GLOVES** are made of synthetic acrylonitrile butadiene rubber (NBR), they are an alternative to latex gloves. They do not contain latex proteins that cause type I allergic reactions, and the content of chemical accelerators is minimized, so they are friendly to the skin. They are characterized by high mechanical strength (also for puncture), chemical and biological resistance (constitute a barrier for viruses and bacteria), flexibility, thanks to which they adjust to the hand, providing the work comfort to the user. Thin walls and textured finger surface enhance grip and guarantee excellent tactile sensitivity. LATEX GLOVES made of natural rubber are recognized as the best barrier against blood-borne pathogens and part of hazardous chemical substances. They are also characterized by high flexibility, providing the work comfort to the user by fitting to the hand like a "second skin". Latex gloves, due to the natural raw material, are biodegradable. However, they can cause type I allergic reactions at people sensitive to proteins of natural rubber. Powdered gloves usually contain modified maize flour to facilitate glove donning. In the case of powder-free gloves, the inner surface is usually chlorinated or polymer-coated.

VINYL GLOVES made of PVC (polyvinyl chloride) are intended for low-risk procedures due to limited mechanical, chemical and biological resistance. The glove material does not cause Type I allergic reactions, only plasticizers can provoke a reaction in sensitized people - phthalates added during production to increase flexibility. ZARYS vinyl gloves are free of harmful phthalates.





# **EXAMINATION GLOVES**

























DIAGNOSTICS NURSING CARE

EMERGENCY MEDICAL SERVICES

STERILISATION UNITS

HYGIENIC AND SANITARY FACTORS

PATIENT HYGIENE

ZARYS



# mediCARE nitrile lime green








# mediCARE nitrile pink



packaging unit/ bulk packaging REF sizes







# mediCARE nitrile AMG

12





# mediCARE nitrile Hydra+



# Nitrile gloves with moisturizing and nourishing coating



# moistuirizes, improves skin appearance, smoothens and soothes irritations

### VITAMIN E

natural antioxidant, vitamin of youth - fights off free radicals and prevents wrinkles, moisturizes for a long time

### **ALMOND OIL**

regenerates, moisturizes, slows the aging process of the skin, smoothens, firms and helps removing dead cells, brighten the skin, soothens irritations and accelerates healing skin lesions, strengthens the hydro-lipid layer of the skin

## DERMATOLOGY SOLUTION FOR SKIN PROBLEMS

- gloves that protect and nourish the hands at the same time, while providing long-lasting hydration



### Anti-wrinkle effects proven in clinical tests\*

\*mediCARE Hydra+ glove presents a better anti-wrinkle effect compared to a traditional nitrile glove with visible reduction of wrinkles from the first day of use with progressive increase in improvment until the end of the tests

13

# easyCARE nitrile black

type

cuff

colour AQL



external surface

14

classification

internal surface

features indication

packaging unit/ bulk packaging REF sizes

examination, **nitrile**, non-sterile textured fingertips chlorinated beaded black 1.0 • class I medical device • personal protective equipment category III masking colour, increased endurance Y. Y ٦ 100 pcs/ 10x100 pcs RNBC03xx (xx:XS-XL)






# easyCARE nitrile pink




# easyCARE nitrile





packaging unit/ bulk packaging REF

sizes



etical mounting option

easyGARE



WWW.ZARYS.PL

Nitrile examination gloves

S

# **UNIVERSAL GLOVE DISPENSER**

for examination gloves

## > TYPES OF DISPENSERS:

- SINGLE • DOUBLE
- TRIPLE
- QUADRUPLE





# **HYGIENE**

- reduces the contamination risk
- of the packaging surface (no need to touch) coated with antibacterial paint to reduce the bacteria (Escherichia coli, MRSA)
- multiplication on the dispenser\* \*on reauest



### **BETTER ORGANIZATION OF THE WORKPLACE** space saving

improvement of putting on gloves procedure •

İij

# **CUSTOMIZED TO INDIVIDUAL** CUSTOMER NEEDS

- wall mounting on a flat surface vertically or horizontally (equipped with a mounting kit)
- vertically or horizontally mounted on a medical trolley, patient's bed or anaesthetic column (assembly kit on request)



# UNIVERSALITY

- vertical or horizontal assembly •
- fit to almost every glove box available on the medical market




# 16 easyCARE latex



type external surface cuff colour AQL classification features indication

packaging unit/ bulk packaging REF sizes



S

XL



# easyCARE latex PF

exa	type
text	external surface
pov	internal surface
bea	cuff
nati	colour
1.5	AQL
• cla • pe	classification
eas	features
ر ا	indication
100 10x RLE	packaging unit/ bulk packaging REF

sizes


examination, latex, non-sterile
textured
powder-free, polymer-coated
beaded
natural latex
1.5
<ul><li>class I medical device</li><li>personal protective equipment category III</li></ul>
easier donning due to inner polymer coating
100 pcs/ 10x100 pcs
RLB03xx (xx:XS-XL)







# easyCARE vinyl



indication

packaging unit/ bulk packaging REF sizes



100 pcs/ 10x100 pcs RWP07xx (xx:XS-XL)





# easyCARE vinyl PF



internal surface cuff colour

AQL

classification

indication

packaging unit/ bulk packaging REF sizes

examination, vinyl, non-sterile smooth powder-free beaded transparent 1.5

• class I medical device • personal protective equipment category III

Ϋ́

100 pcs/ 10x100 pcs



















# **Basic medical devices**

ZARYS International Group prides itself with a long lasting experience gained in the medical industry being a leading manufacturer and supplier of single use medical devices.

A wide range of products and complex solutions make it possible to successfully respond to our customers' needs and meet their expectations.

High standards of design, production, storing and transportation allow to maintain quality, safety and comfort of use to both the medical personnel and patients. A wide range of single use products include basic medical devices and those dedicated to specialized procedures.

ZARYS branded products find a wide application in anesthesiology, urology, gynecology, oxygen and fluid therapy and in many other medical specialities. They make it possible to perform medical procedures including patient diagnostics, prophylaxis, treatment and care.





# BASIC DISPOSABLE EQUIPMENT



W W W . Z A R Y S . C O M

# dispoFINE

# injection needle, sterile

#### features

- · needle tube is made of smooth and durable stainless steel
- sheath and hub made of polypropylene
- thin walls assure adequate stream pace
- tip sharpened in three surfaces (back cut)bevel type: LB/BL standard long bevel
- surface of needle tube covered with sliding agent (silicone)
- glued needle elements to prevent their separation (glue epoxy resin)
- hubs coloured in accordance with ISO code allowing fast recognition
- of needle size • latex-free, phthalate-free
- non-toxic, non-pyrogenic
- EO sterilized
- packaging: 1 pc./paper-foil

### indication

 device designed for intravenous, intramuscular and hypodermic injections

REF	size	colour	G	intermediate packaging	REF	size		colour
IN-03X13	0,3 mm x 13 mm	yellow	30	100 pcs	IN-08X25	0,8 m	m x 25 mm	deep green
IN-04X19	0,4 mm x 19 mm	- medium grey	27	100 pcs	IN-08X40	0,8 m	m x 40 mm	deep greer
IN-045X12	0,45 mm x 12 mm	brown	26	100 pcs	*IN-08X50	0,8 m	m x 50 mm	deep greer
IN-045X16	0,45 mm x 16 mm	brown	26	100 pcs	IN-09X25	0,9 m	m x 25 mm	yellow
IN-045X22	0,45 mm x 22 mm	brown	26	100 pcs	IN-09X40	0,9 m	m x 40 mm	yellow
IN-05X16	0,5 mm x 16 mm	orange	25	100 pcs	IN-11X40	1,1 m	m x 40 mm	cream
IN-05X20	0,5 mm x 20 mm	orange	25	100 pcs	IN-12X40	1,2 m	m x 40 mm	pink
N-05X25	0,5 mm x 25 mm	orange	25	100 pcs	IN-12X50	1,2 m	m x 50 mm	pink
IN-05X40	0,5 mm x 40 mm	orange	25	100 pcs	IN-16X40	1,6 m	m x 40 mm	white
IN-06X25	0,6 mm x 25 mm	deep blue	23	100 pcs	IN-18X40	1,8 m	m x 40 mm	blue-grey
IN-06X30	0,6 mm x 30 mm	deep blue	23	100 pcs	IN-21X40	2,1 m	m x 40 mm	pale green
IN-06X40	0,6 mm x 40 mm	deep blue	23	100 pcs	~~			0
IN-07X30	0,7 mm x 30 mm	black	22	100 pcs			*<	S I
IN-07X40	0,7 mm x 40 mm	black	22	100 pcs		in nu		~
IN-07X50	0,7 mm x 50 mm	black	22	100 pcs	bulk packag 50 x 100 p	0	availa on reg	

# dispoGUARD

# safety injection needle, sterile

#### features

- safety cover can be activated with one hand and it shields needle immediately after injection
- cover integrated with needle
- tip sharpened in three surfaces (back cut)
- colour coding system for size, marked on hub
- activation of safety mechanism signalized by audible click
- needle is securely closed inside activated safety cover
- latex-free, phthalate-free
- non-toxic, non-pyrogenic
- EO sterilized
- packaging: 1 pc./paper-foil

### indication

 device designed for intravenous, intramuscular and hypodermic injections

REF	size	colour	G	intermediate packaging
INB-03X13	0,3 mm x 13 mm	yellow	30	50 pcs
INB-045X12	0,45 mm x 12 mm	brown	26	50 pcs
INB-05X25	0,5 mm x 25 mm	orange	25	50 pcs
INB-05X40	0,5 mm x 40 mm	orange	25	50 pcs
INB-06X25	0,6 mm x 25 mm	deep blue	23	50 pcs
INB-06X30	0,6 mm x 30 mm	deep blue	23	50 pcs
INB-07X30	0,7 mm x 30 mm	black	22	50 pcs
INB-07X40	0,7 mm x 40 mm	black	22	50 pcs
INB-08X40	0,8 mm x 40 mm	deep green	21	50 pcs
INB-09X40	0,9 mm x 40 mm	yellow	20	50 pcs
INB-11X40	1,1 mm x 40 mm	cream	19	50 pcs
INB-12X40	1,2 mm x 40 mm	pink	18	50 pcs



G

21

21

21 20

20

19

18 18

16 15

14

intermediate packaging

100 pcs

100 pcs 100 pcs

100 pcs

100 pcs

100 pcs 100 pcs

100 pcs 100 pcs


100 pcs 100 pcs






bulk packaging: 50 x 50 pcs

# construction of the syringe



# GET TO KNOW OUR SYRINGES

# nozzle types:



# luer slip, centric

tip of the syringe is pushed into the hub of the needle or other device (e.g. catheters) and twisted. It allows to remove the connection immediately if necessary




# luer slip, eccentric


eccentric location of the tip is used in procedures involving the administration of drugs and fluids where closer proximity to skin is necessary



# luer lock

designed to secure connection between the syringe and the device (needle, tube, catheter, stopcock, etc.). The tip is equipped with thread to screw in the syringe to prevent disconnection







# catheter tip 90°

type of connection used for rinsing (cleaning) catheters, stomach tubes and feeding purposes

catheter tip 45° nozzle cut at 45 degrees increases the lumen surface

. . . . . . . . . . . . .

# insulin and tuberculin syringes

special design of the plunger seal reduces the drug/liquid residues

# duoNEX



# syringe type:



2-part - without plunger sealing



3-part - with plunger sealing

# graduation type:

<u>\_\_\_\_\_</u>

extended graduation



regular graduation

catheter syringes



extended graduation

# package colour and shape code:







# **duoNEX**

single-use syringe, 2-part, luer, sterile

### features

- white line
- transparent barrel allows full visual control of content
- white contrasting plunger
- luer nozzle centrally located in 2 ml syringe, eccentric in 5, 10, 20 ml syringes
- smooth plunger sliding
- · limitation flange securing against accidental removal of plunger • clear, legible and permanent black graduation facilitating dosage,
- adequately:
- every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml
- every 0,5 ml for capacity of 10 ml
- every 1 ml for capacity of 20 ml
- extended graduation that includes a 10% extension of the nominal capacity
- name of syringe and manufacturer's logo placed on barrel
- single-use
- non-pyrogenic, non-toxic
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

• 2-piece syringe for subcutaneous, intramuscular and intravenous infusions intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

REF	type	capacity	position of nozzle	intermediate packaging	bulk packaging
002ML-2CZ-B	luer	2 ml	centric	100 pcs	30 x 100 pcs
005ML-2CZ-B	luer	5 ml	eccentric	100 pcs	24 x 100 pcs
010ML-2CZ-B	luer	10 ml	eccentric	100 pcs	16 x 100 pcs
020ML-2CZ-B	luer	20 ml	eccentric	50 pcs	16 x 50 pcs

# dicoNEX SN

single-use syringe, 3-part, luer, with needle alongside the syringe, sterile

#### features

- in set with needle alongside the syringe
- transparent barrel allows full visual control of content, green contrasting plunger
- luer nozzle placed centrally for 2 ml capacity, eccentric for 5, 10, 20, 50 ml capacity
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- limitation flange securing against accidental removal of plunger .
- name of syringe and manufacturer's logo placed on barrel
- clear, legible and permanent black graduation facilitating dosage, adequately:
- every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml (extended to 6 ml)
- every 0,2 ml for capacity of 10 ml
- every 0,5 ml for capacity of 20 ml
- every 1 ml for capacity of 50 ml
- needle tube is made of smooth and durable stainless steel
- needle tip sharpened in three surfaces
- bevel type: LB/BL standard long bevel ٠
- hubs coloured in accordance with ISO code
- non-toxic, non-pyrogeniclatex-free, phthalate-free

### indication

• 3-piece syringe for subcutaneous, intramuscular and intravenous infusions

REF	type	capacity	position of nozzle	needle size	needle colour	intermediate packaging	bulk packaging
002ML-3CZ-SN	luer	2 ml	centric	0,6mm x 30mm	deep blue	100 pcs	30 x 100 pcs
005ML-3CZ-SN	luer	5 ml	eccentric	0,7mm x 30mm	black	100 pcs	24 x 100 pcs
010ML-3CZ-SN	luer	10 ml	eccentric	0,8mm x 40mm	deep green	100 pcs	16 x 100 pcs
020ML-3CZ-SN	luer	20 ml	eccentric	0,8mm x 40mm	deep green	50 pcs	16 x 50 pcs
050ML-3CZ-SN	luer	50 ml	eccentric	0,8mm x 40mm	deep green	25 pcs	16 x 25 pcs





# dicoNEX

# single-use syringe, 3-part, luer, sterile

#### features

- transparent barrel allows full visual control of content
- green contrasting plunger
- luer nozzle

24

- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger
- special design of the 1 ml syringe plunger stopper reduces dead space
- name of syringe and manufacturer's logo placed on barrel
- clear, legible and permanent black graduation facilitating
- dosage, adequately:
- every 0,01ml for capacity of 1ml - every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml
- every 0,2 ml for capacity of 10 ml
- every 0,5 ml for capacity of 20 ml
- every 1 ml for capacity of 50/60 ml
- single-use, non-pyrogenic, non-toxic
- · latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

### indication

· 3-piece syringe intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

REF	type	capacity	position of nozzle	intermediate packaging	bulk packaging
001ML-3CZ-BL	luer	1 ml	centric	100 pcs	32 x 100 pcs
002ML-3CZ-BL	luer	2 ml	centric	100 pcs	30 x 100 pcs
005ML-3CZ-BL	luer	5 ml	centric	100 pcs	24 x 100 pcs
010ML-3CZ-BL	luer	10 ml	centric	100 pcs	16 x 100 pcs
020ML-3CZ-BL	luer	20 ml	eccentric	50 pcs	16 x 50 pcs
050ML-3CZ-BL	luer	50/60 ml	eccentric	25 pcs	16 x 25 pcs

# dicoNEX single-use syringe, 3-part, luer lock, sterile

### features

- transparent barrel allows full visual control of content
- green contrasting plunger (amber-coloured syringe white contrasting plunger)
- luer lock nozzle
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger
- name of syringe and manufacturer's logo placed on barrel
- clear, legible and permanent black graduation facilitating dosage, adequately: - every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml and of 10 ml
- every 0,5 ml for capacity of 20 ml every 1 ml for capacity of 30 ml and of 50/60 ml
- single-use
- non-pyrogenic, non-toxic
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

- · 3-piece syringe intended for drawing up and administering of medicines/liquids
- to patient and collecting body fluid samples from patient
- · amber-coloured syringe is intended for drawing up and administering of photosensitive medicines

						••••••
REF	type	capacity	position of nozzle	colour of plunger	intermediate packaging	bulk packaging
003ML-3CZ-LL-BL	luer lock	3 ml	centric	green	100 pcs	30 x 100 pcs
005ML-3CZ-LL-BL	luer lock	5 ml	centric	green	100 pcs	24 x 100 pcs
010ML-3CZ-LL-BL	luer lock	10 ml	centric	green	100 pcs	16 x 100 pcs
020ML-3CZ-LL-BL	luer lock	20 ml	centric	green	50 pcs	16 x 50 pcs
030ML-3CZ-LL-BL	luer lock	30 ml	centric	green	50 pcs	16 x 50 pcs
050ML-3CZ-LL-BL	luer lock	50/60 ml	centric	green	25 pcs	16 x 25 pcs
amber type						
050ML-3CZ-LL-B-BL	luer lock	50/60 ml	centric	white	25 pcs	16 x 25 pcs







# dicoNEX

# single-use catheter syringe, 3-part, sterile

#### features

- · transparent barrel allows full visual control of content
- green contrasting plunger
- two luer connectors attached
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger
- nozzle cut at 90°
- catheter nozzle secured with cap
- name of syringe and manufacturer's logo placed on barrel
  extended graduation
- clear, legible and permanent black graduation facilitating dosage
- graduated every 1 ml
- single-use
- non-pyrogenic
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

### • 3-piece syringe with catheter nozzle

• intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

REF	ty	ре			intermediate packaging
050ML-3C	Z-CEW-BL ca	theter tip	50/60 ml	centric	25 pcs



# dicoNEX

# single-use catheter syringe, 3-part, sterile

#### features

- 3-piece syringe with centrally located nozzle for catheter
- two luer connectors attached
- made of polypropylene
- transparent barrel allows full visual control of content
- white contrasting plunger
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- limitation flange securing against accidental removal of plunger
- nozzle cut at 45°
- name of syringe and manufacturer's logo placed on barrel
- clear, legible and permanent black graduation facilitating dosage
- graduated every 2 ml
- latex-free, phthalate-free
- single-useEO sterilized
- pockaging: 1 pc./paper-foil

#### indication

- 3-piece syringe with catheter nozzle
- intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

REF	type	capacity	position of nozzle	intermediate packaging
100ML-3CZ-CEW-BL	catheter tip	100 ml	centric	25 pcs

bulk packaging: 4 x 25 pcs





# dicoTUBER

## tuberculin syringe, sterile

#### features

- with bypacked needle 0,45 x 13 mm
- made entirely of polypropylene
- white contrasting plunger
- transparent barrel allows full visual control of content
- luer nozzle

26

- smooth plunger sliding thanks to rubber sealing
- limitation flange securing against accidental removal of plunger
- clear, legible and permanent black graduation every 0,01 ml that allows
   precise dosing of medicine/liquid
- special design of plunger stopper reduces the residual volume of medicine/liquid
- name of syringe and manufacturer's logo placed on barrel
- single-use
- non-pyrogenic, non-toxic
- latex-free, phthalate-free
- EO sterilizedpackaging: 1 pc./paper-foil

#### indication

 3-piece syringe intended for precise drawing up and administering of medicines/liquids

REF	no. of units	capacity	intermediate packaging
TU-U20-1-BL	20	1 ml	100 pcs





# dicoSULIN

# insulin syringe, sterile

### features

- with bypacked needle 0,4 x 13 mm
- transparent barrel made of polypropylene allows full visual
- control of content
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger
- special design of plunger stopper reduces the residual volume
- of medicine/liquid
- clear, legible and permanent black graduation facilitating dosage
- name of syringe and manufacturer's logo placed on barrel
- single-use
- non-pyrogenic, non-toxiclatex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

- 3-piece syringe for subcutaneous, intramuscular and intravenous infusions
- for insulin shots

REF	no. of units	capacity	colour of plunger	intermediate packaging
IN-U40-1-BL	40	1 ml	red	100 pcs
IN-U100-1-BL	100	1 ml	orange	100 pcs





# safeLANCE

# pressure-activated safety lancet, sterile

#### features

- easy to use, minimizes unpleasant feelings
- special construction allows lancet tip to be automatically inserted into finger pad and then to be retracted so that needle is not visible before or after puncturing
- correctly matched size of lancet allows you to get optimal amount of blood with minimal patient's pain
- design ensures safety for medical personnel and patient by eliminating accidental puncture and infection
- · perfect guidance of lancet needle inside housing reduces vibration, thus decreases patient's pain and prevents damage to skin tissue
- available in two sizes:
- 23 G, puncture depth 1,8 mm green, dedicated for women and children 21 G, puncture depth 2,4 mm - pink, dedicated for men
- latex-free, phthalate-free
- single-use
- sterilized by radiation
- packaging: 100 pcs/box

#### indication

 capillary blood sampling from fingertip or ear flap for testing, for example, glucose level

REF	size / puncture depth	colour	intermediate packaging	
NA-18	23 G / 1,8 mm	green	100 pcs	bulk packaging
NA-24	21 G / 2,4 mm	pink	100 pcs	20 x 100 pcs

packaging:

> P bulk packaging: 50 x 100 pcs



# dispoSULIN insulin pen needle, sterile

#### features

- made of polypropylene, polyethylene and stainless steel
- triple-point needle blade allows easy and atraumatic skin penetration
- · covering the needle with silicone lubricant reduces resistance during injection
- · thin walls of needle ensure easy flow of medicine and the right flow rate
- coloured cap and paper cover of packaging make it easy to quickly recognize size of needle
- available in various sizes, allowing you to adapt to individual needs
- universal, compatible with all pen models
- single-use
- non-pyrogenic, non-toxic
- latex-free, phthalate-free
- EO sterilized
- packaging: 100 pcs/box

- · for subcutaneous injection
- · dedicated for use with pre-filled pens (insulin pens)

REF	size	colour	G	intermediate packaging
INP-033X12	0,33 mm x 12 mm	pink	29	100 pcs
INP-030X08	0,30 mm x 8 mm	yellow	30	100 pcs
INP-025X06	0,25 mm x 6 mm	purple	31	100 pcs
INP-025X08	0,25 mm x 8 mm	blue	31	100 pcs



# side hole needle

insulin pen needle, sterile

### features

- made of stainless steel
- the needle hub ideally suited for Luer and Luer-Lock (to ensure full tightnessof the connection with syringes)
- CE sign comply with PN-EN-ISO 7864
- individual packaging comply with PN-EN 556 and PN-EN 868, type blister- pack coloured cap and paper cover of packaging make it easy to quickly recognize
- single-use

28

- EO sterilized
- shelf life: 5 years

#### indication

- made of high quality materials, using the latest technology, which guarantees its high quality and failure-free, intended usage
- special design needle with a pencil cut with a side hole, designed to
- withdraw medication from a vial, including vials with a rubber stopper
- needle construction prevents from accidental, unintended production of rubber scrubs, in order to avoid contamination of the drug
- side hole with special cut, placed on the top of the needle enables free dissolution of the drug and limits foaming by directing the fluid directly onto the walls of the vial

REF	size	colour	G	intermediate packaging
HPP-1830	1,2 x 30 mm	pink	18	100 pcs







bulk packaging: 10 x 100 pcs

# stomach tube/duodenal tube

### sterile

### features

- made of PVC
- made of material that is resistant to bending and kinking due to the appropriate degree of hardness
  atraumatic, soft rounded closed tip
  four lateral eyes with smooth edges

- colour-coded connectors for size recognition, numerical size identification on packaging
  stomach tube is available in 3 lengths: 80 cm, 100 cm, 125 cm and duodenal tube is available in 150 cm length
- latex-free
- single-use
  EO sterilized
  packaging: 1 pc./paper-foil

- inserted through nose or mouth into stomach or duodenum
- drainage of stomach/intestinal content
- decompression of stomach (air or food evacuation)
- gastric lavage
- taking samples of stomach content
- control of bleeding from the upper section of digestive tract
   evaluation of secretory function
- · administration of medicines, fluids and/or food content into stomach or duodenum

REF	size	length	colour	intermediate packaging	bulk packaging
* ZZ-80/CH05	CH 5	800 mm	grey	35 pcs	20 x 35 pcs
* ZZ-80/CH06	CH 6	800 mm	light green	35 pcs	20 x 35 pcs
ZZ-80/CH08	CH 8	800 mm	blue	35 pcs	20 x 35 pcs
ZZ-80/CH10	CH 10	800 mm	black	35 pcs	20 x 35 pcs
ZZ-80/CH12	CH 12	800 mm	white	35 pcs	20 x 35 pcs
ZZ-80/CH14	CH 14	800 mm	green	30 pcs	20 x 30 pcs
ZZ-80/CH16	CH 16	800 mm	orange	30 pcs	20 x 30 pcs
ZZ-80/CH18	CH 18	800 mm	red	25 pcs	20 x 25 pcs
ZZ-80/CH20	CH 20	800 mm	yellow	15 pcs	20 x 15 pcs
ZZ-80/CH22	CH 22	800 mm	purple	15 pcs	20 x 15 pcs
ZZ-80/CH24	CH 24	800 mm	dark blue	15 pcs	20 x 15 pcs
* ZZ-100/CH05	CH 5	1000 mm	grey	30 pcs	20 x 30 pcs
* ZZ-100/CH06	CH 6	1000 mm	light green	30 pcs	20 x 30 pcs
ZZ-100/CH08	CH 8	1000 mm	blue	30 pcs	20 x 30 pcs
ZZ-100/CH10	CH 10	1000 mm	black	30 pcs	20 x 30 pcs
ZZ-100/CH12	CH 12	1000 mm	white	30 pcs	20 x 30 pcs
ZZ-100/CH14	CH 14	1000 mm	green	25 pcs	20 x 25 pcs
ZZ-100/CH16	CH 16	1000 mm	orange	20 pcs	20 x 20 pcs
ZZ-100/CH18	CH 18	1000 mm	red	15 pcs	20 x 15 pcs
ZZ-100/CH20	CH 20	1000 mm	yellow	15 pcs	20 x 15 pcs
ZZ-100/CH22	CH 22	1000 mm	purple	10 pcs	20 x 10 pcs
ZZ-100/CH24	CH 24	1000 mm	dark blue	10 pcs	20 x 10 pcs
* ZZ-125/CH05	CH 5	1250 mm	grey	30 pcs	20 x 30 pcs
* ZZ-125/CH06	CH 6	1250 mm	light green	30 pcs	20 x 30 pcs
ZZ-125/CH08	CH 8	1250 mm	blue	30 pcs	20 x 30 pcs
ZZ-125/CH10	CH 10	1250 mm	black	30 pcs	20 x 30 pcs
ZZ-125/CH12	CH 12	1250 mm	white	20 pcs	20 x 20 pcs
ZZ-125/CH14	CH 14	1250 mm	green	25 pcs	20 x 25 pcs
ZZ-125/CH16	CH 16	1250 mm	orange	20 pcs	20 x 20 pcs
ZZ-125/CH18	CH 18	1250 mm	red	20 pcs	20 x 20 pcs
ZZ-125/CH20	CH 20	1250 mm	yellow	15 pcs	20 x 15 pcs
ZZ-125/CH22	CH 22	1250 mm	purple	10 pcs	20 x 10 pcs
ZZ-125/CH24	CH 24	1250 mm	dark blue	10 pcs	20 x 10 pcs
••••••					
ZD-150/CH5	CH 5	1500 mm	grey	30 pcs	20 x 30 pcs
ZD-150/CH6	CH 6	1500 mm	light green	30 pcs	20 x 30 pcs
ZD-150/CH8	CH 8	1500 mm	blue	30 pcs	20 x 30 pcs
* ZD-150/CH10	CH 10	1500 mm	black	30 pcs	20 x 30 pcs
* ZD-150/CH12	CH 12	1500 mm	white	30 pcs	20 x 30 pcs
* ZD-150/CH14	CH 14	1500 mm	green	20 pcs	20 x 20 pcs
* ZD-150/CH16	CH 16	1500 mm	orange	15 pcs	20 x 15 pcs
* ZD-150/CH18	CH 18	1500 mm	red	15 pcs	20 x 15 pcs
ZD-150/CH20	CH 20	1500 mm	yellow	10 pcs	20 x 10 pcs
ZD-150/CH22	CH 22	1500 mm	purple	10 pcs	20 x 10 pcs
ZD-150/CH24	CH 24	1500 mm	dark blue	10 pcs	20 x 10 pcs







# feeding tube graduated, with X-ray contrast

### sterile

#### features

30

- made of soft PVC with a frosted surface, resistant to bending and kinking due to the appropriate degree of hardness
- · atraumatic, soft rounded, closed tip facilitates precise insertion of catheter
- two alternate end holes with gently stretched edges
- connector compatible with the conical tip of luer syringe
- integrated stoper allows the tube to be closed tightly
- colour-coded connectors for size recognition, numerical size identification on packaging
- radiopaque line along the entire length of the catheter
- digital depth graduation every 1 cm (from 5 cm to 40 cm)
- length 50 cm
- name of the manufacturer placed on the catheter
- individually packed straight
- latex-free
- phthalate-free
- single-useEO sterilized
- shelf life: 5 years
- packaging: 1 pc. / paper-foil

#### indication

· device inserted through patient's nose or mouth to perform enteral nutrition

REF	size	length	colour	intermediate packaging
CKS-04-50	CH 4	500 mm	red	100 pcs
CKS-05-50	CH 5	500 mm	grey	100 pcs
CKS-06-50	CH 6	500 mm	pink	100 pcs
CKS-08-50	CH 8	500 mm	pale green	100 pcs
CKS-10-50	CH 10	500 mm	black	100 pcs



# suction catheter


### sterile

#### features

- made of soft and flexible PVC, that is resistant to bending and kinking
  frozen surface prevents the catheter from sticking to the wall of the
- endotracheal tube
- atraumatic, soft rounded, open tip
- tip with two alternating lateral eyes with smooth edgescolour-coded connectors for size recognition, numerical size
- identification on the packaging
- latex-free
- single-use
- EO sterilizedpackaging: 1 pc./paper-foil

- device designed for suction of secretions from the upper respiratory tract
- recommended for source of suction with its own suction force regulation

REF	size	longth	colour	intermediate
KEF	size	length	colour	packaging
CO-B06-40	CH 6	400 mm	light green	50 pcs
CO-B08-40	CH 8	400 mm	blue	50 pcs
CO-B10-40	CH 10	400 mm	black	50 pcs
CO-B12-50	CH 12	500 mm	white	50 pcs
CO-B12-60	CH 12	600 mm	white	50 pcs
CO-B14-50	CH 14	500 mm	green	50 pcs
CO-B14-60	CH 14	600 mm	green	50 pcs
CO-B16-50	CH 16	500 mm	orange	50 pcs
CO-B16-60	CH 16	600 mm	orange	50 pcs
CO-B18-50	CH 18	500 mm	red	50 pcs
CO-B18-60	CH 18	600 mm	red	50 pcs
CO-B20-50	CH 20	500 mm	yellow	50 pcs
CO-B20-60	CH 20	600 mm	yellow	50 pcs
CO-B22-50	CH 22	500 mm	purple	50 pcs
CO-B22-60	CH 22	600 mm	purple	50 pcs











# suction catheter with vacuum control

sterile


#### features

- made of soft and flexible PVC, that is resistant to bending and kinking
- frozen surface prevents the catheter from sticking to the wall of the
- endotracheal tube • full control over the suction
- atraumatic, soft rounded, open tip
- tip with two alternating lateral eyes with smooth edges
- colour-coded connectors for size recognition, numerical size identification
- on the packaging
- latex-free single-use
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

device designed for suction of secretions from the upper respiratory tract

REF	size	length	colour	intermediate packaging
CO-A06-60	CH 6	600 mm	light green	50 pcs
CO-A08-60	CH 8	600 mm	blue	50 pcs
CO-A10-60	CH 10	600 mm	black	50 pcs
CO-A12-60	CH 12	600 mm	white	50 pcs
CO-A14-60	CH 14	600 mm	green	50 pcs
CO-A16-60	CH 16	600 mm	orange	50 pcs
CO-A18-60	CH 18	600 mm	red	50 pcs





bulk packaging: 10 x 50 pcs

# suction catheter with frozen surface

phthalate-free, sterile

#### features

- made of soft and flexible PVC, that is resistant to bending and kinking
- frozen surface enabling easy insertion
- atraumatic, soft rounded, open tip
- tip with two opposite lateral eyes with smooth edges
- colour-coded connectors for size recognition, numerical size identification on the packaging
   single-use
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

- device designed for suction of secretions from the upper respiratory tract
- recommended for source of suction with its own suction force regulation

REF	size	length	colour	intermediate packaging
CO-BF-ZP-B06-40	CH 6	400 mm	light green	50 pcs
CO-BF-ZP-B06-60	CH 6	600 mm	light green	50 pcs
CO-BF-ZP-B08-40	CH 8	400 mm	blue	50 pcs
CO-BF-ZP-B08-60	CH 8	600 mm	blue	50 pcs
CO-BF-ZP-B10-40	CH 10	400 mm	black	50 pcs
CO-BF-ZP-B10-60	CH 10	600 mm	black	50 pcs
CO-BF-ZP-B12-60	CH 12	600 mm	white	50 pcs
CO-BF-ZP-B14-60	CH 14	600 mm	green	50 pcs
CO-BF-ZP-B16-60	CH 16	600 mm	orange	50 pcs
CO-BF-ZP-B18-60	CH 18	600 mm	red	50 pcs
CO-BF-ZP-B20-60	CH 20	600 mm	yellow	50 pcs
CO-BF-ZP-B22-60	CH 22	600 mm	purple	50 pcs
CO-BF-ZP-B24-60	CH 24	600 mm	dark blue	50 pcs







# vaginal speculum

cusco type, sterile

### features

- both blades made of transparent polystyrene
- locking bolt made of polyethylene
- blades position blocked by a quarter turn of the locking bolt (without screwing in the bolt)
- perfectly smooth and streamlined edges
- colour coded bolts depending on the size
  single-use
- sin

32

- latex-free, phthalate-free EO sterilized
- EO sternized
  packaging: 1 pc./foil
- packaging: 1 pc./ro

### indication

• device designed for gynecological examination of the cervix and vagina

REF	width	colour	size
WG-XXS	16 mm	yellow	XXS
WG-XS	20 mm	white	XS
WG-S	24 mm	blue	S
WG-M	26 mm	red	М
WG-L	30 mm	green	L
WG-XL	34 mm	black	XL





100 pcs

# dicoSPIKE

# withdrawal cannula, sterile

#### features

- easy to use, ergonomic shape
- body with knurled surface for a secure grip
- sharp receiving spike enabling tight and safe connection with containers of different capacity
- receiving spike groove guarantees complete emptying of the container
- hermetic self-locking cap securing the connector against contamination
- available in version with air filter and in version with air and particle filter
  air filter protects against microbiological contamination and against harmful aerosols
- particle filter is an effective barrier to protect drugs and liquids from contamination
- female luer lock connector compatible with products with either male or female luer lock
- free of latex, PVC and phthalates
- non-pyrogenic
- single-use
- EO sterilized
- packaging: 1 pc./ paper-foil

- device for needle-free dissolving and / or withdrawing drugs from multi-dose packages (bottles, vials, etc.) - in both directions
- dicoSPIKE Chemo model is designed for medicines releasing dangerous aerosols such as cytostatics

REF	cap colour	air filter	particle filter	purpose	intermediate packaging
PPLZ	green	0,1 µm	-	standard medicines	100 pcs
PPLN	blue	0,1 μm	5 µm	standard medicines	100 pcs
PPLC	red	0,1 µm	5 µm	cytostatics	100 pcs



# cervical brush

# standard, sterile

#### features

- straight shape enables obtaining samples from the cervical canal
- soft bristle tips to prevent irritation
- optimal outreach thanks to the long handle
- single-use
- latex-free, phthalate-free
- EO sterilized packaging: 1 pc./paper-foil

### indication

• device designed for obtaining cell samples from the cervix, cervical canal and transformation zone in smear tests

REF	type	intermediate packaging
SC-001	standard	100 pcs



# cervical brush

### special, sterile

#### features

- · broom type enables obtaining samples from the cervix, cervical canal and transformation zone
- soft bristle tips enables precise collection of material for examination
- made entirely of polyethylene
- optimal outreach thanks to the long handle
- single-use
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

### indication

• device designed for obtaining cell samples from the cervix, cervical canal and transformation zone in smear tests

REF	type	intermediate packaging
SC-002	special	100 pcs



# umbilical cord clamp

sterile

### features

- made of polypropylene
- safe clamping teeth with rounded edges do not cut the umbilical cord
- serrated clamping surface of the arms stabilize the clamp on the umbilical cord, after closing it is not possible to open the clamp again
- · double lock catch prevents accidental opening
- universal size
- single-use
- non-pyrogenic
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

· device is designed for clamping the umbilical cord of a newborn in a labour ward in order to close the lumen of the umbilical blood vessels after cutting it off of the placenta

REF	size	colour	intermediate packaging	•••••
ZDPB	universal	white	50 pcs	











# **INFUSION THERAPY**



W W W . Z A R Y S . C O M

# SANVIflon

I.V. cannula, sterile

### features

- needle (stylet) made of stainless steel with back-cut blade
- cannula catheter made of Teflon (FEP), equipped with 4 radiopaque lines
- cannula equipped with wings and hand-closed cap of the upper port colour coded according to size
- handle that allows you to insert with one hand
- valve that prevents blood leakage to outside after inserting the cannula into the vessel
- flashback chamber secured with a luer lock cap
  available in various sizes from 14 G to 26 G
- available in var
  single-use
- Ingle-use
  latex-free
- phthalate-free
- EO sterilized
- packaging: 1 pc. hard blister (PVC + lacquered paper) with a colour print in accordance with the colour identification of the size

#### indication

- for cannulation of intravenous, peripheral vessels to:
- multiple injections of intravenous drugs
- intravenous infusion of solutions and parenteral nutrition suitable for peripheral veins
- maintain the patient's hydration and / or reduce dehydration if the patient is unable to take sufficient fluid orally
- blood transfusions and blood-like preparations
- ensuring vascular access in the event of the need for rapid drug delivery,
- especially during diagnostic and therapeutic procedures

REF	size [G]	size [mm]	colour	flow rate [ml/min]	intermediate packaging
KDF-06-19	26 x <sup>3</sup> /4''	0.6 x 19	violet	17	100 pcs
101-00-13	20 X /4	0,0 × 15	violet	11	100 pc3
KDF-07-19	24 x ¾"	0,7 x 19	yellow	23	100 pcs
KDF-09-25	22 x 1"	0,9 x 25	blue	36	100 pcs
KDF-11-32	20 x 1¼''	1,1 x 32	pink	56	100 pcs
KDF-13-45	18 x 1¾"	1,3 x 45	green	90	100 pcs
KDF-15-45	17 x 1¾"	1,5 x 45	white	142	100 pcs
KDF-17-45	16 x 1¾"	1,7 x 45	grey	200	100 pcs
KDF-20-45	14 x 1¾"	2,0 x 45	orange	305	100 pcs



# SANVIIIon S

### I.V. cannula without port, sterile

### features

- needle (stylet) made of stainless steel with back-cut blade
- cannula catheter made of Teflon (FEP), equipped with 4 radiopaque lines
- · cannula equipped with coloured fixing wings, without upper port
- handle that allows you to insert with one hand
- valve that prevents blood leakage to outside after inserting the cannula into the vessel
- flashback chamber secured with a luer lock cap
- available in two sizes: 24 G and 26 G
- single-use
- latex-free
- phthalate-free
- EO sterilized
- packaging: 1 pc. hard blister (PVC + lacquered paper) with a colour print in accordance with the colour identification of the size

#### indication

for cannulation of intravenous, peripheral vessels to:

- multiple injections of intravenous drugs
- intravenous infusion of solutions
- maintain the patient's hydration and / or reduce dehydration if the patient is unable to take sufficient fluid orally
- ensuring vascular access in the event of the need for rapid drug delivery, especially during diagnostic and therapeutic procedures



REF	size [G]	size [mm]	colour	flow rate [ml/min]	intermediate packaging
KDFS-06-19	26 x ¾"	0,6 x 19	violet	17	100 pcs
KDFS-07-19	24 x ¾"	0,7 x 19	yellow	23	100 pcs



# easyFLOW IS

# infusion set, sterile

#### features

36

- 2-way sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with fluid containers/bags
- · air inlet with anti-bacterial filter and coloured (blue) flap
- pliable drip chamber, size 6 cm, equipped with additional pressure wings to facilitate insertion into containers with liquids
- dropper: 20 drops = 1 ml +/- 0,1 ml
- fluid filter protecting against larger particles with pore size 15  $\mu m$
- min. 150 cm-long, soft, flexible tubing with luer lock conical fitting precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
- name of the manufacturer placed on the flow regulator
- · both tips of the device additionally secured with protective caps
- device available in a version containing phthalates and without phthalates single-use
- non-pyrogenic, non-toxic
- latex-free
- EO sterilized
- packaging: 1 pc./paper-foil with a coloured code (blue) identifying the device

#### indication

· device for fluid transfer in infusion therapy

REF	version
IS	with phthalates
IS-BF	without phthalates





# easyFLOW IS ECO

# infusion set, sterile

#### features

- made of high quality materials
- 2-way sharp closure-piercing device of the drip chamber that guarantees
- tight and reliable connection with fluid containers
- air inlet with anti-bacterial filter and colored (blue) flap
- pliable drip chamber, size 4cm
- dropper: 20 drops = 1 ml +/- 0.1 ml
- fluid filter protecting against larger particles with pore size  $~15\,\mu\text{m}$
- min. 150 cm-long, soft, flexible tubing
- universal luer lock conical fitting
- precise safe flow regulator
- both tips of the device additionally secured with protective caps
- latex-free
- · device available in a version containing phthalates and without phthalates
- packaging: 1pc./foil bag
- EO sterilized
- single-use

### indication

• device for fluid transfer in infusion therapy

REF	version
IS-E	with phthalates
IS-BF-E	without phthalates





### easyFLOW LINE

#### extension tube for infusion pump, sterile

#### features

- tubing made of PVC
- · device available in a transparent and amber version for photosensitive medicines
- diameter of the tubing outer: 2,4 mm; inner: 1,24 mm ٠
- luer lock fitting male and female
- single-use
- latex-free, phthalate-free non-pyrogenic, non-toxic
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

· device used to extend infusion line

REF	colour of tube	length	intermediate packaging
P450-BF	transparent	45 cm	50 pcs
P1500-BF	transparent	150 cm	50 pcs
PB1500-BF	amber	150 cm	50 pcs







### easyFLOW TS

#### transfusion set, sterile

#### features

- 2-way sharp closure-piercing device of the drip chamber with anti-slip
- surface that guarantees tight and reliable connection with containers/bags
- air inlet with anti-bacterial filter and a coloured (red) flap
- pliable drip chamber
- dropper: 20 drops = 1 ml +/- 0,1 ml
- special blood filter with large area, pore size 200 μm
- min. 150 cm-long, soft, flexible tubing
- universal luer lock fitting at the end of tubing
  precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
- name of the manufacturer placed on the flow regulator
- · both tips of the device additionally secured with protective caps
- · device available in a version containing phthalates and without phthalates
- single-use
- non-pyrogenic, non-toxic
- latex-free
- EO sterilized
- packaging: 1 pc./paper-foil with a coloured code (red) identifying the device

#### indication

· device used for the administration of blood and blood products

#### REF version

TS	with phthalates	
TS-BF	without phthalates	





### 3-way stopcock

#### sterile

#### features

38

- made of polycarbonate
- transparent housing
- · clear optical identifier of open/close position (marked flow
- directions) operating in 360 degrees
- consistent, straight flow channels and an optimal internal diameter ensure correct pressure measurement
- rotating male luer lock
- all inlets protected with caps
- resistant to lipids
- latex-free, phthalate-free
- single-use
- non-pyrogenic, non-toxic • EO sterilized
- packaging: 1 pc./paper-foil

#### indication

- device designed for managing the flow of infusion fluids
- during administration
- enables switching the infusion flow path
- for blood pressure monitoring
- for direct medicine/fluid supply

REF	colour	intermediate packaging	
KT-002	blue	50 pcs	bulk packaging: 20 x 50 pcs



### 3-way stopcock with extension tube

#### sterile

#### features

- · three-way stopcock with an extension tube of various lengths
- (from 7 cm to 100 cm)
- transparent housing
- clear optical identifier of open/close position (marked flow directions)
- three-way stopcock made of polycarbonate, extension tube made of PVC
- all inlets protected with caps
- resistant to lipids latex-free, phthalate-free
- single-use
- non-pyrogenic, non-toxic
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

- device designed for managing the flow of infusion fluids
- during administration enables switching the infusion flow path
- for blood pressure monitoring
- · for direct medicine/fluid supply

REF	length of tube	intermediate packaging
KTP-7	7 cm	50 pcs
KTP-10	10 cm	50 pcs
KTP-25	25 cm	50 pcs
KTP-50	50 cm	50 pcs
KTP-75	75 cm	50 pcs
KTP-100	100 cm	50 pcs





### needle-free valve

#### sterile

#### features

- body made of transparent material polycarbonate
- high-quality silicone membrane allows multiple port injections
- flat injection surface smooth and evenly ended port membrane ٠ that facilitates disinfection
- metal-free, suitable for MRI examination
- male connector of the transparent port is secured with a cap allowing for aseptic removal of the item from the packaging
- may be used for 7 days or up to 350 activations
- priming capacity: 0,09 ml
- high flow rate: 350 ml/min
- withstands injection fluid pressure of 3 bars (44 PSI)
- withstands back pressure of 2 bars (29 PSI)
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

- · intended for closing infusion ports
- device intended for multiple, needle-free injections
- (liquid drug administration, blood sample taking) protects infusion lines from contamination and infection during infusion therapy
- does not require using safety caps
- · dedicated for standard luer and luer lock connectors

REF	colour	intermediate packaging	
PB-N	blue	100 pcs	bulk packaging:
PB-T	transparent	100 pcs	20 x 100 pcs



39

# cannula stopper

#### sterile

#### features

- precisely and safely closes standard ports
- luer lock stopper has only one male tip
- · combi stopper with male and female tip especially dedicated for securing luer lock fitting in infusion and transfusion sets
- single-use
- non-pyrogenic
- latex-free, phthalate-free • EO sterilized
- packaging: 1 pc./paper-foil

- as a stopper of intravenous cannulas, ports and catheters
- intended to be used in combination with other devices during infusion therapy, protecting against the outflow of blood through the cannula and protecting against contamination and infection

REF	type	intermediate packaging	bulk packaging
KDKLL	luer lock	250 pcs	10 x 250 pcs
KDKC	combi	100 pcs	50 x 100 pcs







# UROLOGICAL CATHETERS AND URINE BAGS



W W W . Z A R Y S . C O M

### foley catheter

### 2-way, silicone-coated latex, sterile

#### features

- made of soft and flexible latex, that is resistant to bending and kinking • atraumatic, soft rounded, special tip facilitating insertion
- of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- universal urine drainage port of the catheter fits into standard connectors in urine bags
- silicone-coated surface that allows atraumatic insertion of the catheter
- · easy to fill and to empty, tear-resistant balloon of various capacities plastic or rubber valve to ensure balloon tightness ٠
- CH 6, CH 8, CH 10 size has a guide to facilitate catheterization
- contains latex, does not contain phthalates
- single-use
- EO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil
- indication
- to drain urine from the bladder

#### with plastic valve

REF	length	balloon	size	colour	intermediate packaging
CFNL-2D-06-03-P	270 mm	3 ml	CH 6	light green	10 pcs
CFNL-2D-08-05-P	270 mm	3-5 ml	CH 8	blue	10 pcs
CFNL-2D-10-05-P	270 mm	3-5 ml	CH 10	black	10 pcs
CFNL-2D-12-10-P	400 mm	5-10 ml	CH 12	white	10 pcs
CFNL-2D-14-10-P	400 mm	5-10 ml	CH 14	green	10 pcs
CFNL-2D-16-10-P	400 mm	5-10 ml	CH 16	orange	10 pcs
CFNL-2D-18-10-P	400 mm	5-10 ml	CH 18	red	10 pcs
CFNL-2D-20-10-P	400 mm	5-10 ml	CH 20	yellow	10 pcs
CFNL-2D-22-10-P	400 mm	5-10 ml	CH 22	purple	10 pcs
CFNL-2D-24-10-P	400 mm	5-10 ml	CH 24	dark blue	10 pcs
CFNL-2D-12-30-P	400 mm	30 ml	CH 12	white	10 pcs
CFNL-2D-26-30-P	400 mm	30 ml	CH 26	pink	10 pcs
CFNL-2D-28-30-P	400 mm	30 ml	CH 28	brown	10 pcs
CFNL-2D-30-30-P	400 mm	30 ml	CH 30	grey	10 pcs
CFNL-2D-14-50-P	400 mm	30-50 ml	CH 14	green	10 pcs
CFNL-2D-16-50-P	400 mm	30-50 ml	CH 16	orange	10 pcs
CFNL-2D-18-50-P	400 mm	30-50 ml	CH 18	red	10 pcs
CFNL-2D-20-50-P	400 mm	30-50 ml	CH 20	yellow	10 szt.
CFNL-2D-22-50-P	400 mm	30-50 ml	CH 22	purple	10 pcs
CFNL-2D-24-50-P	400 mm	30-50 ml	CH 24	dark blue	10 pcs
CFNL-2D-26-50-P	400 mm	30-50 ml	CH 26	pink	10 pcs
CFNL-2D-28-50-P	400 mm	30-50 ml	CH 28	brown	10 pcs
CFNL-2D-30-50-P	400 mm	30-50 ml	CH 30	grey	10 pcs

#### with rubber valve

REF	length	balloon	size	colour	intermediate packaging
CFNL-2D-06-03-L	270 mm	3 ml	CH 6	light green	10 pcs
CFNL-2D-08-05-L	270 mm	3-5 ml	CH 8	blue	10 pcs
CFNL-2D-10-05-L	270 mm	3-5 ml	CH 10	black	10 pcs
CFNL-2D-12-10-L	400 mm	5-10 ml	CH 12	white	10 pcs
CFNL-2D-14-10-L	400 mm	5-10 ml	CH 14	green	10 pcs
CFNL-2D-16-10-L	400 mm	5-10 ml	CH 16	orange	10 pcs
CFNL-2D-18-10-L	400 mm	5-10 ml	CH 18	red	10 pcs
CFNL-2D-20-10-L	400 mm	5-10 ml	CH 20	yellow	10 pcs
CFNL-2D-22-10-L	400 mm	5-10 ml	CH 22	purple	10 pcs
CFNL-2D-24-10-L	400 mm	5-10 ml	CH 24	dark blue	10 pcs
CFNL-2D-26-10-L	400 mm	5-10 ml	CH 26	pink	10 pcs
CFNL-2D-16-30-L	400 mm	30 ml	CH 16	orange	10 pcs
CFNL-2D-18-30-L	400 mm	30 ml	CH 18	red	10 pcs
CFNL-2D-20-30-L	400 mm	30 ml	CH 20	yellow	10 pcs
CFNL-2D-22-30-L	400 mm	30 ml	CH 22	purple	10 pcs
CFNL-2D-24-30-L	400 mm	30 ml	CH 24	dark blue	10 pcs
CFNL-2D-26-30-L	400 mm	30 ml	CH 26	pink	10 pcs









### foley catheter

#### 3-way, silicone-coated latex, sterile

#### features

42

- made of soft and flexible latex, that is resistant to bending and kinking
- universal urine drainage port of the catheter fits into standard connectors in urine bags
- additional inlet for flushing out the bladder
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- · silicone-coated surface that allows atraumatic insertion of the catheter
- easy to fill and to empty, tear-resistant balloon
- plastic valve to ensure balloon tightness
- contains latex, does not contain phthalates
- single-use
- EO sterilized • packaging: 1 pc./double packaging: internal- foil, external- paper-foil
- indication
- to drain urine from the bladder

REF	length	balloon	size	colour	intermediate packaging
CFNL-3D-16-30-P	400 mm	30 ml	CH 16	orange	10 pcs
CFNL-3D-18-30-P	400 mm	30 ml	CH 18	red	10 pcs
CFNL-3D-20-30-P	400 mm	30 ml	CH 20	yellow	10 pcs
CFNL-3D-22-30-P	400 mm	30 ml	CH 22	purple	10 pcs
CFNL-3D-24-30-P	400 mm	30 ml	CH 24	dark blue	10 pcs
CFNL-3D-26-30-P	400 mm	30 ml	CH 26	pink	10 pcs





bulk packaging; 40 x 10 pcs

### foley catheter

#### 2-way, 100% silicone, sterile

#### features

- made of soft and flexible silicone, that reduces the allergic
- reaction of the mucous membranes and increases the patient's comfort
- · atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- universal urine drainage port of the catheter fits into standard connectors in urine bags
- radiopaque line along the entire length of the catheter
- · easy to fill and to empty, tear-resistant balloon of various capacities
- plastic valve to ensure balloon tightness ٠
- CH 6, CH 8, CH 10 size has a nylon guide to facilitate catheterization
- available in version with 10 ml sterile water pre-filled syringe with
- 10% glycerin to seal the balloon
- separately packed Catheter Spigot inside the packaging
- latex-free, phthalate-free single-use ٠
- EO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil
   intermediate packaging: 10 pcs for REF: CFNS-2D-XX-XX-P, 5 pcs for REF: CFNSG-2D-XX-XX-P
- bulk packaging: 40 x 10 pcs for REF: CFNS-2D-XX-XX-P, 50 x 5 pcs for REF: CFNSG-2D-XX-XX-P

#### indication

• to drain urine from the bladder

· recommended for patients with blood clots or sediment in their urine

REF	intermediate packaging	REF with pre-filled syringe	intermediate packaging	length	balloon	size	colour
* CFNS-2D-06-03-P	10 pcs	CFNSG-2D-06-03-P	5 pcs	310 mm	3 ml	CH 6	pink
* CFNS-2D-08-05-P	10 pcs	CFNSG-2D-08-05-P	5 pcs	310 mm	3-5 ml	CH 8	black
* CFNS-2D-10-05-P	10 pcs	CFNSG-2D-10-05-P	5 pcs	310 mm	3-5 ml	CH 10	grey
* CFNS-2D-12-10-P	10 pcs	CFNSG-2D-12-10-P	5 pcs	400 mm	5-10 ml	CH 12	white
* CFNS-2D-14-10-P	10 pcs	CFNSG-2D-14-10-P	5 pcs	400 mm	5-10 ml	CH 14	green
* CFNS-2D-16-10-P	10 pcs	CFNSG-2D-16-10-P	5 pcs	400 mm	5-10 ml	CH 16	orange
* CFNS-2D-18-10-P	10 pcs	CFNSG-2D-18-10-P	5 pcs	400 mm	5-10 ml	CH 18	red
* CFNS-2D-20-10-P	10 pcs	CFNSG-2D-20-10-P	5 pcs	400 mm	5-10 ml	CH 20	yellow
* CFNS-2D-22-10-P	10 pcs	CFNSG-2D-22-10-P	5 pcs	400 mm	5-10 ml	CH 22	purple
* CFNS-2D-24-10-P	10 pcs	CFNSG-2D-24-10-P	5 pcs	400 mm	5-10 ml	CH 24	dark blue









### foley catheter

#### 3-way, 100% silicone, sterile

#### features

- made of soft and flexible silicone, that reduces the allergic reaction of the mucous membranes and increases the patient's comfort
- universal urine drainage port of the catheter fits into standard
- connectors in urine bags • additional inlet for flushing out the bladder
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- radiopaque line along the entire length of the catheter
- easy to fill and to empty, tear-resistant balloon
- plastic valve to ensure balloon tightness
- separately packed Catheter Spigot inside the packaging
  latex-free, phthalate-free
- iatex-free, phthala
  single-use
- EO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

#### indication

- to drain urine from the bladder
- · recommended for patients with blood clots or sediment in their urine

REF	length	balloon	size	colour	intermediate packaging
CFNS-3D-16-30-P	400 mm	30 ml	CH 16	orange	10 pcs
CFNS-3D-18-30-P	400 mm	30 ml	CH 18	red	10 pcs
CFNS-3D-20-30-P	400 mm	30 ml	CH 20	yellow	10 pcs
CFNS-3D-22-30-P	400 mm	30 ml	CH 22	purple	10 pcs
CFNS-3D-24-30-P	400 mm	30 ml	CH 24	dark blue	10 pcs





### foley catheter

### 2-way with tiemann tip,

silicone-coated latex, sterile


#### features

- made of soft and flexible latex, that is resistant to bending and kinking
  universal urine drainage port of the catheter fits into standard connectors
- in urine bags
- angled tip makes it easier to overcome obstacles in the urinary tract
  tip with two lateral eyes with smooth edges
- silicone-coated surface that allows atraumatic insertion of the catheter
- easy to fill and to empty, tear-resistant balloon
- plastic valve to ensure balloon tightness
- contains latex, does not contain phthalates
- single-use
- EO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

#### indication

· to drain urine from the bladder in the case of a narrowed urethra

REF	length	balloon	size	colour	intermediate packaging
CFTL-2D-12-10-P	400 mm	5-10 ml	CH 12	white	10 pcs
CFTL-2D-14-10-P	400 mm	5-10 ml	CH 14	green	10 pcs
CFTL-2D-16-10-P	400 mm	5-10 ml	CH 16	orange	10 pcs
CFTL-2D-18-10-P	400 mm	5-10 ml	CH 18	red	10 pcs
CFTL-2D-20-10-P	400 mm	5-10 ml	CH 20	yellow	10 pcs
CFTL-2D-22-10-P	400 mm	5-10 ml	CH 22	purple	10 pcs
CFTL-2D-24-10-P	400 mm	5-10 ml	CH 24	dark blue	10 pcs
CFTL-2D-12-30-P	400 mm	30 ml	CH 12	white	10 pcs
CFTL-2D-14-30-P	400 mm	30 ml	CH 14	green	10 pcs
CFTL-2D-16-30-P	400 mm	30 ml	CH 16	orange	10 pcs
CFTL-2D-18-30-P	400 mm	30 ml	CH 18	red	10 pcs
CFTL-2D-20-30-P	400 mm	30 ml	CH 20	yellow	10 pcs







bulk packaging 40 x 10 pcs

### pezzer catheter

#### silicone-coated latex, sterile

#### features

- made of soft and flexible silicone-coated latex, resistant to bending and kinking
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with three lateral eyes with smooth edges contains latex, does not contain phthalates
- single-use
- EO sterilized

44

• packaging: 1 pc./double packaging: internal- foil, external- paper-foil

#### indication

- to secure proper urinary drainage directly from the bladder or kidney
- in patients after surgical procedures introduced by the surgical wound, it is used for permanent urinary drainage, facilitating wound healing

 REF	length	size	intermediate packaging
CP-12-40	400 mm	CH 12	10 pcs
CP-14-40	400 mm	CH 14	10 pcs
CP-16-40	400 mm	CH 16	10 pcs
CP-18-40	400 mm	CH 18	10 pcs
CP-20-40	400 mm	CH 20	10 pcs
CP-22-40	400 mm	CH 22	10 pcs
CP-24-40	400 mm	CH 24	10 pcs
CP-26-40	400 mm	CH 26	10 pcs
CP-28-40	400 mm	CH 28	10 pcs
CP-30-40	400 mm	CH 30	10 pcs
CP-32-40	400 mm	CH 32	10 pcs
CP-34-40	400 mm	CH 34	10 pcs
CP-36-40	400 mm	CH 36	10 pcs
CP-38-40	400 mm	CH 38	10 pcs
CP-40-40	400 mm	CH 40	10 pcs







### nelaton catheter

#### sterile

#### features

- made of soft and flexible PVC, that is resistant to bending and kinking
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- + transparent tube allowing visual inspection of excreted urine
- colour-coded connectors for size recognition, numerical size identification on the packaging
- single-uselatex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

· to drain urine from the bladder

REF	length	size	colour	intermediate packaging
•••••				
CN-06-40	400 mm	CH 6	light green	50 pcs
CN-08-40	400 mm	CH 8	blue	50 pcs
CN-10-40	400 mm	CH 10	black	50 pcs
CN-12-40	400 mm	CH 12	white	50 pcs
CN-14-40	400 mm	CH 14	green	50 pcs
CN-16-40	400 mm	CH 16	orange	50 pcs
CN-18-40	400 mm	CH 18	red	50 pcs
CN-20-40	400 mm	CH 20	yellow	50 pcs
CN-22-40	400 mm	CH 22	purple	50 pcs
CN-24-40	400 mm	CH 24	dark blue	50 pcs







### tiemann catheter

#### sterile

#### features

- made of soft and flexible PVC, that is resistant to bending and kinking
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
  cone-shaped angled tip, bent at an angle of 45°
- transparent tube allowing visual inspection of excreted urine
- colour-coded connectors for size recognition, numerical size
- identification on the packaging
- single-use
- latex-free
- EO sterilized • packaging: 1 pc./paper-foil

#### indication

• to drain urine from the bladder in the case of a narrowed urethra

REF	length	size	colour	intermediate packaging	bulk packaging
CT-08-40	400 mm	CH 8	blue	60 pcs	10 x 60 pcs
CT-10-40	400 mm	CH 10	black	60 pcs	10 x 60 pcs
CT-12-40	400 mm	CH 12	white	50 pcs	10 x 50 pcs
CT-14-40	400 mm	CH 14	green	50 pcs	10 x 50 pcs
CT-16-40	400 mm	CH 16	orange	50 pcs	10 x 50 pcs
CT-18-40	400 mm	CH 18	red	40 pcs	10 x 40 pcs
CT-20-40	400 mm	CH 20	yellow	40 pcs	10 x 40 pcs
CT-22-40	400 mm	CH 22	purple	40 pcs	10 x 40 pcs
CT-24-40	400 mm	CH 24	dark blue	40 pcs	10 x 40 pcs





### urine bag

#### sterile

#### features

- latex-free 2000 ml bag made of 0,22 mm thick PVC (double wall)
  double sealed PVC for maximum protection against tearing the bag
- effective non-return valve ٠
- tight T-valve, that can be easily operated with one hand
- connecting tubing with universal stepped connector, standard tubing length 90 cm and 150 cm
- · tubing made of material that is resistant to bending and kinking and which ensures a free and effective urine outflow
- precise, easy-to-read scale (every 100 ml), adequate for diuresis measurement
- strengthened holes for hanging, matching standard hangers
- white back wall making urine easy to notice
- single-use
- latex-free
- EO sterilized • packaging: 1 pc./foil

- for 24-hour extracorporeal collection of urine and other body fluids
- maximum time of use 24 hours

REF	capacity	length of tubing	intermediate packaging
WMD2000-90	2000 ml	90 cm	10 pcs
WMD2000-150	2000 ml	150 cm	10 pcs





### urine bag

#### with sample port, sterile

#### features

46

- latex-free 2000 ml bag made of 0,40 mm thick PVC (double wall)
- effective non-return valve
- tight T-valve, that can be easily operated with one hand
- connecting tubing with universal stepped connector, standard tubing length 120 cm (other lengths available: 90 cm, 150 cm)
- tubing made of material that is resistant to bending and kinking and which ensures a free and effective urine outflow
- self-tightening port for sample taking
  precise, easy-to-read scale, adequate for diuresis measurement, with high measurement accuracy from 25 ml, every 50 ml from 50 ml to 200 ml and every 100 ml to 2000 ml
- strengthened holes for hanging, matching standard hangers
- white back wall making urine easy to notice
- single-use
- latex-free • EO sterilized
- packaging: 1 pc./paper-foil

#### indication

- · for extracorporeal collection of urine and other body fluids
- maximum time of use 7 days

REF	capacity	length of tubing	intermediate packaging
WMT2000-90	2000 ml	90 cm	10 pcs
WMT2000-120	2000 ml	120 cm	10 pcs
WMT2000-150	2000 ml	150 cm	10 pcs





### urine bag

#### with sample port 2W, sterile

#### features

- latex-free 2000 ml bag made of medical grade PVC
- white back wall making urine easy to notice
- effective non-return valve
- tight T-valve, that can be easily operated with one hand and which can be folded and tucked in a hard sheath
- Pasteur type droplet chamber providing a "dry" air barrier
- two hydrophobic air filters (in a bag and chamber), which prevent sucking, balancing the pressure inside the system
- 120 cm long, large diameter tubing, ensuring a free and effective urinary outflow with a universal stepped connector
- self-tightening port for sample taking at the end of tubing
- tubing equipped with a sliding clamp and a fastening pin to stabilize
- the tubing
- at the place where the tubing connects to the chamber, an additional reinforcing element prevents bending and kinking the tubing
- · precise, easy-to-read scale, adequate for diuresis measurement, with high measurement accuracy, every 25 ml from 25 ml to 100 ml and every 100 ml from 100 ml to 2000 ml
- tag to describe the patient's details and date
- integrated universal hanger system matching the round and square bed rails
- additional fixing with a string
- single-use
- latex-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

· device is intended for 14-days urine collection in a closed system

REF	capacity	length of tubing
WM2T2000-120	2000 ml	120 cm





### samples taking urine bag

sterile

#### features

- latex-free bag with a capacity of 100 ml, size 19 cm x 7 cm
- made of delicate, durable PE foil
- easy-to-read scale, adequate for diuresis measurement, highly accurate every 10 ml
- · adhesive part of the bag coated with anti-allergic, high-quality glue
- additional sealing made of anti-allergic foam that protects against the wall sticking
- of the bag, thanks to which free flow of urine is possible
- single-use
- latex-free, phthalate-freeEO sterilized
- packaging: 1 pc./foil

#### indication

· for urine sample collection in infants and small children

REF	capacity	purpose	version	intermediate packaging
WMCH-001	100 ml	for boys	with sponge	100 pcs
	100 1	6 1	111 A	100
WMCH-002	100 ml	for boys	without sponge	100 pcs
WMDZ-001	100 ml	for girls	with sponge	100 pcs
WWDZ-001	100 111	ioi giits	with sponge	100 hcs
WMDZ-002	100 ml	for girls	without sponge	100 pcs
WHIDE 002	100 IIII	IOI BIILD	without sponge	100 pc3



### urine bag hanger

#### non-sterile

#### features

- made of strong and durable plastic
- special clamp that prevents tube bending
- suitable for round and square bed rails
- two handles on each side for hanging urine bags of various sizes
- manufacturer's name placed on the device
- latex-free, phthalate-free
- non-sterile

#### indication

· for hanging standard urine collection bags

REF	intermediate packaging
WNM01	50 pcs





47





### enema bag

#### sterile

#### features

- bag with a capacity of 1750 ml with two holes for hanging
- made of medical PVC
- scaled from 50 ml to 1750 ml
- linear scale every 100 ml, digital scale every 200 ml
- tube with a length of 150 cm with a sliding clamp, ended with an atraumatic opening and one lateral opening
- tube end is secured with a tube cap
- upper part of the bag is reinforced with an additional flange stabilizing the inlet during filling with liquid
- flange is equipped with a non-return valve
- single-use
- latex-free
- EO sterilized
- packaging: 1 pc./foil

#### indication

 bag that allows introduction of fluid into the large bowel through the anus to remove fecal material and stimulate peristaltic bowel movement

REF

WLS



### enema bag

#### non-sterile

#### features

- bag with a capacity of 1500 ml with a hole for hanging
- made of medical PVC
- digitally scaled every 250 ml
- tube with a length of 115 cm with a sliding clamp, ended
- with an atraumatic opening and one lateral opening
- greased end of the tube secured with a cap
  single-use
- 48 single-use
   latex-free
  - latex-free
     non-sterile

#### non otern

#### indication

 bag that allows introduction of fluid into the large bowel through the anus to remove fecal material and stimulate peristaltic bowel movement

•••••	•••	• • •	•••	• • • •	•••	••••	•
1	RE	F.					

.....

### bulk packaging: 100 pcs

### enema set

#### non-sterile

#### features

WLN

- bag with a capacity of 1500 ml with a hole for hanging
- made of medical PVCdigitally scaled every 250 ml
- tube with a length of 115 cm with a sliding clamp, ended with
- an atraumatic opening and one lateral opening • greased end of the tube secured with a cap
- greased end of the tu
  drape 42 cm x 45 cm
- liquid soap (9 ml)
- PE gloves 1 pair
- PE gloves 1
  single-use
- Ingle-use
   latex-free
- non-sterile
- packaging: 1 set/foil

#### indication

REF

ZLN

 set that allows introduction of fluid into the large bowel through the anus to remove fecal material and stimulate peristaltic

bowel movement

bulk packaging 50 x 1 set

### vomit bag

#### non-sterile

#### features

- bag is made of durable, transparent blue foil, enabling the observation of secretions.
- total bag capacity: 2000 ml
- exact measuring graduation (from 10 ml to 100 ml every 10 ml and from 100 to 2000 ml every 100 ml) placed on the bag allows accurate estimation of fluid volume
- direction of graduation reading for small volumes is marked with the tip
- wide bag inlet protected by a plastic collar
   circle-shaped collar ensures a firm grip, makes it easier to maneuver the bag, reducing
- the risk of content contamination
  collar equipped with a special incision to close the bag and hygienic disposal of the vomit "twist and hook"
- device name, manufacturer's name and pictorial instructions for use printed directly on the bag
- single-use
- latex-free
- phthalate-free
- class I non-sterile medical device
- packaging: 20 pcs/foil bag

- recommended in uncomfortable situations, in cases of gastroesophageal reflux requiring immediate action
- used in hospitals, outpatient departments, ambulances, etc.









REF	colour	intermediate packaging
WW-01	blue	20 pcs



# **OXYGEN THERAPY**



W W W . Z A R Y S . C O M

### oxygen mask

with tubing, sterile

#### features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
  equipped with a 210 cm (+/- 5%) long tubing with universal connectors
  tubing with star cross-section that is resistant to bending
  revolving connector that allows adjusting to patient's position

- 50 • single-use
  - latex-free
    - EO sterilized
    - packaging: 1 pc./paper-foil

#### indication

• designed to deliver oxygen with precise regulation of concentration, for medium oxygen concentration

REF	length of tubing	size
MTS-S	210 cm +/- 5%	S
MTS-M	210 cm +/- 5%	М
MTS-L	210 cm +/- 5%	L
MTS-XL	210 cm +/- 5%	XL





### non-rebreather mask

#### with tubing, sterile

#### features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
  equipped with a 210 cm (+/- 5%) long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
- 1000 ml oxygen reservoir
- revolving connector that allows adjusting to patient's position
- · silicone valves on the connector and ventilation openings
- single-use
- latex-free
- EO sterilized packaging: 1 pc./paper-foil

#### indication

- designed to deliver oxygen with precise regulation of concentration, for high oxygen concentration

REF	length of tubing	size
MTW-S	210 cm +/- 5%	S
MTW-M	210 cm +/- 5%	М
MTW-L	210 cm +/- 5%	L
MTW-XL	210 cm +/- 5%	XL





### nebulizer mask

with tubing, sterile

#### features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
- equipped with a 210 cm (+/- 5%) long tubing with universal connectors
  tubing with star cross-section that is resistant to bending
  8 ml nebulizer (graduated every 1 ml)

- .... resource (graduated every 1 ml) • average size of the molecule at a flow of 6l / min 3,23  $\mu m$  single-use
- latex-free • EO sterilized
- packaging: 1 pc./paper-foil

#### indication

• for producing aerosol containing a comminuted medicine, that can be delivered directly to the respiratory system as the patient breathes

REF	length of tubing	size
MTN-S	210 cm +/- 5%	s
MTN-M	210 cm +/- 5%	M
MTN-L	210 cm +/- 5%	L
MTN-XL	210 cm +/- 5%	XL





# multi-vent mask

#### with tubing, sterile

#### features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
- ٠ equipped with a 210 cm (+/- 5%) long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
- equipped with 2 jets white to set mid flow concentration (35%, 40% and 50%) and green to set low concentration (24%, 26%, 28% and 30%)
- single-use
- latex-free • EO sterilized
- packaging: 1 pc./foil

#### indication

• designed to deliver oxygen with precise regulation of concentration, regardless of way of patient's breathing

REF	length of tubing	size
MTDV-L	210 cm +/- 5%	L
MTDV-XL	210 cm +/- 5%	XL





### nasal oxygen cannula

#### sterile

#### features

- made of medical PVC
- · very soft tips with smooth edges
- fastening under the chin with a soft, sliding adjustment ring
  various spacing of the nasal part of the cannula depending on the size
- universal connector compatible with any source of oxygen
- tubing with star cross-section that is resistant to bending
  available in 3 sizes: for adults, for children and for infants with 140 cm,
- 200 cm, 300 cm and 500 cm long tubing single-use
- latex-free

52

- EO sterilized
- packaging: 1 pc./foil

#### inication

• designed to deliver oxygen through the nose

REF	length of tubing	purpose	intermediate packaging	bulk packaging
CTNN-200	200 cm	for infants	25 pcs	10 x 25 pcs
CTNZ-200	200 cm	for children	25 pcs	10 x 25 pcs
CTND-140	140 cm	for adults	30 pcs	10 x 30 pcs
CTND-200	200 cm	for adults	25 pcs	10 x 25 pcs
CTND-300	300 cm	for adults	15 pcs	10 x 15 pcs
CTND-500	500 cm	for adults	10 pcs	10 x 10 pcs

### oxygen tubing

#### sterile

#### features

- made of transparent, non-toxic PVC
- tubing with star cross-section that is resistant to bending, with reinforcing straps along the entire length
- standardized connectors at the end of the tubing allowing
- connection to oxygen devices
- single-use
- latex-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

• for transporting oxygen from an oxygen source to another medical device for oxygen supply, e.g. oxygen mask

REF	length of tubing	bulk packaging
DT2100	210 cm +/- 2%	100 x 1 pc.
DT4260	426 cm +/- 2%	50 x 1 pc.
DT7620	762 cm +/- 2%	25 x 1 pc.

### nebulizer

#### with mouthpiece and tubing, sterile

#### features

- made of transparent, non-toxic PVC
- tubing with star cross-section that is resistant to bending
- 15 ml nebulizer (graduated at 2, 4, 6, 9, 12, 15 ml)
- set components: nebulizer, tubing, mouthpiece
- single-use
- latex-free
- EO sterilized • packaging: 1 pc./foil

#### indication

 designed for producing aerosol containing a comminuted medicine, that can be delivered directly to the respiratory system as the patient breathes

REF length of tubing		
NUD	210 cm +/- 5%	bulk packaging: 50 x 1 pc.







### nebulizer

with T adaptor, mouthpiece, corrugated tube and tubing, sterile

#### features

- made of transparent, non-toxic PVC
- · tubing with star cross-section that is resistant to bending
- 20 ml nebulizer (graduated every 2 ml)
- set components: nebulizer, tubing, T adaptor, mouthpiece and 15 cm long corrugated tube
- diameter of T adaptor: 22F/22F
- single-uselatex-free
- EO sterilized
- packaging: 1 pc./foil

#### indication

 designed for producing aerosol containing a comminuted medicine, that can be delivered directly to the respiratory system as the patient breathes

REF length of tubing		
NUT	210 cm +/- 5%	bulk packaging: 100 x 1 pc.



### nebulizer

# with T F/M adaptor, mouthpiece, corrugated tube and tubing, sterile

#### features

- made of non-toxic materials
- 210 cm long, bend-resistant tubing with star-shaped cross-section
- nebuliser with a 6ml capacity (graduated every 1ml), producing particles with
- an average size of MMAD =  $3.91 \,\mu\text{m}$  at a flow rate of 6 l/min.
- nebuliser, tubing, T-adaptor, mouthpiece, one-way valve and corrugated tube included in a kit
- universal T-adaptor with a diameter of 22M/15F 22F enables the kit to be connected to endotracheal, tracheostomy tube or breathing circuit and perform nebulisation on mechanically ventilated patients
- the universal (detachable) mouthpiece allows nebulisation to be given to patients who are breathing on their own
- silicone one-way valve (detachable) reduces medication waste during nebulisation
- single-use
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./foil

- · device is designed to produce an aerosol containing fragmented medication, which
- can be delivered deep into the respiratory system with the patient's inhalation
- enables nebulisation to be carried out in mechanically ventilated patients by connecting directly to the endotracheal tube, tracheostomy tube or inspiratory arm of the ventilator breathing system
- allows you to perform nebulisation in patient "on one's own breath"

REF length of tubing		
••••••		bulk packaging:
NUTR	210 cm +/- 5%	100 x 1 pc.







# ANAESTHESIA PRODUCTS



W W W . Z A R Y S . C O M

### laryngeal mask

PVC, sterile

#### features

- made of medical PVC
- elastic mask body with description (size, cuff volume, suggested patient's weight) .
- elliptic shape, flexible sealing cuff and universal, transparent connector: 15 mm
- colour-coded valve protection makes it easy to quickly recognize size .
- single-use
- latex-free phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication

- device inserted through the mouth, up to the back of the throat, forms a seal at the top of the glottis in order to:
- direct protecting the patient's airway patency
   artificial ventilation also during cardiopulmonary resuscitation
- maintaining airway patency in: short-term surgical procedures, without the need of muscles relaxation, when endotracheal intubation is not necessary and under anesthesia for diagnostic procedures

REF	size	colour of valve protection	weight of patient	max. inflation	intermediate packaging
MKPJ-10	1	pink	≤ 5 kg	4 ml	5 pcs
MKPJ-15	1.5	orange	5-10 kg	7 ml	5 pcs
MKPJ-20	2	dark green	10-20 kg	10 ml	5 pcs
MKPJ-25	2.5	brown	20-30 kg	14 ml	5 pcs
MKPJ-30	3	green	30-50 kg	20 ml	5 pcs
MKPJ-40	4	purple	50-70 kg	30 ml	5 pcs
MKPJ-50	5	dark blue	70-100 kg	40 ml	5 pcs



### laryngeal mask

#### silicone, sterile

#### features

- made of biocompatible silicone, which ensures excellent flexibility of the device
- elastic mask body with description (size, cuff volume, suggested patient's weight)
- elliptic shape, flexible sealing cuff and universal, transparent connector: 15 mm
- single-use ٠
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

- device inserted through the mouth, up to the back of the throat, forms
- a seal at the top of the glottis in order to:
- direct protecting the patient's airway patency
- artificial ventilation also during cardiopulmonary resuscitation
- maintaining airway patency in: short-term surgical procedures,
- without the need of muscles relaxation, when endotracheal intubation is not necessary and under anesthesia for diagnostic procedures

REF	size	weight of patient	max. inflation	intermediate packaging
MKSJ-10	1	≤ 5 kg	4 ml	5 pcs
MKSJ-15	1.5	5-10 kg	7 ml	5 pcs
MKSJ-20	2	10-20 kg	10 ml	5 pcs
MKSJ-25	2.5	20-30 kg	14 ml	5 pcs
MKSJ-30	3	30-50 kg	20 ml	5 pcs
MKSJ-40	4	50-70 kg	30 ml	5 pcs
MKSJ-50	5	70-100 kg	40 ml	5 pcs







### intubation stylet

#### sterile

#### features

- made of PVC-coated aluminium
- special material for both stiffness and pliability
- intended for use during difficult intubation
- size and diameter markings
- soft distal tip

56

- single-use
- phthalate-free latex-free
- EO sterilized
- packaging: 1 pc./paper-foil
- indication
- allows you to give the endotracheal tube the right curvature • introduced to the endotracheal tube facilitates intubation in routine
- and difficult cases, when the laryngeal inlet is not visible





REF	size	diameter	length	tube size	intermediate packaging	bulk packaging
PRI-06	6 Fr	2,0 mm	310 mm (±10)	from 2,0 to 4,0	20 pcs	20 x 20 pcs
PRI-10	10 Fr	3,3 mm	390 mm (±10)	from 4,5 to 6,5	20 pcs	20 x 20 pcs
PRI-14	14 Fr	4,7 mm	390 mm (±10)	from 7,0 to 10,0	20 pcs	20 x 20 pcs
		,	. ,	, ,		

### tracheal tube introducer

#### curved, sterile

#### features

- made of medical grade low density polyethylene
- curved and rounded tip enables safe intubation without damaging the soft tissues
- optimum flexibility and rigidity facilitate correct insertion and safe removal of the introducer after intubation
- · imprint of the size, diameter and manufacturer's name on the device
- bright yellow body for increased visibility of the introducer in the tracheal tube canal
- clear, easy to read graduation every 1 cm
- available in three sizes and lengths
- single-use
- PVC free
- latex-free
- phthalate-free .
- EO sterilized
- packaging: 1 pc./paper-foil

#### indication • for use during difficult intubation pr

REF	size	diameter O.D.	total length	tube size	graduation	intermediate packaging
PTIZ-06-535	6 Fr	2,0 mm	535 mm	from 2,0 to 4,0	10 – 45 mm	10 pcs
PTIZ-10-600	10 Fr	3,3 mm	600 mm	from 4,5 to 6,0	10 – 50 mm	10 pcs
PTIZ-10-800	10 Fr	3,3 mm	800 mm	from 4,5 to 6,0	10 – 60 mm	10 pcs
PTIZ-15-600	15 Fr	5,0 mm	600 mm	from 6,5 to 10,0	10 – 50 mm	10 pcs
PTIZ-15-800	15 Fr	5,0 mm	800 mm	from 6,5 to 10,0	10 – 60 mm	10 pcs
	REF PTIZ-06-535 PTIZ-10-600 PTIZ-10-800 PTIZ-15-600	REF         size           PTIZ-06-535         6 Fr           PTIZ-10-600         10 Fr           PTIZ-10-800         10 Fr           PTIZ-10-500         15 Fr	REF         size         diameter O.D.           PTIZ-06-535         6 Fr         2,0 mm           PTIZ-10-600         10 Fr         3,3 mm           PTIZ-10-800         10 Fr         3,3 mm           PTIZ-15-600         15 Fr         5,0 mm	REF         size         diameter O.D.         total length           PTIZ-06-535         6 Fr         2,0 mm         535 mm           PTIZ-10-600         10 Fr         3,3 mm         600 mm           PTIZ-10-800         10 Fr         3,3 mm         800 mm           PTIZ-15-600         15 Fr         5,0 mm         600 mm	REF         size         diameter O.D.         total length         tube size           PTIZ-06-535         6 Fr         2,0 mm         535 mm         from 2,0 to 4,0           PTIZ-10-600         10 Fr         3,3 mm         600 mm         from 4,5 to 6,0           PTIZ-10-800         10 Fr         3,3 mm         800 mm         from 4,5 to 6,0           PTIZ-15-600         15 Fr         5,0 mm         600 mm         from 6,5 to 10,0	O.D.         length         size           PTIZ-06-535         6 Fr         2,0 mm         535 mm         from 2,0 to 4,0         10 - 45 mm           PTIZ-10-600         10 Fr         3,3 mm         600 mm         from 4,5 to 6,0         10 - 50 mm           PTIZ-10-800         10 Fr         3,3 mm         800 mm         from 4,5 to 6,0         10 - 60 mm           PTIZ-15-600         15 Fr         5,0 mm         600 mm         from 6,5 to 10,0         10 - 50 mm

intermediate

### endotracheal tube

		packaging	
uncuffed, sterile	RIB-20	I.D. (mm) 2.0	20 pcs
	RIB-25	I.D. (mm) 2.5	20 pcs
features	RIB-30	I.D. (mm) 3.0	20 pcs
made of flexible PVC	RIB-35	I.D. (mm) 3.5	20 pcs
<ul> <li>radiopaque line over the entire tube length</li> <li>equipped with lateral Murphy eye</li> </ul>	RIB-40	I.D. (mm) 4.0	20 pcs
<ul> <li>graduated every 1 cm</li> </ul>	RIB-45	I.D. (mm) 4.5	20 pcs
double depth marker to control the position of the tube	RIB-50	I.D. (mm) 5.0	20 pcs
<ul> <li>tube size marked in 3 places (on the connector and in two</li> </ul>	RIB-55	I.D. (mm) 5.5	20 pcs
places on the tube body)	RIB-60	I.D. (mm) 6.0	20 pcs
<ul> <li>detachable connector with an outer diameter of 15 mm</li> <li>single-use</li> </ul>	RIB-65	I.D. (mm) 6.5	20 pcs
<ul> <li>Ingle-use</li> <li>latex-free</li> </ul>	RIB-70	I.D. (mm) 7.0	20 pcs
• EO sterilized	RIB-75	I.D. (mm) 7.5	20 pcs
<ul> <li>packaging: 1 pc./paper-foil</li> </ul>	RIB-80	I.D. (mm) 8.0	20 pcs
indication	RIB-85	I.D. (mm) 8.5	20 pcs
for both oral and nasal intubation	RIB-90	I.D. (mm) 9.0	20 pcs
<ul> <li>inserted into the trachea to prevent airway obstruction and</li> </ul>	RIB-95	I.D. (mm) 9.5	20 pcs
respiratory failure	RIB-100	I.D. (mm) 10.0	20 pcs

REF

size



Ø bulk packaging: 10 x 20 pcs

### endotracheal tube

cuffed, sterile

#### features

#### • made of flexible PVC

- equipped with lateral Murphy eye and radiopaque line over the entire tube length
   graduated every 1 cm, double depth marker to control the position of the tube

- low-pressure, high-volume cuff
   tube size marked in 4 places (on the connector, on the pilot balloon and in two places on the tube body)
- detachable connector with an outer diameter of 15 mm
- pilot balloon with LOT and size indication
- single-use
- latex-free • EO sterilized
- packaging: 1 pc./paper-foil

#### indication

- for both oral and nasal intubation
- · inserted into the trachea to prevent airway obstruction and respiratory failure

REF	size	intermediate packaging
RIM-30	I.D. (mm) 3.0	20 pcs
RIM-35	I.D. (mm) 3.5	20 pcs
RIM-40	I.D. (mm) 4.0	20 pcs
RIM-45	I.D. (mm) 4.5	20 pcs
RIM-50	I.D. (mm) 5.0	20 pcs
RIM-55	I.D. (mm) 5.5	20 pcs
RIM-60	I.D. (mm) 6.0	20 pcs
RIM-65	I.D. (mm) 6.5	20 pcs
RIM-70	I.D. (mm) 7.0	20 pcs
RIM-75	I.D. (mm) 7.5	20 pcs
RIM-80	I.D. (mm) 8.0	20 pcs
RIM-85	I.D. (mm) 8.5	20 pcs
RIM-90	I.D. (mm) 9.0	20 pcs
RIM-95	I.D. (mm) 9.5	20 pcs
RIM-100	I.D. (mm) 10.0	20 pcs



### reinforced endotracheal tube

#### cuffed, sterile

#### features

- made of flexible PVC
- reinforcing within the wall of the entire tube
- equipped with lateral Murphy eye and radiopaque line over the entire tube length ٠
- graduated every 1 cm, double depth marker to control the position of the tube
- tube size marked in 4 places (on the connector, on the pilot balloon and in two ٠
- places on the tube body) • indissiociable connector with an outer diameter of 15 mm
- single-use
- latex-free
- EO sterilized
- packaging: 1 pc./paper-foil

- for both oral and nasal intubation
- inserted into the trachea to prevent airway obstruction and respiratory failure • for intubation of the patient in an unusual position, with a high risk of kinking the endotracheal tube

 REF	size	intermediate packaging
RIMZP-30	I.D. (mm) 3.0	10 pcs
RIMZP-35	I.D. (mm) 3.5	10 pcs
RIMZP-40	I.D. (mm) 4.0	10 pcs
RIMZP-45	I.D. (mm) 4.5	10 pcs
RIMZP-50	I.D. (mm) 5.0	10 pcs
RIMZP-55	I.D. (mm) 5.5	10 pcs
RIMZP-60	I.D. (mm) 6.0	10 pcs
RIMZP-65	I.D. (mm) 6.5	10 pcs
RIMZP-70	I.D. (mm) 7.0	10 pcs
RIMZP-75	I.D. (mm) 7.5	10 pcs
RIMZP-80	I.D. (mm) 8.0	10 pcs
RIMZP-85	I.D. (mm) 8.5	10 pcs
RIMZP-90	I.D. (mm) 9.0	10 pcs
RIMZP-95	I.D. (mm) 9.5	10 pcs
RIMZP-100	I.D. (mm) 10.0	10 pcs











### endotracheal tube holder

#### sterile

#### features

- made of high-quality materials
- universal Velcro fasteners enabling fast and effective fixing of the holder around the head
- adjustable handle allows for easy and fast fixation of every size of endotracheal tube
- special design of the bite-piece prevents patient trauma or damaging the tube (by biting)
- additional opening allows for chyme, saliva and blood suction from the oral cavity
- single-use

58

- latex-free, phthalate-free
- EO sterilizedpackaging: 1 pc./paper-foil
- packaging, 1 pc./paper

#### indication

 to stabilize the position of the endotracheal tube after its insertion into the trachea

# REF version colour SRI-N horizontal fixation blue-red SRI-A vertical fixation green

bulk packaging: 100 x 1 pc.





### tracheostomy tube

#### cuffed, sterile

#### features

- made of medical, flexible PVC
- swivel 15 mm connector
- flexible neckplates with slits for strap insertion to secure tube position
- soft, low-pressure, thin-walled cuff
- pilot balloon in contrasting colour clearly indicates the level of the cuff's inflation
- 2 fixation straps protected with foil
- clear size marking on the neckplates atraumatic, soft rounded tip, curved cannula
- radiopaque line
- equipped with obturator
- single-use
- latex-free
- EO sterilized
- packaging: 1 pc./rigid blister

- to control respiratory tracts in people after tracheostomy
- introduced to the airways by a surgical incision of the tracheal's
- front wall to allow ventilation • provides free breathing bypassing the nose, throat and larynx

REF	size	intermediate packaging
RTM-40	4.0	10 pcs
RTM-45	4.5	10 pcs
RTM-50	5.0	10 pcs
RTM-60	6.0	10 pcs
RTM-65	6.5	10 pcs
RTM-70	7.0	10 pcs
RTM-75	7.5	10 pcs
RTM-80	8.0	10 pcs
RTM-85	8.5	10 pcs
RTM-90	9.0	10 pcs
RTM-95	9.5	10 pcs
RTM-100	10.0	10 pcs







### dual-limb breathing circuit

sterile

#### features

- made of high quality EVA material, without PVC
- 15F/22M and 22F universal connector fits compatible products
- and devices for mechanical ventilation or anesthesia machine
- expandable arms with maximum length of 160 cm adjustable length of circuit enabling adjustment of length of respiratory system and optimization of dead space
- equipped with latex-free breathing bag: 3l (for adults) or 1l (for children)
  additional expandable tube, length: 100 cm
- elbow with luer lock port
- single-use
- latex-free, phthalate-free
- EO sterilized • packaging: 1 pc./foil
- indication

#### • for patient's mechanical ventilation

- it enables transporting breathing mixture from equipment (ventilator, anesthesia machine) to respiratory protection device (e.g. tracheal tube, tracheostomy mask) and vice versa

REF	purpose	capacity of bag	
OBWD-1601L	for children	1 liter	bulk packaging:
OBWD-1603L	for adults	3 liters	40 x 1 pc.



### catheter mount with double swivel elbow connector sterile

- features
- smooth bore type made of medical PVC
- expandable type made of high quality PE, EVA
- universal patient connection: 22M/15F
- ventilation device side connector: 15M or 22F ٠ expandable version characterized by shape memorizing
- and by the ability to control the position of the connector double swivel elbow connector .
- . port for suctioning and bronchoscope secured with a stopper
- · available in versions: smooth bore and expandable
- single-use
- latex-free
- EO sterilized
- packaging: 1 pc./paper-foil

- used to connect the tracheostomy tube or endotracheal tube with the breathing circuit of the ventilation device
- provides greater mobility and rotation of the connection, increasing patient comfort and safety









REF	type	connection	length (without connectors)	intermediate packaging
PM-15-G	smooth bore (not expandable)	22M/15F-15M	15 cm	10 pcs
PM-22-G	smooth bore (not expandable)	22M/15F-22F	15 cm	10 pcs
PMR-15-K	expandable	22M/15F-15M	8-15 cm	10 pcs
PMR-22-K	expandable	22M/15F-22F	8-15 cm	10 pcs



### oropharyngeal airway/guedel airway

sterile

#### features

- made of high quality PE (polyethylene), free of PVC
- distinctive anatomical curve of airway tube allows for the opening of the upper respiratory tract and the free flow of anaesthetic gases, air
- colour coded sizes
- numeric marking of the size on the airwaysingle-use
- **60** s
  - latex-free, phthalate-free
    - EO sterilized
    - packaging: 1 pc./foil

#### indication

- used during artificial ventilation procedures to ensure the airway patency
- as a tool for paramedic to assist in the management of a casualty's
- airway which increases the effectiveness of the rescue operation
- facilitates the stabilization of the endotracheal tube and protects it against biting

REF	size	number	colour	intermediate packaging
RG-40	40 mm	000	pink	50 pcs
RG-50	50 mm	00	dark blue	50 pcs
RG-60	60 mm	0	black	50 pcs
RG-70	70 mm	1	white	50 pcs
RG-80	80 mm	2	green	50 pcs
RG-90	90 mm	3	yellow	50 pcs
RG-100	100 mm	4	red	50 pcs
RG-110	110 mm	5	blue	50 pcs
RG-120	120 mm	6	orange	50 pcs



bulk packaging: 10 x 50 pcs

bulk packaging: 100 x 1 pc.

### air cushion anaesthetic mask

#### non-sterile

#### features

- due to the high quality of materials used, ventilation is conducted in an atraumatic and safe way for the patient
- anatomically shaped inflatable cuff, ensures perfect tightness of the mask with minimal pressure
- device, when applied to the face, tightly covers the patient's mouth and nose without compressing the eyes
- clear shell provides continuous visual observation of the patient
- colour-coded retaining rings, easy for size recognition
- retaining ring can be removed when it is not needed
- full size range allows to properly choose the size of the mask depending on the age of the patient
- all sizes are equipped with inflatable check valve to adjust air amount for best cushion fitting auto patient's face
- single-uselatex-free, phthalate-free
- non-sterile
- packaging: 1 pc./foil

- intended for the administration of a gas mixture coming through devices for artificial ventilation
- enables performing medical procedures including artificial ventilation, administration of anesthetic gas and carrying out emergency treatment of an unconscious patient

REF	size	age of patient/ size	colour
MAP-0	0	neonate	-
MAP-1	1	infant	-
MAP-2	2	pediatric	yellow
MAP-3	3	adult S	green
MAP-4	4	adult M	red
MAP-5	5	adult L	blue
MAP-6	6	adult XL	orange





# **GENERAL SURGERY**



W W W . Z A R Y S . C O M

### Yankauer surgical suction cannula

standard CH21, sterile

#### features

62

- completely transparent for full possibility of observing the sucked secretion
- light and perfectly balanced handle ensures high comfort of work and reduces hand fatigue
- perfectly smooth inner surface prevents blocking or retaining the sucked content, ensuring correct and optimal flow
- cannula with double bend of curvature
- four lateral eyes, which prevent the tip from sticking to the tissues and one central opening
- universal stepped connector ensures tight connection with a wide range of suction connecting tubes
- single-use
- latex-free, phthalate-freeEO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

#### indication

• evacuating secretions and body fluids from the operating field in general surgery

REF	features	packaging	intermediate packaging
KO-00204	without suction control	1 pc.	50 pcs
KO-00214	with suction control	1 pc.	50 pcs





### Yankauer surgical suction cannula

#### bulb tip CH23, sterile

#### features

- cannula with bulb tip with four lateral eyes, which prevent the tip from sticking to the tissues and one central opening
- completely transparent for full possibility of observing the sucked secretion
   light and perfectly balanced handle ensures high comfort of work and
- reduces hand fatigue
- perfectly smooth inner surface prevents blocking or retaining the sucked content, ensuring correct and optimal flow
- cannula with double bend of curvature
- universal stepped connector ensures tight connection with a wide range of suction connecting tubes
- single-use
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

#### indication

• evacuating secretions and body fluids from the operating field in general surgery

REF	features	packaging	intermediate packaging
KO-00104	without suction control	1 pc.	50 pcs
KO-00114	with suction control	1 pc.	50 pcs







### surgical suction tube

#### sterile

#### features


- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content from settling on the walls, ensuring correct and optimal flow
- connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device (for the diameter of 8 mm to 18 mm)
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection
- single-use
- latex-free EO sterilized
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

#### indication

- device designed for connecting the suction cannula to medical suction device
- Kapkon connector designed to connect the catheter for aspiration of upper respiratory tract to medical suction device

REF	size	lenght of tube	connector
DO-24201	CH24	210 cm	funnel-funnel
DO-24301	CH24	300 cm	funnel-funnel
DO-30201	CH30	210 cm	funnel-funnel
DO-30301	CH30	300 cm	funnel-funnel
DO-24204	CH24	210 cm	funnel-funnel cut-to-fit
DO-24202	CH24	210 cm	funnel-Kapkon







63



### surgical suction set

### with Yankauer cannula, standard, sterile

#### CANNULA

- completely transparent for full possibility of observing the sucked secretion
  light and perfectly balanced handle ensures high comfort of work and reduces
- hand fatigue and particle instruction and construction in the factor of the sector of
- perfectly smooth inner surface prevents blocking or retaining the sucked content, ensuring correct and optimal flow
- cannula with double bend of curvature
- four lateral eyes, which prevent the tip from sticking to the tissues and one central opening
- universal stepped connector ensures tight connection with a wide range of suction connecting tubes

#### TUBE

64

- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content from settling on the walls, ensuring correct and optimal flow
  connector type funnel-funnel cut-to-fit with the possibility of cutting in the
- connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device (for the diameter of 8 mm to 18 mm)
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection

#### features

- available with or without suction control
- single-use
- latex-free
- EO sterilized
- packaging: 1 set/double packaging: internal- foil, external- paper-foil

#### indication

· evacuating secretions and body fluids from the operating field in general surgery

REF	size	lenght of tube	connector
ZO-212400104	CH21/CH24	210 cm	without suction control funnel-funnel
ZO-212400114	CH21/CH24	210 cm	without suction control funnel-funnel
ZO-402423104	CH21/CH24	400 cm	without suction control funnel-funnel
ZO-402423114	CH21/CH24	400 cm	without suction control funnel-funnel
ZO-2124001041	CH21/CH24	210 cm	without suction control funnel-funnel cut-to-fit
ZO-2124001141	CH21/CH24	210 cm	without suction control funnel-funnel cut-to-fit



bulk packaging: 50 x 1 set

### surgical suction set

#### with Yankauer cannula, bulb tip, sterile

#### CANNULA

- cannula with bulb tip with four lateral eyes, which prevent the tip from sticking to the tissues and one central opening
- completely transparent for full possibility of observing the sucked secretion light and perfectly balanced handle ensures high comfort of work
- and reduces hand fatigue perfectly smooth inner surface prevents blocking or retaining the sucked content, ensuring correct and optimal flow
- cannula with double bend of curvature
- universal stepped connector ensures tight connection with a wide range of suction connecting tubes

#### TUBE

- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content
- from settling on the walls, ensuring correct and optimal flow
  connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device (for the diameter of 8 mm to 18 mm)
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection

#### features

- available with or without suction control
- single-use
- latex-free
- EO sterilized
- packaging: 1 set/double packaging: internal- foil, external- paper-foil

#### indication

· evacuating secretions and body fluids from the operating field in general surgery

REF	size	lenght of tube	connector
ZO-2124002041	CH23/CH24	210 cm	with suction control, funnel-funnel cut-to-fit
ZO-2124002141	CH23/CH24	210 cm	with suction control, funnel-funnel cut-to-fit
ZO-3024002041	CH23/CH24	300 cm	with suction control, funnel-funnel cut-to-fit
ZO-3024002141	CH23/CH24	300 cm	with suction control, funnel-funnel cut-to-fit







# **OTHER DISPOSABLE PRODUCTS**



W W W . Z A R Y S . C O M

### latex-free tourniquet

#### non-sterile

#### features

- made of a wide stretch TPE strap (thermoplastic elastomer)
  latex-free material prevents allergic reactions and skin irritation
- high tensile strength
- single-use bands reduce the risk of infection
  total roller length 11,25 m
- single-use
- latex-free, phthalate-free individual packaging contains: 25 tourniquets
  non-sterile
- packaging: 1 roll/paper dispenser that allows convenient sharing of perforated tourniquets

#### indication

• used to induce venous stasis and to make a vein selected for puncture more prominent

REF	size	colour	
S001	25 mm x 450 mm	blue	bulk packaging:
S002	25 mm x 450 mm	pink	100 rolls



# automatic tourniquet

#### non-sterile

#### features

- components: blue stretchy strap and white ABS plastic buckle
- length 480 mm +/- 10 mm, width 25 mm +/- 10 mm
- wide, stretchable band put above the venipuncture site reduces . the pressure intensity
- · protects against the congestion of the arm, minimizes unpleasant sensations, improves the quality of the blood sample
- reusable
- contains latex, phthalate-free
- non-sterile
- packaging: 1 pc./double packaging: internal- foil, external- box

#### indication

 used to induce venous stasis and to make a vein selected for puncture more prominent

REF	intermediate packaging	
SA-001	21	bulk packaging:
SA-001	21 pcs	20 x 21 pcs





### dental brush

#### with suction, non-sterile

#### features

- made of polypropylene
- soft bristles on one side, foam on the other side
- total length 18 cm
- length of the cleaning part 2,5 cm
- suction opening from both the bristle side and the space between
- the foam and the bristles • controlled suction connector at 45° angle for convenient manipulation
- single-use

68

- latex-free, phthalate-free
- non-sterile
- packaging: 1 pc./foil

#### indication

- hygiene of the patient's oral cavity
- function for simultaneous brushing, rinsing and suction

REF	intermediate packaging	
HJM-SO	50 pcs	bu



oral swab

non-sterile

#### features

- overall length 15 cm (+/- 2 mm)
- length of the foam part 2,5 cm
- handle made of polyester
- swab made of polypropylene
- single-use
- latex-free, phthalate-free
- non-sterilepackaging: 1 pc./foil

#### indication

- intended for oral hygiene in patients partially or completely dependent on medical personnel
- allows the oral mucosa to be kept in proper humidity
- allows you to perform a procedure to prevent pathological changes in the mouth

#### 

REF	packaging	
HJM-GB	50 pcs	bulk packaging: 10 x 50 pcs

### tongue depressor

#### wooden

#### features

- made of birch wood
- rounded edges preventing injuries
- natural colour
- long outreach thanks to appropriate length
- durable
- single-use
- available in sterile or non-sterile version
- packaging: sterile 1 pc./paper, non-sterile 100 pcs/box


#### indication

diagnostic examination of the oral cavity

REF	version	dimensions	bulk packaging	
SLN-001	non-sterile	150 mm x 17 mm	50 x 100 pcs	
SLS-001	sterile	150 mm x 17 mm	20 x 100 pcs	conallhealth 袋 Your Partner In Healthcare







### Wash glove

molton

#### features

- wash glove in the form of a glove to put on a handmade of MOLTON material

- ultrasonic sewing
  soft, smooth, gentle and absorbent
  dry, without cleaning agents
  tear resistant

- single-use

14			
		-	-

REF	grammage	dimensions	unit packaging	bulk packaging
MB11	75 g/m²	160 mm x 220 mm	50 pcs	20 x 50 pcs

# Wash glove

#### spunlace

#### features

- wash glove in the form of a glove to put on a hand. Made of MOLTON material
- made of SPUNLACE material
- ultrasonic sewingvery soft, gentle and absorbent
- dry, without cleaning agents tear resistant
- single-use

 REF	grammage	dimensions	unit packaging	bulk packaging
MB12	80 g/m <sup>2</sup>	160 mm x 220 mm	50 pcs	20 x 50 pcs
MDIZ	60 g/111	100 11111 X 220 11111	50 pcs	20 x 30 pcs



### Wash glove

#### with film lining

#### features

- wash glove with film lining in the form of a glove to put on a hand
- made of Airlaid material + PE
- ultrasonic sewing
- gentle and absorbent • dry, without cleaning agents
- tear resistant
- single-use
- REF dimensions unit packaging grammage bulk  $\bigotimes$ packaging MB13PE 80 g/m<sup>2</sup> 160 mm x 220 mm 20 x 50 pcs 50 pcs



### **ID** wristband

#### non-sterile

#### features

- durable wristband coating characterized by high humidity resistance and resistance to the activity of soaps, hand disinfectants and other liquids used in hospitalsmade of soft and delicate vinyl
- flexible, delicate, very light and easy to put on
  rounded edges prevent from injuring

70

- Hogh tensile strength
  colour marking depending on the purpose
  single-use fastening resistant to accidental break
  hypoallergenic, latex-free, phthalate-free
- available in a version with a description field directly on the band and with an insert card
- non-sterile
- packaging: 100 pcs/box

#### indication

· for identification of patients in hospital wards

REF	purpose	colour	version	intermediate packaging
OPD-B-01	for children and adults	white	insert card	100 pcs
OPN-R-01	for infants	pink	insert card	100 pcs
OPN-N-01	for infants	blue	insert card	100 pcs

bulk packaging: 50 x 100 pcs

### entrance mat

non-sterile

#### features

- blue mat consisting of 30 layers thin sticky sheets
  each of the sticky sheets made of polyethylene resin
- does not contain heavy metals
- excellent adhesive properties
- the first and last sheets secured with a non-adhesive foil
- low mat profile: app. 2 mm thick
- the bottom surface adheres very well to the floor over the entire surface, preventing adverse movement of the mat, allowing trouble-free passage of trolleys with equipment or beds, without the risk of curling up the mat or tearing • no need to clean the mat - as soon as the surface of the mat is contaminated,

st bulk packaging: 5 x 10 pcs

- it should be torn off to uncover a new, clean sheet there is numbering of sheets on a white non-adhesive material
- in the mat corner facilitate quick identification and easy removal of subsequent sheets
- light, easy to spread
- non-sterile

#### indication

• preventing spreading of contamination brought on shoe soles and wheelchairs to operating rooms

REF	dimensions	intermediate packaging
M001	45 cm x 115 cm	10 pcs





### **PE gloves**

#### non-sterile

#### features

- made of polyethylenecoarse surface and tear-resistant weld
- ambidextrous
- elongation min. 130%, tensile strength min. 10 MPa
  approved for contact with food
  single-use

- non-sterile
- packaging: 100 pcs/foil bag enabling individual removal of gloves

#### indication

• intended for use as a contribution to other gloves (double-gloving) and for minimally invasive procedures such as drug preparation, feeding the patient

REF	size	length	
BS-RF-M	M (women)	280 mm	bulk pack
BS-RF-L	L (men)	300 mm	100 × 10



### electronic thermometer

#### non-sterile

#### features

- digital thermometer
- fast, safe and reliable temperature readings
- waterproof
- measurement in 10 seconds range 32.0°C-42.9°C (90.0°F-109.9°F)
- accuracy ±0.1°C(±0.2°F)
- has a soft, flexible measuring tip
- sound signal
- automatic turning off
- memory of the last measurement
- low battery indicator
- battery life: approximately 200 hours • package: 1 piece / plastic cover and paper box

#### indication

· for measuring body temperature in the mouth, armpit and rectal area



vitahed TERMOMETR













# Dressings

Dressing materials have for years been the core of patient care. Constant development in this field causes more advanced products to appear next to the basic ones only. Dressing does not translate into using only the basic compresses and tapes any more. Naturally the function of each dressing depends on the requirements, overall patient health and the conditions in which the product is to be used.

The offer of ZARYS Inernational Group falls into a few categories and encompasses a complex range of solutions such as:

 gauze and non-woven dressings, compresses and drapes whose main function is to absorb the exhudate and dress the wound in hospital and home care conditions,

- adhesive tapes made of various materials used for dressing and drain fixation as well as dressings with absorbent pads that enable fast dressing of minor injuries,
- special purpose dressings composing of solutions dedicated to specific conditions such as the dressing of post-operative wounds, cannula and drain fixation,
- supporting dressings- including products such as knitted and elastic bandages and tapes,
- stiffening dressings- including padding and P.O.P. bandages.





# TRADITIONAL DRESSINGS MATERIALS



W W W . Z A R Y S . C O M

# GAUZE AND NON-WOVEN DRESSINGS

# DRESSING GAUZE

### features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- soft and very absorbent
- increase in absorption capacity with increasing number of gauze threads
- homogeneous structure through a strong twist of fibers eliminates the risk of loose fibers on the surface

### indication

• general wound dressing in out-patient clinics, surgical sites, procedure department



# dressing gauze, zig-zag folded

# **GAZA lux**

non-	-sterile <sup>in</sup>	Before use on open wounds, the device must be st in accordance with a validated sterilization method sterilized with steam, ethylene oxide or radiation.		
REF	dimensions	no. of threads	packaging unit foil bag	bulk packaging carton
GR1309100	0,9 m x 100 m	13	1 pc.	5 x 1 pc.
GR1709100	0,9 m x 100 m	17	1 pc.	5 x 1 pc.
n				



edges ES-type folded \* 🔇 available on request

# cutting gauze

### **GAZA lux**

can be cut and formed

non-sterile		Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with steam, ethylene oxide or radiation.				
REF	size	no. of threads	packaging unit paper-foil	bulk packaging carton		
*CG1325SQM	1/4 m <sup>2</sup>	13	1 pc.	100 x 1 pc.		
*CG1305SQM	1/2 m²	13	1 pc.	50 x 1 pc.		
CG131SQM	1 m²	13	1 pc.	25 x 1 pc.		
* CG1725SQM	1/4 m <sup>2</sup>	17	1 pc.	100 x 1 pc.		
* CG1705SQM	1/2 m <sup>2</sup>	17	1 pc.	50 x 1 pc.		
CG171SQM	1 m²	17	1 pc.	25 x 1 pc.		





REF	size	no. of threads	packaging unit paper-foil	bulk packaging carton
CG1325SQM-SS	1/4 m²	13	1 pc.	100 x 1 pc.
CG1305SQM-SS	1/2m <sup>2</sup>	13	1 pc.	50 x 1 pc.
CG131SQM-SS	1 m²	13	1 pc.	25 x 1 pc.
CG1725SQM-SS	1/4 m²	17	1 pc.	100 x 1 pc.
CG1705SQM-SS	1/2 m²	17	1 pc.	50 x 1 pc.
CG171SQM-SS	1 m²	17	1 pc.	25 x 1 pc.

75

Sterilization method: • -S – EO sterilized product - on request • -SS – steam sterilized product

# **GAUZE SWABS**

### features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- edges ES-type folded, which prevents fraying of loose threads
- available with a sewed in X-ray detectable thread, which enables fast and non-invasive identification of the swab left in the operating field

### indication

• general wound dressing in out-patient clinics, surgical sites, procedure departments



Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with steam or ethylene oxide.



# gauze swabs

# **KOMPRI lux**



### without X-ray detectable thread

REF	dimensions	no. of threads	no. of plies	packaging unit paper bag	bulk packaging carton
GS1305081	5 cm x 5 cm	13	8	100 pcs	250 x 100 pcs
GS1375081	7,5 cm x 7,5 cm	13	8	100 pcs	60 x 100 pcs
GS1310081	10 cm x 10 cm	13	8	100 pcs	40 x 100 pcs
GS1305121	5 cm x 5 cm	13	12	100 pcs	80 x 100 pcs
GS1305125	5 cm x 5 cm	13	12	500 pcs	4 x 500 pcs
GS1375121	7,5 cm x 7,5 cm	13	12	100 pcs	40 x 100 pcs
GS1375125	7,5 cm x 7,5 cm	13	12	500 pcs	4 x 500 pcs
GS1310121	10 cm x 10 cm	13	12	100 pcs	20 x 100 pcs
GS1310125	10 cm x 10 cm	13	12	500 pcs	4 x 500 pcs
GS1705081	5 cm x 5 cm	17	8	100 pcs	80 x 100 pcs
GS1775081	7,5 cm x 7,5 cm	17	8	100 pcs	90 x 100 pcs
GS1710081	10 cm x 10 cm	17	8	100 pcs	40 x 100 pcs
GS171020081	10 cm x 20 cm	17	8	100 pcs	40 x 100 pcs
GS1705121	5 cm x 5 cm	17	12	100 pcs	96 x 100 pcs
GS1705125	5 cm x 5 cm	17	12	500 pcs	4 x 500 pcs
GS1775121	7,5 cm x 7,5 cm	17	12	100 pcs	60 x 100 pcs
GS1775125	7,5 cm x 7,5 cm	17	12	500 pcs	4 x 500 pcs
GS1710121	10 cm x 10 cm	17	12	100 pcs	20 x 100 pcs
GS1710125	10 cm x 10 cm	17	12	500 pcs	4 x 500 pcs
GS171020121	10 cm x 20 cm	17	12	100 pcs	20 x 100 pcs
GS1705161	5 cm x 5 cm	17	16	100 pcs	96 x 100 pcs
GS1775161	7,5 cm x 7,5 cm	17	16	100 pcs	60 x 100 pcs
GS1710161	10 cm x 10 cm	17	16	100 pcs	40 x 100 pcs
GS171020161	10 cm x 20 cm	17	16	100 pcs	20 x 100 pcs

### with X-ray detectable thread

REF	dimensions	no. of threads	no. of plies	packaging unit paper bag	bulk packaging carton
• GS1705081X	5 cm x 5 cm	17	8	100 pcs	80 x 100 pcs
GS1775081X	7,5 cm x 7,5 cm	17	8	100 pcs	90 x 100 pcs
GS1710081X	10 cm x 10 cm	17	8	100 pcs	40 x 100 pcs
• GS1705121X	5 cm x 5 cm	17	12	100 pcs	96 x 100 pcs
GS1775121X	7,5 cm x 7,5 cm	17	12	100 pcs	60 x 100 pcs
GS1710121X	10 cm x 10 cm	17	12	100 pcs	20 x 100 pcs
GS171020121X	10 cm x 20 cm	17	12	100 pcs	20 x 100 pcs
GS1705161X	5 cm x 5 cm	17	16	100 pcs	96 x 100 pcs
GS1775161X	7,5 cm x 7,5 cm	17	16	100 pcs	60 x 100 pcs
GS1710161X	10 cm x 10 cm	17	16	100 pcs	40 x 100 pcs
GS171020161X	10 cm x 20 cm	17	16	100 pcs	20 x 100 pcs





KOMPRI Iu	vitaheath	- h
I Internet	KOMPRI Iu	x S yearsouth
I Million		KOMPRE lux S
	10 T	
¥ ===	9=	
((	C(	Company and

### without X-ray detectable thread

REF	dimensions	no. of threads	no. of plies	packaging unit paper-foil	bulk packaging box
GS1705081 <mark>-S</mark>	5 cm x 5 cm	17	8	1 pc.	50 x 1 pc.
GS1705082-SS	5 cm x 5 cm	17	8	2 pcs	50 x 2 pcs
GS1705083-SS	5 cm x 5 cm	17	8	3 pcs	50 x 3 pcs
GS1705085-SS	5 cm x 5 cm	17	8	5 pcs	20 x 5 pcs
GS17050810-SS	5 cm x 5 cm	17	8	10 pcs	40 x 10 pcs
GS17050820-SS	5 cm x 5 cm	17	8	20 pcs	20 x 20 pcs
GS1775081 <mark>-S</mark>	7,5 cm x 7,5 cm	17	8	1 pc.	50 x 1 pc.
GS1775082-SS	7,5 cm x 7,5 cm	17	8	2 pcs	50 x 2 pcs
GS1775083-SS	7,5 cm x 7,5 cm	17	8	3 pcs	50 x 3 pcs
GS1775085-SS	7,5 cm x 7,5 cm	17	8	5 pcs	20 x 5 pcs
GS17750810-SS	7,5 cm x 7,5 cm	17	8	10 pcs	40 x 10 pcs
GS17750820-SS	7,5 cm x 7,5 cm	17	8	20 pcs	20 x 20 pcs
GS1710081 <mark>-S</mark>	10 cm x 10 cm	17	8	1 pc.	50 x 1 pc.
GS1710082-SS	10 cm x 10 cm	17	8	2 pcs	50 x 2 pcs
GS1710083-SS	10 cm x 10 cm	17	8	3 pcs	50 x 3 pcs
GS1710085-SS	10 cm x 10 cm	17	8	5 pcs	50 x 5 pcs
GS17100810-SS	10 cm x 10 cm	17	8	10 pcs	40 x 10 pcs
GS17100820-SS	10 cm x 10 cm	17	8	20 pcs	10 x 20 pcs
GS1705122-SS	5 cm x 5 cm	17	12	2 pcs	50 x 2 pcs
GS1705123-SS	5 cm x 5 cm	17	12	3 pcs	50 x 3 pcs
GS1705123 <mark>-S</mark> -IP	5 cm x 5 cm	17	12	3x1 pcs	20 x 3 pcs
GS1705125-SS	5 cm x 5 cm	17	12	5 pcs	20 x 5 pcs
GS17051210-SS	5 cm x 5 cm	17	12	10 pcs	20 x 10 pcs
GS1775122-SS	7,5 cm x 7,5 cm	17	12	2 pcs	50 x 2 pcs
GS1775123-SS	7,5 cm x 7,5 cm	17	12	3 pcs	50 x 3 pcs
GS1775123 <mark>-S</mark> -IP	7,5 cm x 7,5 cm	17	12	3x1 pcs	20 x 3 pcs
GS1775125-SS	7,5 cm x 7,5 cm	17	12	5 pcs	20 x 5 pcs
GS17751210-SS	7,5 cm x 7,5 cm	17	12	10 pcs	40 x 10 pcs
GS17751220-SS	7,5 cm x 7,5 cm	17	12	20 pcs	20 x 20 pcs
GS1710123-SS	10 cm x 10 cm	17	12	3 pcs	50 x 3 pcs
GS1710123-S-IP	10 cm x 10 cm	17	12	3x1 pcs	20 x 3 pcs
GS1710125-SS	10 cm x 10 cm	17	12	5 pcs	20 x 5 pcs
GS17101210-SS	10 cm x 10 cm	17	12	10 pcs	20 x 10 pcs
GS17101220-SS	10 cm x 10 cm	17	12	20 pcs	20 x 20 pcs
GS171020123-SS	10 cm x 20 cm	17	12	3 pcs	50 x 3 pcs
GS1710201210-SS	10 cm x 20 cm	17	12	10 pcs	10 x 10 pcs
GS171020125-SS	10 cm x 20 cm	17	12	5 pcs	20 x 5 pcs
GS1705163-SS	5 cm x 5 cm	17	16	3 pcs	50 x 3 pcs
GS1705165-SS	5 cm x 5 cm	17	16	5 pcs	30 x 5 pcs
GS1775163-SS	7,5 cm x 7,5 cm	17	16	3 pcs	50 x 3 pcs
GS1775165-SS	7,5 cm x 7,5 cm	17	16	5 pcs	20 x 5 pcs
GS17751610-SS	7,5 cm x 7,5 cm	17	16	10 pcs	16 x 10 pcs
GS1710163-SS	10 cm x 10 cm	17	16	3 pcs	50 x 3 pcs
GS1710165-SS	10 cm x 10 cm	17	16	5 pcs	30 x 5 pcs
GS17101610-SS	10 cm x 10 cm	17	16	10 pcs	20 x 10 pcs
GS171020163-SS	10 cm x 20 cm	17	16	3 pcs	50 x 3 pcs
GS1710201610-SS	10 cm x 20 cm	17	16	10 pcs	10 x 10 pcs

Sterilization method:

-S – EO sterilized product - on request
-SS – steam sterilized product

### with X-ray detectable thread

REF	dimensions	no. of threads	no. of plies	packaging unit paper-foil	bulk packaging box
GS17750810-SSXR	7,5 cm x 7,5 cm	17	8	10 pcs	40 x 10 pcs
GS17750820-SSXR	7,5 cm x 7,5 cm	17	8	20 pcs	20 x 20 pcs
✔ GS17100810-SSXR	10 cm x 10 cm	17	8	10 pcs	40 x 10 pcs
* GS1705125-SSXR	5 cm x 5 cm	17	12	5 pcs	20 x 5 pcs
* GS17051210-SSXR	5 cm x 5 cm	17	12	10 pcs	40 x 10 pcs
✔ GS1775125-SSXR	7,5cm x 7,5cm	17	12	5 pcs	20 x 5 pcs
✔ GS17751210-SSXR	7,5cm x 7,5cm	17	12	10 pcs	10 x 10 pcs
✔ GS17101210-SSXR	10cm x 10cm	17	12	10 pcs	20 x 10 pcs
✔ GS17101220-SSXR	10cm x 10cm	17	12	20 pcs	20 x 20 pcs
✔ GS17101240-SSXR	10cm x 10cm	17	12	40 pcs	10 x 40 pcs
✔ GS1710201210-SSXR	10cm x 20cm	17	12	10 pcs	10 x 10 pcs
✔ GS1710201220-SSXR	10cm x 20cm	17	12	20 pcs	6 x 20 pcs
* GS17051610-SSXR	5 cm x 5 cm	17	16	10 pcs	30 x 10 pcs
✔ GS1775165-SSXR	7,5 cm x 7,5 cm	17	16	5 pcs	20 x 5 pcs
✔ GS17751610-SSXR	7,5 cm x 7,5 cm	17	16	10 pcs	16 x 10 pcs
* GS1710165-SSXR	10 cm x 10 cm	17	16	5 pcs	20 x 5 pcs
✔ GS17101610-SSXR	10 cm x 10 cm	17	16	10 pcs	20 x 10 pcs



 $\checkmark$ self-adhesive labels for medical documentation

\*< available on request

# **GAUZE LAP SPONGES**

### features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- with a sewed in X-ray detectable thread or chip, which enables fast and non-invasive identification of the lap sponge left in the operating field
- with a belt to fix the lap sponge beyond the operating field

### indication

• as a surgical invasive dressing material for intraoperative absorption of large amounts of blood and body fluids from surgical wounds



Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with steam or ethylene oxide.

# gauze lap sponges





### with X-ray thread and belt

REF	dimensions	no. of threads	no. of plies	packaging unit foil bag	bulk packaging carton
LS1745454-XB	45 cm x 45 cm	17	4	20 pcs	10 x 20 pcs
LS1745704-XB	45 cm x 70 cm	17	4	20 pcs	10 x 20 pcs
LS1745456-XB	45 cm x 45 cm	17	6	20 pcs	10 x 20 pcs
LS1745706-XB	45 cm x 70 cm	17	6	20 pcs	10 x 20 pcs
LS1745458-XB	45 cm x 45 cm	17	8	20 pcs	10 x 20 pcs
LS2045454-XB	45 cm x 45 cm	20	4	20 pcs	10 x 20 pcs

### pre-washed with X-ray chip and belt

REF	dimensions	no. of threads	no. of plies	packaging unit foil bag	bulk packaging carton
LS2045454-XB-PW	45 cm x 45 cm	20	4	20 pcs	10 x 20 pcs



The advantages of self-adhesive labels for medical documentation:

- fast data transport without the need to rewrite,
- to rewrite, safe identification of the device -elimination of the possibility of errors, improvement of the preparation of medical documentation.

# **SERVI lux S**



### with X-ray thread and belt

REF	dimensions	no. of threads	no. of plies	packaging unit paper-foil	bulk packaging box
✓ LS17454542-SSXB	45 cm x 45 cm	17	4	2 pcs	20 x 2 pcs
✓ LS17454545-SSXB	45 cm x 45 cm	17	4	5 pcs	8 x 5 pcs
✓ LS17457042-SSXB	45 cm x 70 cm	17	4	2 pcs	20 x 2 pcs
✓ LS17454563-SSXB	45 cm x 45 cm	17	6	3 pcs	10 x 3 pcs
✓ LS17454582-SSXB	45 cm x 45 cm	17	8	2 pcs	10 x 2 pcs
✓ LS17454585-SSXB	45 cm x 45 cm	17	8	5 pcs	5 x 5 pcs
✓ LS17457043-SSXB	45 cm x 70 cm	17	4	3 pcs	10 x 3 pcs
✓ LS17759041-SSXB	75 cm x 90 cm	17	4	1 pc.	15 x 1 pcs





gauze balls

# **TUPFER lux**

non-sterile

### features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- available in a version with or without X-ray detectable thread, which enables fast and non-invasive identification of the gauze ball left in the operating field

### indication

- absorption of blood and body fluids from surgical wounds and hard to reach places

### without X-ray thread

REF	size of gauze	shape	packaging unit foil bag	bulk packaging carton
TK1515250	15 cm x 15 cm	ball	250 pcs	24 x 250 pcs
TK2020250	20 cm x 20 cm	ball	250 pcs	24 x 250 pcs
TK3030100	30 cm x 30 cm	ball	100 pcs	10 x 100 pcs

### with X-ray thread

REF	size of gauze	shape	packaging unit foil bag	bulk packaging carton
TK1515250X	15 cm x 15 cm	ball	250 pcs	24 x 250 pcs
TK2020250X	20 cm x 20 cm	ball	250 pcs	24 x 250 pcs
TK2525250X	25 cm x 25 cm	ball	250 pcs	8 x 250 pcs
TK3030100X	30 cm x 30 cm	ball	100 pcs	10 x 100 pcs
TK3040100X	30 cm x 40 cm	ball	100 pcs	10 x 100 pcs
TK5050100X	50 cm x 50 cm	ball	100 pcs	10 x 100 pcs



# **TUPFER lux S**

sterile

TERIL



REF	size of gauze	shape	packaging unit paper-foil	bulk packaging carton
✓ TF959505-SSXR	9,5 cm x 9,5 cm	bean	5 pcs	50 x 5 pcs
✓ TF959510-SSXR	9,5 cm x 9,5 cm	bean	10 pcs	25 x 10 pcs
✓ TF151505-SSXR	15 cm x 15 cm	bean	5 pcs	25 x 5 pcs
✓ TF151510-SSXR	15 cm x 15 cm	bean	10 pcs	20 x 10 pcs
✓ TR121205-SSXR	12 cm x 12 cm	triangular	5 pcs	15 x 5 pcs
✓ TR181810-SSXR	18 cm x 18 cm	triangular	10 pcs	10 x 10 pcs
✓ TK151510-SSXR	15 cm x 15 cm	ball	10 pcs	20 x 10 pcs
✓ TK202010-SSXR	20 cm x 20 cm	ball	10 pcs	20 x 10 pcs
✓ TK303010-SSXR	30 cm x 30 cm	ball	10 pcs	10 x 10 pcs



# non-woven swabs

### **NONVI lux**

non-sterile

Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with steam, ethylene oxide or radiation.

REF	dimensions	grammage	packaging unit paper bag	bulk packaging carton
NL50-100	5 cm x 5 cm	30 g/m <sup>2</sup>	100 pcs	50 x 100 pcs
NL50-100-40	5 cm x 5 cm	40 g/m <sup>2</sup>	100 pcs	50 x 100 pcs
NL75-100	7,5 cm x 7,5 cm	30 g/m <sup>2</sup>	100 pcs	40 x 100 pcs
NL75-100-40	7,5 cm x 7,5 cm	40 g/m <sup>2</sup>	100 pcs	40 x 100 pcs
NL100-100	10 cm x 10 cm	30 g/m <sup>2</sup>	100 pcs	40 x 100 pcs
NL100-100-40	10 cm x 10 cm	40 g/m <sup>2</sup>	100 pcs	20 x 100 pcs
NL1020-100	10 cm x 20 cm	30 g/m <sup>2</sup>	100 pcs	20 x 100 pcs

# NONVI lux S



REF	dimensions	grammage	packaging unit paper-foil	bulk packaging box
NL052- <mark>S</mark>	5 cm x 5 cm	30 g/m <sup>2</sup>	2 pcs	25 x 2 pcs
NL053- <mark>S</mark>	5 cm x 5 cm	30 g/m <sup>2</sup>	3 pcs	25 x 3 pcs
NL055- <mark>S</mark>	5 cm x 5 cm	30 g/m <sup>2</sup>	5 pcs	25 x 5 pcs
NL0510- <mark>S</mark>	5 cm x 5 cm	30 g/m <sup>2</sup>	10 pcs	10 x 10 pcs
NL052- <mark>S</mark> -40	5 cm x 5 cm	40 g/m <sup>2</sup>	2 pcs	25 x 2 pcs
NL053- <mark>S</mark> -40	5 cm x 5 cm	40 g/m <sup>2</sup>	3 pcs	25 x 3 pcs
NL055- <mark>S</mark> -40	5 cm x 5 cm	40 g/m <sup>2</sup>	5 pcs	25 x 5 pcs
NL0510- <mark>S</mark> -40	5 cm x 5 cm	40 g/m <sup>2</sup>	10 pcs	10 x 10 pcs
NL752- <mark>S</mark>	7,5 cm x 7,5 cm	30 g/m <sup>2</sup>	2 pcs	25 x 2 pcs
NL753- <mark>S</mark>	7,5 cm x 7,5 cm	30 g/m <sup>2</sup>	3 pcs	25 x 3 pcs
NL755- <mark>S</mark>	7,5 cm x 7,5 cm	30 g/m <sup>2</sup>	5 pcs	25 x 5 pcs
NL7510- <mark>S</mark>	7,5 cm x 7,5 cm	30 g/m <sup>2</sup>	10 pcs	10 x 10 pcs
NL752- <mark>S</mark> -40	7,5 cm x 7,5 cm	40 g/m <sup>2</sup>	2 pcs	25 x 2 pcs
NL753- <mark>S</mark> -40	7,5 cm x 7,5 cm	40 g/m <sup>2</sup>	3 pcs	25 x 3 pcs
NL755- <mark>S</mark> -40	7,5 cm x 7,5 cm	40 g/m <sup>2</sup>	5 pcs	25 x 5 pcs
NL7510- <mark>S</mark> -40	7,5 cm x 7,5 cm	40 g/m <sup>2</sup>	10 pcs	10 x 10 pcs
NL7520- <mark>S</mark> -40	7,5 cm x 7,5 cm	40 g/m <sup>2</sup>	20 pcs	5 x 20 pcs
NL102- <mark>S</mark>	10 cm x 10 cm	30 g/m <sup>2</sup>	2 pcs	25 x 2 pcs
NL103- <mark>5</mark>	10 cm x 10 cm	30 g/m <sup>2</sup>	3 pcs	25 x 3 pcs
NL105- <mark>S</mark>	10 cm x 10 cm	30 g/m <sup>2</sup>	5 pcs	20 x 5 pcs
NL1010- <mark>S</mark>	10 cm x 10 cm	30 g/m <sup>2</sup>	10 pcs	10 x 10 pcs
NL102- <mark>S</mark> -40	10 cm x 10 cm	40 g/m <sup>2</sup>	2 pcs	25 x 2 pcs
NL103- <mark>S</mark> -40	10 cm x 10 cm	40 g/m <sup>2</sup>	3 pcs	25 x 3 pcs
NL105- <mark>S</mark> -40	10 cm x 10 cm	40 g/m <sup>2</sup>	5 pcs	20 x 5 pcs
NL1010- <mark>S</mark> -40	10 cm x 10 cm	40 g/m <sup>2</sup>	10 pcs	10 x 10 pcs
NL1020- <mark>S</mark> -40	10 cm x 10 cm	40 g/m <sup>2</sup>	20 pcs	5 x 20 pcs

in out-patient clinics, surgical sites,

**NON-WOVEN SWABS** 

• hydrophilic medical non-woven,

• soft, non-fraying, durable

general wound dressing

procedure department



features

30 or 40 g/m<sup>2</sup> • 4-plies

indication

Sterilization method:

• -S - EO sterilized product

-SS – steam sterilized product





# **NONVI lux S with Y-incision**



### features

- hydrophilic medical non-woven 40 g/m<sup>2</sup>
- 4-plies
  Y-incision for precise, aseptic protection of the outlet of various types of medical drains

### indication

• securing the outlet of drains, catheters, tracheostomy tubes etc.

REF sterilized by EO	REF sterilized by steam	dimensions	packaging unit paper-foil	bulk packaging box
NL105-S-40Y	* NL105-SS-40Y	10 cm x 10 cm	5 pcs	20 x 5 pcs


# **ABSORBENT WOUND DRESSINGS**

### features

- layered structure: non-woven cellulose cellulose pulp
- increased absorption capacity due to the high absorbent cellulose layer
- bottom layer of the dressing covered with waterproof blue non-woven

### indication

- removing the exudate from strongly oozing woundscovering wounds after surgical procedures





# ABD pads, multi-layer

# **MULTIabsorb**

STERILE NON	-sterile <sup>ir</sup>		inds, the device must be sterilized idated sterilization method. It can ne oxide or radiation.
REF	dimensions	packaging unit foil bag	bulk packaging carton
ABD1010	10 cm x 10 cm	50 pcs	40 x 50 pcs
ABD1020	10 cm x 20 cm	50 pcs	20 x 50 pcs
ABD2025	20 cm x 25 cm	50 pcs	12 x 50 pcs



# **MULTIabsorb S**

$\langle \rangle$	
(STERILE)	sterile
$\langle \rangle$	Sterne

REF	dimensions	packaging unit paper-foil	intermediate packaging box	bulk packaging carton
ABD1010-S	10 cm x 10 cm	1 pc.	25 x 1 pc.	24 x 25 pcs
ABD1020-S	10 cm x 20 cm	1 pc.	25 x 1 pc.	18 x 25 pcs
ABD2025-S	20 cm x 25 cm	1 pc.	15 x 1 pc.	18 x 15 pcs
ABD2040-S	20 cm x 40 cm	1 pc.	8 x 1 pc.	12 x 8 pcs

# **EYE DRESSINGS**

### features

- layered structure: non-woven + cellulose pulp
  needs to be fixed with a secondary dressing e.g. with adhesive tape or tubular net bandage
  hypoallergenic

### indication

- covering an eye and its area after ophthalmic procedures and surgical procedures
  sight correction in strabismus therapy



# eye pads

# elastoLUMENAL

STERLE N	on-sterile i	Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can sterilized with ethylene oxide or radiation.		
REF	dimensions	packaging unit paper-paper	intermediate packaging box	bulk packaging carton
803311	5,5 cm x 7,5 cm	1 pc.	50 x 1 pc.	24 x 50 pcs



Conforms well to facial contours thanks to the oval shape.

# elastoLUMENAL S

sterile

REF	dimensions	packaging unit paper-paper	intermediate packaging box	bulk packaging carton
803310	5,5 cm x 7,5 cm	1 pc.	50 x 1 pc.	24 x 50 pcs



# ADHESIVE TAPES AND PLASTERS WITH DRESSING



W W W . Z A R Y S . C O M



# **ADHESIVE TAPES**

### features

- flexible, conforms to the body's contours, guarantees freedom of movement
- · characterized by adequate air and water vapuor permeability
- good adhesion
- · can be divided lengthwise and widthwise without the use of scissors
- hypoallergenic

### indication


- stabilization, protection and fixation of primary, non-adhesive dressings, applied directly to wounds and skin damages
- stabilization, protection and fixation of drains, catheters, medical tubes that do not require specialized fixing in the place of their insertion





# adhesive fabric tape, white

# SENSIplast







# adhesive fabric tape, beige colour

### **PLASTIplast**

STERILE NON-S	sterile		
REF	dimensions	packaging unit box	bulk packaging carton
SZO12505S	1,25 cm x 5 m	24 pcs	30 x 24 pcs
SZO25005S	2,5 cm x 5 m	12 pcs	30 x 12 pcs
SZO50005S	5 cm x 5 m	6 pcs	30 x 6 pcs
$\begin{array}{c} I + I + I \\ I + I + I \\ I + I + I \\ I + I +$	serrated edge	on a roll	beige colour





# adhesive non-woven tape

# SOFTplast




REF	dimensions	packaging unit box	bulk packaging carton
SNW12505W	1,25 cm x 5 m	24 pcs	30 x 24 pcs
SNW25005W	2,5 cm x 5 m	12 pcs	30 x 12 pcs
SNW50005W	5 cm x 5 m	6 pcs	30 x 6 pcs
SNW12509W	1,25 cm x 9,14 m	24 pcs	30 x 24 pcs
SNW25009W	2,5 cm x 9,14 m	12 pcs	30 x 12 pcs
SNW50009W	5 cm x 9,14 m	6 pcs	30 x 6 pcs
$\langle \rangle$			$\diamond$
made of non-woven fabric	straight edge	on a roll	white colour





# adhesive silk tape

# SILKplast















# adhesive PE tape

# FILMplast

EF	dimensions	packaging unit box	bulk packaging carton
PE12505W	1,25 cm x 5 m	24 pcs	30 x 24 pcs
PE25005W	2,5 cm x 5 m	12 pcs	30 x 12 pcs
PE50005W	5 cm x 5 m	6 pcs	30 x 6 pcs
PE12509W	1,25 cm x 9,14 m	24 pcs	30 x 24 pcs
PE25009W	2,5 cm x 9,14 m	12 pcs	30 x 12 pcs
PE50009W	5 cm x 9,14 m	6 pcs	30 x 6 pcs
made of poly- ethylene foil	straight edge	on a roll	kransparent
perforation for easy dividing	<u>↑↑↑↑</u> i : : i air permeabl structure	e exchange b	us structure permea vapor, allowing prop oetween the foil and ents maceration of t





# adhesive perforated non-woven tape

	dimensions	packaging unit box	bulk packaging carton
PW25009W	2,5 cm x 9,14 m	12 pcs	30 x 12 pcs
made of perforated non-woven fabric	straight edge	on a roll	white



# PLASTERS WITH DRESSING

### features

- self-adhesive plasters with absorbent pads
- centrally situated absorbent pad absorbs a small and medium amount of exudate and protects from external factors
- · absorbent pad coated with a polyethylene mesh that prevents sticking to the wound
- characterized by adequate air and water vapuor permeability
- flexible, conforms to the body's contours, guarantees freedom of movement
- silicone backing paper facilitates precise, painless and efficient application
- hypoallergenic


### indication

• protection of small skin injuries, cuts, abrasions and lacerations



# non-woven plaster with dressing

# NONVIplast



REF	dimensions	packaging unit box	bulk packaging carton
NP6100	6 cm x 1 m	1 pc.	200 x 1 pc.
NP6500	6 cm x 5 m	1 pc.	100 x 1 pc.
NP8100	8 cm x 1 m	1 pc.	200 x 1 pc.
NP8500	8 cm x 5 m	1 pc.	100 x 1 pc.





# fabric plaster with dressing

# FABRIplast

STERILE NON-	sterile				
REF	dimensions	pack box	aging unit		k packaging ton
FP6100	6 cm x 1 m	1 pc.		20	0 x 1 pc.
FP6500	6 cm x 5 m	1 pc.		10	0 x 1 pc.
FP8100	8 cm x 1 m	1 pc.		20	0 x 1 pc.
FP8500	8 cm x 5 m	1 pc.		10	0 x 1 pc.
	<u> </u>		of		$\diamond$
made of cotton fabric	air permeabl structure	e	possibility of cutting		beige colour



# SPECIAL PURPOSE DRESSINGS



W W W . Z A R Y S . C O M

# SELF-ADHESIVE NON-WOVEN DRESSINGS

# DRESSING TAPE

### features

- made of hydrophobic non-woven
- covered with hypoallergenic acrylic glue
- micropores in the structure of the non-woven provide adequate vapor permeability, thanks to which
  a proper gas exchange between the dressing and the skin takes place
- flexible, conforms to the body's contours, guarantees freedom of movement
- backing layer made of silicone paper
- · backing paper with metric scale facilitates the division of the tape and precise, painless and efficient application
- packed in a cardboard box a dispenser that allows dosing without removing the tape from the packaging

### indication

- stabilization, protection and fixation of primary, non-adhesive dressings, applied directly to wounds and skin damages
- stabilization, protection and fixation of drains, catheters, medical tubes









colour

# adhesive non-woven dressing tape


## elastopor

non-sterile

REF         dimensions         packaging unit box         bulk packagin carton           806009         2,5 cm x 10 m         1 pc.         200 pcs           806010         5 cm x 10 m         1 pc.         100 pcs           806011         10 cm x 10 m         1 pc.         50 pcs	
806010         5 cm x 10 m         1 pc.         100 pcs           806011         10 cm x 10 m         1 pc.         50 pcs	g
806011 10 cm x 10 m 1 pc. 50 pcs	
806012 15 cm x 10 m 1 pc. 50 pcs	
806013 20 cm x 10 m 1 pc. 25 pcs	
806014 25 cm x 10 m 1 pc. 25 pcs	
806015 30 cm x 10 m 1 pc. 25 pcs	



of backing paper





of

possibility of

cutting

<u>↑↑↑</u> i i i air permeable structure

white

# elastopor E

STERILE NON	-sterile		
REF	dimensions	packaging unit box	bulk packaging carton
806009E	2,5 cm x 10 m	1 pc.	200 pcs
806010E	5 cm x 10 m	1 pc.	100 pcs
806011E	10 cm x 10 m	1 pc.	50 pcs
806012E	15 cm x 10 m	1 pc.	50 pcs
806013E	20 cm x 10 m	1 pc.	25 pcs
806014E	25 cm x 10 m	1 pc.	25 pcs
806015E	30 cm x 10 m	1 pc.	25 pcs





# SELF-ADHESIVE DRESSINGS WITH ABSORBENT PAD

### features

90

- self-adhesive dressings with centrally situated absorbent pad, made of hydrophobic non-woven, covered with acrylic glue • absorbent pad coated with a polyethylene mesh that prevents sticking to the wound, absorbs a small and medium amount
- of exudate and protects against the influence of external factors
- rounded edges prevent accidental removal of the dressing
- silicone backing paper facilitates the precise, painless and efficient application
- micropores in the structure of the non-woven provide adequate vapor permeability, thanks to which a proper gas exchange
- between the dressing and the skin takes place
- hypoallergenic



### made of hydrophobic non-woven





 $\uparrow \uparrow \uparrow$ 

air permeable

structure

**BESTSELLER** 

# non-woven dressing with absorbent pad

## elastopor STERIL



### indication

especially recommended for covering surgical wounds

· for treating wounds with a small and medium amount of exudate

REF	dimensions	packaging unit paper-paper	bulk packaging box
801001	5 cm x 7,2 cm	1 pc.	100 x 1 pc.
801012	6 cm x 10 cm	1 pc.	50 x 1 pc.
801002	8 cm x 10 cm	1 pc.	30 x 1 pc.
801013	10 cm x 10 cm	1 pc.	30 x 1 pc.
801021	10 cm x 12 cm	1 pc.	30 x 1 pc.
801014	10 cm x 15 cm	1 pc.	30 x 1 pc.
801003	6 cm x 15 cm	1 pc.	50 x 1 pc.
801004	8 cm x 15 cm	1 pc.	30 x 1 pc.
801005	8 cm x 20 cm	1 pc.	30 x 1 pc.
801015	10 cm x 20 cm	1 pc.	25 x 1 pc.
801016	10 cm x 25 cm	1 pc.	25 x 1 pc.
801017	10 cm x 30 cm	1 pc.	25 x 1 pc.
801018	10 cm x 35 cm	1 pc.	25 x 1 pc.





BESTSELLER

# non-woven IV cannula dressing

### elastopor IV

### sterile

### features

- · centrally situated absorbent pad protects the non-woven material from sticking to the injection site
- non-woven lining should be placed under cannula's wings, which prevents it from causing damage to the skin

### indication

- dressing application around the inserted cannula
- · stabilization, protection and fixation of intravenous cannulas

REF	dimensions	packaging unit paper-paper	bulk packaging box
811013	6 cm x 8 cm	1 pc.	100 x 1 pc.
811012	5,1 cm x 7,6 cm	1 pc.	100 x 1 pc.
811023	6 cm x 8 cm	1 pc.	50 x 1 pc.
811022	5,1 cm x 7,6 cm	1 pc.	50 x 1 pc.





ተ ተ ተ

air permeable structure





REF	dimensions	packaging unit paper-paper	bulk packaging box
801031	9 cm x 10 cm	1 pc.	30 x 1 pc.
801030	12 cm x 14 cm	1 pc.	25 x 1 pc.

· for treating wounds with a small and medium amount of exudate

non-woven dressing with absorbent

pad, incision and O-hole

 dressing application around the inserted drain securing the outlet of drains and catheters

elastopor STERIL D

sterile

indication

STERILE



### Conforms well to facial contours thanks to the oval shape



# non-woven eye dressing with absorbent pad

# elastopor EYE

### sterile

### indication

- covering an eye and its area after ophthalmic procedures and surgical procedures
- sight correction in strabismus therapy

REF	dimensions	shape	colour	packaging unit paper-paper	bulk packaging box
803010	5,8 cm x 8,3 cm		absorbent pad - white	1 pc.	50 x 1 pc.
803011	6,5 cm x 9,4 cm		non-woven material - white	1 pc.	50 x 1 pc.
803212	5 cm x 7,5 cm		absorbent pad - white	1 pc.	50 x 1 pc.
			non-woven material - beige colour		

# **HOW TO CHOOSE A DRESSING PROPERLY?**

When choosing the proper type of dressing, pay special attention to:

- absorption rate and absorption capacity,
- vapor permeability,
- possibility of an atraumatic dressing change,
- safety of use.

# SELF-ADHESIVE DRESSINGS FOR KIDS

### features

- self-adhesive dressing with absorbent pad, made of hydrophobic non-woven material with color overprint, covered by hypoallergenic glue
- absorbent pad coated with polyethylene mesh to prevent sticking to the wound, absorbing small and medium amounts of exudate and protecting against external influences
- rounded edges to prevent accidental removal of the dressing
- protective layer of silicone paper for precise, painless and effective application
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed



# non-woven IV cannula dressing

# elastoKIDS IV for state features entrally situated absorbent pad protects the non-woven material from sticking to the injection site non-woven lining should be placed under cannula's wings, which prevents it from causing damage to the skin indication application of the dressing around the cannula stabilization, protection and fixation of intravenous cannulas

REF	dimensions	packaging unit paper-paper	bulk packaging box	
811001A	5,1 cm x 7,6 cm	1 pc.	100 x 1 pc.	



# non-woven dressing with absorbent pad

# elastoKIDS STERIL



### indication

especially recommended for covering surgical wounds

• treating wounds with a small and medium amount of exudate

REF	dimensions	packaging unit paper-paper	bulk packaging box
801012A	6 cm x 10 cm	1 pc.	50 x 1 pc.





<u>↑↑↑</u> iiii air permeable structure

beige colour

# non-woven nasal cannula fixation plaster

## elastoNASAL

non-sterile

### features

- flexible fixation plaster made of non-woven
- covered with acrylic glue
  adequate vapor permeability, thanks to which a proper gas exchange between
- the dressing and the skin takes placethin and flexible plaster is comfortable for the patient, conforms to the body's contours
- "two-stage" application type facilitates precise and efficient application
- self-adhesive
- hypoallergenic

### indication

• stabilization of catheters and stomach tubes fixed in the nasal area

REF	dimensions	purpose	packaging unit box	
811014	7 cm x 7,1 cm	for adults	50 x 1 pc.	
811015	5,9 cm x 6,2 cm	for children	50 x 1 pc.	
811016	3,6 cm x 4 cm	for infants	50 x 1 pc.	



blue-white colour

air permeable

structure

# adhesive tape for catheter tubing fixation

# elastoFIXAL

non-sterile

### features

- two-part fixing tape consisting of a part attached to the skin and a part fastening the drain
- part attached to the patient's skin made of non-woven fabric
- drain fixing part is integrated with a part adhered to the patient's skin and has an additional adhesive tape - Velcro and adhesive field for better drain stabilization
- possibility of reposition and correction of the drain position thanks to Velcro
- covered with hypoallergenic glue

### indication

stabilization of various types of drains and catheters

REF	dimensions	packaging unit paper-foil	bulk packaging box	
816101	9 cm x 4 cm	1 pc.	50 x 1 pc.	







# white

# adhesive tape for endotracheal tube

# elastoFIXAL Set

non-sterile

### features

- set consisting of 2 self-adhesive non-woven pads and a fastening strap
- self-adhesive non-woven pads attached to the patient's cheeks with an adhesive part with Velcro on the other side
- a strap for stabilization of the endotracheal or oropharyngeal tube made of nonwoven fabric which should be attached to the velcro layer with which the pads are covered
- adhesive area in the middle of the strip for better stabilization of the medical tube
- the possibility of repositioning and correcting the position of the tube thanks to the Velcro
- pads and strap covered with hypoallergenic glue
- has excellent adhesive properties and high mechanical strength which prevents premature and accidental detachment from the skin

### indication

long-term stabilization of the endotracheal or oropharyngeal tube

REF	kit composition	packaging unit paper-foil	bulk packaging
816102	2 self-adhesive non-woven pads 1 fastening strap	1 set	100x1 set







# band for tracheostomy tube fixation

# elastoFIXAL TUBE

non-sterile

### features

 stabilising band consisting of two parts made of gentle hydrophobic non-woven fabric in flesh-coloured

- on each part there are fastening strips with Velcro on one side enabling to attac
- it to the fastening wings of the tracheostomy tube

 shorter part of the stabilising band is additionally equipped with a second Velcro strap attached to an elastic white rubber band which makes it possible to adjust and determine the total length of the band stabilising the tracheostomy tube

- hypoallergenic
- single use
- non-sterile
- shelf life: 5 years

REF	dimensions	packaging unit paper-foil	bulk packaging
EFT2520	25cm x 2,0 cm	1 set	30x30 pcs
EFT4325	43cm x 2,5 cm	1 set	30x20 pcs



# STABILIZING DRESSINGS



W W W . Z A R Y S . C O M



# SUPPORTING BANDAGES

### features

- delicate, light, air and X-ray permeable
- conforms good to the body's contours
- hypoallergenic

### indication

• fixation of all types of non-adhesive dressings



# knitted supporting bandage elastoBAND BASIC

STERILE NON-	sterile			
features <ul> <li>contains</li> <li>single-us</li> </ul>				
REF	dimensions	packaging unit paper	bulk packaging carton	
OD05400	5 cm x 4 m	1 pc.	250 x 1 pc.	
OD10400	10 cm x 4 m	1 pc.	150 x 1 pc.	
OD15400	15 cm x 4 m	1 pc.	100 x 1 pc.	
0010100	10 0	2 000	100 x 1 pc.	



# elastic supporting bandage, woven elastoBAND FLEX

### non-sterile

### features

- supporting bandage, elastic, woven
- high cotton content makes the bandage safe for the skin and well-tolerated by people with allergic tendency
- equipped with a clip for fixing

### indication

- fixation of all types of non-adhesive dressings, especially where
- a non-elastic bandages can slide downin compression therapy to facilitate correct blood circulation
- stabilization at bone-articular system damage, e.g. dislocation, joint twisting,
- weakness of muscles, tendons and ligaments
- as a sports bandage
- securing and strengthening integument after surgery

REF - cardboard box	REF - paper pouch	dimensions	packaging unit box/paper	bulk packaging carton
OE08400-E	OE08400P-E	8 cm x 4 m	1 pc.	124 x 1 pc.
OE10400-E	OE10400P-E	10 cm x 4 m	1 pc.	108 x 1 pc.
OE12400-E	OE12400P-E	12 cm x 4 m	1 pc.	96 x 1 pc.
OE15400-E	OE15400P-E	15 cm x 4 m	1 pc.	72 x 1 pc.
OE08500-E	OE08500P-E	8 cm x 5 m	1 pc.	124 x 1 pc.
OE10500-E	OE10500P-E	10 cm x 5 m	1 pc.	108 x 1 pc.
OE12500-E	OE12500P-E	12 cm x 5 m	1 pc.	96 x 1 pc.
OE15500-E	OE15500P-E	15 cm x 5 m	1 pc.	72 x 1 pc.





# elastic tubular net bandage

# elastoNET

(STERILE) non-sterile

### features

- elastic net bandage in the form of a tube
- contains polyamide and elastomer
  high stretchability and strength of returning to pre-stretch state
- elastic and very soft
- easy and quick to use
- length before stretching:
- 1 m, stretched: 2,5 m
- 10 m, stretched: 25 m

### indication

- fixation of all dressings on the patient's body, replacing the traditional bandage
- the use of tubular net bandage reduces the time of wound dressing and makes
  - it easier to fix the dressing to the patient's body

REF	lenght before stretching	size	indication	packaging unit box - dispenser
SO-0101	1 m	1	finger	1 pc.
SO-0110	10 m	1	finger	1 pc.
SO-0201	1 m	2	hand, toe	1 pc.
SO-0210	10m	2	hand, toe	1 pc.
SO-0301	1 m	3	hand, foot	1 pc.
SO-0310	10 m	3	hand, foot	1 pc.
SO-0401	1 m	4	lower leg, knee, arm, foot, elbow	1 pc.
SO-0410	10 m	4	lower leg, knee, arm, foot, elbow	1 pc.
SO-0601	1 m	6	head, shoulder, lower leg, knee	1 pc.
SO-0610	10 m	6	head, shoulder, lower leg, knee	1 pc.
SO-0801	1 m	8	thigh, head, hips	1 pc.
SO-0810	10 m	8	thigh, head, hips	1 pc.
SO-1001	1 m	10	hips, trunk	1 pc.
SO-1010	10 m	10	hips, trunk	1 pc.
SO-1401	1 m	14	chest, trunk	1 pc.
SO-1410	10 m	14	chest, trunk	1 pc.


### FIRST AID FOR DRESSING WOUNDS USING BANDAGES

Hippocrates' achievements in treating head wounds were so widespread and advanced that they are still used today. The Hippocrates cap also called "mitra is a headgear that allows you to hold the dressing and stop bleeding of the head wound. It is made of two bandage tapes. The main purpose of using the Hippocrates cap is to hold the dressing directly protecting the bleeding head wound, seal and warm the head, hold the compress applied to the skull area, compression, strengthening and immobilization after surgery and additional protection against the external environment. This solution is used in the absence of access to dressings such as elastic dressings and scarves.

# How to make a HIPOCRATES CAP?

### **RULES FOR BANDAGING**

- 1. Prepare a two-head bandage about 5-6 cm wide by sewing together two bands or by twisting the free end to the middle.
- 2. Put the two-headed bandage on the back of the head and drive bandage to the forehead a circular move.
- 3. Cross the bandages on the forehead and drive one through the top of the head towards the occiput by a longitudinal move.
- 4. Drive the other bandage above the ear to the occiput also (half-round).
- 5. Scroll the longitudinal course through the circular course bandage.
- 6. Drive back the longitudinal course bandage across the top of the head to the forehead.
- 7. Drive the circular course bandage around the other ear to the forehead
- 8. Continue longitudinal and circular moves, the circular course bandage wraps around the scalp
- 9. Little by little, move the longitudinal course bandage towards the temples
- 10. After covering the whole area with a two-head bandage, fix the ends of bandage





# **ORTHOPEDIC PADDINGS**

### features

- very soft and fluffy
- adhere well and adjusts to the shape of the body
- coherent structure resistant to stratification makes it easier placement of the immobilizing dressing
- minimize the risk of skin abrasions and irritations under the P.O.P. bandage
- ${\ensuremath{\cdot}}$  do not absorb moisture and thus do not allow for its contact with the skin
- hypoallergenic

### indication


- as a protective layer put under plaster of Paris bandages
- in compression therapy as a padding under the bandages of various elasticity



 $\bigcirc$ 

white

colour



# orthoNATURAL

non-sterile

REF	dimensions	packaging unit paper cover	bulk packaging foil bag	
* ON-63	6 cm x 3 m	1 pc.	12 x 1 pc.	
* ON-83	8 cm x 3 m	1 pc.	12 x 1 pc.	
ON-103	10 cm x 3 m	1 pc.	12 x 1 pc.	
ON-123	12 cm x 3 m	1 pc.	6 x 1 pc.	
ON-153	15 cm x 3 m	1 pc.	6 x 1 pc.	
* ON-203	20 cm x 3 m	1 pc.	6 x 1 pc.	
* ON-253	25 cm x 3 m	1 pc.	6 x 1 pc.	



available on request

Š





White colour

# orthopedic padding, synthetic

# orthoSYNTHETIC

non-sterile

REF	dimensions	packaging unit paper cover	bulk packaging foil bag	
OS-63	6 cm x 3 m	1 pc.	12 x 1 pc.	
OS-83	8 cm x 3 m	1 pc.	12 x 1 pc.	
OS-103	10 cm x 3 m	1 pc.	12 x 1 pc.	
OS-123	12 cm x 3 m	1 pc.	6 x 1 pc.	
OS-153	15 cm x 3 m	1 pc.	6 x 1 pc.	
OS-203	20 cm x 3 m	1 pc.	6 x 1 pc.	
OS-253	25 cm x 3 m	1 pc.	6 x 1 pc.	





White colour

# triangular sling

# elastoSLING

### non-sterile

### features

- triangle shaped sling made of non-woven or cotton
- thin, light and elastic comfortable for the patient, adjusts to the shape of the body, provides freedom of movement
- possibility of connecting several slings to achieve bigger immobilization
  hypoallergenic

### indication

• short-term immobilization and support of the upper limb in case of various injuries e.g. dislocation or fracture of collar bone, humerus, forearm or hand bone

REF	dimensions	description	packaging unit zipper pouch	bulk packaging foil bag
CHTW9696136	96 cm x 96 cm x 136 cm	non-woven	1 pc.	12 x 1 pc.
CHTB9696130	96 cm x 96 cm x 130 cm	cotton	1 pc.	12 x 1 pc.



### HOW TO PROPERLY PUT ON A TRIANGULAR SLING?

- 1. Place the triangular sling on the patient's body so that one of the corners of the sling is at the level of the elbow of the immobilized hand.
- 2. Place the second corner of the sling on the opposite arm in relation to the immobilized hand.
- 3. Carefully place the sick / injured limb on the sling.
- 4. Wrap the sick/injured limb with dangling, third corner of the sling.
- 5. Make sure that the hand is in the right position, and then tie the corners of the sling around the neck.

# P.O.P. BANDAGES

### features

- made of cotton gauze coated with natural plaster of Paris on both sides
- on a plastic roll

indication

- enable fast and durable immobilization
- before fixing need to be soaked in a container with water at 20°C temperature
- do not absorb X-rays

• large and small immobilizing dressings in orthopedic, ambulatory and surgical procedures



# P.O.P. bandage



### orthoBAND FAST

 $(\mathbb{C})$ fast-setting setting time approx. **3 minutes** fully dry after approx. **24 hours** 

lust setting				
REF	dimensions	type of edge	packaging unit foil bag	bulk packaging carton
OB06300-F	6 cm x 3 m	straight	1 pc.	48 x 1 pc.
OB08300-F	8 cm x 3 m	straight	1 pc.	48 x 1 pc.
OB10300-F	10 cm x 3 m	serrated	1 pc.	48 x 1 pc.
OB12300-F	12 cm x 3 m	serrated	1 pc.	48 x 1 pc.
OB14300-F	14 cm x 3 m	serrated	1 pc.	36 x 1 pc.
OB15300-F	15 cm x 3 m	serrated	1 pc.	36 x 1 pc.
OB20300-F	20 cm x 3 m	serrated	1 pc.	24 x 1 pc.

# orthoBAND NORMAL

fast-setting		me approx. <b>5-6 r</b> fter approx. <b>24 r</b>		
REF	dimensions	type of edge	packaging unit foil bag	bulk packaging carton
OB06300-N	6 cm x 3 m	straight	1 pc.	48 x 1 pc.
OB08300-N	8 cm x 3 m	straight	1 pc.	48 x 1 pc.
OB10300-N	10 cm x 3 m	serrated	1 pc.	48 x 1 pc.
OB12300-N	12 cm x 3 m	serrated	1 pc.	48 x 1 pc.
OB14300-N	14 cm x 3 m	serrated	1 pc.	36 x 1 pc.
OB15300-N	15 cm x 3 m	serrated	1 pc.	36 x 1 pc.
OB20300-N	20 cm x 3 m	serrated	1 pc.	24 x 1 pc.

# orthoBAND LONG



setting time approx. **8-12 minutes** fully dry after approx. **24 hours** 

### slow-setting

REF	dimensions	type of edge	packaging unit foil bag	bulk packaging carton
OB10400-L	10 cm x 4 m	serrated	1 pc.	48 x 1 pc.
OB12400-L	12 cm x 4 m	serrated	1 pc.	36 x 1 pc.
OB15400-L	15 cm x 4 m	serrated	1 pc.	24 x 1 pc.





# **RULES FOR BANDAGING**





Select a bandage of proper type and size for the part of the body to be bandaged.



With your right hand, grasp the beginning of the band, and with your left hand hold the rolled part so called "head".



Put the patient in a comfortable position supporting the bandaged lower or upper limb.



For proper stabilization of dressing, start bandaging with a circular wrap, each subsequent layer covers exactly the previous one in a durable and aesthetic way. Little by little unroll the band during bandaging. Finish the bandaging with a circular wrap, end of the dressing secure with a clip or tie it.



The bandage should adhere correctly to the patient's body and be rolled up with adequate power – too much pressure can lead to impaired circulation, too loose it will slide and fall.

Remove the bandage by transferring it from one hand to other to minimize patient movement.



# SPECIALIST DRESSINGS



W W W . Z A R Y S . C O M



soaked with paraffin

possibility

# paraffin gauze dressing

# paraffiNET

### sterile) sterile

### features

- wide-mesh dressing made of cotton gauze covered with paraffin
- structure of the carrier material allows exudate to be freely drained to the covering secondary dressing
- covering with a layer of paraffin reduces risk of adhesion to the wound and allows a painless change of dressing
- neutral properties of paraffin allow for using the dressing in combination with other preparations containing medicinal substances or antiseptics
- can be shaped and cut according to the needs aseptically
   paged to be covered with a dragging with absorbent properties a
- needs to be covered with a dressing with absorbent properties, e.g. compress, high-absorption dressing, polyurethane foam
- non-adhesive, necessary fixing with a tape, bandage or net tubular dressing
  hypoallergenic

### indication

- keeping proper moisture in the wound
- preventing the wound edges from dying and maceration
- recommended for uninfected wounds
- treatment of chronic wounds: various ulcers, bedsores
- treatment of surgical and burn wounds, skin donor and grafting sites,
- skin damage, cuts
- protection against radiation skin damage

REF	dimensions	packaging unit paper-aluminum foil	bulk packaging box	
PAR-01	5 cm x 5 cm	1 pc.	10 x 1 pc.	
PAR-02	10 cm x 10 cm	1 pc.	10 x 1 pc.	
PAR-06	10 cm x 40 cm	1 pc.	10 x 1 pc.	
PAR-07	15 cm x 20 cm	1 pc.	10 x 1 pc.	

### A NEW APPROACH TO WOUND HEALING

The most important factor affecting the effectiveness of wound healing is maintaining constant humidity within the wound. George D. Winter is a precursor of wet wound therapy, who in 1962 published his observations after experiment of creating a moist environment under polyurethane film and observing wound healing. Experience has shown that the wound heals under these conditions twice as fast as in unprotected conditions with a special film dressing.

Maintaining a moist environment in the wound stimulates the self-cleaning process, facilitates the autolysis of dead tissue, granulation growth, reduce the healing time.

It also prevents the dressing from sticking to the wound, allowing it to be changed without pain. While removing dressing there is no secondary tissue damage.

# FOIL DRESSINGS

# FOIL DRESSINGS

### features

- made of polyurethane film
- covered with hypoallergenic acrylic glue
- transparent they allow observation of the dressed place without the need to remove the dressing
- vapour permeability guarantees correct gas exchange between the dressing and the skin, and prevents it from
  excessive moisture accumulation and maceration
- barrier against microorganisms and bacteria, minimizing the risk of infection
- very thin and elastic adjust to the shape of the body, guarantee freedom of movement
- waterproof
- self-adhesive
- hypoallergenic



1. Remove the back paper protection.



2. Apply the dressing to the desired place.

# foil dressing

### sterile) sterile

### features

- "frame" or three-stage" application system makes dressing placement precise and accurate
- available in several sizes and options:
- with a central absorbent pad
- with notch enabling the dressing to fit the cannulas
- without notch and absorbent pad
- without noten and absorbent pad

### elastoDERM

dressing without absorbent pad, with "three-stage" application type

REF	dimensions	shape	indication	packaging unit paper-paper	bulk packaging box
811210	6 cm x 7 cm	17 17 17	<ul> <li>fixing and securing intravenous cannulas</li> <li>secondary dressing</li> <li>bedsore prevention, works as so-called "second skin"</li> <li>for small epidermal abrasions and damages of the skin surface</li> </ul>	1 pc.	100 x 1 pc.

### elastoDERM F

dressing without absorbent pad, with "frame" application type

REF	dimensions	shape	indication	packaging unit paper-paper	bulk packaging box
812003	4,4 cm x 4,4 cm		- fixing and securing intravenous cannulas - secondary dressing	1 pc.	100 x 1 pc.
812005	6 cm x 7 cm		<ul> <li>bedsore prevention, works as so-called "second skin"</li> </ul>	1 pc.	100 x 1 pc.
with the tape the dressing	ape to describe and damages of		<ul> <li>for small epidermal abrasions and damages of the skin surface</li> </ul>		
812008	10 cm x 12 cm		- fixing and securing central venous catheters	1 pc.	100 x 1 pc.
with the tape the dressing	e to describe		- fixing epidural catheters - secondary dressing - bedsore prevention, works		
812011	10 cm x 15 cm		as so-called "second skin"	1 pc.	100 x 1 pc.
with the tape the dressing	- for vith the tape to describe and		<ul> <li>for small epidermal abrasions and damages of the skin surface</li> </ul>		

11



3. Remove the paper protection on the outer side of the dressing starting from the cut.



4. Put the notes on the tape to describe the dressing and apply it on the dressing without covering the wound/ the place of needle insertion.

# elastoDERM F-IV

IV cannula dressing, with "frame" application type

REF	dimensions	shape	indication	packaging unit paper-paper	bulk packaging box
812006	6 cm x 7 cm		- fixing and securing peripheral cannulas in adults	1 pc.	100 x 1 pc.
with the tape t the dressing	with the tape to describe		- securing cannulas inserted into central veins		

# elastoDERM IV

IV cannula dressing, with "three-stage" application type

REF	dimensions	shape	indication	packaging unit paper-paper	bulk packaging box
811301	5,8 cm x 8 cm	<b>S</b>	<ul> <li>fixing and securing peripheral cannulas in adults</li> </ul>	1 pc.	100 x 1 pc.
with the tape f the dressing	to describe				

# elastoDERM PAD

dressing with central absorbent pad, with "three-stage" application type

REF	dimensions	shape	indication	packaging unit paper-paper	bulk packaging box		
802012	6 cm x 7 cm			1 pc.	100 x 1 pc.		
802014	9 cm x 10 cm			n x 10 cm		1 pc.	50 x 1 pc.
802015	9 cm x 15 cm		- protection of post-operative	1 pc.	50 x 1 pc.		
802016	10 cm x 20 cm			1 pc.	50 x 1 pc.		
802018	10 cm x 25 cm			1 pc.	50 x 1 pc.		





# ADDITIONAL PRODUCTS



WWW.ZARYS.COM



non-sterile

non-sterile

\*< available

on request

conallhealth 郄

# gauze pad moistened with alcohol alkoCLEAN

### features

- made of non-woven
- soaked with isopropyl alcohol (70%)
- two pouches connected in one-piece packaging can be separated thanks to perforation

### indication

• disinfection of small surfaces

REF	dimensions	no. of plies	packaging unit paper-aluminum foil	bulk packaging box
L-02	30 mm x 65 mm (unfolded) 30 mm x 32,5 mm (folded)	2	2 x 1 pc.	100 x 1 pc. (50 x 2 pcs)
*L-03	56 mm x 65 mm (unfolded) 28 mm x 32,5 mm (folded)	4	2 x 1 pc.	100 x 1 pc. (50 x 2 pcs)

gauze pad moistened with alcohol ALKOPAD

### features

- made of non-woven
- each gauze pad is packed separately, in aluminum foil paper
- soaked with isopropyl alcohol (70%)
- two pouches connected in one-piece packaging can be separated thanks to perforation (without L-06)

### indication

disinfection low surface

REF	dimensions	packaging unit paper-aluminum foil	bulk packaging box
L-04	M 32,5 mm x 27,5 mm (folded) 65 mm x 55 mm (unfolded)	2 x 1 pc.	100 x 1 pc. (50 x 2 pcs.)
L-05	XL 55 mm x 45 mm (folded) 110 mm x 90 mm (unfolded)	2 x 1 pc.	100 x 1 pc. (50 x 2 pcs.)
L-06	XXL 190 mm x 140 mm (unfolded)	2 x 1 pc.	50 x 1 pc. (25 x 2 pcs.)

# gauze pad moistened with alcohol **ALKOPAD DUO**

### features

- made of non-woven
- each gauze pad is packed separately, in aluminum foil paper
- gauze pad wet soaked with isopropyl alcohol (70%) and gauze pad dry
- two pouches connected in one-piece packaging can be separated thanks to perforation

### indication

disinfection low surface

REF	dimensions	packaging unit paper-aluminum foil	bulk packaging box	
L-05DUO	125 mm x 125 mm	2 x 1 pc.	50 x 2 pcs	

# rescue blanket, foil, termic thermCARE

### features

- made of metallized PE foil
- protects against adverse weather conditions cold, rain, heat, wind, depending on the side, which is applied towards the victim's body, prevents the life-threatening heat disorders and reduces the risk of shock
- resistant to damage

### indication

- protection against hypothermia and overheating of the victim's body
  in emergency cases such as car accidents, ski accidents, in mountain tourism, in sport

dimensions	purpose	packaging unit zipper pouch	bulk packaging foil bag
160 cm x 210 cm	for adults	1 nc	12 x 1 pc.
	dimensions		zipper pouch





non-sterile















# **Medical clothing**

It is extremely important for us to offer our clients products that ensure high standards of working conditions and comfort of use.

In order to achieve this goal, we have created a line of high quality products made of different types of non-woven materials:

- polypropylene non-woven
- SMS non-woven
- Spunlace non-woven
- viscose non-woven
- Spunlace Mesh non-woven

A combination of several factors such as optimal design, perfect fit, breathability and structure makes BETAtex products comfortable and reliable. A wide range of solutions will certainly meet our clients' individual preferences and secure both medical personnel and patients' comfort and dignity.





# MEDICAL GOWNS



W W W . Z A R Y S . C O M

# **MEDICAL GOWN**

with elastic cuffs, non-woven, non-sterile

#### features:

- made of polypropylene non-woven 20 g/m², 25 g/m², 30 g/m² or 40 g/m²  $\,$ 

- ties at waist and neck
- long sleeves with elastic cuffs
- breathable
- available sizes: L, XL • single-use

REF	grammage	size	colour	packaging unit
BT-019-20B3-L	20 g/m²	L	blue	10 pcs
BT-019-20G3-L	20 g/m <sup>2</sup>	L	green	10 pcs
BT-019-20G4-XL	20 g/m <sup>2</sup>	XL	green	10 pcs
BT-019-20W3-L	20 g/m <sup>2</sup>	L	white	10 pcs
BT-019-25G3-L	25 g/m <sup>2</sup>	L	green	10 pcs
BT-019-25G4-XL	25 g/m <sup>2</sup>	XL	green	10 pcs
BT-019-25B3-L	25 g/m <sup>2</sup>	L	blue	10 pcs
BT-019-30G3-L	30 g/m <sup>2</sup>	L	green	10 pcs
BT-019-30G4-XL	30 g/m <sup>2</sup>	XL	green	10 pcs
BT-019-40G3-L	40 g/m <sup>2</sup>	L	green	10 pcs
BT-019-40B3-L	40 g/m <sup>2</sup>	L	blue	10 pcs







# **MEDICAL GOWN**

with knitted cuffs, non-woven, non-sterile

#### features:

- made of polypropylene non-woven 20 g/m², 25 g/m², 30 g/m² or 40 g/m²  $\,$ 

- ties at waist and neck long sleeves with knitted cuffs
- breathable
- available sizes: L, XL, XXL

• single-use

REF	grammage	size	colour	packaging unit
BT-015-20G3-L	20 g/m <sup>2</sup>	L	green	10 pcs
BT-015-20G4-XL	20 g/m <sup>2</sup>	XL	green	10 pcs
BT-015-25G3-L	25 g/m²	L	green	10 pcs
BT-015-25B3-L	25 g/m²	L	blue	10 pcs
BT-015-25G4-XL	25 g/m²	XL	green	10 pcs
BT-015-30G3-L	30 g/m²	L	green	10 pcs
BT-015-30G4-XL	30 g/m <sup>2</sup>	XL	green	10 pcs
BT-015-30G5-XXL	30 g/m²	XXL	green	10 pcs
BT-015-40G3-L	40 g/m <sup>2</sup>	L	green	10 pcs
BT-015-40B3-L	40 g/m <sup>2</sup>	L	blue	10 pcs
BT-015-40G4-XL	40 g/m <sup>2</sup>	XL	green	10 pcs









# SURGICAL CLOTHING



W W W . Z A R Y S . C O M

## SURGICAL SUIT

short-sleeved shirt+trousers, non-woven, non-sterile

#### features:

 $\bullet$  made of SMS non-woven 35 g/m  $^{2}$ 

- short-sleeved shirt with V-neck and 3 pockets
- trousers with guiding ties
- blue shirt and trousers available also in individual packaging
- single-use

#### SET (SHIRT + TROUSERS) REF size colour packaging unit BT-017-S35B1-S-SMS S blue 1 set BT-017-S35B2-M-SMS М blue 1 set BT-017-S35B3-L-SMS L blue 1 set BT-017-S35B4-XL-SMS XL blue 1 set BT-017-S35B5-XXL-SMS XXL blue 1 set BT-017-S35B6-XXXL-SMS XXXL blue 1 set BT-017-S35B7-XXXXL-SMS XXXXL blue 1 set BT-017-S35G1-S-SMS S 1 set green BT-017-S35G2-M-SMS М 1 set green BT-017-S35G3-L-SMS L 1 set green BT-017-S35G4-XL-SMS XL 1 set green BT-017-S35G5-XXL-SMS XXL 1 set green BT-017-S35V1-S-SMS S purple 1 set BT-017-S35V2-M-SMS М 1 set purple BT-017-S35V3-L-SMS L 1 set purple BT-017-S35V4-XL-SMS XL purple 1 set BT-017-S35V5-XXL-SMS XXL purple 1 set





SHIRT			
REF	size	colour	packaging unit
BT-018-S35B1-S-SMS	S	blue	1 pc.
BT-018-S35B2-M-SMS	М	blue	1 pc.
BT-018-S35B3-L-SMS	L	blue	1 pc.
BT-018-S35B4-XL-SMS	XL	blue	1 pc.
BT-018-S35B5-XXL-SMS	XXL	blue	1 pc.

Ø bulk packaging: 50 x 1 pc.



### TROUSERS

REF		size	colour	packaging unit
BT-018-T35B1-5	S-SMS	S	blue	1 pc.
BT-018-T35B2-N	M-SMS	М	blue	1 pc.
BT-018-T35B3-L	-SMS	L	blue	1 pc.
BT-018-T35B4->	(L-SMS	XL	blue	1 pc.
BT-018-T35B5->	(XL-SMS	XXL	blue	1 pc.

 $\bigotimes$ bulk packaging: 50 x 1 pc.



113

### SURGICAL SHIRT SHORT-SLEEVED

non-woven, non-sterile

#### features:

- ${\rm \cdot}$  made of soft polypropylene non-woven 47  $g/m^2$
- short sleeves
- 3 big pockets
- snap fastener at the neck
- single-use

114

REF	size	packaging unit
BT-018-S47B0-XS-SOFT2	XS	1 pc.
BT-018-S47B1-S-SOFT2	S	1 pc.
BT-018-S47B2-M-SOFT2	М	1 pc.
BT-018-S47B3-L-SOFT2	L	1 pc.
BT-018-S47B4-XL-SOFT2	XL	1 pc.
BT-018-S47B5-XXL-SOFT2	XXL	1 pc.
BT-018-S47B6-XXXL-SOFT2	XXXL	1 pc.
BT-018-S47B7-XXXXL-SOFT2	XXXXL	1 pc.







### SURGICAL SHIRT LONG-SLEEVED

non-woven, non-sterile

#### features:

- ${\rm \cdot}$  made of soft polypropylene non-woven 47  $g/m^2$
- long sleeves with knitted cuffs
- 2 big pockets
- snap fasteners at the front
- single-use

REF	size	packaging unit
BT-018-047B1-S-SOFT3	S	1 pc.
BT-018-047B2-M-SOFT3	М	1 pc.
BT-018-047B3-L-SOFT3	L	1 pc.
BT-018-047B4-XL-SOFT3	XL	1 pc.
BT-018-O47B5-XXL-SOFT3	XXL	1 pc.







# SURGICAL TROUSERS, WITH TIES

non-woven, non-sterile

### features:

made of soft polypropylene non-woven 47 g/m<sup>2</sup>
 trousers with guiding ties

- 1 big pocket with snap fastener single-use

	REF	size	packaging unit
	BT-018-T47B0-XS-SOFT1	XS	1 pc.
	BT-018-T47B1-S-SOFT1	S	1 pc.
	BT-018-T47B2-M-SOFT1	Μ	1 pc.
	BT-018-T47B3-L-SOFT1	L	1 pc.
	BT-018-T47B4-XL-SOFT1	XL	1 pc.
	BT-018-T47B5-XXL-SOFT1	XXL	1 pc.
	BT-018-T47B6-XXXL-SOFT1	XXXL	1 pc.
	BT-018-T47B7-XXXXL-SOFT1	XXXXL	1 pc.









# CAPS AND PROTECIVE MASKS

ZARYS


W W W . Z A R Y S . C O M

### MEDICAL BOUFFANT CAP

non-woven

#### features:

- cap with an elastic rubber band
- diameter: approx. 53 cm
- made of airy polypropylene non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile





#### features

- cap with an elastic rubber band
- diameter: approx. 47 cm
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces

• single-use, non-sterile

REF	colour	packaging unit
BT-008-25G-B1-V	green	100 pcs
BT-008-25B-B1-V	blue	100 pcs
BT-008-25G-B1-V-XL	green	100 pcs





### LILI medical cap with rubber

#### features:

• cap with an elastic rubber band in the back part

- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile

REF	colour	packaging unit
BT-001-25G-R1-V	green	100 pcs
BT-001-25B-R1-V	blue	100 pcs
BT-001-25P-R1-V	pink	100 pcs
BT-001-25Y-R1-V	yellow	100 pcs











### POLA medical cap with rubber

#### features:

118

- cap with an elastic rubber band in the back part
- made of durable Spunlace and polypropylene non-woven
- the used non-woven is characterized by good absorption properties and durability and additionally it is pleasant to the touch, which increases comfort while using
   excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile



### SOL medical cap with ties

#### features:

- cap with ties at the back
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces

• single-use, non-sterile

REF	colour	packaging unit
BT-007-25G-T1-V	green	100 pcs
BT-007-25B-T1-V	blue	100 pcs

bulk packaging: 6 x 100 pcs







### ABI medical cap with ties

#### features:

- deep cap with ties at the back
- being folded out, the front part creates an additional layer protecting against sweat
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile

•	REF	colour	packaging unit
•	BT-007-25G-T2-V	green	100 pcs
	BT-007-25B-T2-V	blue	100 pcs
	BT-007-25P-T2-V	pink	100 pcs
	BT-007-25Y-T2-V	yellow	100 pcs

bulk packaging: 6 x 100 pcs






### NINO medical cap

#### feartures:

- cap without rubber or ties, designed for short-haired people
- two layers of non-woven fabric in the side part
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile



### KOKO long medical cap with snood

#### feartures:

- cap with snood and elastic rubber band in the back part
- optimal protection for long-haired people
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile

REF	colour	packaging unit
BT-001-25G-S1-V	green	100 pcs
BT-001-25B-S1-V	blue	100 pcs
bulk packaging: 6 x 100 pcs	for long hair	available colours:





# DORA long

medical cap with snood

#### feartures:

- cap with anti-sweat insert in the forehead area
- special design allows the cap to adjust to the shape of the head
- and keep one's hair in an excellent manner
- optimal protection for long-haired people
- made of Spunlace and polypropylene non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile





### BEN medical cap with ties

#### features:

- forage cap with ties at the back
- made of perforated Spunlace Mesh non-woven
- excellent breathability and comfort of use
- single-use, non-sterile





### TOM medical cap with ties

#### features:

- cap with ties at the back
- anti-sweat insert on the whole circumference
- made of polypropylene non-woven 20 g/m<sup>2</sup>
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile

REF	packaging unit	
BT-007-P20G	100 pcs	
bulk packaging: 10 x 100 pcs	anti-sweat insert	available colours:



### ROS medical cap with ties

#### features:

- cap with ties at the back
- anti-sweat insert in the forehead area
- made of Spunlace and polypropylene non-woven
- packaging: ergonomic cardboard box to remove individual pieces
- latex-free
- single-use, non-sterile






# MEDICAL MASK

with ties, non-sterile

#### features:

- made of three layers of dust-free and hypoallergenic non-woven
- modelling nose-piece
- protects operative field against exfoliated epidermis, hair and droplet infections
- provides effective protection in patients against infections transmitted from personnel's mouths and noses
- packaging: ergonomic cardboard box to remove individual pieces
- single-use



three layers









# MEDICAL MASK

with earloops, non-sterile

#### features:

• made of three layers of dust-free and hypoallergenic non-woven

modelling nose-piece

- splash resistant
- protects operative field against exfoliated epidermis, hair and droplet infections
- provides effective protection in patients against infections transmitted from personnel's mouths and noses

÷

available colours:

• packaging: ergonomic cardboard box to remove individual pieces

≫

three layers

• single-use


SP

bulk packaging:

20 x 50 pcs

 REF	colour	type in accordance with the standard EN 14683	packaging unit
BT-046-31G	green	type II	50 pcs
BT-046-31B	blue	type II	50 pcs
BT-046-31P	pink	type II	50 pcs
BT-046-31W	white	type II	50 pcs
BT-046-R31B	blue	type II R	50 pcs





# PATIENT CLOTHING



W W W . Z A R Y S . C O M

# PATIENT'S GOWN

tied-up, non-woven, non-sterile

#### features:

- made of SMS non-woven
- short sleeves
- white ties at waist and neck
- easy putting on and taking off in bed-ridden patients
- universal size
- single-use





# **PATIENT'S GOWN**

Y-neck, non-woven, non-sterile

### features:

- made of polypropylene or SMS non-woven
- short sleeves and Y-neck
- slipover
- single-use





### V-neck, non-woven, non-sterile

#### features:

- made of Spunlace or SMS non-woven
- ${\scriptstyle \bullet}$  the used non-woven is characterized by good absorption properties
- and durability and additionally it is pleasant to the touch, which increases

available colours:

- comfort while wearing
- short sleeves and V-neck
- slipover
- single-use

 $\bigotimes$ 

bulk packaging:

10 x 10 pcs

REF	material	size	grammage	colour	packagir unit	ng
BT-0022-45W-SPN	Spunlace non-woven	universal	45 g/m²	white	10 pcs	n pleasant
BT-0022-45B-SPN	Spunlace non-woven	universal	45 g/m²	blue	10 pcs	for the skin





# PATIENT'S GOWN

### for delivery, non-woven, non-sterile

#### features:

- made of SMS non-woven
- short sleeves, ties at waist and neck (other models available on request)

available colours:

- available in 5 sizes: S, M, L, XL, XXL
- single-use

bulk packaging:

10 x 10 pcs

### 124

REF	material	size	grammage	packaging unit
BT-0022-35B-SMS-DM-001-S	SMS non-woven	S	35 g/m²	10 pcs
BT-0022-35B-SMS-DM-001-M	SMS non-woven	Μ	35 g/m²	10 pcs
BT-0022-35B-SMS-DM-001-L	SMS non-woven	L	35 g/m²	10 pcs
BT-0022-35B-SMS-DM-001-XL	SMS non-woven	XL	35 g/m²	10 pcs
BT-0022-35B-SMS-DM-001-XXL	SMS non-woven	XXL	35 g/m²	10 pcs



### **COLONOSCOPY BOXERS**

SMS non-woven, non-sterile

#### features:

 $\bullet$  made of SMS non-woven 35 g/m  $^2$ 

- boxers with a slit in the back part
- completely opaque

bulk packaging: 10 x 10 pcs

single-use

REF	material	size	grammage	packaging unit
BT-0023-35B-SMS	SMS non-woven	universal	35 g/m²	10 pcs
BT-0023-35B-SMS-XL	SMS non-woven	XL	35 g/m²	10 pcs



# COLONOSCOPY BOXERS

available colours:

PP non-woven, non-sterile

#### features:

made of polypropylene non-woven 40 g/m<sup>2</sup>
boxers with a slit in the back part

single-use

REF     material     size     grammage     packaging unit       BT-0023-40DB     PP non-woven     universal     40 g/m <sup>2</sup> 10 pcs       Solution     Solution     Solution     Solution     Solution       bulk packaging: 10 x 10 pcs     available colours:     Solution     Solution						
bulk packaging: available colours:	REF	material	size	grammage		
bulk packaging: available colours:	BT-0023-40DB	PP non-woven	universal	40 g/m <sup>2</sup>	10 pcs	
bulk packaging: available colours:		æ				
		-				



## **GYNECOLOGICAL SKIRT**

non-woven, non-sterile

#### features:

- made of polypropylene non-woven 40 g/m<sup>2</sup>
- dimensions: 160 cm (at the waist after stretching)
- single-use

REF	material	size	grammage	packaging unit	
BT-0025-40DB	PP non-woven	universal	40 g/m <sup>2</sup>	10 pcs	
bulk packaging: 10 x 10 pcs	available colours:				



# **EXAMINATION UNDERPANTS**

single-use, non-sterile

### features:

- ${\scriptstyle \bullet}$  made of polypropylene non-woven 40 g/m²
- ensure patient's intimacy during medical procedures
- single-use





# **EXAMINATION UNDERPANTS**

single-use, non-sterile

#### features:

- made of Spunlace non-woven 45 g/m<sup>2</sup>
- the used non-woven is characterized by good absorption properties and durability and additionally it is pleasant to the touch, which increases comfort while wearing
   ensure patient's intimacy during medical procedures

single-use

	REF	material	size	grammage	packaging unit	
	BT-0026-45W-SPN	Spunlace non-woven	universal	45 g/m²	25 pcs	
	$\langle \rangle$	02				
		(all all all all all all all all all all				
1.	and a set of a set	er in Nederlanden en Leis inner				
C	oulk packaging:	available colours:				
	20 x 25 pcs	$\bigcirc$				





# SHOE COVERS

ZARYS

W W W . Z A R Y S . C O M

# **SHOE COVERS**

non-sterile

### features:

- made of polyethylene or polypropylene non-woven
- $30 \text{ g/m}^2 \text{ or } 40 \text{ g/m}^2$  (for shoe covers with anti-skid layer)
- elastic rubber band at the opening
  protection against contaminating rooms where higher cleanliness standards are required
- increased durability
- single-use

available colours:

REF	material	grammage	size	colour	packaging unit	bulk packaging	
BT-010-G	non-woven	30 g/m²	38 cm x 17 cm	green	100 pcs	10 x 100 pcs	
BT-010-B	non-woven	30 g/m <sup>2</sup>	38 cm x 17 cm	blue	100 pcs	10 x 100 pcs	
BT-011-B	non-woven with anti-skid layer	40 g/m <sup>2</sup>	38 cm x 17 cm	blue	100 pcs	10 x 100 pcs	
BT-002-PEB	polyethylene		41 cm x 15 cm	blue	100 pcs	20 x 100 pcs	<







# **SHOE COVERS**

### calf-high, non-sterile

#### features:

- ${\scriptstyle \bullet}$  made of polypropylene non-woven 30 g/m  ${\scriptstyle ^2}$
- elastic rubber band at the opening
   protection against contaminating rooms where higher cleanliness
- standards are required
- increased durability
- single-use

REF	material	grammage	size	packaging unit
BT-013-B	PP non-woven	30 g/m²	47 cm x 38,5 cm	100 pcs
BT-013-G	PP non-woven	30 g/m <sup>2</sup>	47 cm x 38,5 cm	100 pcs











# **BEDCLOTHES AND ABSORBENT UNDERPADS**



W W W . Z A R Y S . C O M

### **DISPOSABLE MEDICAL BED SHEET**

non-sterile

#### features:

- made of polypropylene non-woven, laminated polypropylene non-woven or SMS non-woven
- dimensions: according to customer's specification
- unit packaging: foil, bulk packaging: cardboard box



REF	material	grammage	dimensions*	colour
BT-030- PP-17-"X" x "Y" - G ( X - width , Y - length)	PP non-woven	17 g/m²		green
BT-030- PP-25-"X" x "Y" - G ( X - width, Y - length)	PP non-woven	25 g/m²	X cm x Y cm	green
BT-030- PP-35-"X" x "Y"- G ( X - width, Y - length)	PP non-woven	35 g/m²	X - width (cm) Y - length (cm)	green
BT-030-PF-25-"X" x "Y" - G ( X - width, Y - length)	laminated PP non-woven	25 g/m²		green
BT-030- PF-45-"X" x "Y" - G ( X - width, Y - length)	laminated PP non-woven	45 g/m²		green
BT-030- SMS-35-"X" x "Y" - B ( X - width, Y - length)	SMS non-woven	35 g/m²		blue

\*according to customer's specification



# DISPOSABLE MEDICAL PILLOW CASE

non-sterile

#### features:

• made of polypropylene non-woven

	REF	grammage	size	packaging unit
	BT-PP-25G1	25 g/m²	60cm x 80 cm	50 pcs
	BT-PP-25G6	25 g/m <sup>2</sup>	80cm x 90 cm	50 pcs
		et e		
b	oulk packaging: 10 x 50 pcs	available colours:		



# DISPOSABLE MEDICAL QUILT COVER

non-sterile

#### features:

made of polypropylene non-woven

REF	grammage	size	packaging unit
BT-PZ-25G1	25 g/m²	160cm x 210 cm	10 pcs
bulk packaging: 10 x 10 pcs	available colo	burs:	



### **DISPOSABLE MEDICAL BEDDING SET**

non-sterile

### features:

130

- made of polypropylene non-woven
- components:
- bed sheet: 150 cm x 210 cm
- quilt cover: 160 cm x 210 cm
- pillow case: 70 cm x 80 cm

REF	grammage	colour	packaging unit
BT-ZES-25G3	25 g/m <sup>2</sup>	green	1 set
BT-ZES-40G3	40 g/m <sup>2</sup>	green	1 set
BT-ZES-40W3	40 g/m <sup>2</sup>	white	1 set



# **DISPOSABLE ABSORBENT UNDERPAD**

non-sterile

bulk packaging: 25 x 1 set

#### features:

#### • composed of five layers:

laminate + cellulose wadding + cellulose pulp + cellulose wadding + PP nonwoven • absorbent insert made of cellulose pulp

available colours:

- impermeable foiled bottom
- single use

REF	dimensions	packaging unit	
BT-PC4060	40 cm x 60 cm	25 pcs	
BT-PC6060	60 cm x 60 cm	25 pcs	
BT-PC6090	60 cm x 90 cm	25 pcs	



4 x 25 pcs

# **DISPOSABLE ABSORBENT UNDERPAD**

for operating table, non-sterile

#### features:

- quilted absorbent layer with addition of SAP (superabsorbent)
- composed of five layers: polypropylene non-woven, cellulose wadding, cellulose pulp
- with superabsorbent SAP, cellulose wadding, PE foil
- impermeable foiled bottom
- single use

REF	dimensions	packaging unit	•••••
BT-PC100225	100 cm x 225 cm	10 pcs	







# DISPOSABLE MEDICAL BED SHEET, REINFORCED

non-sterile

### features:

- made of foil laminated paper
- reinforced with 8 or 48 threads
- grammage 40 g/m<sup>2</sup>
- soft, waterproof, absorbent, tough
- white with blue threads
- product put on mattresses, couches, armchairs and other elements that are the equipment of health care units, ie hospitals, outpatient clinics, clinics, etc.
- the products are designed to protect the patient during medical procedures, to prevent the transmission of contaminants and potentially pathogenic microorganisms located in medical rooms, treatment rooms and equipment that furnish these rooms and that may lead to complications of the patient's treatment process. They are a barrier to contamination
- non-sterile
- single use

REF	type	dimensions	packaging unit
BT-PW101	8 threads	80 cm x 210 cm	25 pcs
BT-PW401	48 threads	80 cm x 210 cm	25 pcs

# THERMAL PATIENT COVER SOFT

non-sterile

#### features:

- made of Spunlace non-woven with polyester filling
- ultrasonic seams over the entire surface
- can be heated up in the heating device
- maintaining patient temperature at a stable level
- non-sterile
- single use

REF	dimensions	outer material	packaging unit	bulk packaging
BT-TOP-21S	110 cm x 120 cm	Spunlace non-woven	1 pc.	50 x 1 pcs
BT-TOP-22S	110 cm x 210 cm	Spunlace non-woven	1 pc.	25 x 1 pcs

## THERMAL PATIENT COVER

non-sterile

#### features:

- made of polypropylene non-woven with polyester filling
- ultrasonic seams over the entire surface
- can be heated up in the heating device
- maintaining patient temperature at a stable level
- non-sterile
- single use

REF	dimensions	outer material	packaging unit	bulk packaging
BT-TOP-11P	110 cm x 120 cm	polypropylene non-woven	1 pc.	50 x 1 pcs
BT-TOP-12P	110 cm x 210 cm	polypropylene non-woven	1 pc.	25 x 1 pcs
BT-TOP-13P	150 cm x 210 cm	polypropylene non-woven	1 pc.	20 x 1 pcs











131













# PROTECTIVE COVERALL SAFE+ AGAINST CHEMICAL AND INFECTIVE AGENTS TYPE 4b - 5b - 6b

Personal protective equipments (PPE) are products that protect against one or more risks to the user's health or safety. These products must comply with Regulation (EU) 2016/425 of the European Parliament and of the Council. PPE is classified into three categories according to the degree of risk.

# The coverall is a PPE category III in accordance with the Regulation (EU) 2016/425 of the European Parliament and of the Council



The category III covers products of complex construction, intended to protect health or life against risks that may cause death or irreversible damage to health.

#### Standards:

requirements

and PB [4])



Coverall is personal protective equipment of category III in conformity with Regulation (EU) 2016/425 of the European Parliament and of the Council



EN 14605:2005+A1:2009 Protective clothing against liquid chemicals. Performance requirements for clothing with liquid-tight (Type 3) or spray-tight (Type 4) connections, including items providing protection to parts of the body only (Types PB [3]

EN ISO 13688:2013 Protective clothing. General



EN ISO 13982-1:2004+A1:2010 Protective clothing for use against solid particulates. Part 1: Performance requirements for chemical protective clothing providing protection to the full body against airborne solid particulates (type 5 clothing)



EN 13034:2005+A1:2009 Protective clothing against liquid chemicals. Performance requirements for chemical protective clothing offering limited protective performance against liquid chemicals (Type 6 and Type PB [6] equipment)



EN 1073-2:2002 Protective clothing against radioactive contamination. Part 2: Requirements and test methods for non-ventilated protective clothing against particulate radioactive contamination



EN 14126:2003 + AC:2004 Protective clothing. Performance requirements and tests methods for protective clothing against infective agents

class 6 Resistance to penetration of contaminated liquids under the influence of hydrostatic pressure class 6 Resistance to penetration by infective agents due to mechanical contact with substances class 3 Resistance to penetration of contaminated liquid aerosols class 3 Resistance to penetration of contaminated



solid particles


EN 1149-5:2018 Protective clothing. Electrostatic properties. Part 5: Material performance and design requirements

EN 14325:2004 Protective clothing against chemicals. Test methods and performance classification of chemical protective clothing materials, seams, joins and assemblages

#### features:

- made of a laminate (non-woven polypropylene and polyethylene) with a weight of 63  $\rm g/m^2$
- polyester seams additionally covered with tape
- equipped with a hood with elastic trim
- equipped with an elastic waist band, wrists and ankles
- zipper covered with a self-adhesive flap
- available in version with boots





REF	size	version	packaging unit
BS-KCB02-S	S	without boots	1 pc.
BS-KCB02-M	М	without boots	1 pc.
BS-KCB02-L	L	without boots	1 pc.
BS-KCB02-XL	XL	without boots	1 pc.
BS-KCB02-XXL	XXL	without boots	1 pc.
BS-KCB02-XXXL	XXXL	without boots	1 pc.
BS-KCB02-B-S	S	with boots	1 pc.
BS-KCB02-B-M	М	with boots	1 pc.
BS-KCB02-B-L	L	with boots	1 pc.
BS-KCB02-B-XL	XL	with boots	1 pc.
BS-KCB02-B-XXL	XXL	with boots	1 pc.
BS-KCB02-B-XXXL	XXXL	with boots	1 pc.

133











# Surgical sets, drapes and surgical gowns

Preventing infections in the operating field is a priority for all medicals. Thanks to our complex offer of single-use protective solutions of ALPHAtex STANDARD, STANDARD+, COMFORT, the post-procedural complications and risk of infection in the operating field are minimized while increasing the protection barrier against the potentially pathogenic germs residing in the human skin layers. We offer:

- patient and equipment drapes
- sets of surgical drapes
- incise films
- surgical covers
- adhesive tapes
- accessories such as fluid collection pouches
- surgical gowns

Properly selected features of surgical drapes provide safety and comfort of use to medical staff and patients. A variety of multi-layer non-woven and laminates reinforced in critical areas in lines STANDARD+ and COMFORT allow to adjust our solutions to particular surgical procedures with maximal cost efficiency.

ALPHAtex products are characterized by parameters frequently exceeding the requirements of EN 13795 for 'Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment'.

Moreover, they provide a successful protection against the risk of infection in the operating field. Practical solutions and composition of our sets enable their quick and simple application while performing medical procedures.



# Prevention of surgical site infection



### Surgical site infection (SSI)

is defined as infection occurring up to 30 days after surgery, when there is a possibility of microorganisms penetration into the open surgical wound. Microorganisms could have various sources EXOGENOUS (carried by an external factor as medical staff or medical devices) or ENDOGENOUS (patient's bacterial flora from the closest area surrounding the surgical postoperative wound). The most important routes of infections are airway or contact path. Various type of barriers are used to minimalize the surgical site infections: surgical drapes, surgical gowns and medical equipment covers.



### **Surgical drapes**

(the names are based on the EN 13795 standard)

Surgical drapes are used to ensure clean microbiological surgical wound surrounding, minimalize the microorganisms migration from patient skin to surgical wound and to supervise the spreading of body fluids from the wound area. Surgical drapes meet the requirements of medical device in accordance with the Directive 93/42/EEC.

### Surgical gowns (the names are based on the EN 13795 standard)

Surgical gowns are used to prevent direct microorganisms transmission from medical staff to patient. Moreover, surgical gowns prevent spreading of exfoliated epidermis in the air of the operating room. The wide range of ALPHAtex products it's an effect of over 30 years of ZARYS International Group existence on medical market. Sharing our experience and knowledge we are convinced that using ALPHAtex products will provide comfort, safety and economic effectiveness during daily work in the operating room.




### **ZONES OF SURGICAL DRAPES ON PATIENT** on the operating room

**ZONE 4 ZONE 2** NOT HAVING ANY CONTACT **DIRECT SURROUNDINGS** WITH PATIENT'S BODY OF OPERATING FIELD performs the function in this zone surgical drape has of anaesthetic cover got contact with not disinfected and operating table cover patient's body **ZONE1** ZONE 3 PARTIAL CONTACT WITH **OPERATING FIELD** NOT OPERATED PARTS disinfected patient's skin in direct surroundings of surgical wound, OF PATIENT'S BODY mostly protected by incise film occasionally exposed to body fluids

# Surgical sets ALPHAtex

contains different disposable medical devices



**ALPHAtex label** 



### Features of disposable surgical drapes ALPHAtex used in the operating room









# **SURGICAL DRAPES**

### BASIC TYPES OF MATERIALS USED IN SURGICAL DRAPES:

PF DRAPE - two layer material of the drape containing absorbency layer of celulose on top and PE foil layer of the drape on back

NF, NFA, NFF, NFFA, NFFAR DRAPE - two layer material of the drape containing absorbency layer of non-woven material on top and PE foil layer on back of the drape


NFN, NFNA, NFNF, NFNFA DRAPE - three layer material on the drape containing absorbency non-woven layer on top, PE foil layer in the middle and non-woven absorbency layer on the back of the drape






### **STERILE DRAPES WITH CELLULOSE LAYER - 2-LAYERED**

ref.	dimensions	packaging unit	bulk packaging
AT-PF-S 1	50cm x 50cm	1 pc.	500 pcs
AT-PF-S 3	50cm x 70cm	1 pc.	400 pcs
AT-PF-S 5	80cm x 80cm	1 pc.	200 pcs
AT-PF-S 6	75cm x 90cm	1 pc.	200 pcs
AT-PF-S 8	90cm x 150cm	1 pc.	100 pcs




### **STERILE SURGICAL DRAPES - 2-LAYERED**

ref.	dimensions	packaging unit	bulk packaging
AT-NF-S1E	45cm x 45cm	100 pcs	4 x 100 pcs
AT-NF-S3E	45cm x 60cm	100 pcs	4 x 100 pcs
AT-NF-S4E	45cm x 75cm	100 pcs	4 x 100 pcs
AT-NF-S 5 E	75cm x 75cm	50 pcs	4 x 50 pcs
AT-NF-S9E	75cm x 90cm	50 pcs	4 x 50 pcs
AT-NF-S 13 E	100cm x 150cm	50 pcs	4 x 50 pcs
AT-NF-S 15 E	150cm x 180cm	20 pcs	4 x 20 pcs
AT-NF-S 21 E	37cm x 45cm	100 pcs	4 x 100 pcs




### **STERILE SURGICAL DRAPES - 2-LAYERED**

ref.	dimensions	packaging unit	bulk packaging
AT-NF-S 1	45cm x 45cm	1 pc.	2 x 240 pcs
AT-NF-S 2	50cm x 50cm	1 pc.	2 x 200 pcs
AT-NF-S 3	50cm x 60cm	1 pc.	2 x 280 pcs
AT-NF-S4	50cm x 70cm	1 pc.	2 x 160 pcs
AT-NF-S 5	75cm x 75cm	1 pc.	2 x 150 pcs
AT-NF-S8	80cm x 90cm	1 pc.	2 x 150 pcs
AT-NF-S 9	75cm x 90cm	1 pc.	2 x 150 pcs
AT-NF-S 10	90cm x 120cm	1 pc.	2 x 80 pcs
AT-NF-S 13	100cm x 150cm	1 pc.	2 x 50 pcs
AT-NF-S 14	120cm x 150cm	1 pc.	2 x 50 pcs
AT-NF-S 15	150cm x 180cm	1 pc.	2 x 20 pcs
AT-NF-S 16	150cm x 200cm	1 pc.	2 x 20 pcs
AT-NF-S 17	180cm x 180cm	1 pc.	2 x 20 pcs
AT-NF-S 19	160cm x 220cm	1 pc.	2 x 35 pcs
AT-NF-S 20	170cm x 300cm	1 pc.	2 x 20 pcs



### STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 2-LAYERED

ref.	dimensions	packaging unit	bulk packaging
AT-NFA-S1E	45cm x 45cm	100 pcs	4 x 100 pcs
AT-NFA-S4E	45cm x 75cm	100 pcs	4 x 100 pcs
AT-NFA-S9E	75cm x 90cm	50 pcs	4 x 50 pcs
AT-NFA-S 10 E	90cm x 120cm	50 pcs	4 x 50 pcs
AT-NFA-S 13 E	100cm x 150cm	50 pcs	4 x 50 pcs
AT-NFA-S 15 E	150cm x 175cm	20 pcs	4 x 20 pcs
AT-NFA-S 16 E	150cm x 200cm	20 pcs	4 x 20 pcs
AT-NFA-S 20 E	150cm x 240cm	20 pcs	4 x 20 pcs






economic version

141


### STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 2-LAYERED

ref.	dimensions	packaging unit	bulk packaging
AT-NFA-S 1	45cm x 45cm	1 pc.	2 x 230 pcs
AT-NFA-S 2	50cm x 50cm	1 pc.	2 x 200 pcs
AT-NFA-S 3	50cm x 60cm	1 pc.	2 x 100 pcs
AT-NFA-S 4	50cm x 70cm	1 pc.	2 x 210 pcs
AT-NFA-S 5	75cm x 75cm	1 pc.	2 x 135 pcs
AT-NFA-S 9	75cm x 90cm	1 pc.	2 x 125 pcs
AT-NFA-S 10	90cm x 120cm	1 pc.	2 x 80 pcs
AT-NFA-S 13	100cm x 150cm	1 pc.	2 x 50 pcs
AT-NFA-S 15	150cm x 180cm	1 pc.	2 x 20 pcs
AT-NFA-S 16	150cm x 200cm	1 pc.	2 x 20 pcs
AT-NFA-S 20	150cm x 240cm	1 pc.	2 x 20 pcs
AT-NFA-S 22	170cm x 300cm	1 pc.	2 x 20 pcs



### STERILE SURGICAL DRAPES WITH CENTRAL FENESTRATION - 2-LAYERED

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFF-S 2 E	45cm x 75cm	6cm x 8cm	100 pcs	4 x 100 pcs
AT-NFF-S3E	45cm x 60cm	ø 6cm	100 pcs	4 x 100 pcs
AT-NFF-S4E	75cm x 90cm	ø7cm	50 pcs	4 x 50 pcs
AT-NFF-S 5 E	75cm x 90cm	6cm x 8cm	50 pcs	4 x 50 pcs




### STERILE SURGICAL DRAPES WITH CENTRAL FENESTRATION - 2-LAYERED

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFF-S 1	45cm x 45cm	ø 6cm	1 pc.	2 x 200 pcs
AT-NFF-S 2	50cm x 75cm	6cm x 8cm	1 pc.	2 x 120 pcs
AT-NFF-S 3	50cm x 75cm	ø7cm	1 pc.	2 x 120 pcs
AT-NFF-S4	75cm x 90cm	ø7cm	1 pc.	2 x 120 pcs
AT-NFF-S 5	75cm x 90cm	6cm x 8cm	1 pc.	2 x 120 pcs
AT-NFF-S 6	90cm x 120cm	ø 10cm	1 pc.	2 x 120 pcs
AT-NFF-S 21	50cm x 60cm	ø 7cm	1 pc.	2 x 300 pcs




### STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL FENESTRATION - 2-LAYERED

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFFA-S1E	45cm x 45cm	ø 5cm	100 pcs	4 x 100 pcs
AT-NFFA-S 2 E	45cm x 75cm	6cm x 8cm	100 pcs	4 x 100 pcs
AT-NFFA-S 3 E	50cm x 60cm	ø 6cm	100 pcs	4 x 100 pcs
AT-NFFA-S4E	75cm x 90cm	ø 7cm	50 pcs	4 x 50 pcs
AT-NFFA-S 5 E	75cm x 90cm	6cm x 8cm	50 pcs	4 x 50 pcs
AT-NFFA-S 6 E	90cm x 120cm	ø8cm	50 pcs	4 x 50 pcs



### STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL FENESTRATION - 2-LAYERED

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFFA-S 1	45cm x 45cm	ø 6cm	1 pc.	2 x 200 pcs
AT-NFFA-S 2	50cm x 75cm	6cm x 8cm	1 pc.	2 x 120 pcs
AT-NFFA-S 3	50cm x 75cm	ø7cm	1 pc.	2 x 120 pcs
AT-NFFA-S 4	75cm x 90cm	ø7cm	1 pc.	2 x 100 pcs
AT-NFFA-S 5	75cm x 90cm	6cm x 8cm	1 pc.	2 x 100 pcs
AT-NFFA-S 6	90cm x 120cm	ø 8cm	1 pc.	2 x 50 pcs
AT-NFFA-S 7	120cm x 150cm	ø7cm	1 pc.	2 x 50 pcs
AT-NFFA-S 8	120cm x 150cm	ø 10cm	1 pc.	2 x 50 pcs
AT-NFFA-S 11	50cm x 60cm	6cm x 8cm	1 pc.	2 x 200 pcs





### STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL **ADJUSTABLE FENESTRATION - 2-LAYERED**

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFFAR-S 2 E	45cm x 75cm	6cm x 8cm	100 pcs	4 x 100 pcs
AT-NFFAR-S4E	75cm x 90cm	ø7cm	50 pcs	4 x 50 pcs





version

4ja

60 g/m²


### STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL **ADJUSTABLE FENESTRATION - 2-LAYERED**


ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFFAR-S 2	45cm x 75cm	6cm x 8cm	1 pc.	2 x 150 pcs
AT-NFFAR-S 4	75cm x 90cm	ø7cm	1 pc.	2 x 90 pcs

### **STERILE ADHESIVE EPIDURAL DRAPES - 2-LAYERED**

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFZRS-S 1	75cm x 80cm	6cm x 15cm	1 pc.	2x60 pcs
AT-NFZRS-S 11	75cm x 90cm	10cm x 15cm	1 pc.	2x60 pcs

### indication: epidural anesthesia






adhesive tape in the corners of drape

### STERILE INSTRUMENT TABLE COVERS


ref.	dimensions	packaging unit	bulk packaging
AT-SD-S 3	100cm x 150cm	1 pc.	100 pcs
AT-SD-S 4	150cm x 190cm	1 pc.	60 pcs
AT-SD-S 5	150cm x 260cm	1 pc.	45 pcs





#### 144 **STERILE SURGICAL DRAPES - 3-LAYERED**

ref.	dimensions	packaging unit	bulk packaging
*AT-NFN-S1	45cm x 45cm	1 pc.	300 pcs
AT-NFN-S 2	50cm x 50cm	1 pc.	300 pcs
AT-NFN-S 3	50cm x 60cm	1 pc.	150 pcs
AT-NFN-S 4	50cm x 70cm	1 pc.	160 pcs
AT-NFN-S 5	75cm x 75cm	1 pc.	100 pcs
AT-NFN-S 9	75cm x 90cm	1 pc.	80 pcs
*AT-NFN-S 10	90cm x 120cm	1 pc.	75 pcs
AT-NFN-S 13	100cm x 150cm	1 pc.	50 pcs
AT-NFN-S 15	150cm x 180cm	1 pc.	50 pcs
AT-NFN-S 16	150cm x 210cm	1 pc.	25 pcs
* AT-NFN-S 20	150cm x 240cm	1 pc.	25 pcs






4<u>1</u>4 75 g/m²

### **STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 3-LAYERED**

ref.	dimensions	packaging unit	bulk packaging
AT-NFNA-S 1	45cm x 45cm	1 pc.	300 pcs
AT-NFNA-S 2	50cm x 50cm	1 pc.	256 pcs
AT-NFNA-S 3	50cm x 60cm	1 pc.	200 pcs
AT-NFNA-S 4	50cm x 70cm	1 pc.	160 pcs
AT-NFNA-S 5	75cm x 75cm	1 pc.	100 pcs
AT-NFNA-S 9	75cm x 90cm	1 pc.	90 pcs
AT-NFNA-S 10	90cm x 120cm	1 pc.	80 pcs
AT-NFNA-S 13	100cm x 150cm	1 pc.	80 pcs
AT-NFNA-S 15	150cm x 180cm	1 pc.	25 pcs
AT-NFNA-S 16	150cm x 210cm	1 pc.	50 pcs
AT-NFNA-S 20	150cm x 240cm	1 pc.	25 pcs



adhesive tape

4<sup>°</sup>

on the longer side of drape

75 g/m²

### STERILE SURGICAL DRAPES WITH CENTRAL FENESTRATION - 3-LAYERED

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFNF-S 1	45cm x 45cm	ø 6cm	1 pc.	300 pcs
AT-NFNF-S 2	50cm x 75cm	6cm x 8cm	1 pc.	160 pcs
AT-NFNF-S 3	50cm x 75cm	ø 7cm	1 pc.	200 pcs
AT-NFNF-S4	75cm x 90cm	ø 7cm	1 pc.	80 pcs
AT-NFNF-S 5	75cm x 90cm	6cm x 8cm	1 pc.	80 pcs
AT-NFNF-S 6	90cm x 120cm	ø 10cm	1 pc.	100 pcs





#### STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL FENESTRATION - 3-LAYERED

ref.	dimensions	fenestration	packaging unit	bulk packaging
*AT-NFNFA-S 1	45cm x 45cm	ø 6cm	1 pc.	300 pcs
AT-NFNFA-S 2	50cm x 70cm	6cm x 8cm	1 pc.	200 pcs
AT-NFNFA-S 3	50cm x 75cm	ø7cm	1 pc.	100 pcs
AT-NFNFA-S 4	75cm x 90cm	ø7cm	1 pc.	100 pcs
AT-NFNFA-S 5	75cm x 90cm	6cm x 8cm	1 pc.	100 pcs
AT-NFNFA-S 6	90cm x 120cm	ø 10cm	1 pc.	63 pcs



لم 75 g/m² 0


#### STERILE SURGICAL DRAPES WITH ADHESIVE FENESTRATION, REINFORCED

ref.	dimensions	fenestration	packaging unit	bulk packaging
AT-NFZRF-S 1	120cm x 170cm	ø 12cm	1 pc.	50 pcs
AT-NFZRF-S 2	200cm x 280cm	ø 10cm	1 pc.	20 pcs
AT-NFZRF-S 3	200cm x 320cm	ø 15cm	1 pc.	20 pcs



hydrophobic SMS material








#### **GENERAL SURGERY**

AT-SD-ABD1-IF-C	ABDOMIN	AL DRAPE
ref.	dimensions	quantity

f.	dimensions	quantity	component
F-SD-ABD1-IF-C	200cm x 320cm	1 pc.	drape with adhesive fenestration 25cm x 30cm filled with an incise film, reinforced
r-sd-s4	150cm x 190cm	1 pc.	instrument table cover, reinforced

1-IF






AT-SD-ABPER1-S ABDO-PERINEAL DRAPE 1

ref.	dimensions	quantity	component
AT-SD-ABPER1-S	180/280cm x 220cm	1 pc.	drape integrated with limb covers, adhesive trapezium-shaped fenestration 30cm x 20cm x 20cm and perineal fenestration 10cm x 15cm with covering flap
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced






150x190cm

180/280x220cm

#### AT-SD-LAP3-SP LAPAROSCOPY DRAPE 3

ref.	dimensions	quantity	component
AT-SD-LAP3-SP	200/270cm x 300cm	1 pc.	drape integrated with limb covers, adhesive fenestration 25cm x 30cm, reinforced, with 8 cable holders Velcro type
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced






150x190cm

200/270x300cm



AT-SD-ORT1-SP	ORTHOPAED	ORTHOPAEDIC DRAPE 1					
ref.	dimensions	quantity	component				
AT-SD-ORT1-SP	150cm x 250cm	1 pc.	drape with adhesive U-split 7cm x 70cm, reinforced				
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced				



bulk packaging: 20 pcs

AT-SD-ORT2-SP

bulk packaging: 16 pcs

AT-SD-ORT3-C

AT-SD-ORT3-C

AT-SD-S4


 $\bigotimes$ bulk packaging: 20 pcs

ref.

Ŵ



Ŵ



150x190cm

230x260cm

ref.	dimensions	quantity	component
AT-SD-ORT2-SP	230cm x 260cm	1 pc.	drape with adhesive U-split 20cm x 100cm, reinforced
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced

**ORTHOPAEDIC DRAPE 2** 

**ORTHOPAEDIC DRAPE 3** 

quantity

1 pc.

1 pc.

component

drape with adhesive U-split 6cm x 65cm, reinforced


instrument table cover, reinforced

dimensions

150cm x 200cm

150cm x 190cm







150x190cm

150x200cm

AT-SD-ORT4-C **ORTHOPAEDIC DRAPE 4** ref. dimensions quantity component








147


AT-SD-EXT1-SP	EXTREMITY I	DRAPE 1		V		<b>.</b>
ref.	dimensions	quantity	component	Ĩ		
AT-SD-EXT1-SP	150cm x 180cm	1 pc.	drape with elastic fenestration ø 3,5cm, reinforced			
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced			
					<b>150</b> x <b>190</b> cm	<b>150</b> x <b>180</b> cm
bulk packaging: 24 pcs						
AT-SD-EXT2-SP	EXTREMITY	DRAPE 2		¥.		
ref.	dimensions	quantity	component			
AT-SD-EXT2-SP	230cm x 300cm	1 pc.	drape with elastic fenestration ø 3,5cm, reinforced			
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced			
bulk packaging:						
18 pcs						
					<b>150</b> x <b>190</b> cm	230x300cm
AT-SD-EXT3-SP	EXTREMITY I					
				۲.		
ref. AT-SD-EXT3-SP	dimensions 230cm x 300cm	quantity 1 pc.	component drape with elastic fenestration ø 7cm, reinfo	read		
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover. reinforced	in Ceu		
	1500117 170011	1 pc.				
bulk packaging: 18 pcs						
ia pos						
					150x190cm	<b>230</b> x <b>300</b> cm
				- Mb		
AT-SD-EXT4-C	EXTREMITY I	DRAPE 4				



**150**x**190**cm

200x320cm

**148** 



#### AT-SD-GYN2-S **GYNAECOLOGY DRAPE 2**


160/240cm x 180cm

150cm x 190cm

dimensions

<b>K</b>	
----------	--





150x190cm

160/240x180cm

8 pcs

AT-SD-DEL1-SP

AT-SD-S 4

Ø


bulk packaging:

ref.

AT-SD-GYN2-S

**DELIVERY DRAPE 1** 





quantity

1 pc.

1 pc.

component

drape integrated with limb covers, adhesive perineal fenestration 10cm x 15cm, integrated fluid collection pouch 35cm x 50cm with formable molder, filter and tubing port

instrument table cover, reinforced





Ø bulk packaging: 45 pcs







#### AT-SD-OPHT1-IF-S



#### **OPHTHALMIC DRAPE 1-IF**








۲

#### Ø bulk packaging:

50 sets

AT-SD-OPHT2-IF-C


#### **OPHTHALMIC DRAPE 2-IF**



#### AT-SD-OPHT3-IF-C

#### **OPHTHALMIC DRAPE 3-IF**

ref.	dimensions	quantity	component
AT-SD-OPHT3-IF-C	120cm x 170cm	1 pc.	ophthalmic drape with fenestration 10cm x 10cm filled with an incise film, integrated with 2 fluid collection pouches and with 1 formable molder, reinforced
AT-SD-S 3	100cm x 150cm	1 pc.	instrument table cover, reinforced



bulk packaging: 40 sets



ĥ

hydrophobic SMS material




formable molder to fit the patient's face

#### AT-SD-OPHT4-IF-C

#### **OPHTHALMIC DRAPE 4-IF**

ref.	dimensions	quantity	component
AT-SD-OPHT4-IF-C	140cm x 240cm	1 pc.	renforced ophthalmic drape with a fenestration of 10cm x 10cm filled with an incise film, integrated with 2 fluid collection pouches and 1 formable molder
AT-SD-S 4	150cm x 190cm	1 pc.	instrument table cover, reinforced



150x190cm

140x240cm



Ø

bulk packaging:

25 sets




# SURGICAL COVERS AND ACCESSORIES



W W W . Z A R Y S . P L



#### STERILE MAYO STAND COVER, REINFORCED - RED







80x140cm

#### STERILE SET OF LIMB COVERS

152






7x120cm

STERILE LIMB COVERS

ref.	dimensions	packaging unit	bulk packaging
AT-SD-S 6	30cm x 60cm	1 pc. + OP tape 10cm x 50cm	200 pcs
AT-SD-S 7	30cm x 80cm	1 pc. + OP tape 10cm x 50cm	120 pcs
AT-SD-S 8	35cm x 120cm	1 pc. + OP tape 10cm x 50cm	100 pcs



#### STERILE CAMERA CABLE COVERS - PERFORATED, WITH 2 ADHESIVE TAPES

ref.	dimensions	packaging unit	bulk packaging
AT-SP-S 3	14cm x 250cm	1 pc.	160 pcs
AT-SP-S 4	16cm x 200cm	1 pc.	200 pcs
AT-SP-S 5	16cm x 250cm	1 pc.	160 pcs
AT-SP-S 6	16cm x 300cm	1 pc.	200 pcs






sterile camera cover

#### STERILE CAMERA CABLE COVER - UNPERFORATED, WITHOUT ADHESIVE TAPE


ref.	dimensions	packaging unit
AT-SP-S 7	7cm x 120cm	1 pc.
bulk packaging: 200 pcs		

#### STERILE CAMERA CABLE COVER - UNPERFORATED, WITH AN ADHESIVE TAPE




#### STERILE CIRCULAR BANDED COVERS FOR MEDICAL DEVICES

ref.	dimensions in stretched state	maximum size of covered device	packaging unit	bulk packaging
AT-SP-S 9	65cm x 85cm	35cm x 45cm	1 pc.	200 pcs
AT-SP-S 10	85cm x 85cm	50cm x 50cm	1 pc.	100 pcs
AT-SP-S 11	105cm x 105cm	60cm x 60cm	1 pc.	100 pcs
AT-SP-S 12	120cm x 120cm	70cm x 70cm	1 pc.	100 pcs
AT-SP-S 13	140cm x 140cm	80cm x 80cm	1 pc.	50 pcs



#### STERILE RECTANGULAR BANDED COVERS FOR MEDICAL DEVICES

ref.	dimensions	packaging unit	bulk packaging
AT-SP-S 14	70cm x 80cm	1 pc.	100 pcs
AT-SP-S 15	90cm x 100cm	1 pc.	50 pcs
AT-SP-S 16	80cm x 150cm	1 pc.	50 pcs






# STERILE C-ARM COVER SET










#### STERILE SURGICAL TAPES





OP tape

154

Velcro surgical tape

#### **STERILE ADHESIVE POUCHES**

ref.	dimensions	description	packaging unit	bulk packaging	
AT-SP-S 20	30cm x 40cm	adhesive pouch, one-chamber	1 pc.	180 pcs	
AT-SP-S 21	50cm x 60cm	adhesive pouch, one-chamber	1 pc.	180 pcs	
AT-SP-S 22	40cm x 30cm	adhesive pouch, two-chamber 15cm x 30cm left chamber 25cm x 30cm right chamber		250 pcs	

#### STERILE ADHESIVE FLUID COLLECTION POUCH



#### STERILE ABSORBENT DRAPE





500 pcs

#### STERILE ABSORBENT DRAPE





1//

#### 155

#### STERILE INCISE FILM elastoFILM

ref.	total dimensions	adhesive filed dimensions	packaging unit	bulk packaging
808067	56cm x 84cm	56cm x <b>78cm</b>	1 pc.	6 pcs
808066	40cm x 42cm	40cm x <b>36cm</b>	1 pc.	20 pcs
808065	45cm x 55cm	45cm x <b>49cm</b>	1 pc.	20 pcs
808064	45cm x 28cm	45cm x <b>22cm</b>	1 pc.	20 pcs
808063	30cm x 28cm	30cm x <b>22cm</b>	1 pc.	20 pcs
808062	15cm x 28cm	15cm x <b>22cm</b>	1 pc.	40 pcs
808061	14cm x 25cm	14cm x <b>21cm</b>	1 pc.	40 pcs
805018	55cm x 55cm	55cm x <b>51cm</b>	1 pc.	20 pcs
805025	30cm x 40cm	30cm x <b>36cm</b>	1 pc.	20 pcs
805016	20cm x 30cm	20cm x <b>26cm</b>	1 pc.	30 pcs





- thin incise film made of polyurethane
- flexible
- perfectly fits the patient's body
- double packed: externally paper-foil, internally paper



#### CELLULOSE

Cellulose is a solid with typically fiber structure. Two-step process is used for industrial manufacturing of cellulose fibers, which was developed in 1892. Cellulose raw material is used in different industries as textiles and paper as well as in production of explosives.



# SURGICAL SETS



W W W . Z A R Y S . C O M



#### AT-S-BASIC1-S


#### **BASIC SET 1**

quantity	dimensions	component
1 pc.	75cm x 90cm	drape with central adhesive fenestration
1 pc.	10cm x 17cm	tray
10 pcs	10cm x 10cm	gauze swab
5 pcs	30cm x 30cm	gauze ball
1 pc.	100cm x 150cm	instrument table cover, reinforced

# $\bigotimes$


bulk packaging:

28 sets



#### AT-S-UNI1-S **UNIVERSAL SET 1**

quantity	dimensions	component
4 pcs	75cm x 90cm	drape with adhesive edge
1 pc.	150cm x 190cm	instrument table cover, reinforced
bulk packag 50 sets	ing:	



AT-S-UNI2-S

**UNIVERSAL SET 2** 

quantity	dimensions	component					
2 pcs	150cm x 180cm	drape with adhesive edge					
2 pcs	75cm x 90cm	drape with adhesive edge					
4 pcs	30cm x 30cm	absorbent towel					
1 pc.	150cm x 190cm	instrument table cover, reinforced					
bulk packagir 20 sets	ıg:		2x	<b>150x180</b> cm	2x <b>75x90</b> cm	4x 30x30cm	<b>150</b> x <b>190</b> cm

#### AT-S-UNI4-C **UNIVERSAL SET 4**

quantity	dimensions	component
1 pc.	170cm x 240cm	drape with adhesive edge, reinforced
1 pc.	170cm x 180cm	drape with adhesive edge, reinforced
2 pcs	75cm x 90cm	highly absorbent drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced



MA

**2**x


**75**x**90**cm


**4**x

30x30cm

Г


**10**x**50**cm





170x240cm

170x180cm



#### AT-S-UNI5-S

**UNIVERSAL SET 5** 


quantity	dimensions	component
1 pc.	150cm x 240cm	drape with adhesive edge
1 pc.	180cm x 180cm	drape with adhesive edge
2 pcs	75cm x 90cm	drape with adhesive edge
4 pcs	30 cmx30 cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced



150x240cm





180x180cm




#### AT-S-UNI6-SP **UNIVERSAL SET 6**

quantity	dimensions	component
1 pc.	150cm x 240cm	drape with adhesive edge, reinforced
1 pc.	180cm x 180cm	drape with adhesive edge, reinforced
2 pcs	75cm x 90cm	drape with adhesive edge, reinforced
4 pcs	30,5cm x 34cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced





80x145cm



150x190cm

#### AT-S-UNI7-C **UNIVERSAL SET 7**

ij


highly absorbent

material

quantity	dimensions	component
1 pc.	150cm x 240cm	drape with adhesive edge
1 pc.	180cm x 180cm	drape with adhesive edge
2 pcs	75cm x 90cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced

**2**x

**75**x**90**cm



**4**x

30,5x34cm

**10**x**50**cm

158

Ø

bulk packaging:


9 sets

#### AT-S-UNI8-C **UNIVERSAL SET 8**

quantity	dimensions	component				E	
1 pc.	150cm x 240cm	drape with adhesive edg	e, reinforced				
1 pc.	180cm x 180cm	drape with adhesive edg	e, reinforced				L L
2 pcs	75cm x 90cm	drape with adhesive edg	e, reinforced				
4 pcs	30cm x 30cm	absorbent towel				E	
1 pc.	10cm x 50cm	OP tape				•	L L
1 pc.	80cm x 145cm	Mayo stand cover, reinfo	orced				
1 pc.	150cm x 190cm	instrument table cover,	reinforced				
					150×24	<b>0</b> cm	<b>180</b> × <b>180</b> cm
		¢.				_	
bulk packagir 7 sets	ng: 3-layer dra PREMIU		nt				
			2x	4x			
			<b>75</b> x <b>90</b> cm	<b>30</b> x <b>30</b> cm	<b>10</b> x <b>50</b> cm	<b>80</b> x <b>145</b> cm	<b>150</b> x <b>190</b> cm

# AT-S-UNI9-SP PEDIATRIC UNIVERSAL SET 9

quantity	/ dimensions	component
1 pc.	150cm x 240cm	drape with adhesive edge
1 pc.	180cm x 180cm	drape with adhesive edge
2 pcs	75cm x 75cm	drape with adhesive edge, reinforced
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced






٢ 150x240cm 180x180cm bulk packaging: 10 sets pediatric set **2**x **75**x**75**cm **30**x**30**cm **10**x**50**cm 80x145cm 150x190cm thinner tape in drapes - 2,5cm

### AT-S-UNI10-S UNIVERSAL SET 10

quantity	dimensions	component
1 pc.	150cm x 230cm	drape with adhesive edge
1 pc.	150cm x 180cm	drape with adhesive edge
2 pcs	75cm x 90cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced








150x230cm






80x145cm


150x190cm

159

#### AT-S-UNI5-G-S UNIVERSAL SET 5-G


quantity	dimensions	component
1 pc.	150cm x 240cm	drape with adhesive edge
1 pc.	180cm x 180cm	drape with adhesive edge
2 pcs	75cm x 90cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced
1 pc.	L	surgical gown STANDARD
2 pcs	XL	surgical gown STANDARD






150x240cm

180x180cm




150x190cm




# AT-S-UNI5-GP-S UNIVERSAL SET 5-GP

quantity	dimensions	component
1 pc.	150cm x 240cm	drape with adhesive edge
1 pc.	180cm x 180cm	drape with adhesive edge
2 pcs	75cm x 90cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced
1 pc.	XL	surgical gown STANDARD
1 pc.	XXL	surgical gown, reinforced STANDARD PLUS
1 pc.	30cm x 40cm	adhesive pouch, one-chamber



150x240cm

180x180cm



**160** 

bulk packaging: 5 sets


# AT-S-ABD1-IF-C ABDOMINAL SET 1-IF

 $\widehat{}$ 

quantity	dimensions	component
1 pc.	200cm x 320cm	drape with adhesive fenestration 25cm x 30cm filled with an incise film, reinforced
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced


Ø bulk packaging: 10 sets





10x50cm


161





# AT-S-ABPER1-S ABDO-PERINEAL SET 1

quantity	dimensions	component
1 pc.	180/270cm x 220cm	drape integrated with limb covers, adhesive trapezium-shaped fenestration 30cm x 20cm x 20cm and perineal fenestration 10cm x 15cm with covering flap
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	2cm x 22cm	surgical tape, Velcro type
1 pc.	14cm x 250cm	camera cable cover
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced



180/270x220cm

**30**x**30**cm

 $\bigotimes$ 

bulk packaging: 10 sets





**4**x

30x30cm

10:550cm








80x140cm

150x190cm

flap to secure the fenestration

#### AT-S-LAP1-IF-C LAPAROSCOPY SET 1-IF




#### AT-S-LAP2-SP LAPAROSCOPY SET 2

quantity	dimensions	component
1 pc.	200/270cm x 300cm	drape integrated with limb covers, adhesive fenestration 25cm x 30cm, reinforced, with 8 cable holders Velcro type
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	2cm x 22cm	surgical tape, Velcro type
1 pc.	14cm x 250cm	camera cable cover
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced




highly absorbent material



cable holders, Velcro type








200/270x300cm


30x30cm



AT-S-LAP3-C


#### LAPAROSCOPY SET 3

quantity	dimensions	component
1 pc.	200/270cm x 300cm	drape integrated with limb covers, adhesive fenestration 25cm x 30cm, reinforced, with 8 cable holders Velcro type
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	2cm x 22cm	surgical tape, Velcro type
1 pc.	14cm x 250cm	camera cable cover
1 pc.	30cm x 40cm	adhesive pouch, one-chamber
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced



163

200/270x300cm



÷ highly absorbent

material

bulk packaging: 10 sets

 $\bigotimes$ 

bulk packaging: 8 sets

AT-S-VAR1-S VARICOSE VEIN SET 1

 $\widehat{}$ 


hydrophobic SMS material

quantity	dimensions	component
1 pc.	200cm x 260cm	drape with adhesive U-split 8,5cm x 85cm
1 pc.	160cm x 180cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

**10**x**50**cm

2x22cm

14x250cm




**4**x

**30**x**30**cm


200x260cm

**4**x




150x190cm


160x180cm



200x260cm







30x30cm

**4**x



80x140cm

150x190cm

IJ

drape with

U-split

#### **VARICOSE VEIN SET 2**

quantity	dimensions	component
1 pc.	200cm x 260cm	drape with adhesive U-split 8,5cm x 95cm, integrated genital area cover, reinforced
1 pc.	170cm x 180cm	drape with adhesive edge, reinforced
4 pcs	30cm x 30cm	absorbent towel
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

integrated genital area cover

 $\bigotimes$ bulk packaging: 10 sets

> U drape with

U-split

AT-S-VAR2-C



÷ĵ highly absorbent material

 $\widehat{}$ 



164


#### AT-S-EXT1-C EXTREMITY SET 1

quantity	dimensions	component	
1 pc.	200cm x 300cm	drape with elastic fenestration ø 3,5cm, reinforced	
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red	
1 pc.	150cm x 190cm	instrument table cover, reinforced	


	ç Î
bulk packaging:	highly absorbent
12 sets	material



<i>.</i>
limb
surgery



200x300cm



# AT-S-EXT2-C EXTREMITY SET 2


quantity	dimensions	component
1 pc.	200cm x 300cm	drape with elastic fenestration ø 3,5cm, reinforced
1 pc.	160cm x 180cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced



ċĵ highly absorbent material









**4**×

30x30cm

160x180cm




200x300cm



80x140cm

10x50cm

150x190cm



#### AT-S-ART2-C ARTHROSCOPY SET 2

 $\widehat{}$ 


quantity	dimensions	component
1 pc.	200cm x 320cm	drape with elastic fenestration 6cm x 8cm, reinforced
2 pcs	30cm x 30cm	absorbent towel
2 pcs	10cm x 50cm	OP tape
1 pc.	14cm x 250cm	camera cable cover
1 pc.	30cm x 60cm	elastic limb cover
1 pc.	60cm x 100cm	fluid collection pouch with elastic fenestration, formable molder and tubing port
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

 $\bigotimes$ bulk packaging: 14 sets






200x320cm



#### AT-S-ART2-IP-C ARTHROSCOPY SET 2-IP



10x50cm

14x250cm

AT-S-ART3-C ARTHROSCOPY SET 3

 $\widehat{}$ 

hydrophobic SMS material

quantity	dimensions	component
1 pc.	200cm x 320cm	drape with elastic fenestration 6cm x 8cm, reinforced
2 pcs	30cm x 30cm	absorbent towel
3 pcs	10cm x 50cm	OP tape
1 pc.	14cm x 250cm	camera cable cover
1 pc.	30cm x 60cm	elastic limb cover
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

÷

highly absorbent

material


**30**x**30**cm

200x320cm

80x140cm

30x60cm

150x190cm



bulk packaging: 12 sets AT-S-ART4-S ARTHROSCOPY SET 4


quantity	dimensions	component
1 pc.	230cm x 320cm	drape with elastic fenestration ø 7cm, integrated fluid collection pouch
4 pcs	30cm x 30cm	absorbent towel
3 pcs	10cm x 50cm	OP tape
1 pc.	22cm x 75cm	elastic limb cover
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
2 pcs	150cm x 190cm	instrument table cover, reinforced

 $\bigotimes$ bulk packaging: 7 sets





230x320cm




#### AT-S-ART5-SP ARTHROSCOPY SET 5

quantity	dimensions	component
1 pc.	230cm x 300cm	drape with elastic fenestration ø 7cm, reinforced
4 pcs	30cm x 30cm	absorbent towel
3 pcs	10cm x 50cm	OP tape
1 pc.	22cm x 75cm	elastic limb cover
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
2 pcs	150cm x 190cm	instrument table cover, reinforced




ŝ highly absorbent material





230x300cm



AT-S-HIP1-C HIP SET 1

**168** 

quantity	dimensions	component	
1 pc.	200cm x 260cm	drape with adhesive U-split 8,5cm x 85cm, reinforced	
1 pc.	170cm x 240cm	drape with adhesive edge, reinforced	
4 pcs	30cm x 30cm	absorbent towel	
2 pcs	10cm x 50cm	OP tape	
1 pc.	35cm x 120cm	elastic limb cover	
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red	
1 pc.	150cm x 190cm	instrument table cover, reinforced	
ulk packagir 8 sets	ng: hydrophobi SMS materia		<b>200</b> x <b>260</b> cm
		E	

30x30cm

U

**2**x [

10x50cm

AT-S-HIP2-C HIP SET 2

 $\widehat{}$ 

hydrophobic

SMS material

component


170x240cm

quantity	dimensions	component
1 pc.	200cm x 260cm	drape with adhesive U-split 8,5cm x 85cm, reinforced
1 pc.	170cm x 300cm	drape with adhesive edge, reinforced
1 pc.	160cm x 180cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
2 pcs	10cm x 50cm	OP tape
1 pc.	35cm x 120cm	elastic limb cover
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

÷


highly absorbent material

160x180cm



35x120cm

80x140cm



150x190cm

drape with U-split **4**x **2**x [ **30**x**30**cm **10**x**50**cm 35x120cm 150x190cm 80x140cm

www.zarys.com

Ø

bulk packaging:

8 sets


#### AT-S-HIP3-SP HIP SET 3

quantity	dimensions	component	
1 pc.	230cm x 260cm	drape with adhesive U-split 20cm x 100cm, reinforced	
1 pc.	150cm x 240cm	drape with adhesive edge, reinforced	
1 pc.	180cm x 180cm	drape with adhesive edge, reinforced	
2 pcs	75cm x 90cm	drape with adhesive edge, reinforced	
4 pcs	30cm x 30cm	absorbent towel	
3 pcs	10cm x 50cm	OP tape	
1 pc.	30cm x 120cm	elastic limb cover	
1 pc.	10cm x 450cm	elastic bandage	
1 pc.	80cm x 145cm	Mayo stand cover, reinforced	
1 pc.	150cm x 190cm	instrument table cover, reinforced	

U

drape with

. U-split



230x260cm

150x240cm



÷ĵ

highly absorbent


material

elastic bandage

Ø

bulk packaging:

5 sets



AT-S-SHLR1-C SHOULD

**SHOULDER SET 1** 




bulk packaging: 8 sets hydrophobic SMS material material

150×180cm



30x30cm

170x180cm



200x260cm

**10**x**50**cm



169


150x190cm

AT-S-SHLR2-S SHOULDER SET 2

Ø bulk packaging:

8 sets

quantity	dimensions	component
1 pc.	263cm x 406cm	drape with elastic fenestration 13cm x 14cm, integrated fluid collection pouch
4 pcs	30cm x 30cm	absorbent towel
2 pcs	10cm x 50cm	OP tape
1 pc.	22cm x 75cm	elastic limb cover
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced



#### **TRUE OR FALSE?**


Ħ

chair position

"I am too young / too old for hip replacement."


#### NOTHING COULD BE MORE WRONG!

Hip endoprosthesis is performed in patients from a wide age range. Very important is patient condition, however this is only one factor which classifies patient for procedure.



22x75cm





150x190cm






#### AT-S-GYN1-S GYNAECOLOGY SET 1

quantity	/ dimensions	component
1 pc.	160/240cm x 180cm	drape integrated with limb covers, adhesive perineal fenestration 10cm x 15cm, integrated fluid collection pouch 35cm x 50cm with formable molder, filter and tubing port
1 pc.	75cm x 90cm	drape without adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	150cm x 190cm	instrument table cover, reinforced




bulk packaging: 6 sets




160/240x180cm

171









150x190cm

fluid collection pouch with filter and tubing port

#### AT-S-GYN2-C **GYNAECOLOGY SET 2**

quantity	dimensions	component
1 pc.	160/240cm x 180cm	drape integrated with limb covers, adhesive perineal fenestration 10cm x 15cm, integrated fluid collection pouch 35cm x 50cm with formable molder, filter and tubing port, reinforced
1 pc.	75cm x 90cm	drape without adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced





150x190cm

fluid collection pouch with filter and tubing port

#### AT-S-GYN3-C **GYNAECOLOGY SET 3**

75x90cm



30x30cm


10x50cm



#### **CYSTOSCOPY SET 1**

100×150cm

quantity	dimensions	component
1 pc.	75cm x 90cm	drape with adhesive fenestration 6cm x 8cm
2 pcs	70cm x 120cm	limb cover
1 pc.	100cm x 150cm	instrument table cover, reinforced



150x190cm

80x140cm

70x120cm



hydrophobic SMS material


ċĵ

highly absorbent material

#### AT-S-UROG2-S URO/GYNAECOLOGY SET 2

quantity	dimensions	component	
1 pc.	160cm x 180cm drape with adhesive edge		
1 pc.	75cm x 100cm	drape with adhesive fenestration 6cm x 14cm	
1 pc.	10cm x 50cm	OP tape	
2 pcs	70cm x 120cm	limb cover	
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red	
1 pc.	150cm x 190cm	instrument table cover, reinforced	






70x120cm

80x140cm

150x190cm

# AT-S-UROG3-S URO/GYNAECOLOGY SET 3

quantity	dimensions component	
1 pc.	75cm x 200cm	drape with adhesive diamond-shaped fenestration 8cm x 12cm
4 pcs	30cm x 30cm	absorbent towel
2 pcs	70cm x 120cm limb cover	
1 pc.	100cm x 150cm	instrument table cover, reinforced




bulk packaging: 20 sets

 $\bigotimes$ 

AT-S-TUR1-S TUR SET 1

quantity	dimensions	component
1 pc.	190cm x 230cm	drape integrated with limb covers, adhesive abdominal fenestration ø 6cm and adhesive perineal fenestration ø 6cm, finger coat, fluid collection pouch with formable molder, filter and tubing port
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	150cm x 190cm	instrument table cover, reinforced










 $\widehat{}$ 

hydrophobic SMS material




made of two different materials

AT-S-TUR1-S-P1


**TUR SET 1** 

quantity	dimensions	component
1 pc.	175cm x 240cm	drape with integrated leggings, abdominal fenestration, perineal fenestration with finger coat and fluid collection pouch. Drape made of hydrophobic 2 layer non-woven material PP and PE foil
2 pcs	30cm x 30cm	absorbent towel
1 pc.	2cm x 22cm	surgical tape, Velcro type
1 pc.	150cm x 190cm	instrument table cover, reinforced






175x240cm



173

# AT-S-PCNL1-IP-S PERCUTANEOUS LITHOTRIPSY SET 1-IP

quantity	dim	ensions	component
1 pc.	180	cm x 300cm	drape with adhesive fenestration 21,5cm x 24cm, integrated with long fluid collection pouch
1 pc.	150	cm x 190cm	instrument table cover, reinforced
bulk packag 12 sets	jing:	မြုံ့ kidney stone removal set	





#### AT-S-CSEC1-C C-SECTION SET 1


quantity	dimensions	component
1 pc.	200cm x 320cm	drape with adhesive fenestration 25cm x 30cm, integrated with fluid collection pouch 75cm x 85cm with tubing ports
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	75cm x 80cm	absorbent baby towel
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

bulk packaging: 10 sets




 $\widehat{}$ 

hydrophobic SMS material











**75**x**80**cm

174


# AT-S-CSEC1-IF-C C-SECTION SET 1-IF

 $\widehat{}$ 

		component
1 pc.	200cm x 320cm	drape with adhesive fenestration 25cm x 30cm filled with an incise film, integrated with fluid collection pouch 75cm x 85cm with tubing ports
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	75cm x 80cm	absorbent baby towel
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced

bulk packaging: 10 sets





200x320cm

**4**x **10**x**50**cm **75%800**cm 80x140cm 150x190cm


AT-S-CSEC2-IF-C C-SECTION SET 2-IF

quantity	dimensions	component		
1 pc.	200cm x 320cm	drape with adhesive fenestration 25cm x 30cm filled with an incise film, integrated with fluid collection pouch 75cm x 85cm with tubing ports		
4 pcs	30cm x 30cm	absorbent towel		
1 pc.	10cm x 50cm	OP tape		
1 pc.	2cm x 22cm	surgical tape, Velcro type		
1 pc.	75cm x 80cm	absorbent baby towel	111	
10 pcs	40cm x 40cm	gauze lap sponge with X-ray thread and belt		
20 pcs	30cm x 30cm	gauze ball		
1 pc.	L	surgical gown STANDARD		
3 pcs	XL	surgical gown STANDARD		
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red		
1 pc.	150cm x 190cm	instrument table cover, reinforced		<b>4</b> x
~			<b>200</b> x <b>320</b> cm	<b>30</b> x <b>30</b> cm



#### AT-S-CSEC3-IF-SP C-SECTION SET 3-IF

quantity	dimensions	component
1 pc.	180/240cm x 215cm	drape integrated with limb covers, adhesive fenestration 23,5cm x 23,5cm x 22,5cm x 33cm filled with an incise film, integrated with fluid collection pouch 56cm x 80cm with tubing port, reinforced
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	75cm x 80cm	absorbent baby towel
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced






176

ċ highly absorbent material





#### AT-S-DEL1-S


 $\bigotimes$ 

bulk packaging:

20 sets

#### **DELIVERY SET 1**

quantity	dimensions	component
1 pc.	90cm x 120cm	drape without adhesive edge
2 pcs	30cm x 30cm	absorbent towel
3 pcs	75cm x 80cm	absorbent baby towel
1 pc.	60cm x 90cm	absorbent underpad
1 pc.	80cm x 130cm	foil apron
2 pcs	-	umbilical cord clamp




80x130cm

AT-S-DEL2-S **DELIVERY SET 2** 

#### quantity dimensions component 90cm x 120cm drape without adhesive edge 6 pcs 30cm x 30cm absorbent towel 75cm x 80cm absorbent baby towel 2 pcs 60cm x 90cm absorbent underpad



bulk packaging: 20 sets



2x

**75**x**80**cm

**60**x**90**cm


# AT-S-DEL3-S DELIVERY SET 3

quantity	dimensions	component
1 pc.	90cm x 120cm	drape without adhesive edge
2 pcs	30cm x 30cm	absorbent towel
3 pcs	75cm x 80cm	absorbent baby towel
1 pc.	60cm x 90cm	absorbent underpad
1 pc.	80cm x 130cm	foil apron
2 pcs	-	umbilical cord clamp
10 pcs	30cm x 30cm	gauze ball

# 90x120cm 3x **75**x**80**cm

 $\bigotimes$ 


bulk packaging: 15 sets



177

# AT-S-DEL4-S DELIVERY SET 4

quantity	dimensions	component
2 pcs	90cm x 120cm	drape without adhesive edge
6 pcs	30cm x 30cm	absorbent towel
1 pc.	75cm x 80cm	absorbent baby towel
1 pc.	60cm x 90cm	absorbent underpad
1 pc.	80cm x 130cm	foil apron
2 pcs	-	umbilical cord clamp







#### AT-S-DEL5-SP DELIVERY SET 5


quantity	dimensions	component
1 pc.	114cm x 150cm	delivery drape integrated with fluid collection pouch 78cm x 80cm with formable molder, filter and tubing port, reinforced
2 pcs	30cm x 30cm	absorbent towel
1 pc.	75cm x 80cm	absorbent baby towel
1 pc.	150cm x 190cm	instrument table cover, reinforced


 $\bigotimes$ 

bulk packaging: 18 sets



pouch with a very large scale





75x85cm


150x190cm



178

#### AT-S-OPHT2-IF-C **OPHTHALMIC SET 2-IF**


quantity	dimensions	component		
1 pc.	120cm x 170cm	ophthalmic drape with fenestration 10cm x 10cm filled with an incise film, integrated with 2 fluid collection pouches and with 1 formable molder, reinforced		
2 pcs	30cm x 30cm	absorbent towel	-	2×
1 pc.	40cm x 65cm	Phaco tray cover	-	
2 pcs	35cm x 75cm	armboard cover	- <b>120</b> x <b>170</b> cm	<b>30</b> x <b>30</b> cm
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red	-	
1 pc.	150cm x 190cm	instrument table cover, reinforced		



# AT-S-OPHT3-IF-C OPHTHALMIC SET 3-IF

hydrophobic SMS material

quantity	dimensions	component			
1 pc.	140cm x 240cm	ophthalmic drape with fenestration 10cm x 10cm filled with an incise film, integrated with 2 fluid collection pouches and with 1 formable molder, reinforced			
4 pcs	30cm x 30cm	absorbent towel			
1 pc.	40cm x 65cm	Phaco tray cover			
2 pcs	35cm x 75cm	armboard cover			
1 pc.	L	surgical gown STANDARD			
1 pc.	XL	surgical gown STANDARD		4x	
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red	140x240cm	<b>30</b> x <b>30</b> cm	<b>40</b> x <b>65</b> cm
1 pc.	150cm x 190cm	instrument table cover, reinforced			



bulk packaging: 12 sets



#### AT-S-OTOL1-C OTOLARYNGOLOGY SET 1

hydrophobic

SMS material

 $\frac{1}{2}$ 

hydrophobic

SMS material

41 head surgery

quantity	dimensions	component
1 pc.	200cm x 280cm	drape with adhesive fenestration ø 10cm, reinforced
2 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	150cm x 190cm	instrument table cover, reinforced

2x

30x30cm

10x50cm

**4**x

30x30cm

10x


10x50cm


10x50cm

**20**x

80x145cm

**75**x**90**cm





÷ highly absorbent

bulk packaging:

20 sets

material

Ø

bulk packaging:

12 sets

ŝ

highly absorbent material


Ø


bulk packaging: 12 sets

÷ highly absorpent material

#### AT-S-OTOL2-C OTOLARYNGOLOGY SET 2


quantity	dimensions	component
1 pc.	150cm x 200cm	drape with adhesive U-split 6cm x 65cm, reinforced
1 pc.	150cm x 180cm	drape with adhesive edge
1 pc.	75cm x 90cm	drape with adhesive edge
4 pcs	30cm x 30cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	80cm x 140cm	Mayo stand cover, reinforced - red
1 pc.	150cm x 190cm	instrument table cover, reinforced





150x200cm







150x190cm

#### AT-S-CRA1-IF-SP CRANIOTOMY SET 1-IF

quantity	dimensions	component
1 pc.	230cm x 300cm	drape with adhesive fenestration 20cm x 30cm filled with anincise film, with fluid collection pouch, reinforced
4 pcs	50cm x 50cm	drape with adhesive edge
2 pcs	30cm x 34cm	absorbent towel
1 pc.	10cm x 50cm	OP tape
1 pc.	2cm x 25cm	surgical tape, velcro type
1 pc.	80cm x 145cm	Mayo stand cover, reinforced
1 pc.	150cm x 190cm	instrument table cover, reinforced


2


\_

2x25cm

2x

30x34cm





179

# CARDIOSURGERY AND ANGIOGRAPHY


180

# AT-S-ABL1-IF-SP ABLATION SET 1-IF

quantity	dimensions	component
1 pc.	180cm x 305cm	drape with 2 adhesive fenestrations 12cm x 15cm filled with an incise film, reinforced
2 pcs	30cm x 30cm	absorbent towel
1 pc.	ø 80cm	circular banded cover for medical equipment
1 pc.	XL	surgical gown STANDARD
1 pc.	150cm x 190cm	instrument table cover, reinforced




ċĵ highly absorbent material





180x305cm

**30**x**30**cm



#### AT-S-ANG2-C ANGIOGRAPHY SET 2


quantity	dimensions	component
1 pc.	220cm x 320cm	drape with 2 adhesive fenestrations ø 10cm, integrated with transparent foil inserts protecting the control panel reinforced
1 pc.	10cm x 50cm	OP tape
1 pc.	XL	surgical gown STANDARD
1 pc.	100cm x 150cm	drape without adhesive edge
1 pc.	100cm x 150cm	instrument table cover, reinforced



 $\widehat{}$ hydrophobic SMS material



highly absorbent material





# SURGICAL GOWNS



W W W . Z A R Y S . P L

#### **STERILE GOWN NORMAL**





#### STERILE SURGICAL GOWN STANDARD

ref.	size	trimming	packaging unit
AT-SGS-M 1	М		1 pc.
AT-SGS-L 1	L		1 pc.
AT-SGS-XL 1	XL		1 pc.
AT-SGS-XXL 1	XXL		1 pc.
AT-SGS-XXXL 1	XXXL		1 pc.

bulk packaging: 25 pcs



BESTSELLER



#### STERILE SURGICAL GOWN STANDARD

ref.	size	trimming	packaging unit
AT-SGS-M 1E	М		1 pc.
AT-SGS-L 1E	L		1 pc.
AT-SGS-XL 1E	XL		1 pc.
AT-SGS-XXL 1E	XXL		1 pc.

bulk packaging: 25 pcs 2





economic

version

#### STERILE SURGICAL GOWN STANDARD

ref.	size	trimming	packaging unit
AT-SGS-M 1A	М		1 pc.
AT-SGS-L 1A	L		1 pc.
AT-SGS-XL 1A	XL		1 pc.
AT-SGS-XXL 1A	XXL		1 pc.






**30**x**40**cm

30×30cm

#### STERILE SURGICAL GOWN STANDARD

ref.	size	trimming	packaging unit	
AT-SGS-M	М		1 pc.	
AT-SGS-L	L		1 pc.	
AT-SGS-XL	XL		1 pc.	
AT-SGS-XXL	XXL		1 pc.	
bulk packaging: 25 pcs	includ 2 absorben 30cm x 3	t towels	ک <u>آ</u> 45 g/㎡	



#### STERILE SURGICAL GOWN SOFT





#### STERILE SURGICAL GOWN STANDARD PLUS





#### STERILE SURGICAL GOWN STANDARD PLUS

ref.	size	trimming	packaging unit	
AT-SGSP-M 1E	М		1 pc.	
AT-SGSP-L 1E	L		1 pc.	
AT-SGSP-XL 1E	XL		1 pc.	
AT-SGSP-XXL 1E	XXL		1 pc.	
	Ġ	]		4 <sup>1</sup>
bulk packaging: 25 pcs	witho 2 absorben 30cm x 3	t towels	economic version	35 g/m²



#### 183

#### STERILE SURGICAL GOWN STANDARD PLUS

Я

includes

2 absorbent towels

30cm x 40cm

ref.	size	trimming	packaging unit
AT-SGSP-M 1A	М		1 pc.
AT-SGSP-L 1A	L		1 pc.
AT-SGSP-XL 1A	XL		1 pc.
AT-SGSP-XXL 1A	XXL		1 pc.

without

colorful trimming

4<u>1</u>4

35 g/m²



#### STERILE SURGICAL GOWN STANDARD PLUS

ref.	size	trimming	packaging unit	
AT-SGSP-M	М		1 pc.	
AT-SGSP-L	L		1 pc.	
AT-SGSP-XL	XL		1 pc.	
AT-SGSP-XXL	XXL		1 pc.	
AI-SGSP-XAL	AAL	4	I þc.	
		ปี 👘	5 <sup>1</sup>	

45 g/m²



184

bulk packaging:

50 pcs



#### STERILE SURGICAL GOWN EXTRA SAFE





#### STERILE SURGICAL GOWN COMFORT

AT-SGC-XL XL 1 pc. AT-SGC-XXL 1 pc. bulk packaging: includes 68 g/m <sup>2</sup>	AT-SGC-M	М	1 pc.
AT-SGC-XXL 1 pc.	AT-SGC-L		1 pc.
bulk packaging: includes 68 g/m <sup>2</sup>	AT-SGC-XL	XL	1 pc.
	AT-SGC-XXL	XXL	1 pc.
	bulk packaging: 25 pcs	includes 2 absorbent towels	 68 g/㎡



www.zarys.com

#### STERILE SURGICAL GOWN COMFORT PLUS

AT-SGCP-M	M	1 pc.
AT-SGCP-L	L L	1 pc.
AT-SGCP-XL	XL	1 pc.
AT-SGCP-XXL	XXL	1 pc.





**30**x**30**cm

### NON-STERILE SURGICAL GOWNS

#### NON-STERILE SURGICAL GOWN STANDARD

ref.	size	trimming	packaging unit
AT-NS-SGS-M 1E	М		1 pc.
AT-NS-SGS-L 1E	L		1 pc.
AT-NS-SGS-XL 1E	XL		1 pc.
AT-NS-SGS-XXL 1E	XXL		1 pc.








#### NON-STERILE SURGICAL GOWN PREMIUM

ref.	size	trimming	packaging unit	
AT-NS-SGP-L 1E	L		10 pcs	
AT-NS-SGP-XL 1E	XL		10 pcs	
bulk packaging: 10 x 10 pcs		inforced H+PE		
meet the standards of EN 14126	35	ſ∆ g/m²		



### **Material legends:**

detailed information and material structures available in **material templet ALPHAtex** 



in the critical area

# STANDARDS GUIDE

### SURGICAL GOWNS

				REQUIF	REMENT	
CHARACTERISTICS	STANDARD	UNIT	STANDARD R	EQUIREMENT	HIGH REQ	UIREMENT
			CRITICAL PRODUCT AREA	LESS CRITICAL PRODUCT AREA	CRITICAL PRODUCT AREA	LESS CRITICAL PRODUCT AREA
microbial penetration - dry	EN ISO 22612	CFU	not required	≤ 300 ª	not required	≤ 300 ª
microbial penetration - wet	EN ISO 22610	I <sub>B</sub>	≥2,8 <sup>b</sup>	not required	6,0 <sup>bc</sup>	not required
cleanliness - microbial bioburden	EN ISO 11737-1	CFU/100 cm <sup>2</sup>	≤ 300	≤ 300	≤ 300	≤ 300
particulate release	EN ISO 9073-10	IPM	≤ 4,0	≤ 4,0	≤ 4,0	≤ 4,0
linting	EN ISO 9073-10	log₁₀ (lint count)	≥40	≥40	≥40	≥40
liquid penetration	EN ISO 811	cm H <sub>2</sub> O	≥30	≥10	≥100	≥10
bursting strength - dry	EN ISO 13938-1	kPa	≥40	≥40	≥40	≥40
bursting streangth - wet	EN ISO 13938-1	kPa	≥40	not required	≥40	not required
tensile strength - dry	EN 29073-3	Ν	≥20	≥20	≥20	≥20
tensile strength - wet	EN 29073-3	Ν	≥20	not required	≥20	not required

### SURGICAL DRAPES

			REQUIREMENT					
CHARACTERISTICS	STANDARD	UNIT	STANDARD R	EQUIREMENT	HIGH REQUIREMENT			
			CRITICAL PRODUCT AREA	LESS CRITICAL PRODUCT AREA	CRITICAL PRODUCT AREA	LESS CRITICAL PRODUCT AREA		
microbial penetration - dry	EN ISO 22612	CFU	not required	≤ 300 ª	not required	≤ 300 ª		
microbial penetration - wet	EN ISO 22610	I <sub>B</sub>	≥2,8 <sup>b</sup>	not required	6,0 <sup>bc</sup>	not required		
cleanliness - microbial bioburden	EN ISO 11737-1	CFU/100 cm <sup>2</sup>	≤ 300	≤ 300	≤ 300	≤ 300		
particulate release	EN ISO 9073-10	IPM	≤ 4,0	≤ 4,0	≤ 4,0	≤ 4,0		
linting	EN ISO 9073-10	log₁₀ (lint count)	≥40	≥40	≥40	≥40		
liquid penetration	EN ISO 811	cm H <sub>2</sub> O	≥30	≥10	≥100	≥10		
bursting strength - dry	EN ISO 13938-1	kPa	≥40	≥40	≥40	≥40		
bursting streangth - wet	EN ISO 13938-1	kPa	≥40	not required	≥40	not required		
tensile strength - dry	EN 29073-3	Ν	≥15	≥15 ≥15		≥20		
tensile strength - wet	EN 29073-3	Ν	≥15	≥15 not required ≥		not required		

### Standard EN 13795-1:2019

Essential requirements for surgical drapes and surgical gowns results from directive 93/42/EEC and have been included in EN 13795 standard "Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment". The purpose of EN 13795 standard is to maintain stable level of surgical drapes and surgical gowns safety for whole using period as well as minimalization of direct microorganisms transmission from medical staff to patient during surgical procedures.



187











# Cohesive bandages, kinesiotapes and others...

**yellow**SPORT is a brand offering comprehensive solutions in the field of sports medicine - according to the best practices that have been developed during physical activity tests. Based on the knowledge and experience of experts, the product line is dedicated not only to sportswomen and sportsmen but also to people taking their first steps to improve their physical activity (recreation).

**yellow**SPORT is primarily a tribute to, among others, physiotherapists, personal trainers, medical staff, sports clubs and owners of rehabilitation clinics. It creates an excellent base and product base while maintaining the expected ratio between cost and required quality. The offer of products from the **yellow**SPORT family was prepared by specialists in the field of medical devices and physiotherapy. Therefore, from the very first day of its presence on the market, the "yellow brand" is a concrete response to the real needs of users, regardless of the level of advancement in the field of physical activity or just recreation.

don't be shy, be yellow



### yellow**BAND**

#### **COHESIVE BANDAGE**

#### USE:

190

- fixation of all types of non-adhesive dressings, especially where the elastic band slides
- fixation of splinting rails
- in compression therapy
- reducing joints mobility after an injury

#### FEATURES:

- 1-layer, breathable, made of non-woven
- cohesive it coheres to itself, does not adhere to skin, hair or other materials
- latex-free

#### **AVAILABLE SIZES AND COLOURS:**



# O, O, O, O, O, O, O, O, O, Ø, Ø, Ø, Ø, O,

2,5cm x 4,5m				YCB-2545-X		YCB-2545-B						_		
5cm x 4,5m	YCB-545-S	YCB-545-Y	YCB-545-G	YCB-545-X	YCB-545-R	YCB-545-B	YCB-545-W	YCB-545-BC	YCB-545-IO	YCB-545-IG	YCB-545-GP	YCB-545-VP	YCB-545-XP	YCB-545-PH
7,5cm x 4,5m	YCB-7545-S	YCB-7545-Y	YCB-7545-G	YCB-7545-X	YCB-7545-R	YCB-7545-B	YCB-7545-W	YCB-7545-BC	YCB-7545-IO	YCB-7545-IG	YCB-7545-GP	YCB-7545-VP	YCB-7545-XP	YCB-7545-PH
10cm x 4,5m	YCB-1045-S	YCB-1045-Y	YCB-1045-G	YCB-1045-X	YCB-1045-R	YCB-1045-B	YCB-1045-W	YCB-1045-BC	YCB-1045-IO	YCB-1045-IG	YCB-1045-GP	YCB-1045-VP	YCB-1045-XP	YCB-1045-PH
15cm x 4,5m				YCB-1545-X										

# yellowTAPE

#### KINESIOTAPE

#### USE:

- kinesiotaping (dynamic taping)
- ensures nutrition of tissues and accelerates their recovery
- / improves the biomechanics of movement without restraining joint's mobility
- reduces pain by affecting sensory receptors

#### **FEATURES:**

- 🥖 elastic, breathable tape made of cotton
- ø wave-like acrylic coating

#### AVAILABLE SIZES AND COLOURS:





YT-55-4



don't be shy, be yellow CHECK THE PRODUCT



# yellow**RIGID TAPE**

### SPORTS TAPE, LINE GLUE COATING

#### USE:

- / athletic (rigid) taping
- sports injuries prevention
- protecting the joints from excessive mobility after the injury, which allows for activity continuation
- protection of the areas exposed to abrasions

#### FEATURES:

- non-elastic, breathable tape made of cotton
- acrylic coating
- serrated edges no need for scissors
- alternate arrangement of the stripes covered with glue and those without it form grooves that allow free sweat and water flow that reduces the tape's sensivity and increases its resitance to enviromental conditions

### AVAILABLE SIZES AND COLOURS:



# yellowSPORTS TAPE

#### SPORTS TAPE


#### USE:

- 🥖 athletic (rigid) taping
- sports injuries prevention
- protecting the joints from excessive mobility after the injury, which allows for activity continuation
- protection of the areas exposed to abrasions

#### FEATURES:

- non-elastic, breathable tape made of cotton
- acrylic coating
- serrated edges no need for scissors

#### AVAILABLE SIZES AND COLOURS:



yellow WRAP

### FOAM UNDERWRAP

#### USE:

- skin protection during athletic taping
- ensuring correct base for athletic taping in hairy areas
- protecting from abrasions caused by orthoses or sport equipment

#### FEATURES:

- made of polyurethane foam
- porous structure allows for air flow and moisture drainage

#### AVAILABLE SIZES AND COLOURS:









# yellow**FOAM**

#### FOAM TAPE

#### USE:

192

- stabilization, protection and fixation of non-adhesive primary dressings put directly on wounds or injured skin
- ${\star}$  protection of the areas exposed to abrasions when the epidermis
  - is undamaged e.g. while wearing new shoes or another equipment

#### FEATURES:

- / waterproof tape made of polyethylene foam
- acrylic coating

#### **AVAILABLE SIZES AND COLOURS:**





### scissors

#### FOR SPORTS TAPES

#### USE:

for cutting sports tapes or other adhesive tapes

#### **FEATURES:**

- / sharp blades ensure precision of the cut
- covered with teflon preventing the tapes from adhesing to the blade
- ergonomic shape fits perfectly in hand ensuring comfort of using
- Iength of the blade: 21 cm

#### **REF number: YC-X**



# yellow**BOX**

#### DISPENSER FOR KINESIOTAPES

#### USE:

- storing and transporting of tapes
- protects from contamination and damage
- "rolling" core makes it easy to cut the desired length of tape

#### FEATURES:

- made of hard PVC
- suitable for tapes of 5 cm width

#### **REF number: YB**







let us introduce you to the new sport line





ANY QUESTIONS?

**Contact us!** 

ESPORT@ZARYS.PL
www.yellowsport.pl/en/







### **PROCEDURE SETS**

Operating theatre is a specific environment where time is the most valuable factor influencing the success of the operation. It is a place where every minute - or even a second - counts, because it can affect the patient's life. To meet the expectations of medicine, ZARYS introduces operating sets to its offer that meet the highest standards, helping to improve the preparation of the operating theatre before the procedure

### ADVANTAGES of using omegapack operating sets



### **REDUCTION OF ERRORS**

the complete set contains all disposable medical devices necessary for the surgical operation, which eliminates the possibility of errors during the selection of individual products



### WIDE RANGE OF PRODUCTS

a wide range of medical devices selected individually in accordance with the Customer's preferences The sets are assembled on the basis of individual requirements and preferences of Customers



Medical devices under the omegapack brand by ZARYS are operating sets containing complete equipment with all disposable products necessary to perform a given medical procedure to which the set is dedicated. The operating set contains carefully selected disposable components, such as:

- surgical drapes, pads
- accessories, e.g. fluid collection bags, tapes
- surgical gowns
- surgical gloves
- covers for equipment
- shoe covers, masks, caps

- gauze balls
- gauze and non-woven swabs bandages, dressings
- catheters, urine collection bags
- trays, bowls
- taps, three-way stopcocks
- manifolds, connectors
- infusion sets
- suction sets
- disposable surgical instruments
- injection needles and syringes
- disposable electrosurpicel electrodes
   and accessories

and many more..





vitahealth line of adhesive plasters sets are designed to provide optimum protection in every situation and to satisfy the needs of the most demanding users. They are an indispensable element of everyday life.

### Always with you and your loved ones.



### The vitahealth family includes:

- plasters for cutting with scissors
- sets of plasters suited to your needs
- adhesives

- gauze compresses and bandages
- medical masks
- thermometers



# since 1989...

# We create good solutions

W W W . Z A R Y S . C O M

# ZARYS International Group sp. z o.o. sp.k. Please contact us if you require further details:

S	Area	Speaking Language	Direct Line	Mobile
	Baltics, CIS, Central Europe, Finland,	English	+48 32 376 08 00	+48 518 130 314
	Denmark, Spain			
	UK, Benelux, France, Sweden, Norway,	English, French	+48 32 797 05 71	+48 532 085 454
	North America, Australia, Portugal,			
	Central & South America, East Asia			
	Germany, Balcans, South Europe	English, Spanish	+48 32 376 07 02	+48 531 743 734
	Middle East, Africa , Italy	English	+48 32 376 07 06	+48 533 331 695

Pod Borem 18 41-808 Zabrze, Poland

(⊠)

(@) e-mail: export@zarys.pl



### **EC Certificate** Directive 93/42/EEC Annex V **Production Quality Assurance Medical Devices**

Registration No.: DD 60139535 0001

**Report No.:** 

26300232 017

Manufacturer: **ZARYS** International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

**Products:** 

(see attachments for products and sites included) Replaces EC Certificate, Registration No.: DD 60117020 0001

### **Expiry Date:**

2024-05-27

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date:

Date:

2019-06-09

2019-05-27

**Notified Body** TÜVRheini Tifizierung

Rafal Byczkowski

TUV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



Doc. 1/6, Rev. 0

# TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

- Sterile and non-sterile cutting gauze
- Non-sterile dressing gauze
- Sterile and non-sterile gauze swabs (with or without X-ray thread)
- Sterile and non-sterile gauze lap sponges (with X-ray thread/ with X-ray chip)
- Sterile and non-sterile gauze balls (with or without X-ray thread)
- Sterile and non-sterile gauze rolls (with or without X-ray thread)
- Sterile and non-sterile non-woven swabs (with or without X-ray thread)
- Sterile paraffin gauze dressings
- Sterile three-way stopcocks
- Sterile transfusion sets for single use
- Sterile infusion sets for single use
- Sterile extension tubes for infusion pump

**Notified Body** 



Rafal Byczkowski

Date: 2019-05-27

0/020 h 04.08 @ TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approva



Doc. 2/6, Rev. 0

# TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

- Sterile endotracheal tubes
- Sterile tracheostomy tubes
- Sterile breathing circuits
- Sterile catheter mounts
- Non-sterile anaesthetic masks
- Sterile laryngeal masks
- Sterile oxygen masks
- Sterile Multi-Vent masks
- Sterile non-rebreath masks
- Sterile nebulizer masks
- Sterile nasal oxygen cannulas
- Sterile nebulizer sets
- Sterile oxygen tubing
- Sterile suction catheters
- Sterile abdominal drains
- Sterile feeding tubes
- Sterile stomach and duodenal tubes
- Sterile urology catheters
- Sterile surgical suction sets

TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior app

- Sterile surgical suction cannulas
- Sterile syringes for single use

Notified Body



Rafal Byczkowski

Date: 2019-05-27



Doc. 3/6, Rev. 0

# TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

- Sterile insulin syringes
- Sterile tuberculin syringes
- Sterile hypodermic needles
- Sterile insulin pen needles
- Sterile blood lancets
- Sterile IV cannulas
- Sterile needle free valves
- Sterile surgical gloves

TÜVRheinlan **Notified Body** Tinzierungs

Date: 2019-05-27

10/020 h 04.08 (8) TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approval

Rafal Byczkowski



Doc. 4/6, Rev. 0

# **TÜV Rheinland** LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate **Registration No.: Report No.:** 

DD 60139535 0001 26300232 017

Manufacturer:

**ZARYS** International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

For the following medical devices the scope covers only the aspects of manufacture concerned with securing and maintaining sterile conditions:

- Adhesive cannulla fixation dressings
- Adhesive wound dressings
- Eye pads
- Incise films
- Transparent film dressings
- Foam dressings
- Absorbent wound dressings
- Surgical gowns
- Surgical drapes
- Sets of surgical drapes
- Fluid collection pouches
- Nelaton catheters



Rafal Byczkowski

Date: 2019-05-27

10/020 h 04.08 
TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approval



Doc. 5/6, Rev. 0

# TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

For the following medical devices the scope covers only the aspects of manufacture concerned with securing and maintaining sterile conditions:

- Vaginal speculums
- Cervical brushes
- Urine bags
- Tongue depressors
- Guedel airways
- Intubation stylets
- Endotracheal tube holders
- Suction tubes
- Withdrawal cannulas
- Alginate dressings
- Cannula stoppers
- Umbilical cord clamps



Rafal Byczkowski

Date: 2019-05-27

0/020 h 04.08 ®

TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approva



Doc. 6/6, Rev. 0

# TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Gustawa Eiffel'a 15 44-109 Gliwice Poland

Activity: Production

**Notified Body** TÜVRheinla Tifizierung

Date: 2019-05-27

04.08 ®

TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior app

Rafal Byczkowski



# Certificate

The Certification Body of TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

### ZARYS International Group Spolka z ograniczona odpowiedzialnoscia spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Poland

has established and applies a quality management system for medical devices for the following scope:

#### Production and distribution of disposable and reusable medical devices and distribution of in-vitro medical devices (see attachments for sites included)

Proof has been furnished that the requirements specified in

### EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date:

2017-11-07

Certificate Registration No.: SX 60120471 0001

An audit was performed. Report No.: 26300232 006

This Certificate is valid until: 2020-06-08

**Certification Body** 





TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail:cert-validity@de.tuv.com http://www.tuv.com/safety

Date 2017-11-07



Doc. 1/1, Rev. 0

# TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

SX 60120471 0001 26300232 006

**Organization:** 

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Poland

Scope:

Sites included:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia spolka komandytowa ul. Gustawa Eiffel'a 15 44-109 Gliwice Poland

Activity: Production, storage and distribution of medical devices



Date: 2017-11-07

**Certification Body** 

