

LOGIQ Totus™

Probe guide



Featuring XDclear™ Technology

The LOGIQ Totus is GE HealthCare's ultrasound imaging system designed for general imaging applications including abdominal, vascular, obstetric, gynecologic, neonatal, pediatric, urological, transcranial, cardiac, and small parts applications.

	Description	Applications	FOV	Bandwidth	Biopsy Guide	Volume Navigation
Convex						
 C1-6-D C1-6VN-D*	XDclear broad-spectrum convex probe	Abdominal, Obstetrics, Gynecology, Vascular, Musculoskeletal, Pediatric	80°	1 – 6 MHz	Multi-angle disposable with a reusable bracket	Yes * Internal VNav sensor, does not require an external bracket
 C2-7-D C2-7VN-D*	Broad-spectrum convex probe	Abdominal, Pediatric	110°	1 – 6 MHz	Multi-angle disposable with reusable bracket options	Yes * Internal VNav sensor, does not require an external bracket
 C3-10-D	XDclear broad-spectrum convex probe	Neonatal, Pediatrics, Vascular, Small Parts, Neonatal Transcranial, Abdominal	95°	2 – 11 MHz	No	Yes
 IC5-9-D	Broad-spectrum micro-convex intra-cavitary probe	Obstetrics, Gynecology, Urology	179°	3 – 10 MHz	Single-angle disposable or single-angle reusable	Yes
Linear						
 ML6-15-D	Broad-spectrum linear matrix array probe	Vascular, Small Parts, Neonatal, Pediatrics, Neonatal Transcranial, Musculoskeletal, Breast, Abdominal	50 mm	4 – 16 MHz	Multi-angle disposable with a reusable bracket	Yes
 L3-12-D	Broad-spectrum linear probe	Neonatal, Obstetric, Vascular, Musculoskeletal, Small Parts, Pediatric, Neonatal Transcranial, Abdominal, Breast	51 mm	2 – 11 MHz	Multi-angle disposable with a disposable bracket	Yes
 L6-24-D	Broad-spectrum linear probe	Musculoskeletal, Breast, Small Parts, Abdominal, Neonatal Abdominal, Neonatal Transcranial, Peripheral Vascular	26 mm	6 – 20 MHz	No	No
 9L-D	Broad-spectrum linear probe	Vascular, Small Parts, Pediatric, Abdominal, OB/GYN, Neonatal, Neonatal Transcranial, Musculoskeletal, Breast	44 mm	2 – 8 MHz	Multi-angle, disposable with a reusable bracket	Yes



M5Sc-D



6S-D



12S-D



RAB6-D



RIC5-9-D



P2D



P6D



Vscan Air™ CL
convex probe



Vscan Air CL
linear probe

Description	Applications	FOV	Bandwidth	Biopsy Guide	Volume Navigation
Sector					
XDclear broad-spectrum sector probe	Cardiac, Transcranial, Abdominal, Adult Cephalic	120°	1 – 5 MHz	Multi-angle disposable with a reusable bracket	Yes
Broad-spectrum sector probe	Pediatric Cardiac, Pediatric Abdominal	115°	2 – 8 MHz	No	No
Broad-spectrum sector probe	Pediatric, Pediatric Cardiac, Neonatal, Neonatal Cardiac	90°	4 – 12 MHz	No	No
Real-time 4D					
Broad-spectrum real-time 4D probe	Abdominal, Obstetrics, Gynecology, Pediatrics, Neonatal	80°	2 – 8 MHz	Single-angle disposable with a reusable bracket	No
Broad-spectrum real-time 4D micro-convex probe	Obstetrics, Gynecology, Urology	179°	3 – 10 MHz	Single-angle reusable	No
Specialty					
CW split crystal pencil probe	Cardiac, Vascular, Adult Cephalic	N/A	1 – 3 MHz	No	No
CW split crystal pencil probe	Cardiac, Vascular, Adult Cephalic	N/A	5 – 7 MHz	No	No
Wireless					
Broad-spectrum convex probe	Abdominal, Obstetrics, Vascular, Musculoskeletal, Cardiac	60°	2 – 5 MHz	No	No
Broad-spectrum linear probe	Vascular, Thyroid, Scrotal, Breast, Musculoskeletal, Pediatric	38 mm	3– 12 MHz	No	No



GE HealthCare

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

© 2024 GE HealthCare. LOGIQ Totus, XDclear, and Vscan Air are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

March 2024
DOC2976282