


UNDERPAD
REF.NO:800.001C-1.02
PRODUCT DATA SHEET



**PRODUCT DATA SHEET
UNDERPAD-NON STERILE
800.001C-1.02**

1) PRODUCT COMPONENTS



ITEM:1
Underpad
60x60 cm
Absorbent Area
50x50 cm

2) DETAILS OF PRODUCT COMPONENTS

| ITEM | QTY | PRODUCT | PRODUCTION TYPE | DIMENSION | COLOUR |
|------|-----|----------|---------------------------|---------------|--------|
| 1 | 1 | Underpad | Total 26 gr (±2) | 60x60 cm (±3) | White |
| | | | Backsheet-PE Film 15 gsm | | |
| | | | Topsheet 10 gsm | | |
| | | | Pulp-Absorbent Area-16 gr | 50x50 cm (±3) | |

3) TEST RESULT (ISO 11948-1)

| | |
|----------------------------------|-----|
| ABSORPTION CAPACITY (ML)* | 511 |
|----------------------------------|-----|

**PRODUCT DATA SHEET
UNDERPAD-NON STERILE
800.001C-1.02**

4)GENERAL INFORMATION

| | |
|-------------------------------|---|
| Product Description | UNDERPAD Non-Sterile |
| Certificates | EN ISO 13485:2016 |
| Standards | EN ISO 14971:2019 EN ISO 15223-1:2021 |
| Intended Use | Underpad product is used as a mattress protector in health institutions such as hospitals, polyclinics, intensive care units, and home care patients. |
| Instructions For Usage | Cut the products bag open with a scissor from the bottom side of the packaging. Doing so will present you a better place to hold onto the pad as you take it out of the packaging. Take the underpad out and place it in a folded-up state, onto the surface you'll be using it on. Unfold the pad. Take the edge of the product and unfurl it outward, away from you. Allow it to unfold completely. Place the pad down on the surface, with the non plastic side up. The soft side absorbs the urine, while the plastic-looking side will prevent any liquid from going through and onto the surface. The person will lie down on the softer side. |

5) PACKAGING

Packaging

Bag in 30 pcs