

# Sterling™

PTA Balloon Dilatation Catheter



# Demand

from your 0.018" / 0.46 mm Balloon Dilatation Catheter

# More

**The #1 0.018"/ 0.46 mm\*\*  
PTA Balloon Dilatation Catheter**

with the most comprehensive  
size offering available

**2-10 mm diameters  
10-220 mm lengths\***

The most comprehensive  
size offering available

**10 seconds\***

Best-in-class  
deflation time

**0.020"/ 0.51 mm\***

Lowest in-class  
lesion entry profile

The Sterling™ family† is designed to offer best-in-class lesion entry profile and deflation times with the most clinically comprehensive size offering and indications available.



**Demand Exceptional Crossability**

- Exceptionally low 0.020"/ 0.51 mm lesion entry profile on a 0.018"/ 0.46 mm wire
- 7x220 though 5F / 0.65mm



**Demand the Broadest Size Matrix**

- The most comprehensive size offering available in both Over-The-Wire and Monorail™ rapid-exchange
- Indicated for wider acute use including carotid‡ and post stent dilatation

\* Average measurements taken by Boston Scientific (n=3, 5 x 150 mm balloons). Data on file. Bench testing may not be representative of clinical performance.


\*\* US and EU market share data per internal and external analysis.

† The Sterling family includes Sterling and Sterling SL.

‡ Monorail only


# Sterling™ Monorail™

PTA Balloon Dilatation Catheter



# Sterling™ Over-the-Wire

PTA Balloon Dilatation Catheter



Diameter (mm)	Length (mm)												
	Balloon	10	15	20	30	40	60	80	100	120	150	200	220
2.0 mm													
2.5 mm													
3.0 mm													
3.5 mm													
4.0 mm													
4.5 mm													
5.0 mm													
5.5 mm													
6.0 mm													
6.5 mm													
7.0 mm													
8.0 mm													
9.0 mm													
10.0 mm													

Platform	Shaft Length (cm)			
	40	80/90	135	150
OTW				
Rapid Exchange				

**Boston Scientific**  
Advancing science for life™

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome.

PI-195601-AB OCT2014 Printed in Germany by medicalvision.

[www.bostonscientific-international.com](http://www.bostonscientific-international.com)

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved. DINPER4642EB