

# Rotation, tilting, and vibrating

## Versatile and flexible HulaMixer™ Sample Mixer



- Tilt, rotate, and/or vibrate your samples
- Continuous or timed operation
- All settings easily adjustable
- For use at 4°C to 40°C

A lot of labs may not think about it much, but good mixing is important for optimal results. The HulaMixer™ Sample Mixer is perfect for sample preparation with Dynabeads® products and for any other application needing thorough mixing.

### Flexible mixing

The multipurpose HulaMixer™ helps mix practically any sample in any tube. Speed is adjustable from 1 to 100 rpm, and it can rotate in three ways: orbital (end-over-end), reciprocal (tilting), and vibrating (vortex). These motions can be performed individually or in

combinations. The mixer comes with two separate platforms, accommodating tubes ranging from 0.5 to 50 mL.

### User-friendly and portable

The HulaMixer™ Sample Mixer is flexible, yet intuitive and easy to operate (Figure 1). You adjust speed, time, and movement combinations to match your changing needs. The HulaMixer™ comes with a manual, a laminated Quick Guide, and AC/DC adaptors. The motor is extremely quiet. The mixer is small; just carry it to the incubator or coldroom.

### Personalize your HulaMixer™

Every lab needs a flexible, easy-to-use sample mixer. Or maybe you need a couple of them! Want one all for yourself? The HulaMixer™ comes with a selection of durable stickers, allowing you to add your personal touch to your sample mixer.

### Ordering information

Product	Cat. No.
HulaMixer™ Sample Mixer	159-20D

For current prices and further information, visit

[www.invitrogen.com/hulamixer](http://www.invitrogen.com/hulamixer).

**Stop waiting for an available mixer. Order one for yourself or your group today!**



Figure 1. The versatile HulaMixer™ Sample Mixer comes with two tube holder plates to suit all kinds of microtubes and vacutainers (0.5–50 mL).