





Our quality target; to offer a high standard products to our all customers and stakeholders.

To provide employment by opening fields of business areas in the medical sector. To assure qualified and standard services globally.











Karadeniz Medikal started its first activity in the medical sector in 1994. Our company, whose foundations were laid by our founder Mr. Metin Azak in Istanbul Gungoren; started to show itself by making a rapid start in the medical sector, especially with product designs in the first years of its activities.

In the following years, our organization continued to develop its product range with its business partners and experience, which played an important role in its activities. Additionally, It continued to increase this momentum in the following years with the priority of quality products and unconditional customer satisfaction in the medical sector, where it has been operating in a short time with alternative product groups and services offered according to needs. Karadeniz Medikal has managed to become one of the leading brands in Turkey, especially in the "Wheelchair" sector in less than 10 years. In the following years, it has carried the experience gained in this success process beyond the borders of the country to the Asian, Middle Eastern, African and European markets.

The fact that the architects of the adventure, which started with 5 people, insisting of the principle of "Perfect Service" under any circumstances, constituted one of the cornerstones of our success. The motto of "Perfect Service" has led us to catch the "Producing Perfect" perspective. In this direction; it has created an infrastructure for us to create a staff who can compete, is open to innovations and is aware of the concept of innovation, is sensitive to the environment and has a professional perspective. Our organization's perspective of "Right manpower for the right job" in employment indirectly formed the basis of corporate belonging and guided the spontaneous formation of the concept of competent personnel.











We have accelerated and strengthened all our processes by blending our experience in production, import and export processes, which we have gained in Turkey and abroad for many years, with our 4 brands in our new modern production facility, which we established on a 2.000m2 closed area in Kocaeli Dilovası region in 2020. We are proud to be known as the only Turkish brand that produces certified products at international standards in Turkey's own modern facilities.

Today, Karadeniz Medikal is trusted, respected and known in its field among the medical companies in Turkey and aims at this awareness abroad as well. It continues on its way with confident steps by adopting the principle of transferring the learning and development process that has been going on since 1994 to the next generation without compromising its prestige.

## Vision

To be a leading and permanent company among innovative, domestic and foreign medical companies that have production standards and awareness in the global market, where our users who prefer our daily life aid products feel safe in all processes and product expectations are met.

### Mission

To increase the employment areas where our member of workforce work happily and safely in the field of the production of medical products, which we also serve, by offering quality products to our local and international users and stakeholders.

CHAIRMAN OF THE BOARD

Kemal AZAK

V. Anak

## INDEX

# 

## **GOLFI**MANUAL WHEELCHAIR



G099	06
G100	07
G100Y	08
G101	10
G103	11
G130	12
G131	13
G132	14
G133	15
G135	16

G305



G120



G551

G561

PR771

56 JETTY JETTY POWER WHEELCHAIRS

JT099	56	JT200	62
JT100	58	JT311	64
ITW111A	60		

Jetty<sup>®</sup>

SCUBA

SCUBA POWER WHEELCHAIRS & SCOOTERS

S150	67
S190	69
S220	71
S301	73
S550	75
S350	77
REHA2	78

88 LOCO PATIENT AID PRODUCTS



PR891	82	PR444	86	PR104	92
PR449	83	PR445	86	PR105	92
PR888	83	PR447	86	PR106	92
PR881	83	PR 448	86	PR107	93
PR882	83	PR691	88	PR843	93
PR491	84	PR991	89	PR844	93
PR440	85	PR981	90	PR603C	93
PR441	85	PR999	91	PR602	94
PR442	85	PR102	92	PR603	94
PR443	85	PR103	92	PR911	95









Standard Manual Wheelchair

Extra durable structure with double cross bars Foldable frame for easy transport and storage







# GOLFI G100

Standard Manual Wheelchair

Special designed structure for attendant usages

> Foldable frame for easy transport and storage.












Standard Manual Wheelchair



Foldable frame for easy transport and storage  $\, \blacktriangleleft \,$ 











Standard Manual Wheelchair











Functional Manual Wheelchair

- Easy transferable to bed, car etc. with swivel-up armrests
- > Armrest design simplify to use dining table and school desk.
- > Foldable frame for easy transport and storage.
- Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc.











Easy transferable to bed, car etc. with swivel-up armrests Armrest design simplify to use dining table and school desk

Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc

Foldable backrest and quick release axles 
save more space while transported

200\*50 PU front wheels provides ◀ easy driving on rough terrains and street drains



Functional Manual Wheelchair

Easy transferable to bed, car etc. with swivel-up armrests.

Armrest design simplify to use dining table and school desk.

Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc.

The footrests designed to elevate 90° for relaxing with lenghten legs while seated.

Foldable backrest and quick release axles save more space while transported.

200\*50 PU front wheels provides easy driving on rough terrains and street drains.









Functional Manual Wheelchair

Easy transferable to bed,car etc. with swivel-up armrests <a>Armrest design simplify to use dining table and school desk</a> <a></a>

Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc

Foldable backrest and quick release axles <a>save more space while transported</a>

200\*50 PU front wheels provides easy ◀ driving on rough terrains and street drains





## GOLFI G133

Functional Manual Wheelchair

Easy transferable to bed, car etc. with swivel-up armrests

Armrest design simplify to use dining table and school desk

Detachable footrests offering easy storage and extra maneuvera bility into narrow places, elevators etc.

➤ The footrests designed to elevate 90° for relaxing with lenghten legs while seated.

- Foldable backrest and quick release axles save more space while transported
- 200\*50 PU front wheels provides easy driving on rough terrains and street drains











Extraordinary strenghten frame structure and double cross 
brace system offer bariatric users comfortable usage

Easy transferable to bed, car etc. with swivel-up armrests <

Armrest design simplify to use dining table and school desk  $\, \blacktriangleleft \,$ 

Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc

200\*50 PU front wheels provides easy driving ◀ on rough terrains and street drains

Seat and back cushion provides <a>comfort in a long term usages</a>







BariatricManual Wheelchair

Extraordinary strenghten frame structure and double cross brace system offer bariatric users comfortable usage

Easy transferable to bed, car etc. with swivel-up armrests

Armrest design simplify to use dining table and school design

Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc.

- 200\*50 PU front wheels provides easy driving on rough terrains and street drains.
- Seat and back cushion provides comfort in a long term usages.














Multi-functional Manual Wheelchair



- The headrest available positioning up to ergonomics of users <
  - The footrests designed to elevate 90° for relaxing with ◀ lenghten legs while seated
  - The footrest and legrest height designed adjustable up to 
    user leg length
    - 200\*50 PU front wheels provides easy driving ◀ on rough terrains and street drains





# Comfort

It has been produced for the complete comfort of the user, with all the details in mind. It is a class-leading in functional wheelchairs.











Cerebral Palsy Wheelchair



The backrest designed to recline 160° like as sunbed with considering users comfort in a long term usage

Easy to adjust chest holders to fix the user body upright position <

Adjustable positioning front padded support protects 
the users regarding to fall down from the chair

The frame allows tilting between 3°-22° up to user demand <

The footrests designed to elevate 90° for relaxing with lenghten legs while seated

The anti-tipper heel protects the user if fall down while using the chair on the sloping areas and in a sunbed position









# GGG30

Aluminum Wheelchair

Aluminum lightweight foldable framefor easy transport and storage

Armrest design simplify to use dining table and school desk

> Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc

> Foldable backrest and quick release axles save more space while transported.

> Additional back rear wheels allow functional usage indoors.







GGGGG





Aluminum lightweight foldable frame for easy transport and storage <

- Height adjustable armrest design simplify to use dining table and school desk
- Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc
- Curb climb sticks, handles and extra attendant united 
  brakes for convenient attendant usage
  - Foldable backrest and quick release axles save 
    more space while transported
- Additional seat and back cushions provide extra comfortability






GOLFI G637

Aluminum Wheelchair







Height adjustable armrest design simplify to use dining table and school desk

Detachable footrests offering easy storage and extra maneuverability into narrow places, elevators etc

Curb climb sticks, handles and extra attendant united brakes for convenient attendant usage

> Quick release axles save more space while transported

> Additional seat and back cushions provide extra comfortability









Aluminum Patient Transfer Wheelchair

- Special innovative design allows less storage in < public areas(Hospitals, Rehabilition Centers, Airports etc.)
- Rigidstrengthen aluminum durable frame for multiple usages <
- Central self locked brake allows safe and comfort attendant usage 🔌
  - Automatic safety belt, IV Stand and back pocket are more < selective functions
    - Extra durable aesthetical plastic parts for longterm <
      - Custumized color option <



**FEET PUTTING** 





- \* 68% Space saving
- \* Ergonomic and comfortable seat
  - \* Accessory options
  - \* Innovative design
    - \* Colour options
  - \* Noiseless and easy to use
  - \* Footrest and foldable armrest







# PEDIATRIC MANUAL WHEELCHAIR





Multi-functional Pediatric Wheelchair

> The backrest designed to recline 6 level/160° like as sunbed with considering users comfort in a long term usage

> Seat width, depth and backrest height designed for pediatric users

- > The footrests designed to elevate 90° for relaxing with lenghten legs while seated
- The anti-tipper wheel protects the user if fall down while using the chair on the sloping areas and in a sunbed position













Seat width, depth and backrest height designed for pediatric users <

- Easy transferable to bed,car etc. with flip-up armrests <
  - Height adjustable footrests can be adjusted againist to of user leg height
  - Foldable frame for easy transport and storage <
- Additional seat and back cushions provide extra comfortability





GOLFI G333

Aluminum Pediatric Wheelchair

- Seat width, depth and backrest height designed for pediatric users
- > Special designed chest belt stabilize the patient againist to posture defects.
- > Height adjustable footrests can be adjusted againist to of user leg height
- > Adjustable seat width provides longterm usage for pediatric patient growth
- > Additional seat cushion provide extra comfortability









Designed for patients with Cerebral Palsy and Abnormal Posture

> Seat width, depth, footrests and backrest height designed for pediatric users

> Special designed seat for pediatric users comfortability

Dining table offers easy feeding at all areas

> Easy to adjust chest holders to fix the user body upright position

The frame allows tilting between 3°-22° up to user demand

> The anti-tipper heel protects the user if fall down while using the chair on the sloping areas and in a sunbed position









# TRANSFER










Transfer Wheelchair



- Foldable backrest and footrests save more space 
  while transported
  - Curb climb sticks, handles and extra attendant brakes for convenient attendant usage
    - Anti-bacterial upholstery protects the user <a href="againist to bacterial infections">againist to bacterial infections while seated</a>






### Golfi 37



Aluminum Transfer Wheelchair



> Foldable backrest and footrests save more space while transported.

> Curb climb sticks, handles and extra attendant united brakes for convenient attendant usage.













Foldable strengthen aluminum frame for easy transport and storage Foldable backrest and footrests save more space while transported

Curb climb sticks, handles and extra attendant united 
brakes for convenient attendant usage

Anti-bacterial upholstery protects the user <a>againist to bacterial infections while seated</a>




Golfi



G503





Foldable backrest and footrests save more space while transported

Curb climb sticks, handles and extra attendant united brakes for convenient attendant usage

Anti-bacterial upholstery protects the user againist to bacterial infections while seated









GOLFI G605

Aluminum Wheelchair







> Swivel-up footrests offering easy storage and extra maneuverability into narrow places, elevators etc.

> Curb climb sticks, handles and extra attendant united brakes for convenient attendant usage.

> Aluminum forks protect againist to defects in a longterm usages.

> Foldable backrestsave more space while transported.







### Lightweight

Portable, Long Lasting

Fully meet your expectation...



G637
Suits you perfectly!

combination of aesthetics and quality

#mygolfi♥







## SHOWER COMMODE WHEELCHAIR









Easy transferable to bed, car etc. with detachable armrests 

Detachable footrests offering easy storage and extra 

maneuverability into narrow places, elevators etc

The seat designed as a toilet apparatus for <

Anti-bacterial upholstery protects the user <a>againist to bacterial infections while seated</a>





G124

Multi-functional Wheelchair



Dining table offers easy feeding at all areas

➤ The footrests designed to elevate 90° for relaxing with lenghten legs while seated

200\*50 PU front wheels provides easy driving on rough terrains and street drains

The seat designed as a toilet apparatus for eliminate needs with deodizer bucket







G125



Shower Commode Wheelchair



Waterproof upholstery provides a long term usage for bathing

Easy transferable to bed, car etc. with flip-down armrests <

The padded cover of toilet apparatus provides comfort while seated <

Adaptable on closet with excluding the bucket <

The wheels with brake help to immobilize <a> chair for slippery slopes while eliminating</a>







Shower Commode Wheelchair

- Multi-functional structure allows to use indoor and on the closet
- Waterproof upholstery provides a long term usage for bathing
- Easy transferable to bed, car etc. with flip-down armrests
- ➤ The padded cover of toilet apparatus provides comfort while seated. The middle part of seat can slide out to offer easy usage while seated
- > Adaptable on closet with excluding the bucket
- The wheels with brake help to immobilize chair for slippery slopes while eliminating













Transfer Commode Wheelchair

- Easy patient transfer from bed with its special structure <
  - Special designed frame allows long lifetime for usage 
    into bath, toilet etc
    - Easy usage on the closet and indoor with 
      multi-functional structure
    - Polyurethane seat convenient for easy bathing <
      - Hydroulic height adjustability «
  - Easy to adapt on the closet after removing the bucket <









### GOLFI G551

Aluminum Shower Commode Wheelchair

- Easy foldable aluminum structure allows to use indoor and on the closet
- Waterproof PE seat&back provides a long term usage for bathing
- > EVA padded cover of toilet apparatus provides comfort while seated
- > Adaptable on closet with excluding the bucket
- The wheels with brake help to immobilize chair for slippery slopes while eliminating









G561



Aluminum Shower Commode Wheelchair

- Stainless aluminum frame allows long lifetime for <a href="text-usage">usage</a> into bath, toilet etc
  - Easy usage on the closet and indoor with 
    multi-functional structure
    - EVA padded cover of toilet apparatus 

      provides comfort while seated
- Waterproof PE seat&back provides a long term usage for bathing <
- Easy transferable to closet, bed, car etc. with swivel-up armrests <
  - Easy to adapt on the closet after removing the bucket <
- Quick release axles allow to use chair indoor and easy transport <







Adaptable on closet with excluding the bucket

The wheels with brake help to immobilize chair for slippery

slopes while eliminating











# POWER WHEELCHAIRS & SCOOTERS









Easy transferable to bed, car etc. with swivel-up armrests

Foldable frame for easy transport and storage

16" Pneumatic rear wheels& 10" Pneumatic front wheels are protected againist to melt and cracks

Available manual usage

after disengage the motor lever once needed

The anti-tipper wheel protects the user if fall down while using the chair on the sloping areas





Maximum Speed : 6.0 km/h : 15 km

Maximum Range

Turning Radius : 58 cm

: 2xDC 24V 250W Motors

Controller : European Type 45-50 AH

: 2x12 AH Lead Acid Batteries

Charger : 2.8 A Off-Board

Brake : Electromagnetic

Maximum Safe Slope : 12°

Maximum Curb Climb : 15 cm

Backrest Recline :-

Body Type : Folding





JT100





Easy transferable to bed, car etc. with <a href="windstanding-right">swivel-up armrests</a>

Foldable frame for easy transport and storage «

12" PU designed rear wheels are protected 
againist to melt and cracks

200\*50 PU front wheels provides easy ◀ driving on rough terrains and street drains

Available manual usage after <a>disengage the motor lever once needed</a>

The anti-tipper wheel protects the user if 
fall down while using the chair on
the sloping areas





Maximum Speed : 6.0 km/h : 25 km

Maximum Range

Turning Radius : 58 cm

: 2×DC 24V 200W Motors

Controller : European Type 45-50A

Batteries : 2x28 AH Gel

: DC 24V 5A Charger

Brake : Electromagnetic

Maximum Safe Slope :12°

Maximum Curb Climb : 10 cm

Backrest Recline

Body Type : Folding


Product Weight : 60 kg

Carrying Capacity : 136 kg









Maximum Range : 28 km

Turning Radius : 58 cm

: 2×DC 24V 340W Motors

Controller : European Type 45-50A

Batteries : 2x26 AH Gel

: DC 24V 5A Charger

Brake : Electromagnetic

:12° Maximum Safe Slope

Maximum Curb Climb : 10 cm

Backrest Recline :-

Body Type : Folding

Product Weight : 59 kg

Carrying Capacity : 136 kg









Multi-functional Power Wheelchair

- Easy transferable to bed, car etc. with flip away 
  and height adjustable armrests
- Positioning the backrest offers extra comfort to the users <
  - Foldable backrest save more space while transported <
    - 12" PU designed rear wheels are protected 
      againist to melt and cracks
    - 200\*50 PU front wheels provides easy <a href="https://driving.com/driving-nc-up/4">driving on rough terrains and street drains</a>
      - Available manual usage after <a>disengage the motor lever once needed</a>
    - The anti-tipper wheel protects the user if fall down while using the chair on the sloping areas
      - Ergonomic captain seat provides <a>comfort in a long term usages</a>



Maximum Speed : 6.0 km/h : 25 km Maximum Range Turning Radius : 58 cm : 2×DC 24V 250W Motors Controller : European Type 45-50A Batteries : 2x32 AH Gel : DC 24V 5A Charger Brake : Electromagnetic Maximum Safe Slope :12° Maximum Curb Climb : 10 cm Backrest Recline Body Type : Folding Product Weight : 65 kg Carrying Capacity : 136 kg











One button foldability, lithium batteries allows easy <a href="carriage">carriage</a> and transport

Easy transferable to bed, car etc. with swivel-up armrests

Carriage bag for belongings <

12" PU designed rear wheels are protected againist to melt and cracks «

200\*50 PU front wheels provides easy driving on rough terrains and street drains


Available manual usage <a> after disengage the motor lever once needed</a>

The anti-tipper wheel protects the user if fall down while using the chair on the sloping areas

Seat&back cushions provides comfort in a long term usages <







Maximum Range : 30 km

Turning Radius : 50 cm

Motors : 2×DC 24V 250W

Controller : European Type 45-50A

Batteries : 2x7,5 AH Lithium Batteries

Charger : DC 24V 2A

Brake : Electromagnetic

Maximum Safe Slope : 12°

Maximum Curb Climb : 10 cm

Backrest Recline :-

Body Type : Folding









# SCUBA S150 Stair Climbing Power Wheelchair

- Special designed frame with rubber tracks allows movement onto stairs
- Extraordinary strenghten frame structure offers confidential usage
- Foldable frame for easy transport and storage
- Back and front handles allows to carry as stretcher as needed
- Available manual usage after disengage the motor lever once needed
- Seat&back cushion provides comfort in a long term usages



Body Material	: 6063 Aluminum
Battery	: 2x12 AH Lead Acid
Engine	: 120 W - 3000rpm
Maximum Slope	: 5°
Stair Slope	:35±5°
Minimum Ladder Width	: 800mm
Dimensions	: 1.012mm / 512mm / 1.250mm
Maximum Loading Capacity	: 90kg
Product Weight	: 43kg
TSE Standard	: TS EN ISO 7176-28







### SCUBA **5190**

Multi-functional Power Wheelchair

- Extraordinary strenghten frame structure offers confidential usage
- Easy transferable to bed, car etc. with flip away and height adjustable armrests
- Positioning the backrest offers extra comfort to the users
- Foldable backrest save more space while transported
- Back suspension system protects the user with absorbing the shock hits
- ➤ 14"(3.00-8) Pneumatic drive wheels help the user drive the chair on rough terrains, lands and high sidewalked cities easily
- 9"(2.80/2.50-4) Pneumatic front wheels provides easy driving on rough terrains and street drains
- Available manual usage after disengage the motor lever once needed
- Ergonomic captain seat provides comfort in a long term usages



Maximum Speed	: 6.0 km/h
Maximum Range	: 40 km
Turning Radius	: 58 cm
Motors	: 2×DC 24V 450W
Controller	: PG VR2 60AH
Batteries	: 2x55 AH Gel
Charger	: DC 24V 5A
Brake	: Electromagnetic
Maximum Safe Slope	: 23°
Maximum Curb Climb	: 26 cm
Backrest Recline	: 45°
Product Weight	: 90 kg
Maximum Loading Capacity	: 136 kg









### SCUBA **S220**

Multi-functional Power Wheelchair

- Extraordinary strenghten frame structure offers confidential usage.
- Easy transferable to bed, car etc. with flip away and height adjustable armrests.
- Positioning the backrest offers extra comfort to the users.
- Back warning lamb and lighting system offers more visibility.
- Foldable backrest save more space while transported.
- Back suspension system protects the user with absorbing the shock hits.
- ➤ 14"(3.00-8) Pneumatic drive wheels help the user drive the chair on rough terrains, lands and high sidewalked cities easily.
- 9"(2.80/2.50-4) Pneumatic front wheels provides easy driving on rough terrains and street drains.
- Available manual usage after disengage the motor lever once needed.
- Ergonomic captain seat provides comfort in a long term usages.



Maximum Speed	:6KM/H
Maximum Range	: 40 km
Turning Radius	: 82 cm
Motors	: 2xDC 24V 450W
Controller	: PG VR2 70A
Batteries	: 2x55 AH Gel
Charger	: 5 Amp Off-Board
Brake	: Electromagnetic
Maximum Safe Slope	:12°
Maximum Curb Climb	: 10 cm
Backrest Recline	: 45°
Body Type	: Fixed
Product Weight	: 104 kg
Maximum Loading Capacity	: 180 kg













- Easy transferable to bed, car etc. with flip away and height adjustable armrests.
- Positioning the backrest&legrests offers extra comfort to the users.
- Back suspension system protects the user with absorbing the shock hits.
- 14"(3.00-8) Pneumatic drive wheels help the user drive the chair on rough terrains, lands and high sidewalked cities easily.
- 10"(3.00-4) Pneumatic front wheels provides easy driving on rough terrains and street drains.
- Available manual usage after disengage the motor lever once needed.
- The anti-tipper wheel protects the user if fall down while using the chair on the sloping areas.
- Back&Seat cushion provides comfort in a long term usages.
- Lighting system offers more visibility.



SCUBA **5301** 

Standing Power Wheelchair

 Maximum Speed
 : 10 km/h

 Maximum Range
 : 22 km

 Motors
 : 2x DC 24V 340W

 Controller
 : Dynamic 60A


 Batteries
 : 2x55 AH Gel

 Maximum Safe Slope
 : 12°

 Product Weight
 : 143 kg

Maximum Loading Capacity: 140 kg









	Maximum Speed	: 6.4 km/h			
	Max. Range Distance	: 20 km			
	Turning Radius	: 1.35 m (53").			
	Motors	: DC Motor			
	Batteries	: 48V Lithium			
	Charger	: 110-240V AC			
	Brake	: Electromagnetic			
	Maximum Safe Slope	: 6°			
	Body Type	: Folding			















- Extraordinary strenghten frame structure offers confidential usage.
- Back warning lamb and lighting system offers more visibility
- Foldable backrest save more space while transported
- Back suspension system protects the user with absorbing the shock hits
- ➤ 410/350-5" Pneumatic drive&front wheels help the user drive the chair on rough terrains, lands and high sidewalked cities easily
- Available manual usage after disengage the motor lever once needed
- Ergonomic captain seat provides comfort in a long term usages

Maximum Speed	: 12 KM/H
Maximum Range	: 40 km
Turning Radius	: 82 cm
Motors	: 1x DC 24V 600W (MTM)
Controller	: PG S-Drive 120A
Batteries	: 2x40 AH Gel
Charger	: 5 Amp Off-Board
Brake	: Electromagnetic
Maximum Safe Slope	: 12°
Maximum Curb Climb	: 10 cm
Backrest Recline	: 45°
Body Type	: Fixed







Max. Capacity

Carton dimensions

: 136 kg

: 44x40x7 cm



#### AIR CUSHION!

Exclusive for seekers to want comfort with affordable price

it is one step ahead of its peers with 1.5 kg product weight in 44x40x7 cm dimensions







### PATIENT AID PRODUCTS



## PR891 Aluminum Luxury Rollator















Barcode No: 8698785010325





Barcode: 8698785010349

















Barcode No: 8698785010967





Barcode No: 8698785010950















Barcode No: 8698785010271





Barcode No: 8698785010998



LOCO
PR448
Aluminum Gradual Walker



LOCO

#### **PR448**

**Aluminum Gradual** Walker

You are safe with non-slip rubber tips.



























Barcode: 8698785010837



Barcode No: 8698785010851



Barcode No: 8698785010820



Barcode No: 8698785010844











Barcode No: 8698785010936























Barcode No: 8698785010882 (S)

8698785010912 (M)

8698785010318 (L)




# TECHNICAL SPECIFICATIONS




















			Ţ		Ī	Ţ	1	
			PR-102	PR-103	PR-104	PR-105	PR-106	PR-107
d.	Height	:	71-93.5 cm	71-93.5 cm	80-90 cm	80-90 cm	72 cm	116 cm
NG	Product Weight	:	0.34 kg	0.34 kg	0.36 kg	0.34kg	0.88 kg	0.16 kg
≉.	Material	:	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
*	Color of Parts	:	Bronze	Grey	Black	Bronze	Grey	Grey-Bronze
8	Dimensions of Carton	:	65x15x28 cm	65x15x28 cm	63x16x27 cm	63x16x27 cm	86x31x41 cm	39x13x29 cm
8	Piece / Carton	:	20 PCS / Carton	10 PCS / Carton	20 PCS/ Carton			

			Ţ Ţ	<u></u>	7	7	7	F
			PR-843	PR-844	PR-602	PR-603	PR-603C	PR-691
क्ष	Height	:	74-98 cm	74-98 cm	86-116 cm	104-132 cm	60-77 cm	84-120 cm
MS	Product Weight	:	0.74 kg	0.84 kg	0.50 kg	0.62 kg	0.50 kg	0.92 kg
≉.	Material	:	Aluminum	Aluminum	Aluminum	Alüminyum	Aluminum	Aluminum
٤.	Color of Parts	:	Grey	Grey	Black	Grey	Grey	5 Colors
8	Dimensions of Carton	:	81x33x37 cm	80x26x45 xm	86x26x29 cm	100x34x29 cm	26x26x29 cm	79x11x21 cm
8	Piece / Carton	:	10 PCS / Carton	10 PCS / Carton	10 Pair / Carton	10 Pair / Carton	10 Pair / Carton	1 Pair / Carton




			PR-981	PR-991	PR-999	
th.	Height	:	107-130 cm	106-144 cm	88-129 cm	
r r	Product Weight	:	1.02 kg	1.08 kg	1.08 kg	
≉.	Material	3	Aluminum	Aluminum	Aluminum	
£.	Color of Parts	;	Grey	Purple-Blue-Red-Yellow	Grey	
8	Dimensions of Carton	3	20x105x130 cm	106x12x23 cm	81x7x26 cm	
8	Piece / Carton	:	10 Pair / Carton	1 Pair	4 Pair / Carton	



**PR882** 

Rollator



			PR-447	PR-448	PR-449	PR-491
0	Length	:	49 cm	48 cm	48 cm	45 cm
$\longleftrightarrow$	Width	:	64 cm	66 cm	66 cm	63 cm
1	Height		75-94 cm	81-99 cm	75-93 cm	74-93 cm
Ê	Product Weight	:	2.24 kg	2.24 kg	2.64 kg	2.88 kg
҈.	Material	:	Aluminum	Aluminum	Aluminum	Aluminum
<b>*</b>	Color of Parts	:	Bronze	Grey	Grey	Blue
8	Dimensions of Carton	:	55x30x80 cm	55x28x82 cm	57x44x79 cm	60x10x81 cm
8	Piece / Carton	:	4 PCS / Carton	4 PCS / Carton	4 PCS / Carton	1 PCS / Carton





PR888 Rollator



NOTES
***************************************









