




Home » Products » Medical Imaging Products » Angiographic Syringes Accessories



Angiographic Syringes Accessories

· The high pressure angiography connecting tubing is used together with Angiographic Syringe for angiography clinically.

INQUIRY

FEATURES

- VENTED SPIKE.
- "J" Fill Tube
- Connecting Tubing (300psi)



VENTED SPIKE

Product No.	Length(mm)
42.18.20000	20 ,Bule
42.18.20001	34,Red

“J” Fill Tube

Product No.
42.18.10000



PRODUCT QUICK LOOK

Product at a glance -



D A P A M E T E D



42.02.90000	Y tube with two check valves	150	1.8	3.6
42.02.90001	Coiled tube	150	1.6	3
42.02.90002	Y tube with two check valves	250	1.8	3.6
42.02.90003	Y tube with one check valve	150	1.8	3.6
42.02.90004	Y tube with one check valve	250	1.8	3.6
42.02.90005	Double infusion line with two check valves	150	2.8	4
42.02.90006	Straight tube	150	1.6	3
42.02.90007	Straight tube	20	1.6	3
42.02.90008	Straight tube with two male check valves	150	1.8	3.6
42.02.90009	Straight tube with two male check valves	250	1.8	3.6
42.02.90010	Coiled tube	150	1.8	3.6
42.02.90011	Straight tube with one needles adapter and spike	25	1.8	3.6
42.02.90012	Y tube with two check valves	15	1.8	3.6



42.02.90015	Straight tube with one check valve	12	1.8	3.6
42.02.90016	T tube with one coild tube 150cm and 100cm spike tube	150	1.8	3.6

PRODUCT RECOMMENDATION



Angiographic Syringes

What else do you want to know? Contact us right now!



*Verify Code

hpk U

*We respect your privacy. When you submit your contact information, we agree to only contact you in accordance with our [Privacy Policy](#).

Submit

We'd Love to Hear From You

Curious about product features, price quote, registration, delivery time, etc? We're ready to answer any and all questions within one business day.

Contact Us

Contact



NO.4 Baolong 6th Road, Baolong industrial Town,
Longgang District, Shenzhen, P.R. China



Copyright © SCW Medicath Ltd All rights reserved. [Privacy](#) [Site Map](#) [Terms of use](#)