Cordis STABILIZER[™] Support Guidewire and Cordis STABILIZER[™] Extra Support Guidewire

0.014" Peripheral Steerable Guidewires



Providing All the Support You Need for 0.014" Peripheral Interventions



Cordis STABILIZER[™] Support Guidewire and Cordis STABILIZER[™] Extra Support Guidewire

0.014" Peripheral Steerable Guidewires

Versatility

Broad transitions and high support for device deliverability and compatibility with OUTBACK™ Elite Re-Entry Catheter

Enhanced Lubricity

Proximal DURAGLIDE™ PTFE Spray and distal PTFE Sleeve for excellent surface lubricity

Atraumatic Tip

SuperSoft tip flexibility provides atraumatic vessel entry and navigation

Ordering Information

Cordis STABILIZER[™] Support Guidewire

| Catalog Numbers | Diameter (inches) | Tip Flexibility | Length (cm) | Tip Shape |
|--------------------|-----------------------------|--------------------|-----------------------|--------------|
| 507 - 180S | 0.014 | SuperSoft | 180 | Straight |
| 507 - 300S | 0.014 | SuperSoft | 300 | Straight |

Cordis STABILIZER[™] Extra Support Guidewire


| Catalog Numbers | Diameter (inches) | Tip Flexibility | Length (cm) | Tip Shape |
|--------------------|-----------------------------|--------------------|-----------------------|--------------|
| 527 - 180E | 0.014 | SuperSoft | 180 | Straight |
| 527 - 300E | 0.014 | SuperSoft | 300 | Straight |

Compatible with **OUTBACK[™] Elite Re-Entry Catheter**

| Codes | Length (cm) | Sheath (F) | - | - |
|------------|----------------|---------------|---|---|
| OTB 59080A | 80 | 6 | - | - |
| OTB 59120A | 120 | 6 | - | - |


Cordis **STABILIZER™** Support Guidewire

Support for the delivery of advanced peripheral devices in complex and tortuous anatomy



Cordis **STABILIZER™** Extra-Support Guidewire

When a higher level of support is needed for the most challenging peripheral interventions



For Healthcare Professionals Only. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy. Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO, DURAGLIDE, OUTBACK and STABILIZER are trademarks of Cordis and may be registered in the US and/or in other countries.

