




**We do ostomy bag**




桐乡市宝禾医疗器械有限公司

**Our story-since 2012**

**We find eureguard hydrocolloid**

**Soft  
comfortable**



# Skin high-match technology

## Soft & comfortable

Instead of using high adhesive material, we improve the adhesive ability by increasing the adhesive area onto the skin.

The leakage may occur through the wrinkles.



BAO-Health adhesive could fill in the wrinkles to avoid leakage



Adhesive material matches to the micro-surface of the skin, enter into the wrinkles and block the leakage channels.

Less skin irritation will occur by the Skin High-Match technology

No hurt when peel off the pouch from skin  
Anti-leakage



## A Simple test

With a slight pressure, the decorative pattern shape was printed on the surface of the adhesive






Plate cause maceration:  
The skin turn to white




Cureguard plate:  
No maceration



Cureguard has super-absorbent capacity, can absorb 800% of its own weight moisture, and quickly transferred to a material inside, does not form jelly-like gel that can cause skin maceration.





Easy to peel off whenever you want to change the bag. Minimal or NO damage to epithelial cells, prevent skin tearing and prevent irritation.










# Clinic report from Germany market

PUBLISHED ON:

<https://de.surveymonkey.net/results/SM-DZ2TQJ2F/>


Wie würden Sie den Zustand der Umgebung der Haut um das Stoma nach Anwendung des BAO FREE beschreiben?

Beantwortet: 32 Umfrageteilnehmer: 1



Antwortkategorie	Beurteilungstypus	Anteil
unauffällig (keine Probleme)		33
entzündet		0
gereizt		1
irritiert		2
andere		0
Beitrag insgesamt: 32		

How skin situation after using BAO HEALTH?



Antwortkategorie	Erwartungen	Anteil
unauffällig (keine Probleme)		29
entzündet		0
gereizt		1
irritiert		3
andere		0
Beitrag insgesamt: 32		



# Tender item introduction



CODE: 1166000F



1. Flange size: Diameter=9cm  
Thickness 0.75 mm
2. Fasten flange size: Diameter:11.5cm
3. Maximum cutting: up to 60mm
4. Pouch size: 15cm x 28cm
5. With filter
6. Water-resistance nonwoven cover (window design)
7. Pouch color: semi-transparent (skin color)
8. Clamp closure: 1clamp for 5 pouches
9. Volume: ≥550ml







Thank you