



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

---

## STATEMENT

We, Zhejiang Orient Gene Biotech Co., Ltd , having a registered office at 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China assign SRL SANMEDICO having a registered office at A. Corobceanu street 7A, apt. 9, Chişinău MD-2012, Moldova, as exclusive authorized representative for Orient Gene Brand product in correspondence with the conditions of directive 98/79/EEC in Moldova only. The detailed product list is in the Annex 1 in the following pages.

We declare that the company mentioned above is authorized to register, notify, renew or modify the registration of medical devices on the territory of the Republic of Moldova.

This Statement letter will be valid from Mar.10th,2025to Mar.09th, 2027.

Zhejiang Orient Gene Biotech Co., Ltd

General Manager

Date:2025/3/10



---

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号  
Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China  
电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

Annex 1

Numărul de catalog (referință)*	Denumire generică (denumirea dispozitivului)	Denumire comercială (brand)*	Modelul
GCCOV-702a-H1	TEST RAPID	Orient Gene	Rapid COVID-19 Antigen Oral Fluid Self-Test
GCCOV-702a-H5	TEST RAPID	Orient Gene	Rapid COVID-19 Antigen Oral Fluid Self-Test
GCCOV-702a-H20	TEST RAPID	Orient Gene	Rapid COVID-19 Antigen Oral Fluid Self-Test
GCCOV-502a-H1OGE	TEST RAPID	Orient Gene	Rapid COVID-19 (Antigen) Self-test
GCCOV-502a-H5OGE	TEST RAPID	Orient Gene	Rapid COVID-19 (Antigen) Self-test
GCCOV-502a-H20OGE	TEST RAPID	Orient Gene	Rapid COVID-19 (Antigen) Self-test

GIHSA-102a	TEST RAPID	Orient Gene	One Step Microalbumin Test Cassette
GIHSA-101a	TEST RAPID	Orient Gene	One Step Microalbumin Test Strip (Urine)
GCROA-602a	TEST RAPID	Orient Gene	Rotavirus rapid test cassette (feces)
GCMAL(pf/pv)-402a	TEST RAPID	Orient Gene	Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)
GCROA/ADE-602a	TEST RAPID	Orient Gene	Rotavirus and Adenovirus Combo Rapid Test Cassette (Feces)
GASPE-902a	TEST RAPID	Orient Gene	Male Fertility Rapid Test Cassette (Semen)
GAFSH-101a	TEST RAPID	Orient Gene	One Step Menopause Test Strip (Urine)
GAFSH-102a	TEST RAPID	Orient Gene	One Step Menopause Test Cassette (Urine)
GAIGF1-502a	TEST RAPID	Orient Gene	iGFBP-1 Rapid test Cassette (Swab)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GAIGF1-501a	TEST RAPID	Orient Gene	iGFBP-1 Rapid Test Strip (Cervical Secretion)
GALH-101a	TEST RAPID	Orient Gene	One Step Ovulation Test Strip (Urine) (25mlU)
GALH-101b	TEST RAPID	Orient Gene	One Step Ovulation Test Strip (Urine) (40mlU)
GALH-102a	TEST RAPID	Orient Gene	One Step Ovulation Test Cassette (Urine) (25mlU)
GALH-102b	TEST RAPID	Orient Gene	One Step Ovulation Test Cassette (Urine) (40mlU)
GCTYP-302a	TEST RAPID	Orient Gene	Typhoid IgG/IgM Rapid Test Cassette (serum/plasma)
GCMAL(pf/pan)-402a	TEST RAPID	Orient Gene	Malaria P.f./Pan Ag Rapid Test Cassette (Whole Blood)
GCDEN-425a	TEST RAPID	Orient Gene	Dengue NS1+IgM/IgG Combo Rapid Test Cassette (Whole blood/serum/plasma)
GCDEN(NS)-402c	TEST RAPID	Orient Gene	Dengue NS1 Antigen Rapid Test Cassette (Whole blood/serum/plasma)
GCDEN(ab)-402c	TEST RAPID	Orient Gene	Dengue IgM/IgG Rapid Test Cassette (Whole blood/serum/plasma)
GCVCH(O1/O9)-602a	TEST RAPID	Orient Gene	V.cholerae O1/O139 Ag Combo Rapid Test Cassette (Feces)
GCMAL(pf)-402a	TEST RAPID	Orient Gene	Malaria Pf Ag Rapid Test Cassette (Whole blood)
GCSAL(ST)-602a	TEST RAPID	Orient Gene	S. typhi Ag Rapid Test Cassette (Serum/plasma/Feces)
GCCHK(IgM)-402a	TEST RAPID	Orient Gene	Chikungunya IgM Rapid Test Cassette (Whole blood/Serum/Plasma)
GCCOV-502a-NN	TEST RAPID	Orient Gene	Coronavirus Ag Rapid Test Cassette (Swab)
GCCOV (Nab)-402b	TEST RAPID	Orient Gene	SARS-CoV-2 Neutralizing Antibody Rapid Test Cassette (Whole blood/Serum/Plasma)
GCCOV-502a-NA	TEST RAPID	Orient Gene	Coronavirus Ag Rapid Tests Cassette (Swab)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GCCOV (Ag)-PN10	TEST RAPID	Orient Gene	COVID-19 Antigen Control Kit
GCCOV (Ag)-PN20	TEST RAPID	Orient Gene	COVID-19 Antigen Control Kit
GCFERA-545a	TEST RAPID	Orient Gene	Flu, COVID-19,RSV&Adeno Ag Combo Tests Cassette (Swab)
GCTV-502a	TEST RAPID	Orient Gene	Trichomonas Ag Rapid Test Cassette (Swab)
GCVCH(O1)-602a	TEST RAPID	Orient Gene	V.cholerae O1 Ag Rapid Test Cassette (Feces)
GCCHA-402a	TEST RAPID	Orient Gene	Chagas Ab Rapid Test Cassette (Whole blood/serum/plasma)
GCMAL(pf/pv Ab)-302a	TEST RAPID	Orient Gene	Malaria Pf/Pv Ab Rapid Test Cassette (Serum/plasma)
GCMAL(pf/pv Ab)-402a	TEST RAPID	Orient Gene	Malaria Pf/Pv Ab Rapid Test Cassette (Whole blood/Serum/plasma)
GCMKP-502b	TEST RAPID	Orient Gene	Monkeypox Ag Rapid Test Cassette (Swab)
GCCOV(Del)-T502a	TEST RAPID	Orient Gene	SARS-CoV-2 Delta-series Mutant Strain Ag Rapid Test cassette (Swab)
GCCOV-PN10	TEST RAPID	Orient Gene	COVID-19 IgG/IgM Control kit
GCCOV-PN20	TEST RAPID	Orient Gene	COVID-19 IgG/IgM Control kit
GCFER-T525a	TEST RAPID	Orient Gene	COVID-19/Flu A&B/ RSV Ag Combo Rapid Test Cassette (Swab)
GCCOV(B117)-525a	TEST RAPID	Orient Gene	COVID-19 Ag&B.1.1.7 Mutant Strain Combo Test Cassette (Swab)
GCFERA-T525a	TEST RAPID	Orient Gene	COVID-19/Flu A&B/ RSV/Adeno Ag Combo Rapid Test Cassette (Swab)
GCCOV-702a	TEST RAPID	Orient Gene	COVID-19 Ag Rapid Test Cassette (Oral fluid)
GCFE-T502a	TEST RAPID	Orient Gene	COVID-19/Flu A&B Ag Combo Rapid Test Cassette (Swab)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax:+86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GCMKP-402a	TEST RAPID	Orient Gene	Monkeypox IgG/IgM Rapid Test Cassette (Whole blood/serum/plasma)
GCKa1-401a	TEST RAPID	Orient Gene	Leishmania Ab Rapid Test strip (Whole blood/serum/plasma)
GCKa1-T402a	TEST RAPID	Orient Gene	Leishmania IgG/IgM Rapid test cassette (whole blood/serum/plasma)
GCCOV-503a	TEST RAPID	Orient Gene	Smart Rapid COVID-19 Ag Test Device
GCBRU-402a	TEST RAPID	Orient Gene	Brucella Antibody Rapid Test Cassette (Whole blood/serum/plasma)
GCCHA-302a	TEST RAPID	Orient Gene	Chagas Ab Rapid Test Cassette (Serum/plasma)
GCCHK(IgM)-302a	TEST RAPID	Orient Gene	Chikungunya IgM Rapid Test Cassette (Serum/Plasma)
GCCOV-501a	TEST RAPID	Orient Gene	Rapid COVID-19 Antigen Test Strip
GCMON-352a	TEST RAPID	Orient Gene	Mononucleosis IgG/IgM Rapid Test Cassette (Serum/plasma)
GCMON-402a	TEST RAPID	Orient Gene	Mononucleosis Rapid Test Cassette (Whole blood/Serum/plasma)
GCMON-425a	TEST RAPID	Orient Gene	Mononucleosis IgG/IgM Rapid Test Cassette (Whole blood/Serum/plasma)
GCEV71 (IgM)-302a	TEST RAPID	Orient Gene	EV71 IgM Rapid Test Cassette (Serum/plasma)
GCEV71 (IgM)-402a	TEST RAPID	Orient Gene	EV71 IgM Rapid Test Cassette (Whole blood/Serum/plasma)
GENMP22-102a	TEST RAPID	Orient Gene	One Step Nuclear Matrix Protein 22 Test Cassette (Urine)
GEFOB/TF-602a	TEST RAPID	Orient Gene	Fecal Occult Blood and Transferrin Combo Rapid Test Cassette (Feces)
GCHEV-302a	TEST RAPID	Orient Gene	HEV IgM Rapid Test Cassette (Serum/Plasma)
GCMP (IgM)-302a	TEST RAPID	Orient Gene	M.pneumonia IgM Rapid Test Cassette (Serum/plasma)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

FCCOV-502a	TEST RAPID	Orient Gene	SARS-CoV-2 Ag Fluorescence Rapid Test Cassette (Swab)
GCCOV-402Ba	TEST RAPID	Orient Gene	COVID-19 IgG/IgM Rapid test cassette (whole blood/serum/plasma)
GCCOV-402a	TEST RAPID	Orient Gene	COVID-19 IgG/IgM Rapid test cassette (whole blood/serum/plasma)
GCCOV-502a	TEST RAPID	Orient Gene	Coronavirus Ag Rapid Test Cassette (Swab)
GCFC-T503a	TEST RAPID	Orient Gene	Smart Rapid COVID-19 &Flu A/B Ag Test Device
GAHCG-101a	TEST RAPID	Orient Gene	One step pregnancy test strip (urine)
GAHCG-101d	TEST RAPID	Orient Gene	One step pregnancy test strip (urine)
GAHCG-101b	TEST RAPID	Orient Gene	One step pregnancy test strip (urine)
GAHCG-102a	TEST RAPID	Orient Gene	One step pregnancy test cassette (urine)
GAHCG-102d	TEST RAPID	Orient Gene	One step pregnancy test cassette (urine)
GAHCG-102b	TEST RAPID	Orient Gene	One step pregnancy test cassette (urine)

GEFOB-602c	TEST RAPID	Orient Gene	Fecal Occult Blood Rapid Test Cassette (Feces)
GEFOB-602b	TEST RAPID	Orient Gene	Fecal Occult Blood Rapid Test Cassette (Feces)
GEFOB-601c	TEST RAPID	Orient Gene	Fecal Occult Blood Rapid Test Strip (Feces)
GEFOB-601b	TEST RAPID	Orient Gene	Fecal Occult Blood Rapid Test Strip (Feces)
GECEA-402a	TEST RAPID	Orient Gene	Carcinoembryonic Antigen Rapid Test Cassette (Whole blood/serum/plasma)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号  
Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China  
电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GECEA-401a	TEST RAPID	Orient Gene	Carcinoembryonic Antigen Rapid Test Strip (Whole blood/serum/plasma)
GETF-602a	TEST RAPID	Orient Gene	Transferrin Rapid Test Cassette (Feces)
GETF-601a	TEST RAPID	Orient Gene	Transferrin Rapid Test Strip (Feces)
GEAFP-401a	TEST RAPID	Orient Gene	Alpha-Fetoprotein Rapid Test Strip (Whole blood/serum/plasma)
GEAFP-402a	TEST RAPID	Orient Gene	Alpha-Fetoprotein Rapid Test Cassette (Whole blood/serum/plasma)
GIHSA-101a	TEST RAPID	Orient Gene	One step microalbumin test strip (urine)
GIHSA-102a	TEST RAPID	Orient Gene	One step microalbumin test cassette (urine)
GDCAR-335a	TEST RAPID	Orient Gene	Myoglobin/CK-MB/Troponin I Combo Test Cassette (Serum/plasma)
GDCKM-302a	TEST RAPID	Orient Gene	One step CK-MB Test Cassette (Serum/Plasma)
GDCKM-402a	TEST RAPID	Orient Gene	CK-MB Rapid Test Cassette (Whole blood/serum/plasma)
GDCRP-402a	TEST RAPID	Orient Gene	C-Reactive Protein Semi-Quantitative Rapid Test Cassette (Whole Blood/serum/plasma)
GDTRO-302a	TEST RAPID	Orient Gene	Troponin I Rapid Test Cassette (Serum/Plasma)
GDTRO-402a	TEST RAPID	Orient Gene	Troponin I Rapid Test Cassette (Whole blood/Serum/Plasma) <b>(Except the tender No. ocds-b3wdp1-MD-1722410248839 din 05.09.2024, limited to the quantity 28060 pcs only, as per the tender)</b>
GDTRO-402b	TEST RAPID	Orient Gene	Troponin I Rapid Test Cassette (Whole blood/Serum/Plasma)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GDMYO-402a	TEST RAPID	Orient Gene	Myoglobin Rapid Test Cassette (Whole blood/serum/plasma)
GDCAR-W435a	TEST RAPID	Orient Gene	Myoglobin/CK-MB/Troponin I Combo Rapid Test Cassette (whole blood/Serum/plasma)
GDPCT-402a	TEST RAPID	Orient Gene	Procalcitonin Rapid Test Cassette (Whole blood/serum/plasma)
GDPCT-T402a	TEST RAPID	Orient Gene	Procalcitonin Semi-Quantitative Rapid Test Cassette (Whole blood/serum/plasma)
GDPCT-T401a	TEST RAPID	Orient Gene	Procalcitonin Semi-Quantitative Rapid Test Strip (Whole blood/serum/plasma)
FDPCT -302a	TEST RAPID	Orient Gene	Procalcitonin Rapid Test Kit (serum/plasma)
GDDDI-402b	TEST RAPID	Orient Gene	D-Dimer Rapid Test Cassette (Whole blood/plasma)
FDCAR-T302a	TEST RAPID	Orient Gene	Troponin I/CK-MB/Myoglobin Fluorescence Combo Test Kit (Serum/plasma)
FDTRO-302a	TEST RAPID	Orient Gene	Troponin I Fluorescence Rapid Test Kit (Serum/plasma)
FDBNP-302a	TEST RAPID	Orient Gene	NT-ProBNP Fluorescence Rapid Test Kit (Serum/plasma)
FDCRP-402a	TEST RAPID	Orient Gene	C-Reactive Protein Rapid Test Kit (Whole blood/serum/plasma)
GDCKM-402a	TEST RAPID	Orient Gene	CK-MB Rapid Test Cassette (whole blood/serum/plasma)
GAHCG-201a	TEST RAPID	Orient Gene	One step pregnancy test strip (Urine/serum)
GAHCG-202a	TEST RAPID	Orient Gene	One step pregnancy test cassette (Urine/serum)
GAHCG-201b	TEST RAPID	Orient Gene	One step pregnancy test strip (Urine/serum)
GAHCG-202b	TEST RAPID	Orient Gene	One step pregnancy test cassette (Urine/serum)
GCHAV(IgM)-302Ba	TEST RAPID	Orient Gene	HAV IgM Rapid Test Cassette (Serum/plasma)
GCHAV-602a	TEST RAPID	Orient Gene	HAV Ag Rapid Test Cassette (Feces)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GCHAV(IgG/IgM) -302a	TEST RAPID	Orient Gene	HAV IgG/IgM Rapid Test Cassette (Serum/plasma)
GCHSV(IgG)-402a	TEST RAPID	Orient Gene	HSV IgG Rapid Test Cassette (Whole blood/serum/plasma)
GCHSV(IgM)-302a	TEST RAPID	Orient Gene	HSV IgM Rapid test Cassette (serum/plasma)
GCHP-601a	TEST RAPID	Orient Gene	H.pylori Ag Rapid Test Strip (feces)
GCHP-602a	TEST RAPID	Orient Gene	H.pylori Ag Rapid Test Cassette(feces)
GCTB-302a	TEST RAPID	Orient Gene	Tuberculosis IgG/IgM Rapid Test Cassette (serum/plasma)
GCTB-402a	TEST RAPID	Orient Gene	Tuberculosis IgG/IgM Rapid Test Cassette (whole blood/serum/plasma)
GCFLU(A/B)-501a	TEST RAPID	Orient Gene	Influenza A&B Ag Rapid Test Strip (Swab)
GCFLU(A/B)-502a	TEST RAPID	Orient Gene	Influenza A&B Ag Rapid Test Cassette (Swab)
GCFLU(A/B)-502Ca	TEST RAPID	Orient Gene	Influenza A&B Ag Rapid Test Cassette (Swab)
GCFLU(A)-501a	TEST RAPID	Orient Gene	Influenza A Ag Rapid Test Strip (Swab)
GCFLU(A)-502a	TEST RAPID	Orient Gene	Influenza A Ag Rapid Test Cassette (Swab)
GCHP-301a	TEST RAPID	Orient Gene	H.Pylori Ab Rapid Test Strip (serum/plasma)
GCHP-302a	TEST RAPID	Orient Gene	H.pylori Ab Rapid Test Cassette (serum/plasma)
GCHP-401a	TEST RAPID	Orient Gene	H.pylori Ab Rapid Test Strip (Whole blood/serum/plasma)
GCHP-402a	TEST RAPID	Orient Gene	H.pylori Ab Rapid Test Cassette (Whole blood/serum/plasma)
GCCA-502a	TEST RAPID	Orient Gene	Candida albicans Antigen rapid test cassette (swab)
GCGON-502b	TEST RAPID	Orient Gene	Gonorrhea Rapid Test Cassette (Swab)
GCGIA-602a	TEST RAPID	Orient Gene	Giardia lamblia Antigen Rapid tests cassette (feces)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GCSTR-501a	TEST RAPID	Orient Gene	Strep A Rapid Test Strip (Throat swab)
GCSTR-502a	TEST RAPID	Orient Gene	Strep A Rapid Test Cassette (Throat swab)
GCFC-525a	TEST RAPID	Orient Gene	Rapid COVID-19 + Influenza Antigen Test
GCRSV-502a	TEST RAPID	Orient Gene	RSV Antigen Rapid Test Cassette (swab)
GCADE-502a	TEST RAPID	Orient Gene	Adenovirus antigen rapid test cassette (swab)
GCADE-602a	TEST RAPID	Orient Gene	Adenovirus Rapid test cassette (feces)
GCCD(GDH)-602a	TEST RAPID	Orient Gene	Clostridium difficile Antigen GDH Rapid Test cassette (feces)
GCCD (Toxin A/B)-602a	TEST RAPID	Orient Gene	Clostridium difficile Toxin A&B rapid test cassette (feces)
GCCD-602a	TEST RAPID	Orient Gene	Clostridium difficile GDH & Toxin A/B Rapid Test Cassette (Feces)
GCHSV (IgM)-402a	TEST RAPID	Orient Gene	HSV IgM Rapid test Cassette (whole blood/serum/plasma)
GCHSV(IgG)-302a	TEST RAPID	Orient Gene	HSV IgG Rapid Test Cassette (serum/plasma)
GCSYP-301a	TEST RAPID	Orient Gene	Syphilis Ab Rapid Test Strip (serum/plasma)
GCSYP-302a	TEST RAPID	Orient Gene	Syphilis Ab Rapid Test Cassette (serum/plasma)
GCSYP-401a	TEST RAPID	Orient Gene	Syphilis Ab Rapid test strip (whole blood/serum/plasma)
GCSYP-402a	TEST RAPID	Orient Gene	Syphilis Ab Rapid test cassette (whole blood/serum/plasma)
GBBAR-101a	TEST RAPID	Orient Gene	One Step Barbiturates Test Strip (Urine)
GBBAR-102a	TEST RAPID	Orient Gene	One Step Barbiturates Test Cassette (Urine)
GBAMP-101a	TEST RAPID	Orient Gene	One Step Amphetamine Test Strip (Urine)
GBAMP-102a	TEST RAPID	Orient Gene	One Step Amphetamine Test Cassette (Urine)
GBAMP-105a	TEST RAPID	Orient Gene	One Step Amphetamine Dip Card (Urine)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GBPPX-101a	TEST RAPID	Orient Gene	One Step Propoxyphene Test Strip (Urine)
GBPPX-102a	TEST RAPID	Orient Gene	One Step Propoxyphene Test Cassette (Urine)
GBDSA-XXXXFX	TEST RAPID	Orient Gene	Oral Fluid Drug test Cube
GBDSA-XXXXEX	TEST RAPID	Orient Gene	Oral Fluid Drug test
GBDSA-XXXXFSI	TEST RAPID	Orient Gene	Oral Fluid Drug test Cube
GBDSA-XXXXCX	TEST RAPID	Orient Gene	Oral Fluid Drug test cylinder
GBOPI-102a	TEST RAPID	Orient Gene	One Step Opiate Test Cassette (Urine)
GBOPI-101a	TEST RAPID	Orient Gene	One Step Opiate Test Strip (Urine)
GBETG-101b	TEST RAPID	Orient Gene	One Step Ethyl Glucoronide Test Strip (urine)
GBETG-102b	TEST RAPID	Orient Gene	One Step Ethyl Glucoronide Test Cassette (urine)
GBMOP-101a	TEST RAPID	Orient Gene	One step Morphine Test strip (urine)
GBMOP-102a	TEST RAPID	Orient Gene	One step Morphine Test Cassette (urine)
GBMOP-105a	TEST RAPID	Orient Gene	One step Morphine Test dip card (urine)
GBTHC-101a	TEST RAPID	Orient Gene	One Step Marijuana Test Strip (Urine)
GBTHC-102a	TEST RAPID	Orient Gene	One Step Marijuana Test Cassette (Urine)
GBTHC-105a	TEST RAPID	Orient Gene	One Step Marijuana Test Dip Card (Urine)
GBMTD-101a	TEST RAPID	Orient Gene	One step Methadone Test strip (urine)
GBMTD-102a	TEST RAPID	Orient Gene	One step Methadone Test cassette (urine)
GBXXX-101	TEST RAPID	Orient Gene	One Step Drugs of Abuse Test Strip (Urine)
GBXXX-102	TEST RAPID	Orient Gene	One Step Drugs of Abuse Test Cassette (Urine)
GBXXX-105	TEST RAPID	Orient Gene	One Step Drugs of Abuse Test Dip Card (Urine)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GBDSA-XXXXJSI	TEST RAPID	Orient Gene	Oral fluid drug test cylinder
GBDSA-XXXXJX	TEST RAPID	Orient Gene	Oral fluid drug test cylinder
GBDSA-XXXXKX	TEST RAPID	Orient Gene	Oral fluid drug test cylinder
GBDSA-XXXXMX	TEST RAPID	Orient Gene	Oral fluid drug test device
GBDSA-XXXXA/B/G/H/I	TEST RAPID	Orient Gene	Multi-drug rapid screen test cassette (oral fluid)
GBMTC-101a	TEST RAPID	Orient Gene	One Step Methcathinone Test Strip (Urine)
GBMTC-102a	TEST RAPID	Orient Gene	One Step Methcathinone Test Cassette (Urine)
GBKRA-101a	TEST RAPID	Orient Gene	One step kratom test strip (urine)
GBKRA-102a	TEST RAPID	Orient Gene	One step kratom test cassette (urine)
GBLSD-101a	TEST RAPID	Orient Gene	One Step Lysergic Acid Diethylamide Test Strip (Urine)
GBLSD-102a	TEST RAPID	Orient Gene	One Step Lysergic Acid Diethylamide Test Cassette (Urine)
FBXXX-1102	TEST RAPID	Orient Gene	Hair Multi-drug rapid test kit
GBETG-105a	TEST RAPID	Orient Gene	One step ethyl glucuronide test dip card (urine)
GBPGB-102b	TEST RAPID	Orient Gene	One step pregabalin test cassette (urine)
GBTRA-101a	TEST RAPID	Orient Gene	One step tramadol test strip (urine)
GBTRA-102a	TEST RAPID	Orient Gene	One step tramadol test cassette (urine)
GBOXY-101a	TEST RAPID	Orient Gene	One step oxycodone Test strip (urine)
GBOXY-102a	TEST RAPID	Orient Gene	One step oxycodone Test cassette (urine)
GBMDP-101a	TEST RAPID	Orient Gene	One step 3,4-Methylenedioxypropylvalerone Test strip (urine)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel: +86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.: 313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GBMDP-101a	TEST RAPID	Orient Gene	One step 3,4-Methylenedioxypropylone Test cassette (urine)
GBMQL-102a	TEST RAPID	Orient Gene	One step Methaqualone Test cassette (urine)
GBMQL-101a	TEST RAPID	Orient Gene	One step Methaqualone Test strip (urine)
GBMPD-101a	TEST RAPID	Orient Gene	One step Methylphenidate Test strip (urine)
GBMPD-102a	TEST RAPID	Orient Gene	One step Methylphenidate Test cassette (urine)
GBUR-101a	TEST RAPID	Orient Gene	One step UR-144 test strip (urine)
GBUR-102a	TEST RAPID	Orient Gene	One step UR-144 test cassette (urine)
GBBUP-101a	TEST RAPID	Orient Gene	One step buprenorphine test strip (urine)
GBBUP-102a	TEST RAPID	Orient Gene	One step buprenorphine test cassette (urine)
GBPCP-101a	TEST RAPID	Orient Gene	One step Phencyclidine Test strip (urine)
GBPCP-102a	TEST RAPID	Orient Gene	One step Phencyclidine Test cassette (urine)
GBTCA-101a	TEST RAPID	Orient Gene	One step Tricyclic Antidepressants test strip (urine)
GBTCA-102a	TEST RAPID	Orient Gene	One step Tricyclic Antidepressants test cassette (urine)
GBEDD-101a	TEST RAPID	Orient Gene	One step EDDP test strip (urine)
GBEDD-102a	TEST RAPID	Orient Gene	One step EDDP test cassette (urine)
GBFEN-101b	TEST RAPID	Orient Gene	One step Fentanyl Test strip (urine)
GBFEN-102b	TEST RAPID	Orient Gene	One step Fentanyl Test cassette (urine)
GBALC-101a	TEST RAPID	Orient Gene	Urine Alcohol Test Strip
GBMAM-S102	TEST RAPID	Orient Gene	One step 6-MAM Test cassette (urine)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GBMAM-S101	TEST RAPID	Orient Gene	One step 6-MAM Test cassette (urine)
GBHCD-101a	TEST RAPID	Orient Gene	One step Hydrocodone test strip (urine)
GBHCD-102a	TEST RAPID	Orient Gene	One step Hydrocodone test cassette (urine)
GBNFT-101c	TEST RAPID	Orient Gene	One step Norfentanyl test strip (urine)
GBNFT-102c	TEST RAPID	Orient Gene	One step Norfentanyl test cassette (urine)
GBXXX-1102	TEST RAPID	Orient Gene	Hair Multi-Drug Rapid Test Kit (ICA)
GBDSA-XXXXLX	TEST RAPID	Orient Gene	Oral Fluid Drug Test Mini Cube
GBDUA-1X4	TEST RAPID	Orient Gene	One Step Multi-Drug Screen Test Dip Card (urine)
GBDOA-1X5	TEST RAPID	Orient Gene	One Step Multi-Drug Screen Test Cassette (urine)
GBDUA-1X6	TEST RAPID	Orient Gene	One Step Multi-Drug Screen Test Cup (urine)
GBCOT-102a	TEST RAPID	Orient Gene	One step cotinine test cassette (urine)
GBK2-101a	TEST RAPID	Orient Gene	One step K2 Test strip (urine)
GBK2-102a	TEST RAPID	Orient Gene	One step K2 Test cassette (urine)
GBKET-101a	TEST RAPID	Orient Gene	One step Ketamine Test strip (urine)
GBKET-102a	TEST RAPID	Orient Gene	One step Ketamine Test cassette (urine)
GBBZO-101a	TEST RAPID	Orient Gene	One step Benzodiazepines Test Strip (urine)
GBBZO-102a	TEST RAPID	Orient Gene	One step Benzodiazepines Test Cassette (urine)
GBCOC-101a	TEST RAPID	Orient Gene	One step Cocaine Test strip (urine)
GBCOC-102a	TEST RAPID	Orient Gene	One step Cocaine Test cassette (urine)
GBCOC-105a	TEST RAPID	Orient Gene	One step Cocaine Test dip card (urine)
GBMDM-101a	TEST RAPID	Orient Gene	One step ecstasy Test strip (urine)

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax:+86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GBMDM-102a	TEST RAPID	Orient Gene	One step ecstasy Test cassette (urine)
GBMET-101a	TEST RAPID	Orient Gene	One step Methamphetamine test strip (urine)
GBMET-102a	TEST RAPID	Orient Gene	One step Methamphetamine test cassette (urine)
GBMET-105a	TEST RAPID	Orient Gene	One step Methamphetamine test dip card (urine)
GCTOXI(IgG/IgM)-302a	TEST RAPID	Orient Gene	Toxoplasma gondii test cassette (serum/plasma)
GCTOXI(IgG)-302a	TEST RAPID	Orient Gene	Toxoplasma gondii IgG test cassette (serum/plasma)
GCTOXI(IgM)-302a	TEST RAPID	Orient Gene	Toxoplasma gondii IgM test cassette (serum/plasma)
GCCHL-502a	TEST RAPID	Orient Gene	Chlamydia Trachomatis Antigen test cassette (swab/urine)
GEPSA-402a	TEST RAPID	Orient Gene	Prostate specific antigen test cassette (whole blood/serum/plasma)
GEPSA-401a	TEST RAPID	Orient Gene	Prostate specific antigen test strip (whole blood/serum/plasma)
GEPSA-302a	TEST RAPID	Orient Gene	Prostate specific antigen test cassette (serum/plasma)
GEPSA-301a	TEST RAPID	Orient Gene	Prostate specific antigen test strip (serum/plasma)
GALH-101a-1T	TEST RAPID	Orient Gene	LH Ovulation Test Strip
GALH-101a-5T	TEST RAPID	Orient Gene	LH Ovulation Test Strip
GALH-101a-7T	TEST RAPID	Orient Gene	LH Ovulation Test Strip
GALH-101b-1T	TEST RAPID	Orient Gene	LH Ovulation Test Strip
GALH-101b-5T	TEST RAPID	Orient Gene	LH Ovulation Test Strip
GALH-101b-7T	TEST RAPID	Orient Gene	LH Ovulation Test Strip
GALH-102a-1T	TEST RAPID	Orient Gene	LH Ovulation Test Cassette
GALH-102a-5T	TEST RAPID	Orient Gene	LH Ovulation Test Cassette
GALH-102a-7T	TEST RAPID	Orient Gene	LH Ovulation Test Cassette
GALH-102b-5T	TEST RAPID	Orient Gene	LH Ovulation Test Cassette
GALH-102b-1T	TEST RAPID	Orient Gene	LH Ovulation Test Cassette

地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号

Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China

电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

GALH-102b-7T	TEST RAPID	Orient Gene	LH Ovulation Test Cassette
VPH-502a	TEST RAPID	Orient Gene	Vaginal pH test cassette (Vaginal secretions)
URS-1T to 14T with various combination	STRIPURI DE URINA	Orient Gene	LEU/NIT/URO/MA/PRO/PH/BLO/S G/ASC/CRE/KET/BIL/GLU/CA
GCHCV-302a	TEST RAPID	Orient Gene	HCV Hepatitis C Virus Rapid Test (serum/plasma) cassette
GCHCV-402a	TEST RAPID	Orient Gene	HCV Hepatitis C Virus Rapid Test (whole blood/serum/plasma) cassette
GCHIV-302a	TEST RAPID	Orient gene	HIV 1/2 Human Immunodeficiency Virus (Serum/Plasma) cassette
GCHIV-402a	TEST RAPID	Orient gene	HIV 1/2 Human Immunodeficiency Virus (Whole blood/serum/plasma)cassette
GCHBsg-302a	TEST RAPID	Orient gene	HBsAg Hepatitis B Surface Antigen Rapid Test (Serum/Plasma)
GCHBsg-402a	TEST RAPID	Orient gene	HBsAg Hepatitis B Surface Antigen Rapid Test(Whole Blood/Serum/Plasma)

The end.



地址：浙江省湖州市安吉县递铺镇阳光大道东段 3787 号  
Add: 3787#, East Yangguang Avenue, Dipu Street Anji 313300, Huzhou, Zhejiang, China  
电话 Tel:+86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.:313300



# Certificate

No. Q5 092305 0001 Rev. 02

**Holder of Certificate:** **Zhejiang Orient Gene Biotech Co., Ltd.**  
3787#, East Yangguang Avenue, Dipu Street Anji  
313300 Huzhou, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA

**Certification Mark:**



**Scope of Certificate:** **Design and Development, Production and Distribution of In Vitro Diagnostic Reagent and Instrument for the Detection of Drugs of Abuse, Fertility, Infectious Diseases, Oncology, Biochemistry, Cardiac Diseases, Allergic Disease based on Rapid Test, PCR and Liquid Biochip Method.**

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: [www.tuvsud.com/ps-cert?q=cert:Q5 092305 0001 Rev. 02](http://www.tuvsud.com/ps-cert?q=cert:Q5 092305 0001 Rev. 02)

**Report No.:** SH2398804

**Valid from:** 2024-03-17  
**Valid until:** 2027-03-16

**Date,** 2024-03-01



Christoph Dicks  
Head of Certification/Notified Body





浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co., LTD



CE-DOC-OG029  
Version 4.0

# EC Declaration of Conformity

In accordance with Directive 98/79/EC

**Legal Manufacturer:** *Zhejiang Orient Gene Biotech Co., Ltd*

**Legal Manufacturer Address:** *3787#, East Yangguang Avenue, Dipu Street,  
Anji 313300, Huzhou, Zhejiang, China*

Declares, that the products  
Product Name and Model(s)

Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)	GCMAL(pf/pv)-402a
--	-------------------

**Classification:** *Other*  
**Conformity assessment route:** *Annex III (EC DECLARATION OF CONFORMITY)*

We, the Manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

We hereby explicitly appoint

**EC Representative's Name:** QARAD BV

**EC Representative's Address:** Ciplastraat 3, 2440 Geel, BELGIUM

to act as our European Authorized Representative as defined in the aforementioned Directive.

I, the undersigned, hereby declare that the medical devices specified above conform with the directive 98/79/EC on in vitro diagnostic medical devices and pertinent essential requirements

Date Signed: March 4, 2022

Name of authorized signatory: Joyce Pang  
Position held in the company: Vice-President



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co., LTD



CE-DOC-OG060  
Version 1.0

# EC Declaration of Conformity

In accordance with Directive 98/79/EC

**Legal Manufacturer:** *Zhejiang Orient Gene Biotech Co., Ltd*

**Legal Manufacturer Address:** *3787#, East Yangguang Avenue, Dipu Street,  
Anji 313300, Huzhou, Zhejiang, China*

Declares, that the products  
Product Name and Model(s)

Fecal Occult Blood Rapid Test Strip (Feces)	GEFOB-601b
Fecal Occult Blood Rapid Test Cassette (Feces)	GEFOB-602b

**Classification:** *Other*  
**Conformity assessment route:** *Annex III (EC DECLARATION OF CONFORMITY)*

We, the Manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

We hereby explicitly appoint

**EC Representative's Name:** Shanghai International Holding Corp. GmbH (Europe)

**EC Representative's Address:** Eiffestrasse 80, 20537 Hamburg, Germany

to act as our European Authorized Representative as defined in the aforementioned Directive.

I, the undersigned, hereby declare that the medical devices specified above conform with the directive 98/79/EC on in vitro diagnostic medical devices and pertinent essential requirements

Date Signed: November 28, 2017

Name of authorized signatory: Joyce Pang  
Position held in the company: Vice-President



# EC Declaration of Conformity

In accordance with Directive 98/79/EC

**Legal Manufacturer:** Healgen Scientific Limited Liability Company

**Legal Manufacturer Address:** 3818 Fuqua Street, Houston, TX 77047, USA.

Declares, that the products  
Product Name and Model(s)

Orient Gene HCV Hepatitis C Virus Rapid Test (Serum/Plasma) (Cassette)	GCHCV-302a
Orient Gene HCV Hepatitis C Virus Rapid Test (Whole blood/Serum/Plasma)(Cassette)	GCHCV-402a

EDMA Code: 15 70 02 02

**Classification:** Annex II List A  
**Conformity assessment route:** Annex IV (Full Quality Assurance)

Compliance of the designated product with the Directive 98/79/EC has been assessed and certified by the Notified Body

**Notified Body:** TÜV SÜD Product Service GmbH

**Notified Body Address:** Munich Branch Ridlerstraße 65 80339 München Germany

**EC Certificate No.:** V1 092378 0004 Rev. 02 Valid until: 2025-05-26

**EC Design-Examination Certificate No.:** V7 092378 0009 Rev. 00 Valid until: 2025-05-26

It bears the mark

**CE 0123**

We, the Manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.


We hereby explicitly appoint

**EC Representative Name:** QARAD b.v.b.a.

**EC Representative Address:** Cipalstraat 3, B-2440 Geel, Belgium

to act as our European Authorized Representative as defined in the aforementioned Directive.

I, the undersigned, hereby declare that the medical devices specified above conform with the directive 98/79/EC on in vitro diagnostic medical devices and pertinent essential requirements

Signature:   
Name of authorized signatory: Joyce Pang  
Position held in the company: Vice-President  
Date: 2022.4.22



Zhejiang Orient Gene Biotech Co., LTD

CERTIFICATE OF ANALYSIS

Product Name: HBsAg Rapid Test (Whole blood/Serum/Plasma) (Cassette)

Catalog NO.: GCHBsg-402a

Purchase NO.: 2025-IEU224#

Lot NO.: 2512141

Quantity:7000 pcs

Expiration Date: 2027.11

CONTROLS		SPECIFICATION	TEST RESULT	CONCLUSION
Negative Specimens		Negative	Negative	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
Positive Specimens	1ng/ml	Positive	Positive	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
	2ng/ml	Positive	Positive	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
	3ng/ml	Positive	Positive	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
	5ng/ml	Positive	Positive	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail

Conclusion: Pass: All results meet QC standard.  
Fail



Test by:

Date: 2025.12.06

QC Supervisor:

Date: 2025.12.06

# Hepatitis B Surface Antigen Rapid Test Cassette (Whole blood/Serum/Plasma)

## INTENDED USE

The Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) is a lateral flow chromatographic immunoassay for the qualitative detection of Hepatitis B surface antigen (HBsAg) in human whole blood, serum or plasma. It is intended to be used as a screening test and as an aid in the diagnosis of infection with Hepatitis B virus (HBV). Any reactive specimen with the HBsAg Rapid Test Cassette must be confirmed with alternative testing method(s) and clinical findings.

## INTRODUCTION

Viral hepatitis is a systemic disease primarily involving the liver. Most cases of acute viral hepatitis are caused by Hepatitis A virus, Hepatitis B virus (HBV) or Hepatitis C virus. The complex antigen found on the surface of HBV is called HBsAg. The presence of HBsAg in serum or plasma is an indication of an active Hepatitis B infection, either acute or chronic. In a typical Hepatitis B infection, HBsAg will be detected 2 to 4 weeks before the ALT level becomes abnormal and 3 to 5 weeks before symptoms or jaundice develop. HBsAg has four principal subtypes: adw, ayw, adr and ayr. Because of antigenic heterogeneity of the determinant, there are 10 major serotypes of Hepatitis B virus. The HBsAg Test Cassette (Whole Blood/Serum/Plasma) is a rapid test to qualitatively detect the presence of HBsAg in whole blood, serum or plasma specimens. The test utilizes a combination of double monoclonal antibodies to selectively detect elevated levels of HBsAg in whole blood, serum or plasma.

## PRINCIPLE

The Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) is a lateral flow chromatographic immunoassay based on the principle of the double antibody-sandwich technique. The membrane is pre-coated with anti-HBsAg antibodies on the test line region of the test. During testing, Hepatitis B Surface Antigen in the whole blood, serum or plasma specimen reacts with the particle coated with anti-HBsAg antibody. The mixture migrates upward on the membrane, chromatographically by capillary action to react with anti-HBsAg antibodies on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that the proper volume of specimen has been added and membrane wicking has occurred.

## PRODUCT CONTENTS

The Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) containing anti-HBsAg antibodies particles and anti-HBsAg antibodies coated on the membrane.

## MATERIALS SUPPLIED

Test cassette	Dropper	Buffer	Package insert
---------------	---------	--------	----------------

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Specimen collection containers
2. Lancets (for fingerstick whole blood only)
3. Centrifuge (for plasma only)
4. Timer
5. Heparinized capillary tubes and dispensing bulb (for fingerstick whole blood only)

## STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

## WARNINGS AND PRECAUTIONS

1. For professional In Vitro diagnostic use only. Do not use after expiration date.
2. Warning: the reagents in this kit contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.

3. Do not use it if the tube/pouch is damaged or broken.
4. Test is for single use only. Do not re-use under any circumstances.
5. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
6. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
7. Humidity and temperature can adversely affect results.
8. Do not perform the test in a room with strong air flow, i.e. electric fan or strong air-conditioning.

## SPECIMEN COLLECTION AND PREPARATION

1. Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.
2. To collect Fingerstick Whole Blood specimens:
  - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
  - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
  - Puncture the skin with a new sterile lancet for each person. Wipe away the first sign of blood.
  - Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
  - Add the Fingerstick Whole Blood specimen to the test cassette by using a capillary tube:
    - Touch the end of the capillary tube to the blood until filled to approximately 50µL. Avoid air bubbles.
    - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood into the specimen well (S) of the test cassette.
  - Add the Fingerstick Whole Blood specimen to the test cassette by using hanging drops:
    - Position the patient's finger so that the drop of blood is just above the specimen well (S) of the test cassette.
    - Allow 2 hanging drops of fingerstick whole blood to fall into the center of specimen well (S) on the test cassette, or move the patient's finger so that the hanging drop touches the center of the specimen well (S). Avoid touching the finger directly to the specimen well (S).
3. Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
4. Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
5. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
6. If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

## TEST PROCEDURE

**Allow test cassette, specimen, buffer and/or controls to equilibrate to room temperature (15-30°C) prior to testing.**


1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface.

**For Serum or Plasma Specimens:** Hold the dropper vertically and transfer 2-3 drops of serum or plasma (approximately 60-90µL) to the specimen well (S) of the test cassette. See illustration below.

**For Venipuncture Whole Blood Specimens:** Hold the dropper vertically and transfer 2 drops of venipuncture whole blood (approximately 50µL) to the specimen well (S) of the test cassette, then add 1 drop of buffer (approximately 40 µL) and start the timer. See illustration below.

**For Fingerstick Whole Blood Specimens:** Allow 2 hanging drops of fingerstick whole blood (approximately 50 µL) to fall into the center of the specimen well (S) on the test cassette, then add 1 drop of buffer (approximately 40 µL) and start the timer. See illustration below.

3. Wait for the red line(s) to appear. The result should be read in 15 minutes. Do not interpret the result after 15 minutes



### INTERPRETATION OF RESULTS



**Positive Negative Invalid**

**POSITIVE:** Two distinct red lines appear. One line should be in the control region (C) and another line should be in the test region (T).

**NEGATIVE:** One red line appears in the control region (C). No apparent red or pink line appears in the test region (T).

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test Cassette. If the problem persists discontinue using the test kit immediately and contact your local distributor

### QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

### LIMITATIONS

1. Though the Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) is a reliable screening assay, it should not be used as a sole criterion for diagnosis of HBV infection.
2. The HBsAg Rapid Test Cassette is limited to the qualitative detection of HBsAg in human whole blood, serum or plasma. The intensity of the test band does not have linear correlation with HBsAg titer in the specimen.
3. A negative test result does not preclude the possibility of exposure to or infection with HBV. Infection through recent exposure (seroconversion) to HBV may not be detectable.
4. A negative result can occur if the quantity of HBsAg present in the specimen is below the detection limits of the assay (lower than 1 ng/mL), or the HBsAg that are detected are not present during the stage of disease in which a sample is collected.
5. Interference due to heterophile antibodies, Rheumatoid Factors and other nonanalyte substances in patient's serum, capable of binding antibodies multivalently and providing erroneous analyte detection in immunoassays, has been reported in various studies. Both laboratory professionals and clinicians must be vigilant to this possibility of antibody interference. Results that appear to be internally inconsistent or incompatible with the clinical presentation should invoke suspicion of the presence of an endogenous artifact and lead to appropriate in vitro investigative action.
6. This kit is intended ONLY for testing of individual samples. Do not use it for testing of cadaver samples, saliva, urine or other body fluids, or pooled (mixed) blood.
7. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

### PERFORMANCE CHARACTERISTICS

#### Sensitivity:

The Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested with a sensitivity panel ranging from 0 to 300 ng/mL. All 10 HBsAg subtypes produced positive results on the Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma). The test can detect 5ng/mL of HBsAg in 10 minutes, and 1 ng/mL of HBsAg in 15 minutes.

#### Specificity:

Antibodies used for the Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) were developed against whole Hepatitis B antigen isolated from Hepatitis B virus. Specificity of the Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) was also tested with laboratory strains of Hepatitis A and Hepatitis C. They all yielded negative results.

#### Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma) vs. EIA test

Method	EIA		Total Results	
	Results	Positive		Negative
Hepatitis B Surface Antigen Rapid Test Cassette (Whole Blood/Serum/Plasma)	Positive	345	5	350
	Negative	2	980	982
Total Results		347	985	1332

Relative sensitivity: 99.4%

Relative specificity 99.5%

Accuracy: 99.5%

### REFERENCE

1. Blumberg, B. S. The Discovery of Australian Antigen and its relation to viral hepatitis. *Vitro*. 1971; 7: 223

Catalogue No:GCHBsg-402a

Effective Date: 2023-08-22

B20137-05

# Fecal Occult Blood Rapid Test Cassette (Feces)



## INTENDED USE

Fecal Occult Blood Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of human occult blood in feces by professional laboratories or physician's offices. It is useful to detect bleeding caused by a number of gastrointestinal disorders, e.g., diverticulitis, colitis, polyps, and colorectal cancer.

Fecal Occult Blood Rapid Test Cassette (Feces) is recommended for use in 1) routine physical examinations, 2) hospital monitoring for bleeding in patients, and 3) screening for colorectal cancer or gastrointestinal bleeding from any source.

## INTRODUCTION

Most of diseases can cause hidden blood in the stool. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood. Traditional guaiac-based method lacks sensitivity and specificity, and has diet-restriction prior to the testing.

Fecal Occult Blood Rapid Test Cassette (Feces) is a rapid test to qualitatively detect low levels of fecal occult blood in feces. The test uses double antibody-sandwich assay to selectively detect as low as 50 ng/mL of hemoglobin or 6 µg hemoglobin/g feces. In addition, unlike the guaiac assays, the accuracy of the test is not affected by the diet of the patients.

## PRINCIPLE

Fecal Occult Blood Rapid Test Cassette (Feces) is a lateral flow chromatographic immunoassay based on the principle of the double antibody-sandwich technique. The membrane is pre-coated with anti-hemoglobin antibodies on the test line region of the device. During testing, the specimen reacts with the colloidal gold coated with anti-hemoglobin antibodies. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-hemoglobin antibodies on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

## MATERIALS PROVIDED

- 20 Test cassettes
- 20 Specimen collection tubes with buffer
- 1 Package insert

## MATERIALS REQUIRED BUT NOT PROVIDED

1. Specimen collection containers
2. Clock or timer

## STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened at 2°C-30°C. If stored at 2°C-8°C, ensure that the test device is brought to room temperature before opening. The test is not stable out of the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

## PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
3. Do not use it if the tube/pouch is damaged or broken.
4. Test is for single use only. Do not re-use under any circumstances.
5. **Do not use specimen with visible blood for the testing.**
6. Handle all specimens as if they contain infectious agents. Observe established standard procedure for proper disposal of specimens.
7. Specimen extraction buffer contains Sodium Azide (0.1%). Avoid contact with skin or eyes. Do not ingest.
8. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
9. Humidity and temperature can adversely affect results.
10. Do not perform the test in a room with strong air flow, ie. electric fan or strong air conditioning.

## PATIENT PREPARATION

1. A specimen should not be collected from a patient with following conditions that may interfere with the test results:

- Menstrual bleeding
  - Bleeding hemorrhoids
  - Constipating bleeding
  - Urinary bleeding.
2. Dietary restrictions are not necessary.
  3. Alcohol and certain medications such as aspirin, indomethacin, phenylbutazone, reserpine, cortocosteroids, and nonsteroidal anti-inflammatory drugs may cause gastrointestinal irritation and subsequent bleeding, thus gives positive reactions. On the advice of the physician, such substances should be discontinued at least 48 hours prior to testing.

## SPECIMEN COLLECTION AND PREPARATION

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.


1. Collect a random sample of feces in a clean, dry receptacle.
2. Unscrew the top of the collection tube and remove the applicator stick.
3. Randomly pierce the fecal specimen in at least five (5) different sites.
4. Remove excess sample off the shaft and outer grooves. Be sure sample remains on inside grooves.
5. Replace the stick in the tube and tighten securely.
6. Shake the specimen collection bottle so that there is proper homogenisation of feces in buffer solution.

**Note:** Specimens prepared in the specimen collection tube may be stored at room temperature (15-30°C) for 3 days maximum, at 2-8°C for 7 days maximum or at -20°C for 3 months maximum if not tested within 1 hour after preparation.

## TEST PROCEDURE

**Allow the test cassette, specimen, and/or controls to reach room temperature (15-30°C) prior to testing.**

1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean, flat surface.
3. Shake the specimen collection tube several times.
4. Hold the specimen collection tube upright and then unscrew and open the upper cap.
5. Squeeze 3 drops (~90 µL) of the sample solution in the sample well of the cassette and start the timer.
6. Wait for the colored line(s) to appear. Read results in 5 minutes. Do not interpret the result after 5 minutes.



## INTERPRETATION OF RESULTS

(Please refer to the illustration above)

**Positive:** Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

**Negative:** One colored line appears in the control line region (C). No line appears in the test line region (T).

**Invalid:** Control line fails to appear. The test should be repeated using a new cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

**NOTE:**

1. The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test region should be considered positive. Note that this is a qualitative test only, and

# Fecal Occult Blood Rapid Test Cassette (Feces)

cannot determine the concentration of analytes in the specimen.

2. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

## QUALITY CONTROL

An internal procedural control is included in the test. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

## LIMITATIONS

1. This test kit is to be used for the qualitative detection of human hemoglobin in fecal samples. A positive result suggests the presence of human hemoglobin in fecal samples. In addition to intestinal bleeding the presence of blood in stools may have other causes such as hemorrhoids, blood in urine etc.
2. Not all colorectal bleedings are due to precancerous or cancerous polyps. The information obtained by this test should be used in conjunction with other clinical findings and testing methods, such as colonoscopy gathered by the physician.
3. Negative results do not exclude bleeding since some polyps and colorectal region cancers can bleed intermittently or not at all. Additionally, blood may not be uniformly distributed in fecal samples. Colorectal polyps at an early stage may not bleed.
4. Urine and excessive dilution of sample with water from toilet bowl may cause erroneous test results. The use of a receptacle is recommended.
5. Feces specimens should not collect during the menstrual period and not three day before or afterwards, at bleeding due to constipation, bleeding haemorrhoids, or at taking rectally administered medication. It could cause false positive results.
6. This test may be less sensitive for detecting upper g.i. Bleeding because blood degrades as it passes through the g.i. Track.
7. The Fecal Occult Blood Rapid Test Cassette (Feces) is to aid in diagnosis and is not intended to replace other diagnostic procedures such as G.I. fibroscope, endoscopy, colonoscopy, or X-ray analysis. Test results should not be deemed conclusive with respect to the presence or absence of gastrointestinal bleeding or pathology. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source for the occult blood in the feces.

## PERFORMANCE CHARACTERISTICS

### 1. Sensitivity: 99.6%

Fecal Occult Blood Rapid Test Cassette (Feces) can detect the levels of human occult blood as low as 50 ng/mL hemoglobin or 6 µg hemoglobin/g feces.

### 2. Prozone Effect:

It is observed that this FOB test can detect 2 mg/mL hemoglobin.

### 3. Specificity: 99.9%

Fecal Occult Blood Rapid Test Cassette (Feces) is specific to human hemoglobin. Specimen containing the following substances at the standard concentration was tested on both positive and negative controls and showed no effects on test results at standard concentration.

Substances	Concentrations (Diluted with the extraction buffer)
Beef hemoglobin	2 mg/mL
Chicken hemoglobin	0.5 mg/mL
Pig hemoglobin	0.5 mg/mL
Goat hemoglobin	0.5 mg/mL
Horse hemoglobin	20 mg/mL
Rabbit hemoglobin	0.06 mg/mL

## REFERENCES

1. Simon J.B. Occult Blood Screening for Colorectal Carcinoma: A Critical Review, Gastroenterology, Vol. 1985;88:820.
2. Blebae J. and Napherson RA. False-Positive Guaiac Testing With Iodine, Arch Pathol Lab Med, 1985;109:437-40.

## INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

 Zhejiang Orient Gene Biotech Co., Ltd  
Address: 3787#, East Yangguang Avenue, Dipu Street,  
Anji 313300, Huzhou, Zhejiang, China  
Tel: +86-572-5226111 Fax: +86-572-5226222  
Website: www.orientgene.com

 Shanghai International Holding Corp. GmbH (Europe)  
Add: Eiffestrasse 80, 20537 Hamburg, Germany

 GEFOB-602b

# HCV Ab Rapid Test Cassette (Whole Blood/Serum/Plasma)

## INTENDED USE

The HCV Ab Rapid Test Cassette (Whole Blood/Serum/Plasma) is a sandwich lateral flow chromatographic immunoassay for the qualitative detection of antibodies (IgG, IgM, and IgA) anti- Hepatitis C virus (HCV) in human whole blood, serum or plasma. It is intended to be used as a screening test and as an aid in the diagnosis of infection with HCV. Any reactive specimen with the HCV Ab Rapid Cassette must be confirmed with alternative testing method(s) and clinical findings.

## INTRODUCTION

Hepatitis C Virus (HCV) is a small, enveloped, positive-sense, single-stranded RNA Virus. Antibody to HCV is found in over 80% of patients with well-documented non-A, non-B hepatitis. Conventional methods fail to isolate the virus in cell culture or visualize it by electron microscope. Cloning the viral genome has made it possible to develop serologic assays that use recombinant antigens<sup>(1, 2)</sup>. Compared to the first generation HCV EIAs using single recombinant antigen, multiple antigens using recombinant protein and/or synthetic peptides have been added in new serologic tests to avoid nonspecific cross-reactivity and to increase the sensitivity of the HCV antibody tests<sup>(3, 4)</sup>.

HCV Ab Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid test to qualitatively detect the presence of antibody to HCV in a whole blood, serum or plasma specimen. The test utilizes a combination of recombinant antigen to selectively detect elevated levels of HCV antibodies in whole blood, serum or plasma.

## PRINCIPLE

The HCV Ab Rapid Test Cassette is a lateral flow chromatographic immunoassay based on the principle of the double antigen-sandwich technique. The test cassette consists of: 1) a burgundy colored conjugate pad containing HCV antigens conjugated with colloidal gold (HCV Ag conjugates) and rabbit IgG-gold conjugates, 2) a nitrocellulose membrane cassette containing a test band (T band) and a control band (C band). The T band is pre-coated with non-conjugated HCV antigens, and the C band is pre-coated with goat anti-rabbit IgG. When an adequate volume of test specimen is dispensed into the sample well of the cassette, the specimen migrates by capillary action across the cassette. The antibodies: either the IgG, the IgM, or the IgA, to HCV if present in the specimen will bind to the HCV Ag conjugates. The immunocomplex is then captured on the membrane by the pre-coated HCV antigens, forming a burgundy colored T band, indicating a HCV Ab positive test result. Absence of the T band suggests a negative result. The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex of goat anti-rabbit IgG/rabbit IgG-gold conjugate regardless the presence of any antibodies to HCV. Otherwise, the test result is invalid and the specimen must be retested with another device.

## PRODUCT CONTENTS

HCV Ab Rapid Test Cassette(Whole Blood/Serum/Plasma) containing HCV antigen coated particles and HCV antigen coated on the membrane.

## MATERIALS SUPPLIED

1. Test Cassette (25 sealed pouches) 2. Disposable pipette for each test 3.Desiccant 4. 1 Buffer (4 ml) 5.Package Insert 6. Sterile disposable lancet for each test 7. Alcohol swab for each test

## MATERIAL REQUIRED BUT NOT PROVIDED

1.Specimen collection containers 2.Centrifuge (for plasma only) 3.Timer

## STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

## WARNINGS AND PRECAUTIONS

- For professional In Vitro diagnostic use only. Do not use after expiration date.
- Warning: the reagents in this kit contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Do not use it if the tube/pouch is damaged or broken.

- Test is for single use only. Do not re-use under any circumstances.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Humidity and temperature can adversely affect results.

## SPECIMEN COLLECTION

1.The HCV Rapid Test Cassette (Whole Blood/Serum/Plasma) can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

2.To collect Fingerstick Whole Blood specimens:

Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.

Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger. Puncture the skin with a new sterile lancet for each person. Wipe away the first sign of blood. Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site. Add the Fingerstick Whole Blood specimen to the test cassette by using a capillary tube:

Touch the end of the capillary tube to the blood until filled to approximately 50 µL. Avoid air bubbles. Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood into the sample pad of the test cassette.

Add the Fingerstick Whole Blood specimen to the test cassette by using hanging drops:

Position the patient's finger so that the drop of blood is just above the sample pad of the test cassette.

Allow 2 hanging drops of fingerstick whole blood to fall into the center of sample pad of the test cassette or, move the patient's finger so that the hanging drop touches the center of the sample pad. Avoid touching the finger directly to the sample pad.

3.Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.

4.Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.

5.Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

6.If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

## TEST PROCEDURE

**Allow test cassette, specimen, buffer and/or controls to equilibrate to room temperature (15-30°C) prior to testing.** 1.Remove the test device from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.

2. Place the test device on a clean and level surface.


**For Serum or Plasma Specimens:** Hold the dropper vertically and transfer 1 drop of serum or plasma (approximately 30 µL) to the sample pad of the test cassette, then add 1 drop of buffer (approximately 30 µL) and start the timer. See illustration below.

**For Venipuncture Whole Blood Specimens:** Hold the dropper vertically and transfer 2 drops of venipuncture whole blood (approximately 50µL) to the sample pad of the cassette, then add 1 drop of buffer (approximately 30 µL) and start the timer. See illustration below.

**For Fingerstick Whole Blood Specimens:** Allow 2 hanging drops of fingerstick whole blood (approximately 50 µL) to fall into the center of the sample pad on the test cassette, then add 1 drop of buffer (approximately 30 µL) and start the timer. See illustration below.

For venipuncture whole blood and plasma: K2EDTA, Sodium Heparin, Sodium citrate Sterile, and Lithium heparin should be used as the anticoagulant. Other anticoagulants have not been tested and may give incorrect results.

3. Wait for the red line(s) to appear. The result should be read in 15 minutes. Do not interpret the result after 15 minutes.



### INTERPRETATION OF RESULTS

(please refer to the illustration above)

**Positive:** Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

**Negative:** One colored line appears in the control line region (C). No line appears in the test line region (T).

**Invalid:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

### QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this test. However, it is recommended that positive and negative controls are sourced from a local competent authority and tested as a good laboratory practice, to confirm the test procedure and verify the test performance.

### LIMITATIONS

1. The HCV Ab Rapid Test Cassette (Whole Blood/ Serum/Plasma) is for in vitro diagnostic use only. This test should be used for the detection of antibodies to HCV in whole blood, serum or plasma specimen.
2. The HCV Ab Rapid Test Cassette (Whole Blood/Serum/Plasma) will only indicate the presence of antibodies to HCV in the specimen and should not be used as the sole criteria for the diagnosis of Hepatitis C viral infection.
3. As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
4. If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is recommended. A negative result at any time does not preclude the possibility of Hepatitis C Virus infection.
5. A negative result can occur if the quantity of the antibodies to HCV present in the specimen is below the detection limits of the assay, or the antibodies that are detected are not present during the stage of disease in which a sample is collected.
6. Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.

### PERFORMANCE CHARACTERISTICS

**Sensitivity:** HCV Ab Rapid Test Cassette (Whole Blood/ Serum/Plasma) has passed a seroconversion panel and compared with leading commercial HCV EIA test using clinical specimens.

**Specificity:** The recombinant antigens used for HCV Ab Rapid Test Cassette (Whole Blood/Serum/Plasma) are encoded by genes for both structural (nucleocapsid) and non-structural proteins. HCV Ab Rapid Test Cassette (Whole Blood/Serum/Plasma) is highly specific for antibodies to Hepatitis C Virus compared with a leading commercial HCV EIA test.

#### The HCV Ab Rapid Test Cassette vs.EIA test

Method	EIA		Total Results
	Positive	Negative	
HCV Ab Rapid Test	105	19	124
	2	1760	1762
Total Results	107	1779	1886

Relative sensitivity: 99.9%

Relative specificity: 99.9%

### REFERENCE

1. Choo, Q.L., G.Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 189;244:359
  2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An assay for circulating antibodies to a major etiologic Virus of human non-A, non-B hepatitis. *Science* 1989; 244:362.
  3. Van der Poel, C.L., H.T.M. Cuyper, H.W. Reesink, and P.N. Lelie. Confirmation of hepatitis C Virus infection by new four- antigen recombinant immunoblot assay. *Lancet* 1991;337:317
- Wilber, J.C. Development and use of laboratory tests for hepatitis C infection: a review. *J. Clin. Immunoassay* 1993;16:204

# Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)



## INTENDED USE

The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) is a rapid lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of Malaria P.falciparum specific histidine rich protein-2 (Pf HRP-II) and Malaria P.vivax specific lactate dehydrogenase (Pv-LDH) in human blood specimen as an aid in the diagnosis of Malaria infection. It is for *In-Vitro* Diagnostic use only.

## INTRODUCTION

Malaria is a serious, sometimes fatal, parasitic disease characterized by fever, chills, and anemia and is caused by a parasite that is transmitted from one human to another by the bite of infected Anopheles mosquitoes. There are four kinds of malaria that can infect humans: Plasmodium falciparum, P. vivax, P. ovale, and P. malariae. In humans, the parasites (called sporozoites) migrate to the liver where they mature and release another form, the merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year. At the present, malaria is diagnosed by looking for the parasites in a drop of blood. Blood will be put onto a microscope slide and stained so that the parasites will be visible under a microscope.

## PRINCIPLE

The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) contains a membrane, which is precoated with mouse monoclonal antibodies specific to HRP-II of P. falciparum on test line Pf region and with mouse monoclonal antibodies specific to lactate dehydrogenase of P.vivax species on test line Pv region respectively. Conjugate pad is dispensed with monoclonal antibodies conjugated to colloidal gold, which are specific to P.falciparum histidine rich protein-2 (Pf HRP-II) and specific to the lactate dehydrogenase of P.vivax.

During the assay, an adequate volume of the blood specimen is dispensed into the sample well (S) of the test cassette, a lysis buffer is added to the buffer well (B). The buffer contains a detergent that lyses the red blood cells and releases various antigens, which migrate by capillary action across the strip held in the cassette. Pv-LDH if presents in the specimen will bind to the Pv-LDH-gold conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-Pv-LDH antibody, forming a burgundy colored Pv band, indicating a Pv positive test result.

Alternatively, pHRP-II if presents in the specimen will bind to the pHRP-II-gold conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-pHRP-II antibodies, forming a burgundy colored Pf band, indicating a Pf positive test result.

Absence of any T bands suggests a negative result. The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex of goat anti- mouse IgG I mouse IgG (anti-Pv-LDH and anti-pHRP-II)-gold conjugates regardless of the color development on any of the T bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

## MATERIALS SUPPLIED

25 Sealed pouches each containing a test cassette, a dropper and a desiccant  
1 Buffer, 7.0 mL  
1 Package insert

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Clock or timer
2. Collection by venipuncture: collection tube (containing EDTA, citrate or heparin)
3. Collection using a lancet: sterile lancet

## STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened, preferably at 2°C-30°C. Do not expose the kit over 30°C. Do not freeze the kit. Ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch if it is stored at 2°C-30°C.

## WARNINGS AND PRECAUTIONS

1. For professional *in vitro* diagnostic use only. Do not use after expiration date.
2. The instruction must be followed exactly to get accurate results. Failure to follow the insert gives inaccurate test results.
3. Do not eat, drink or smoke in the area where the specimens or kits are handled.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.

5. Hemolized blood may be used for the testing, but do not take precipitants.
6. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
7. Humidity and temperature can adversely affect results.
8. Do not perform the test in a room with strong air flow, ie. an electric fan or strong airconditioning.

## SPECIMEN COLLECTION

### Collection by venipuncture:

- 1) Collect whole blood into a collection tube (containing EDTA, citrate or heparin) by venipuncture.
- 2) If specimens are not immediately tested, they should be refrigerated at 2-8°C. For storage periods greater than three days, freezing is recommended. They should be brought to room temperature prior to use. Using the specimen after long-term storage of more than three days can cause non-specific reaction.
- 3) When stored at 2-8°C, the whole blood sample should be used within three days.

### Collection using a lancet:

- 1) Clean the area to be lanced with an alcohol swab.
- 2) Squeeze the end of the fingertip and pierce with a sterile lancet.
- 3) Wipe away the first drop of blood with sterile gauze or cotton.
- 4) Using the dropper provided, while gently squeezing the tube, immerse the open end in the blood drop and then gently release the pressure to draw blood into the dropper.


## TEST PROCEDURE

Allow the test device, specimen, buffer, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface. Be sure to label the device with specimen's ID number.
3. With a 5 µL mini plastic dropper provided, draw whole blood specimen to exceed the specimen line as showed in the following image and then transfer drawn whole blood into the sample well (S). Then add 3 drops (about 120 µL) of Lysis Buffer to the buffer well (B) immediately.

*Note: Practice a few times prior to testing if you are not familiar with the mini dropper. For better precision, transfer specimen by pipette capable to deliver 5 µL of volume.*

4. Set up timer.
- If preferred, after 5 minutes of adding specimen and buffer, you may add one more drop of Lysis Buffer to help the background become clearer.
5. Results can be read in 20 to 30 minutes. It may take more than 20 minutes to have the background become clearer. Don't read results after 30 minutes. To avoid confusion, discard the test cassette after interpreting the result.



# Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)

## INTERPRETATION OF RESULTS

(Please refer to the illustration above)

### POSITIVE:

**P.f Positive:** One line appears in the control region, and one line appears in P.f. line region.

**P.v Positive:** One line appears in the control region and one line appears in Pv line region.

**P.f and P.v Positive:** One line appears in the control region, one line appears in Pv line region and one line appears in P.f. line region.

**NEGATIVE:** Only one colored line appears in the control region.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

## QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

## LIMITATIONS

- The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) is for *in vitro* diagnostic use only. This test should be used for the detection of P.f and P.v antigens in whole blood specimens only. Neither the quantitative value nor the rate of increase in P.f and P.v concentration can be determined by this qualitative test.
- The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) will only indicate the presence of antigens of P.f and / or P.v in the specimen and should not be used as the sole criterion for the diagnosis of malaria infection.
- As known relevant interference, haemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples can lead to impair the test results.
- The test is limited to the detection of antigen to Malaria Plasmodium sp. Although the test is very accurate in detecting HRP-II specific to P.f or pLDH specific to P.v, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of malaria infection.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

## PERFORMANCE CHARACTERISTICS

### 1. Clinical Performance for P.f Ag test:

A total of 352 samples from susceptible subjects were tested by the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) and by thick blood smear test.

Method		Smear Test		Total Results
Malaria Pf/Pv Ag Rapid Test	Results	Positive	Negative	
	Positive	50	4	54
	Negative	0	298	298
Total Results		50	302	352

Relative Sensitivity: 100%

Relative Specificity: 98.7%

Overall Agreement: 98.9%

### 2. Clinical Performance for P.v Ag test:

A total of 289 samples from susceptible subjects were tested by the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) and by thick blood smear test.

Method		Smear Test		Total Results
Malaria Pf/Pv Ag Rapid Test	Results	Positive	Negative	
	Positive	63	3	66
	Negative	0	223	223
Total Results		63	226	289

Relative Sensitivity: 100%

Relative Specificity: 98.7%

Overall Agreement: 99.0%

**3. Precision:** Within-run and between-run have been determined by the testing 10 replicates of four specimens: a negative, a low positive, a medium positive and a strong positive. All values were correctly identified 100% of the time.

**4. Interference:** To evaluate the interference of Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) with known relevant interfering specimens, the haemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples were investigated. In these studies, those specimens did not interfere with the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood).

## REFERENCE

- Leonard K. Basco, Frederique Marquet, Michael M. Makler, and Jacques Le Bras.: Plasmodium falciparum and Plasmodium vivax: Lactate Dehydrogenase Activity and its Application for *in vitro* Drug Susceptibility Assay. Experimental Parasitology 80, 260-271 (1995)
- David L. Vander Jagt, Lucy A. Hunsaker and John E. Heidrich : Partial Purification and Characterization of Lactate Dehydrogenase from Plasmodium falciparum. Molecular and Biochemical Parasitology, 4 (1981) 255-264
- David J. Bzik, Barbara A. Fox and Kenneth Gonyer : Expression of Plasmodium falciparum lactate dehydrogenase in Escherichia coli Molecular and Biochemical Parasitology, 59(1993) 155-166
- Histidine-Rich Protein II: a Novel Approach to Malaria Drug Sensitivity Testing ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, June 2002, p. 1658©1664 Vol. 46, No. 6

## INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

Zhejiang Orient Gene Biotech Co., Ltd  
Address: 3787#, East Yangguang Avenue, Dipu Street,  
Anji 313300, Huzhou, Zhejiang, China  
Tel: +86-572-5226111 Fax: +86-572-5226222  
Website: www.orientgene.com

QARAD BV  
Cipalstraat 3, 2440 Geel BELGIUM

GCMAL(pf/pv)-402a