

ePTFE (PTFE EXPANDIDO) PRÓTESIS VASCULARES

Aplicación
quirúrgica

Acceso Vascular

Recto



SEALPTFE™
Wrap

injerto de pared estándar

Cónico



taperflo™

Soporte
central



taperflo™
con soporte central

Resistencia
superior al
aplastamiento



SEALPTFE™
Ultrathin

con soporte PTFE sólido

Soporte



SEALPTFE™
Ultrathin

de pared ultradelgada con
construcción **unity™**

Sin soporte



SEALPTFE™
Ultrathin

Periférico

Soporte



COBRAHOOD™
con soporte PTFE sólido

Flujo

hemodinámico


Sin soporte



COBRAHOOD™

SEALPTFE™ Ultrathin


Pared ultradelgada (0,41 mm)
ePTFE sellado con gelatina



Longitud	Diámetro interno	Código de producto	Código de producto	Código de producto
10 cm	5 mm	T1005S	T1005ES	T1005PS
	6 mm	T1006S	T1006ES	T1006PS
	7 mm	T1007S	T1007ES	T1007PS
	8 mm	T1008S	T1008ES	T1008PS
	10 mm	T1010S	T1010ES	T1010PS
20 cm	5 mm	T2005S	T2005ES	T2005PS
	6 mm	T2006S	T2006ES	T2006PS
	7 mm	T2007S	T2007ES	T2007PS
	8 mm	T2008S	T2008ES	T2008PS
	10 mm	T2010S	T2010ES	T2010PS
30 cm	5 mm	T3005S	T3005ES	T3005PS
	6 mm	T3006S	T3006ES	T3006PS
	7 mm	T3007S	T3007ES	T3007PS
	8 mm	T3008S	T3008ES	T3008PS
	10 mm	T3010S	T3010ES	T3010PS
40 cm	5 mm	T4005S	T4005ES	T4005PS
	6 mm	T4006S	T4006ES	T4006PS
	7 mm	T4007S	T4007ES	T4007PS
	8 mm	T4008S	T4008ES	T4008PS
	10 mm	T4010S	T4010ES	T4010PS
50 cm	5 mm	T5005S	T5005ES	T5005PS
	6 mm	T5006S	T5006ES	T5006PS
	7 mm	T5007S	T5007ES	T5007PS
	8 mm	T5008S	T5008ES	T5008PS
	10 mm	T5010S	T5010ES	T5010PS
60 cm	5 mm	T6005S	T6005ES	T6005PS
	6 mm	T6006S	T6006ES	T6006PS
	7 mm	T6007S	T6007ES	T6007PS
	8 mm	T6008S	T6008ES	T6008PS
	10 mm	T6010S	T6010ES	T6010PS
70 cm	5 mm	T7005S	T7005ES	T7005PS
	6 mm	T7006S	T7006ES	T7006PS
	7 mm	T7007S	T7007ES	T7007PS
	8 mm	T7008S	T7008ES	T7008PS
	10 mm	T7010S	T7010ES	T7010PS
80 cm	5 mm	T8005S	T8005ES	T8005PS
	6 mm	T8006S	T8006ES	T8006PS
	7 mm	T8007S	T8007ES	T8007PS
	8 mm	T8008S	T8008ES	T8008PS
	10 mm	T8010S	T8010ES	T8010PS

SEALPTFE™ Wrap

Pared estándar (0,59 mm)
ePTFE sellado con gelatina



Longitud	Diámetro interno	Código de producto	Código de producto	Código de producto
10 cm	5 mm	S1005S	S1005ES	S1005PS
	6 mm	S1006S	S1006ES	S1006PS
	7 mm	S1007S	S1007ES	S1007PS
	8 mm	S1008S	S1008ES	S1008PS
	10 mm	S1010S	S1010ES	S1010PS
20 cm	5 mm	S2005S	S2005ES	S2005PS
	6 mm	S2006S	S2006ES	S2006PS
	7 mm	S2007S	S2007ES	S2007PS
	8 mm	S2008S	S2008ES	S2008PS
	10 mm	S2010S	S2010ES	S2010PS
30 cm	5 mm	S3005S	S3005ES	S3005PS
	6 mm	S3006S	S3006ES	S3006PS
	7 mm	S3007S	S3007ES	S3007PS
	8 mm	S3008S	S3008ES	S3008PS
	10 mm	S3010S	S3010ES	S3010PS
40 cm	5 mm	S4005S	S4005ES	S4005PS
	6 mm	S4006S	S4006ES	S4006PS
	7 mm	S4007S	S4007ES	S4007PS
	8 mm	S4008S	S4008ES	S4008PS
	10 mm	S4010S	S4010ES	S4010PS
50 cm	5 mm	S5005S	S5005ES	S5005PS
	6 mm	S5006S	S5006ES	S5006PS
	7 mm	S5007S	S5007ES	S5007PS
	8 mm	S5008S	S5008ES	S5008PS
	10 mm	S5010S	S5010ES	S5010PS
60 cm	5 mm	S6005S	S6005ES	S6005PS
	6 mm	S6006S	S6006ES	S6006PS
	7 mm	S6007S	S6007ES	S6007PS
	8 mm	S6008S	S6008ES	S6008PS
	10 mm	S6010S	S6010ES	S6010PS
70 cm	5 mm	S7005S	S7005ES	S7005PS
	6 mm	S7006S	S7006ES	S7006PS
	7 mm	S7007S	S7007ES	S7007PS
	8 mm	S7008S	S7008ES	S7008PS
	10 mm	S7010S	S7010ES	S7010PS
80 cm	5 mm	S8005S	S8005ES	S8005PS
	6 mm	S8006S	S8006ES	S8006PS
	7 mm	S8007S	S8007ES	S8007PS
	8 mm	S8008S	S8008ES	S8008PS
	10 mm	S8010S	S8010ES	S8010PS

taperflo™

ePTFE cónico sellado con gelatina
Arterial (0,61 mm)*
Venoso (0,38 mm)*



Longitud	Diámetro interno	Código de producto
20 cm	4 - 6 mm	W20ST0406S
20 cm	4 - 7 mm	W20ST0407S
25 cm	4 - 6 mm	W25ST0406S
25 cm	4 - 7 mm	W25ST0407S
30 cm	4 - 6 mm	W30ST0406S
30 cm	4 - 7 mm	W30ST0407S
35 cm	4 - 6 mm	W35ST0406S
35 cm	4 - 7 mm	W35ST0407S
35 cm	5 - 8 mm	W35ST0508S
40 cm	4 - 6 mm	W40ST0406S
40 cm	4 - 7 mm	W40ST0407S
40 cm	5 - 8 mm	W40ST0508S
45 cm	4 - 6 mm	W45ST0406S
45 cm	4 - 7 mm	W45ST0407S
45 cm	5 - 8 mm	W45ST0508S
50 cm	4 - 6 mm	W50ST0406S
50 cm	4 - 7 mm	W50ST0407S
50 cm	5 - 8 mm	W50ST0508S



Soporte central cónico escalonado

Longitud	Diámetro interno	Código de producto
40 cm	4 - 7 mm	W40ST0407CRS
45 cm	4 - 7 mm	W45ST0407CRS
45 cm	5 - 8 mm	W45ST0508CRS
50 cm	4 - 7 mm	W50ST0407CRS



Cono corto

Longitud	Diámetro interno	Código de producto
30 cm	4 - 7 mm	W30T0407S
40 cm	4 - 6 mm	W40T0406S
40 cm	4 - 7 mm	W40T0407S
45 cm	4 - 6 mm	W45T0406S
45 cm	4 - 7 mm	W45T0407S
50 cm	4 - 7 mm	W50T0407S




Soporte central cónico corto

Longitud	Diámetro interno	Código de producto
40 cm	4 - 6 mm	W40T0406CRS
40 cm	4 - 7 mm	W40T0407CRS
45 cm	4 - 7 mm	W45T0407CRS

SEALPTFE™ Wrap


Pared estándar (0,59 mm)



Longitud	Diámetro interno	Código de producto
40 cm	6 mm	S4006CPS
	7 mm	S4007CPS
	8 mm	S4008CPS
50 cm	6 mm	S5006CPS
	7 mm	S5007CPS
	8 mm	S5008CPS

SEALPTFE™ Wrap

Pared estándar (0,59 mm)




Longitud	Diámetro interno	Código de producto
50 cm	6 mm	S5006EPS



Soporte central cónico corto

Longitud	Diámetro interno	Código de producto
40 cm	4 - 6 mm	W40T0406CRS
40 cm	4 - 7 mm	W40T0407CRS
45 cm	4 - 7 mm	W45T0407CRS

* Datos obtenidos de producto de 4-7 mm



SEALPTFE™

El único ePTFE recubierto de gelatina del mundo

AHORRE TIEMPO


- Reduzca al mínimo el tiempo de espera entre operaciones
 - Reducción de la hemorragia de línea de sutura^{1,2}
- Potencial para reducir al mínimo el trauma del tejido
 - Permite una tunelización más atraumática
- Posicionamiento más fácil del injerto
 - Gracias a la menor fricción que aporta la gelatina
- Reducción del tiempo de operación por paciente
 - 3,4 minutos de media hasta la hemostasia²

AHORRE DINERO

- No se precisan suturas especiales
 - La gelatina sella los agujeros de sutura cuando se utiliza sutura monofilamentosa¹

REDUZCA LAS COMPLICACIONES

- Potencial de reducción del seroma alrededor del injerto
 - El gel actúa como una protección contra la exudación²
- Opción de enlace de antibiótico³




taperflo™



 **MAXIFLO™ Ultrathin**


Pared ultradelgada (0,35 mm)
ePTFE no sellado



Longitud	Diámetro interno	Código de producto Sin soporte	Código de producto Construcción de soporte unity™	Código de producto Soporte total de PTFE
10 cm	5 mm	T1005	T1005E	T1005P
	6 mm	T1006	T1006E	T1006P
	7 mm	T1007	T1007E	T1007P
	8 mm	T1008	T1008E	T1008P
	10 mm	T1010	T1010E	T1010P
20 cm	5 mm	T2005	T2005E	T2005P
	6 mm	T2006	T2006E	T2006P
	7 mm	T2007	T2007E	T2007P
	8 mm	T2008	T2008E	T2008P
	10 mm	T2010	T2010E	T2010P
30 cm	5 mm	T3005	T3005E	T3005P
	6 mm	T3006	T3006E	T3006P
	7 mm	T3007	T3007E	T3007P
	8 mm	T3008	T3008E	T3008P
	10 mm	T3010	T3010E	T3010P
40 cm	5 mm	T4005	T4005E	T4005P
	6 mm	T4006	T4006E	T4006P
	7 mm	T4007	T4007E	T4007P
	8 mm	T4008	T4008E	T4008P
	10 mm	T4010	T4010E	T4010P
50 cm	5 mm	T5005	T5005E	T5005P
	6 mm	T5006	T5006E	T5006P
	7 mm	T5007	T5007E	T5007P
	8 mm	T5008	T5008E	T5008P
	10 mm	T5010	T5010E	T5010P
60 cm	5 mm	T6005	T6005E	T6005P
	6 mm	T6006	T6006E	T6006P
	7 mm	T6007	T6007E	T6007P
	8 mm	T6008	T6008E	T6008P
	10 mm	T6010	T6010E	T6010P
70 cm	5 mm	T7005	T7005E	T7005P
	6 mm	T7006	T7006E	T7006P
	7 mm	T7007	T7007E	T7007P
	8 mm	T7008	T7008E	T7008P
	10 mm	T7010	T7010E	T7010P
80 cm	5 mm	T8005	T8005E	T8005P
	6 mm	T8006	T8006E	T8006P
	7 mm	T8007	T8007E	T8007P
	8 mm	T8008	T8008E	T8008P
	10 mm	T8010	T8010E	T8010P

 **MAXIFLO™ Wrap**

Pared estándar (0,47 mm)
ePTFE no sellado



Longitud	Diámetro interno	Código de producto Sin soporte	Código de producto Construcción de soporte unity™	Código de producto Soporte total de PTFE
10 cm	5 mm	S1005	S1005E	S1005P
	6 mm	S1006	S1006E	S1006P
	7 mm	S1007	S1007E	S1007P
	8 mm	S1008	S1008E	S1008P
	10 mm	S1010	S1010E	S1010P
20 cm	5 mm	S2005	S2005E	S2005P
	6 mm	S2006	S2006E	S2006P
	7 mm	S2007	S2007E	S2007P
	8 mm	S2008	S2008E	S2008P
	10 mm	S2010	S2010E	S2010P
30 cm	5 mm	S3005	S3005E	S3005P
	6 mm	S3006	S3006E	S3006P
	7 mm	S3007	S3007E	S3007P
	8 mm	S3008	S3008E	S3008P
	10 mm	S3010	S3010E	S3010P
40 cm	5 mm	S4005	S4005E	S4005P
	6 mm	S4006	S4006E	S4006P
	7 mm	S4007	S4007E	S4007P
	8 mm	S4008	S4008E	S4008P
	10 mm	S4010	S4010E	S4010P
50 cm	5 mm	S5005	S5005E	S5005P
	6 mm	S5006	S5006E	S5006P
	7 mm	S5007	S5007E	S5007P
	8 mm	S5008	S5008E	S5008P
	10 mm	S5010	S5010E	S5010P
60 cm	5 mm	S6005	S6005E	S6005P
	6 mm	S6006	S6006E	S6006P
	7 mm	S6007	S6007E	S6007P
	8 mm	S6008	S6008E	S6008P
	10 mm	S6010	S6010E	S6010P
70 cm	5 mm	S7005	S7005E	S7005P
	6 mm	S7006	S7006E	S7006P
	7 mm	S7007	S7007E	S7007P
	8 mm	S7008	S7008E	S7008P
	10 mm	S7010	S7010E	S7010P
80 cm	5 mm	S8005	S8005E	S8005P
	6 mm	S8006	S8006E	S8006P
	7 mm	S8007	S8007E	S8007P
	8 mm	S8008	S8008E	S8008P
	10 mm	S8010	S8010E	S8010P

 **MAXIFLO™**

ePTFE cónico no sellado Arterial (0,59 mm)*
Venoso (0,36 mm)*



Longitud	Diámetro interno	Código de producto Cónico escalonado
20 cm	4 - 6 mm	W20ST0406
20 cm	4 - 7 mm	W20ST0407
25 cm	4 - 6 mm	W25ST0406
25 cm	4 - 7 mm	W25ST0407
30 cm	4 - 6 mm	W30ST0406
30 cm	4 - 7 mm	W30ST0407
35 cm	4 - 6 mm	W35ST0406
35 cm	4 - 7 mm	W35ST0407
35 cm	5 - 8 mm	W35ST0508
40 cm	4 - 6 mm	W40ST0406
40 cm	4 - 7 mm	W40ST0407
40 cm	5 - 8 mm	W40ST0508
45 cm	4 - 6 mm	W45ST0406
45 cm	4 - 7 mm	W45ST0407
45 cm	5 - 8 mm	W45ST0508
50 cm	4 - 6 mm	W50ST0406
50 cm	4 - 7 mm	W50ST0407
50 cm	5 - 8 mm	W50ST0508



Soporte central cónico escalonado

Longitud	Diámetro interno	Código de producto
40 cm	4 - 7 mm	W40ST0407CR
45 cm	4 - 7 mm	W45ST0407CR
45 cm	5 - 8 mm	W45ST0508CR
50 cm	4 - 7 mm	W50ST0407CR



Cono corto

Longitud	Diámetro interno	Código de producto
30 cm	4 - 7 mm	W30T0407
40 cm	4 - 6 mm	W40T0406
40 cm	4 - 7 mm	W40T0407
45 cm	4 - 6 mm	W45T0406
45 cm	4 - 7 mm	W45T0407
50 cm	4 - 7 mm	W50T0407




Soporte central cónico corto

Longitud	Diámetro interno	Código de producto
40 cm	4 - 6 mm	W40T0406CR
40 cm	4 - 7 mm	W40T0407CR
45 cm	4 - 7 mm	W45T0407CR

 **MAXIFLO™ Wrap**

Pared estándar (0,47 mm)



Longitud	Diámetro interno	Código de producto Soporte central de PTFE
40 cm	6 mm	S4006CP
	7 mm	S4007CP
	8 mm	S4008CP
50 cm	6 mm	S5006CP
	7 mm	S5007CP
	8 mm	S5008CP

 **MAXIFLO™ Wrap**

Pared estándar (0,47 mm)



Longitud	Diámetro interno	Código de producto Soporte de extremo de PTFE
50 cm	6 mm	S5006EP

* Datos obtenidos de producto de 4-7 mm



MAXIFLO™



La **solución rentable** para todas sus necesidades de **ePTFE**

MAXIFLO™ Ultrathin

- ePTFE de paredes ultradelgadas con características de manejo excelentes²
- Flexibilidad mejorada que ofrece conformabilidad con los vasos²

MAXIFLO™ Wrap



- ePTFE de pared estándar con mayor retención de sutura¹

Gama cónica **MAXIFLO™**

- Diseño cónico creado específicamente para el acceso vascular
- Diseño cónicos cortos y cónicos escalonados con y sin soporte de polipropileno central

MAXIFLO™ con construcción **Unity™**


- El soporte PTFE expandido exclusivo se puede retener en la anastomosis, creando un “sistema de protección anastomótica”¹






Longitud	Diámetro interno	Sin reforzar	Soporte total de PTFE reforzado
		Código de producto	Código de producto
40 cm	6 mm	H4006	H4006P
	7 mm	H4007	H4007P
	8 mm	H4008	H4008P
50 cm	6 mm	H5006	H5006P
	7 mm	H5007	H5007P
	8 mm	H5008	H5008P
60 cm	6 mm	H6006	H6006P
	7 mm	H6007	H6007P
	8 mm	H6008	H6008P
70 cm	6 mm	H7006	H7006P
	7 mm	H7007	H7007P
	8 mm	H7008	H7008P
80 cm	6 mm	H8006	H8006P
	7 mm	H8007	H8007P
	8 mm	H8008	H8008P

DELGADO Y RESISTENTE




TECNOLOGÍA DE DOBLE CAPA DE POROS ABIERTOS



DISEÑO DE CAMPANA DE FLUJO DISTAL



CAMPANA DISTAL TOTALMENTE PERSONALIZABLE





 **COBRAHOOD™**

Prótesis vasculares **ePTFE** con **Campana distal**

DELGADO Y RESISTENTE

- El mismo grosor de pared que la competencia¹
- 51% más de retención de sutura¹

TECNOLOGÍA DE DOBLE CAPA DE POROS ABIERTOS

- La envuelta tiene la misma porosidad que el injerto base, por lo que permite la incorporación de tejido⁴

DISEÑO DE CAMPANA DE FLUJO DISTAL

Comparado con una prótesis ePTFE estándar de 6 mm, **COBRAHOOD™** ofrece:

- 5,6 veces más superficie⁵
- 10 veces más volumen anastomótico⁵

CAMPANA DISTAL TOTALMENTE PERSONALIZABLE

- El diseño de cubierta de pared delgada se puede cortar y conformar fácilmente



Referencias bibliográficas

1. Datos Clínicos del soporte post-marketing en archivo de Vascutek
2. Supporting Post-Market Clinical Data on file at Vascutek
3. BSI Ref EQ# 10020927
4. Moore WS & Malone JM: Vascular Repair. In Hunt TK & Dumphrey JE (eds) **Fundamentals of Wound Management** New York *Appleton-Century-Crofts* 1979
5. Based on direct comparison between 6mm straight ePTFE prosthesis and 6mm Cobrahood™. Datos Clínicos del soporte post-marketing en archivo de Vascutek. 2007

PARA MÁS INFORMACIÓN
(+44) 141 812 5555
ePTFE@vascutek.com

Distribuido por~

CARDIVA S.L.

c/ General Eguía, 48

48013 Bilbao (Vizcaya), España

Tel : 34 94 439 50 98

Fax : 34 94 439 52 68

SEALPTFE™, *taperflo*™, MAXIFLO™ y COBRAHOOD™ son marcas registradas de Vascutek.
unity™ è un prodotto della

 **VASCUTEK**
TERUMO

VASCUTEK, a TERUMO Company
Newmains Avenue, Inchinnan
Renfrewshire PA4 9RR Escocia, Reino Unido

Tel : (+44) 141 812 5555

Fax : (+44) 141 812 7170

ePTFE@vascutek.com

www.vascutek.com

B253/3E

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

No.**CE 00422****Issued To:**


**Vascutek Limited
a Terumo Company
Newmains Avenue
Inchinnan
Renfrewshire
Scotland
PA4 9RR
United Kingdom**

In respect of:

The design, development and manufacture of Sealed and Unsealed, Woven and Knitted Polyester, PTFE, Stented Nitinol Vascular Prostheses and Single Use Vascular Instruments.

on the basis of our examination of the quality assurance system under the requirements of Council Directive 93/42/EEC, Annex II excluding section 4. The quality assurance system meets the requirements of the directive. For the placing on the market of class III products an Annex II section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):



Albert Roossien, Regulatory Lead

First Issued: **1994-12-21**

Date: **2019-02-05**

Expiry Date: **2023-03-02**

...making excellence a habit.™

Page 1 of 1

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 00422**
 Date: **2019-02-05**
 Issued To: **Vasutek Limited
 a Terumo Company
 Newmains Avenue
 Inchinnan
 Renfrewshire
 Scotland
 PA4 9RR
 United Kingdom**

Subcontractor:	Service(s) supplied
Gelita USA, Inc. 2445 Port Neal Industrial Road Sergeant Bluff Iowa 51054 USA	Animal Tissues / Derivatives
Sterigenics Cotes Park Industrial Estate Somercotes, Alfreton Derbyshire, DE55 4NJ United Kingdom	ETO Sterilization
Synergy Health Sterilisation UK Ltd (Synergy Health – AST – Bradford) Roysdale Way Euroway Industrial Estate Bradford BD4 6SE United Kingdom	Gamma Sterilization

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 00422**
 Date: **2019-02-05**
 Issued To: **Vascutek Limited
 a Terumo Company
 Newmains Avenue
 Inchinnan
 Renfrewshire
 Scotland
 PA4 9RR
 United Kingdom**

Subcontractor:	Service(s) supplied
Synergy Health Sterilisation UK Ltd Brunel Close Drayton Fields Industrial Estate Daventry NN11 8RB United Kingdom	Gamma Sterilization
Terumo Vietnam Co., Ltd. Lot 44-B-C Quang Minh Industrial Zone Me Linh District, Hanoi City Vietnam	Manufacture
Tyson Fresh Meat Inc. 1131 Dakota Avenue Dakota City Nebraska 68731 USA	Animal Tissues / Derivatives

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Certificate History

Certificate No: **CE 00422**
Date: **2019-02-05**
Issued To: **Vascutek Limited
a Terumo Company
Newmains Avenue
Inchinnan
Renfrewshire
Scotland
PA4 9RR
United Kingdom**

Date	Reference Number	Action
21 December 1994		First Issue.
19 May 1995		Griffith Micro Science Limited added to the sub-contractor list.
1 February 1996		Tissue Patches added to the scope.
13 September 1996		Isotron PLC (Bedford) added to the sub-contractor list.
21 September 1999		PTFE Vascular Protheses added to the scope and certificate renewal.
12 January 2001		Stented Vascular Protheses added to the scope.
03 March 2003		Nitinol Vascular Protheses added to the scope.

...making excellence a habit.™

Page 1 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Certificate - Full Quality Assurance System

Certificate History

Certificate No: **CE 00422**
 Date: **2019-02-05**
 Issued To: **Vascutek Limited
 a Terumo Company
 Newmains Avenue
 Inchinnan
 Renfrewshire
 Scotland
 PA4 9RR
 United Kingdom**

Date	Reference Number	Action
28 September 2005		Removal of Tissue Patches from the scope, and sterilizer name change from Griffith Micro Science Ltd to IBA and re-formatted company address adding Scotland.
21 September 2006		Heart Valve Conduits added to the scope. Kohler Medical Ltd added as a sub-contract manufacturer and sterilizer.
29 June 2007		Acquisition of subcontractor Koehler Medical Ltd. Change of name from Koehler Medical Ltd to Vascutek Limited, a Terumo Company and update to list of significant subcontractors.
03 March 2008	7005772	Certificate renewal including change of subcontractor name from IBA to Sterigenics. Also change "Disposable" to "Single Use" consistent with current terminology.

...making excellence a habit.™

Page 2 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Certificate - Full Quality Assurance System

Certificate History

Certificate No: **CE 00422**
 Date: **2019-02-05**
 Issued To: **Vascutek Limited
 a Terumo Company
 Newmains Avenue
 Inchinnan
 Renfrewshire
 Scotland
 PA4 9RR
 United Kingdom**

Date	Reference Number	Action
15 July 2010	7546406	Addition of Terumo Vietnam Co., to the list of significant subcontractors.
25 April 2012	7794623	Update the name of Isotron, Ltd to Synergy Health Sterilisation UK Ltd.
27 February 2013	7905922	Certificate renewal.
27 August 2015	8405023	Removal of "Heart Valve Conduits" from the scope and the significant subcontractor Vascutek Ltd (Swillington, Leeds).
23 October 2017	8595440	Addition of animal substance subcontractors Tyson Fresh Meat Inc. and Gelita USA Inc.
02 March 2018	8891380	Certificate renewal and addition of Synergy Health Daventry as a subcontractor for sterilization.

...making excellence a habit.™

Page 3 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Certificate - Full Quality Assurance System

Certificate History

Certificate No: **CE 00422**
Date: **2019-02-05**
Issued To: **Vascutek Limited
a Terumo Company
Newmains Avenue
Inchinnan
Renfrewshire
Scotland
PA4 9RR
United Kingdom**

Date	Reference Number	Action
Current	7778835	Traceable to NB 0086. Administrative Subcontractor Service wording update for: Gelita USA Inc, Iowa 51054 from "Animal Substances" to "Animal Tissues / Derivatives". Sterigenics Derbyshire DE55 4NJ from "Sterilization" to "ETO Sterilization". Synergy Health Sterilisation UK Ltd, Daventry NN11 8RB from "Sterilization" to "Gamma Sterilization". Synergy Health Sterilisation UK Ltd, Bradford BD4 6SE from "Sterilization" to "Gamma Sterilization". Tyson Fresh Meat Inc, Nebraska 68731 from "Animal Substances" to "Animal Tissues / Derivatives".

...making excellence a habit.™

Page 4 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

Vascutek Limited
Newmains Avenue
Inchinnan
PA4 9RR
United Kingdom

Holds Certificate Number:

MD 93897

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

The design and manufacture of sealed and unsealed, woven, knitted, stented and PTFE vascular prostheses and disposable vascular instruments.
Contract sterilization of medical devices in accordance with EN ISO 11135:2014.



For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2005-02-14

Latest Revision Date: 2020-05-14

Effective Date: 2020-05-14

Expiry Date: 2022-12-25

Page: 1 of 1



003

...making excellence a habit.™

Declaration of Conformity Record

Manufacturer's Name: Vascutek Ltd, a TERUMO Company**Manufacturer's Address:** Newmains Avenue, Inchinnan, Renfrewshire PA4 9RR Scotland**Product Family:** Woven Vascular Grafts**Product Class:** Class III**Product Range: Twillweave Vascular Grafts**Straight
Bifurcate**Gelweave Vascular Grafts**Straight
Bifurcate
Trifurcate
Trifurcate with Side Branch With or Without Radiopaque Marker
Double Bifurcates
Curved Graft
Curved Ante-Flo
Precurved 4 Branch Plexus
Ante-Flo
Trifurcate Arch Graft
Thoracic Arch Graft with or Without Radiopaque Marker
Three Branch Plexus
Aortic Arch
Trifurcate Arch Graft with Side Branch
Four Branch Plexus
Plexus (With Extra Black Line)
Thoracoabdominal
Aortic Arch Trifurcate Graft
Branched Arch Graft With or without Radiopaque Markers
Ante-flo with Collar
Four Branch Plexus with Collar
Siena Collared Straight with Radiopaque Markers
Siena Collared Ante-flo with Radiopaque Markers
Siena Collared Four Branch Plexus with Radiopaque Markers
Valsalva
Valsalva with Collar
Valsalva + Ante-Flo**Revision Number:** EUDOC-001W-001

I, the undersigned, hereby declare that the medical device(s) specified above conforms with the Essential Requirements listed in Annex I of EC Directive 93/42/EEC.

This declaration is supported by:

EC Quality System approval statement (Annex II) No.: **CE00422**


Issued by:	BSI Netherlands (NB No: 2797) *	First Issued On:	21/12/94
Current Issue No.:	18	Dated:	05/02/19

EC Design Dossier examination No. **CE00514, CE01230**

Issued by:	BSI Netherlands (NB No: 2797) *	First Issued On:	14/02/95, 06/03/96
Current Issue No.:	23, 10	Dated:	05/02/19, 05/02/19

Declaration of Conformity Record

*Previous issues of these certificates were issued by BSI UK (NB No: 0086). Vascutek completed Notified Body migration from BSI UK to BSI Netherlands in February 2019 as part of preparations for the UK leaving the EU.

Signed: 
Name: Neil McLACHLAN

Date: 23 APRIL 2019

VP Regulatory Affairs

Declaration of Conformity Record

Catalogue Numbers Covered by Attached Woven Vascular Grafts Declaration of Conformity

CE01230 TWILLWEAVE VASCULAR GRAFTS

RANGE: TWILLWEAVE STRAIGHTS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
6 mm	536006	60 cm
7 mm	536007	60 cm
8 mm	536008	60 cm
9 mm	536009	60 cm
10 mm	536010	60 cm
11 mm	536011	60 cm
12 mm	536012	60 cm
14 mm	536014	60 cm
16 mm	536016	60 cm
18 mm	536018	60 cm
20 mm	536020	60 cm
22 mm	536022	60 cm
24 mm	536024	60 cm
26 mm	536026	60 cm
28 mm	536028	60 cm
30 mm	536030	60 cm
32 mm	536032	60 cm
34 mm	536034	60 cm
36 mm	536036	60 cm
38 mm	536038	60 cm
18 mm	535018	50 cm
20 mm	535020	50 cm
22 mm	535022	50 cm
24 mm	535024	50 cm
26 mm	535026	50 cm
28 mm	535028	50 cm
30 mm	534030	40 cm
32 mm	534032	40 cm
34 mm	534034	40 cm
36 mm	534036	40 cm
38 mm	534038	40 cm
6 mm	533006	30 cm
7 mm	533007	30 cm
8 mm	533008	30 cm
9 mm	533009	30 cm
10 mm	533010	30 cm
11 mm	533011	30 cm
12 mm	533012	30 cm
14 mm	533014	30 cm
16 mm	533016	30 cm
18 mm	533018	30 cm
20 mm	533020	30 cm
22 mm	533022	30 cm
24 mm	533024	30 cm
26 mm	533026	30 cm
28 mm	533028	30 cm
30 mm	533030	30 cm
32 mm	533032	30 cm
34 mm	533034	30 cm
36 mm	533036	30 cm
38 mm	533038	30 cm

Declaration of Conformity Record

RANGE: TWILLWEAVE STRAIGHTS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
18 mm	532518	25 cm
20 mm	532520	25 cm
22 mm	532522	25 cm
24 mm	532524	25 cm
26 mm	532526	25 cm
28 mm	532528	25 cm
30 mm	532030	20 cm
32 mm	532032	20 cm
34 mm	532034	20 cm
36 mm	532036	20 cm
38 mm	532038	20 cm
6 mm	531506	15 cm
7 mm	531507	15 cm
8 mm	531508	15 cm
9 mm	531509	15 cm
10 mm	531510	15 cm
11 mm	531511	15 cm
12 mm	531512	15 cm
14 mm	531514	15 cm
16 mm	531516	15 cm
18 mm	531218	12.5 cm
20 mm	531220	12.5 cm
22 mm	531222	12.5 cm
24 mm	531224	12.5 cm
26 mm	531226	12.5 cm
28 mm	531228	12.5 cm
30 mm	531030	10 cm
32 mm	531032	10 cm
34 mm	531034	10 cm
36 mm	531036	10 cm
38 mm	531038	10 cm

RANGE: TWILLWEAVE BIFURCATES		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
12 x 6 mm	531206	45 cm
14 x 7 mm	531407	45 cm
16 x 8 mm	531608	45 cm
18 x 9 mm	531809	45 cm
20 x 10 mm	532010	45 cm
22 x 11 mm	532211	45 cm
24 x 12 mm	532412	45 cm

Declaration of Conformity Record

RANGE: UNCRIMPED TWILLWEAVE STRAIGHTS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
14 mm	553014	30 cm
16 mm	553016	30 cm
18 mm	553018	30 cm
20 mm	553020	30 cm
24 mm	553024	30 cm
26 mm	553026	30 cm
28 mm	553028	30 cm
30 mm	553030	30 cm
32 mm	553032	30 cm
34 mm	553034	30 cm
36 mm	553036	30 cm
38 mm	553038	30 cm
14 mm	555014	50 cm
16 mm	555016	50 cm
18 mm	555018	50 cm
20 mm	555020	50 cm
22 mm	555022	50 cm
24 mm	555024	50 cm
26 mm	555026	50 cm
28 mm	555028	50 cm
30 mm	555030	50 cm
32 mm	555032	50 cm
34 mm	555034	50 cm
36 mm	555036	50 cm
38 mm	555038	50 cm

RANGE: UNCRIMPED TWILLWEAVE BIFURCATES		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
12 x 6 mm	551206	45 cm
14 x 7 mm	551407	45 cm
16 x 8 mm	551608	45 cm
18 x 9 mm	551809	45 cm
20 x 10 mm	552010	45 cm
22 x 11 mm	552211	45 cm
24 x 12 mm	552412	45 cm

SEALED GRAFTS: WOVEN – GELWEAVE

The base fabric for Gelweave is Twillweave

RANGE: GELWEAVE STRAIGHTS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
6 mm	736006	60 cm
7 mm	736007	60 cm
8 mm	736008	60 cm
9 mm	736009	60 cm
10 mm	736010	60 cm
11 mm	736011	60 cm
12 mm	736012	60 cm
14 mm	736014	60 cm
16 mm	736016	60 cm
6 mm	734006	40 cm
8 mm	734008	40 cm
6 mm	733006	30 cm
7 mm	733007	30 cm
8 mm	733008	30 cm
9 mm	733009	30 cm
10 mm	733010	30 cm
11 mm	733011	30 cm
12 mm	733012	30 cm
14 mm	733014	30 cm

Declaration of Conformity Record

16 mm	733016	30 cm
6 mm	731506	15 cm
7 mm	731507	15 cm
8 mm	731508	15 cm
9 mm	731509	15 cm
10 mm	731510	15 cm
11 mm	731511	15 cm
12 mm	731512	15 cm
14 mm	731514	15 cm
16 mm	731516	15 cm
18 mm	735018	50 cm
20 mm	735020	50 cm
22 mm	735022	50 cm
24 mm	735024	50 cm
26 mm	735026	50 cm
28 mm	735028	50 cm
30 mm	735030	50 cm
32 mm	735032	50 cm
34mm	735034	50 cm
26 mm	730526	5 cm
28 mm	730528	5 cm
18 mm	732518	25 cm
20 mm	732520	25 cm
22 mm	732522	25 cm
24 mm	732524	25 cm
26 mm	732526	25 cm
28 mm	732528	25 cm
18 mm	731218	12.5 cm
20 mm	731220	12.5 cm
22 mm	731222	12.5 cm
24 mm	731224	12.5 cm
26 mm	731226	12.5 cm
28 mm	731228	12.5 cm
18 mm	733018	30 cm
20 mm	733020	30 cm
22 mm	733022	30 cm
24 mm	733024	30 cm
26 mm	733026	30 cm
28 mm	733028	30 cm
30 mm	733030	30 cm
32 mm	733032	30 cm
34 mm	733034	30 cm
36 mm	733036	30 cm
38 mm	733038	30 cm
30 mm	731030	10 cm
32 mm	731032	10 cm
34 mm	731034	10 cm
38 mm	731038	10 cm
30 mm	732030	20 cm
32 mm	732032	20 cm
34 mm	732034	20 cm
38 mm	732038	20 cm

Declaration of Conformity Record

30 mm	734030	40 cm
32 mm	734032	40 cm
34 mm	734034	40 cm
36 mm	734036	40 cm
38 mm	734038	40 cm
18 mm	736018	60 cm
20 mm	736020	60 cm
22 mm	736022	60 cm
24 mm	736024	60 cm
26 mm	736026	60 cm
28 mm	736028	60 cm
30 mm	736030	60 cm
32 mm	736032	60 cm
34 mm	736034	60 cm
36 mm	736036	60 cm
38 mm	736038	60 cm
24 mm	7310024	100 cm
26 mm	7310026	100 cm
28 mm	7310028	100 cm

RANGE: GELWEAVE STRAIGHTS WITH 3 BLACK LINES (VS192)

SIZE	CATALOGUE NUMBER	USABLE LENGTH
18 mm	7312183BL	12.5 cm
20 mm	7312203BL	12.5 cm
22 mm	7312223BL	12.5 cm
24 mm	7312243BL	12.5 cm
26 mm	7312263BL	12.5 cm
28 mm	7312283BL	12.5 cm
18 mm	7325183BL	25 cm
20 mm	7325203BL	25 cm
22 mm	7325223BL	25 cm
24 mm	7325243BL	25 cm
26 mm	7325263BL	25 cm
28 mm	7325283BL	25 cm
18 mm	7350183BL	50 cm
20 mm	7350203BL	50 cm
22 mm	7350223BL	50 cm
30 mm	7310303BL	10 cm
30 mm	7320303BL	20 cm
30 mm	7340303BL	40 cm

RANGE: GELWEAVE BIFURCATES

SIZE	CATALOGUE NUMBER	USABLE LENGTH
12 x 6 mm	731206	45 cm
14 x 7 mm	731407	45 cm
16 x 8 mm	731608	45 cm
18 x 9 mm	731809	45 cm
20 x 10 mm	732010	45 cm
22 x 11 mm	732211	45 cm
24 x 12 mm	732412	45 cm

RANGE: GELWEAVE TRIFURCATE

SIZE	CATALOGUE NUMBER	USABLE LENGTH
12 x 6 x 6 x 7 mm	731206/7	45 cm
14 x 7 x 7 x 7 mm	731407/7	45 cm
16 x 8 x 8 x 7 mm	731608/7	45 cm

Declaration of Conformity Record

RANGE: GELWEAVE TRIFURCATE WITH SIDE BRANCH		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
12 x 10 x 6 x 6 x 10 mm	731206/10/10	45 cm + 1 @ 15cm
14 x 10 x 7 x 7 x 10 mm	731407/10/10	45 cm + 1 @ 15cm
16 x 10 x 8 x 8 x 10 mm	731608/10/10	45 cm + 1 @ 15cm
16 x 8 x 8 x 8 x 10 mm	731608/8/10	45 cm + 1 @ 15cm

RANGE: GELWEAVE™ TRIFURCATE WITH SIDE BRANCH WITH RADIOPAQUE MARKER		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
12/10/6/6 x 10mm	731206/10/10RM	45 cm + 1 @ 15cm
14/10/7/7 x 10mm	731407/10/10RM	45 cm + 1 @ 15cm
16/10/8/8 x 10mm	731608/10/10RM	45 cm + 1 @ 15cm
16/8/8/8 x 10mm	731608/8/10RM	45 cm + 1 @ 15cm

RANGE: GELWEAVE DOUBLE BIFURCATES		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
16 x 7 x 8 mm	731607DB8	45 cm
18 x 8 x 9 mm	731808DB9	45 cm
20 x 9 x 10 mm	732009DB10	45 cm

RANGE: GELWEAVE EXTRA LONG BIFURCATES				
CATALOGUE NUMBER	GRAFT DIAMETER	BODY LENGTH	LEG LENGTH	USEABLE LENGTH
731206XL55	12 x 6mm	10cm	45cm	55cm
731407XL55	14 x 7mm	10cm	45cm	55cm
731608XL55	16 x 8mm	10cm	45cm	55cm
731809XL55	18 x 9mm	10cm	45cm	55cm
732010XL55	20 x 10mm	10cm	45cm	55cm
732211XL55	22 x 11mm	10cm	45cm	55cm
732412XL55	24 x 12mm	10cm	45cm	55cm

RANGE: GELWEAVE CURVED GRAFT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22 mm	735022ST	50 cm
24 mm	735024ST	50 cm
26 mm	735026ST	50 cm
28 mm	735028ST	50 cm
30 mm	735030ST	50 cm
32 mm	735032ST	50 cm
34 mm	735034ST	50 cm

RANGE: GELWEAVE CURVED ANTE-FLO GRAFT OFFSET BRANCH		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
26 x 8 mm	735026ST/8	50cm + 1 @ 15 cm
28 x 8 mm	735028ST/8	50cm + 1 @ 15 cm
30 x 8 mm	735030ST/8	50cm + 1 @ 15 cm

RANGE: GELWEAVE™ PRE-CURVED ANTE-FLO		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22 x 8mm	735022ST/8B	50cm + 1 @ 15cm
24 x 8mm	735024ST/8B	50cm + 1 @ 15cm
26 x 8mm	735026ST/8B	50cm + 1 @ 15cm
28 x 8mm	735028ST/8B	50cm + 1 @ 15cm
30 x 8mm	735030ST/8B	50cm + 1 @ 15cm
32 x 8mm	735032ST/8B	50cm + 1 @ 15cm
34 x 8mm	735034ST/8B	50cm + 1 @ 15cm
22 x 10mm	735022ST/10B	50cm + 1 @ 15cm
24 x 10mm	735024ST/10B	50cm + 1 @ 15cm
26 x 10mm	735026ST/10B	50cm + 1 @ 15cm
28 x 10mm	735028ST/10B	50cm + 1 @ 15cm

Declaration of Conformity Record

30 x 10mm	735030ST/10B	50cm + 1 @ 15cm
32 x 10mm	735032ST/10B	50cm + 1 @ 15cm
34 x 10mm	735034ST/10B	50cm + 1 @ 15cm
22 x 8 mm	735022ST/8A	50cm + 1 @ 15cm
24 x 8 mm	735024ST/8A	50cm + 1 @ 15cm
26 x 8 mm	735026ST/8A	50cm + 1 @ 15cm
28 x 8 mm	735028ST/8A	50cm + 1 @ 15cm
30 x 8 mm	735030ST/8A	50cm + 1 @ 15cm
32 x 8 mm	735032ST/8A	50cm + 1 @ 15cm
34 x 8 mm	735034ST/8A	50cm + 1 @ 15cm

RANGE: GELWEAVE™ PRECURVED 4 BRANCH PLEXUS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20/12/8/10 x 10mm	732012810/10STA	40cm + 4 @ 15cm
22/12/8/10 x 10mm	732212810/10STA	40cm + 4 @ 15cm
24/12/8/10 x 10mm	732412810/10STA	40cm + 4 @ 15cm
26/12/8/10 x 10mm	732612810/10STA	40cm + 4 @ 15cm
28/12/8/10 x 10mm	732812810/10STA	40cm + 4 @ 15cm
30/12/8/10 x 10mm	733012810/10STA	40cm + 4 @ 15cm
32/12/8/10 x 10mm	733212810/10STA	40cm + 4 @ 15cm
34/12/8/10 x 10mm	733412810/10STA	40cm + 4 @ 15cm
20/10/8/8 x 10mm	73201088/10STA	40cm + 4 @ 15cm
22/10/8/8 x 10mm	73221088/10STA	40cm + 4 @ 15cm
24/10/8/8 x 10mm	73241088/10STA	40cm + 4 @ 15cm
26/10/8/8 x 10mm	73261088/10STA	40cm + 4 @ 15cm
28/10/8/8 x 10mm	73281088/10STA	40cm + 4 @ 15cm
30/10/8/8 x 10mm	73301088/10STA	40cm + 4 @ 15cm
32/10/8/8 x 10mm	73321088/10STA	40cm + 4 @ 15cm
34/10/8/8 x 10mm	73341088/10STA	40cm + 4 @ 15cm
20/10/8/8 x 8mm	73201088/8SD	40cm + 4 @ 15cm
22/10/8/8 x 8mm	73221088/8SD	40cm + 4 @ 15cm
24/10/8/8 x 8mm	73241088/8SD	40cm + 4 @ 15cm
26/10/8/8 x 8mm	73261088/8SD	40cm + 4 @ 15cm
28/10/8/8 x 8mm	73281088/8SD	40cm + 4 @ 15cm
30/10/8/8 x 8mm	73301088/8SD	40cm + 4 @ 15cm
32/10/8/8 x 8mm	73321088/8SD	40cm + 4 @ 15cm
34/10/8/8 x 8mm	73341088/8SD	40cm + 4 @ 15cm

RANGE: GELWEAVE ANTE-FLO		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
18 x 8mm	734018/8	40cm + 1 @ 15cm
20 x 8 mm	734020/8	40cm + 1 @ 15cm
22 x 8 mm	734022/8	40cm + 1 @ 15cm
24 x 8 mm	734024/8	40cm + 1 @ 15cm
26 x 8 mm	734026/8	40cm + 1 @ 15cm
28 x 8 mm	734028/8	40cm + 1 @ 15cm
30 x 8 mm	734030/8	40cm + 1 @ 15cm
32 x 8 mm	734032/8	40cm + 1 @ 15cm
34 x 8 mm	734034/8	40cm + 1 @ 15cm
22 x 8 mm	732022/8	20cm + 1 @ 15cm
24 x 8 mm	732024/8	20cm + 1 @ 15cm
26 x 8 mm	732026/8	20cm + 1 @ 15cm
28 x 8 mm	732028/8	20cm + 1 @ 15cm
30 x 8 mm	732030/8	20cm + 1 @ 15cm
32 x 8 mm	732032/8	20cm + 1 @ 15cm
34 x 8 mm	732034/8	20cm + 1 @ 15cm
24 x 10 mm	735024/10	50cm + 1 @ 15cm
28 x 10 mm	735028/10	50cm + 1 @ 15cm

Declaration of Conformity Record

32 x 10 mm	735032/10	50cm + 1@15cm
20 x 10 mm	734020/10	40cm + 1@15cm
22 x 10 mm	734022/10	40cm + 1@15cm
22 x 10 mm	734022/10S	40cm + 1@15cm
24 x 10 mm	734024/10	40cm + 1@15cm
24 x 10 mm	734024/10S	40cm + 1@25cm
26 x 10 mm	734026/10	40cm + 1@15cm
26 x 10 mm	734026/10S	40cm + 1@25cm
28 x 10 mm	734028/10	40cm + 1@15cm
28 x 10 mm	734028/10S	40cm + 1@25cm
30 x 10 mm	734030/10	40cm + 1@15cm
20 x 10 mm	734020/10S	40cm + 1@25cm
30 x 10 mm	734030/10S	40cm + 1@25cm
32 x 10 mm	734032/10	40cm + 1@15cm
34 x 10 mm	734034/10	40cm + 1@15cm
22 x 10 mm	732022/10	20cm + 1@15cm
24 x 10 mm	732024/10	20cm + 1@15cm
26 x 10 mm	732026/10	20cm + 1@15cm
28 x 10 mm	732028/10	20cm + 1@15cm
30 x 10 mm	732030/10	20cm + 1@15cm
32 x 10 mm	732032/10	20cm + 1@15cm
34 x 10 mm	732034/10	20cm + 1@15cm

RANGE: GELWEAVE EXTRA LENGTH ANTE-FLO WITH OFFSET BRANCH		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22 x 8 mm	735022/8S	50cm + 1@15cm
24 x 8 mm	735024/8S	50cm + 1@15cm
26 x 8 mm	735026/8S	50cm + 1@15cm
28 x 8 mm	735028/8S	50cm + 1@15cm
30 x 8 mm	735030/8S	50cm + 1@15cm
32 x 8 mm	735032/8S	50cm + 1@15cm

RANGE: GELWEAVE™ LONGER ANTE-FLO WITH OFFSET BRANCH (VS076)		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22 x 8 mm	736022/8S	60 cm + 1@15cm
24 x 8 mm	736024/8S	60 cm + 1@15cm
26 x 8 mm	736026/8S	60 cm + 1@15cm
28 x 8 mm	736028/8S	60 cm + 1@15cm
30 x 8 mm	736030/8S	60 cm + 1@15cm
32 x 8 mm	736032/8S	60 cm + 1@15cm

RANGE: GELWEAVE™ ANTE-FLO WITH OFFSET SIDE BRANCH		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
26 x 10 mm	733026/10S	30cm + 1@15cm
28 x 10 mm	733028/10S	30cm + 1@15cm
30 x 10 mm	733030/10S	30cm + 1@15cm
32 x 10 mm	733032/10S	30cm + 1@15cm
26 x 10 mm	736026/10S	60cm + 1@15cm
28 x 10 mm	736028/10S	60cm + 1@15cm
30 x 10 mm	736030/10S	60cm + 1@15cm
32 x 10 mm	736032/10S	60cm + 1@15cm

RANGE: GELWEAVE TRIFURCATE ARCH GRAFT		
SIZE	CATALOGUE NUMBER	SIZE
12/8/8mm	732012/8X2	20cm+2@10cm
12/10/10mm	732012/10X2	20cm+2@10cm
14/8/8mm	732014/8X2	20cm+2@10cm
14/10/10mm	732014/10X2	20cm+2@10cm

Declaration of Conformity Record

RANGE: GELWEAVE™ THORACIC ARCH GRAFT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
14X10X10mm	732214/10/10A	22cm + 10cm + 13cm
12X10X8mm	732212/10/8A	22cm + 10cm + 13cm

RANGE: GELWEAVE™ THORACIC ARCH GRAFT WITH RADIOPAQUE MARKER		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
14X10X10mm	732214/10/10ARM	22cm + 10cm + 13cm
12X10X8mm	732212/10/8ARM	22cm + 10cm + 13cm

RANGE: GELWEAVE –THREE BRANCH PLEXUS		
SIZE	CATALOGUE NUMBER	SIZE
20/10/8/8 mm	73201088	40 cm + 3 @ 15 cm
22/10/8/8 mm	73221088	40 cm + 3 @ 15 cm
24/10/8/8 mm	73241088	40 cm + 3 @ 15 cm
26/10/8/8 mm	73261088	40 cm + 3 @ 15 cm
28/10/8/8 mm	73281088	40 cm + 3 @ 15 cm
30/10/8/8 mm	73301088	40 cm + 3 @ 15 cm
32/10/8/8 mm	73321088	40 cm + 3 @ 15 cm
34/10/8/8 mm	73341088	40 cm + 3 @ 15 cm

RANGE: GELWEAVE™ AORTIC ARCH		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
14/10/10/10mm	732514/10X3	25cm + 3@10cm

RANGE: GELWEAVE™ TRIFURCATE ARCH GRAFT WITH SIDE BRANCH		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
12/10/10X8mm	73128/10X2	30cm + 3@10cm
10/8/8X8mm	73108/8X2	30cm + 3@10cm
10/10/8X8mm	731010/8X2	30cm + 3@10cm
12/10/8X8mm	731210/8X2	30cm + 3@10cm

RANGE: GELWEAVE – FOUR BRANCH PLEXUS		
SIZE	CATALOGUE NUMBER	SIZE
20/10/8/8 x 8 mm	73201088/8	40 cm + 4 @ 15 cm
22/10/8/8 x 8 mm	73221088/8	40 cm + 4 @ 15 cm
24/10/8/8 x 8 mm	73241088/8	40 cm + 4 @ 15 cm
26/10/8/8 x 8 mm	73261088/8	40 cm + 4 @ 15 cm
28/10/8/8 x 8 mm	73281088/8	40 cm + 4 @ 15 cm
30/10/8/8 x 8 mm	73301088/8	40 cm + 4 @ 15 cm
32/10/8/8 x 8 mm	73321088/8	40 cm + 4 @ 15 cm
34/10/8/8 x 8 mm	73341088/8	40 cm + 4 @ 15 cm
20/10/8/8 x 10 mm	73201088/10	40 cm + 4 @ 15 cm
22/10/8/8 x 10 mm	73221088/10	40 cm + 4 @ 15 cm
24/10/8/8 x 10 mm	73241088/10	40 cm + 4 @ 15 cm
26/10/8/8 x 10 mm	73261088/10	40 cm + 4 @ 15 cm
28/10/8/8 x 10 mm	73281088/10	40 cm + 4 @ 15 cm
30/10/8/8 x 10 mm	73301088/10	40 cm + 4 @ 15 cm
32/10/8/8 x 10 mm	73321088/10	40 cm + 4 @ 15 cm
34/10/8/8 x 10 mm	73341088/10	40 cm + 4 @ 15 cm
28/12/8/8X8 mm	73281288/8	40 cm + 4 @ 15 cm

RANGE: GELWEAVE™ PLEXUS (WITH EXTRA BLACK LINE)		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20/12/8/8X8mm	73201288/8A	40cm + 4@15cm
22/12/8/8X8mm	73221288/8A	40cm + 4@15cm
24/12/8/8X8mm	73241288/8A	40cm + 4@15cm
26/12/8/8X8mm	73261288/8A	40cm + 4@15cm

Declaration of Conformity Record

28/12/8/8X8mm	73281288/8A	40cm + 4@15cm
30/12/8/8X8mm	73301288/8A	40cm + 4@15cm
32/12/8/8X8mm	73321288/8A	40cm + 4@15cm
34/12/8/8X8mm	73341288/8A	40cm + 4@15cm

RANGE: THORACOABDOMINAL GRAFT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
26/12/12/12 x12mm	736026/12 x 4	60cm + 4 @ 12 cm
28/12/12/12 x 12mm	736028/12 x 4	60cm + 4 @ 12 cm
20/6/6/8 x 6mm	7320668/6	60cm + 4@10cm
22/6/6/8 x 6mm	7322668/6	60cm + 4@10cm
24/6/6/8 x 6mm	7324668/6	60cm + 4@10cm
26/6/6/8 x 6mm	7326668/6	60cm + 4@10cm
28/6/6/8 x 6mm	7328668/6	60cm + 4@10cm
30/6/6/8 x 6mm	7330668/6	60cm + 4@10cm
20/8/8/8 x 8mm	7320888/8	60cm + 4@10cm
22/8/8/8 x 8mm	7322888/8	60cm + 4@10cm
24/8/8/8 x 8mm	7324888/8	60cm + 4@10cm
26/8/8/8 x 8mm	7326888/8	60cm + 4@10cm
28/8/8/8 x 8mm	7328888/8	60cm + 4@10cm
30/8/8/8 x 8mm	7330888/8	60cm + 4@10cm

RANGE: COSELLI THORACO ABDOMINAL		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20/10/10/8X8 mm	732010108/8S4	60cm + 4 @ 10 cm
22/10/10/8X8 mm	732210108/8S4	60cm + 4 @ 10 cm
24/10/10/8X8 mm	732410108/8S4	60cm + 4 @ 10 cm
26/10/10/8X8 mm	732610108/8S4	60cm + 4 @ 10 cm
28/10/10/8X8 mm	732810108/8S4	60cm + 4 @ 10 cm
30/10/10/8X8 mm	733010108/8S4	60cm + 4 @ 10 cm

RANGE: GELWEAVE™ AORTIC ARCH TRIFURCATE GRAFT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
24/12/8/8mm	73248812	15cm + 3@30cm

RANGE: GELWEAVE™ THORACIC ARCH GRAFT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22/12/8/8 x 10mm	732022CX4	20cm + 2@15cm + 2@30cm
24/12/8/8 x10mm	732024CX4	20cm + 2@15cm + 2@30cm
26/12/8/8 x 10mm	732026CX4	20cm + 2@15cm + 2@30cm
28/12/8/8 x 10mm	732028CX4	20cm + 2@15cm + 2@30cm
30/12/8/8 x 10mm	732030CX4	20cm + 2@15cm + 2@30cm
32/12/8/8 x 10mm	732032CX4	20cm + 2@15cm + 2@30cm
34/12/8/8 x 10mm	732034CX4	20cm + 2@15cm + 2@30cm

RANGE: GELWEAVE™ BRANCHED ARCH GRAFT WITH RADIOPAQUE MARKERS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22/12/8/8 x 10mm	732022CX4RM	20cm + 2@15cm + 2@30cm
24/12/8/8 x10mm	732024CX4RM	20cm + 2@15cm + 2@30cm
26/12/8/8 x 10mm	732026CX4RM	20cm + 2@15cm + 2@30cm
28/12/8/8 x 10mm	732028CX4RM	20cm + 2@15cm + 2@30cm
30/12/8/8 x 10mm	732030CX4RM	20cm + 2@15cm + 2@30cm
32/12/8/8 x 10mm	732032CX4RM	20cm + 2@15cm + 2@30cm
34/12/8/8 x 10mm	732034CX4RM	20cm + 2@15cm + 2@30cm

RANGE: GELWEAVE LUPIAE BRANCHED GRAFT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
22/10/10/8 x 10mm	734020CX4RMS	40cm + 40cm + 3@30cm
22/10/10/8 x 10mm	734022CX4RMS	40cm + 40cm + 3@30cm

Declaration of Conformity Record

24/10/10/8 x10mm	734024CX4RMS	40cm + 40cm + 3@30cm
26/10/10/8 x 10mm	734026CX4RMS	40cm + 40cm + 3@30cm
28/10/10/8 x 10mm	734028CX4RMS	40cm + 40cm + 3@30cm
30/10/10/8 x 10mm	734030CX4RMS	40cm + 40cm + 3@30cm
32/10/10/8 x 10mm	734032CX4RMS	40cm + 40cm + 3@30cm
34/10/10/8 x 10mm	734034CX4RMS	40cm + 40cm + 3@30cm

RANGE: GELWEAVE SIENA™ COLLARED ANTE-FLO		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20/8mm	734520/8ET	45cm X 1 @ 30cm
22/8mm	734522/8ET	45cm X 1 @ 30cm
24/8mm	734524/8ET	45cm X 1 @ 30cm
26/8mm	734526/8ET	45cm X 1 @ 30cm
28/8mm	734528/8ET	45cm X 1 @ 30cm
30/8mm	734530/8ET	45cm X 1 @ 30cm
32/8mm	734532/8ET	45cm X 1 @ 30cm
34/8mm	734534/8ET	45cm X 1 @ 30cm
20/10mm	734520/10ET	45cm X 1 @ 30cm
22/10mm	734522/10ET	45cm X 1 @ 30cm
24/10mm	734524/10ET	45cm X 1 @ 30cm
26/10mm	734526/10ET	45cm X 1 @ 30cm
28/10mm	734528/10ET	45cm X 1 @ 30cm
30/10mm	734530/10ET	45cm X 1 @ 30cm
32/10mm	734532/10ET	45cm X 1 @ 30cm
34/10mm	734534/10ET	45cm X 1 @ 30cm

RANGE: GELWEAVE SIENA™ COLLARED 4 BRANCH PLEXUS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20/10/8/8/8mm	73201088/8ET	30cm X 4 @ 15cm
22/10/8/8/8mm	73221088/8ET	30cm X 4 @ 15cm
24/10/8/8/8mm	73241088/8ET	30cm X 4 @ 15cm
26/10/8/8/8mm	73261088/8ET	30cm X 4 @ 15cm
28/10/8/8/8mm	73281088/8ET	30cm X 4 @ 15cm
30/10/8/8/8mm	73301088/8ET	30cm X 4 @ 15cm
32/10/8/8/8mm	73321088/8ET	30cm X 4 @ 15cm
34/10/8/8/8mm	73341088/8ET	30cm X 4 @ 15cm
20/10/8/8/10mm	73201088/10ET	30cm X 4 @ 15cm
22/10/8/8/10mm	73221088/10ET	30cm X 4 @ 15cm
24/10/8/8/10mm	73241088/10ET	30cm X 4 @ 15cm
26/10/8/8/10mm	73261088/10ET	30cm X 4 @ 15cm
28/10/8/8/10mm	73281088/10ET	30cm X 4 @ 15cm
30/10/8/8/10mm	73301088/10ET	30cm X 4 @ 15cm
32/10/8/8/10mm	73321088/10ET	30cm X 4 @ 15cm
34/10/8/8/10mm	73341088/10ET	30cm X 4 @ 15cm

RANGE: GELWEAVE SIENA™ COLLARED STRAIGHT WITH RADIOPAQUE MARKERS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20mm	735020RM	50cm
22mm	735022RM	50cm
24mm	735024RM	50cm
26mm	735026RM	50cm
28mm	735028RM	50cm
30mm	735030RM	50cm
32mm	735032RM	50cm
34mm	735034RM	50cm

RANGE: GELWEAVE SIENA™ COLLARED ANTE-FLO WITH RADIOPAQUE MARKERS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20mm	735020/8RM	50cm
22mm	735022/8RM	50cm
24mm	735024/8RM	50cm
26mm	735026/8RM	50cm

Declaration of Conformity Record

28mm	735028/8RM	50cm
30mm	735030/8RM	50cm
32mm	735032/8RM	50cm
34mm	735034/8RM	50cm
20mm	735020/10RM	50cm
22mm	735022/10RM	50cm
24mm	735024/10RM	50cm
26mm	735026/10RM	50cm
28mm	735028/10RM	50cm
30mm	735030/10RM	50cm
32mm	735032/10RM	50cm
34mm	735034/10RM	50cm

RANGE: GELWEAVE SIENA™ COLLARED 4 BRANCH PLEXUS WITH RADIOPAQUE MARKERS		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
20mm	73201088/8RM	50cm
22mm	73221088/8RM	50cm
24mm	73241088/8RM	50cm
26mm	73261088/8RM	50cm
28mm	73281088/8RM	50cm
30mm	73301088/8RM	50cm
32mm	73321088/8RM	50cm
34mm	73341088/8RM	50cm
20mm	73201088/10RM	50cm
22mm	73221088/10RM	50cm
24mm	73241088/10RM	50cm
26mm	73261088/10RM	50cm
28mm	73281088/10RM	50cm
30mm	73301088/10RM	50cm
32mm	73321088/10RM	50cm
34mm	73341088/10RM	50cm

RANGE: GELWEAVE VALSALVA		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
16mm	730016ADP	15cm
18mm	730018ADP	15cm
20mm	730020ADP	15cm
22mm	730022ADP	15cm
24 mm	730024ADP	15cm
26 mm	730026ADP	15cm
28 mm	730028ADP	15cm
30 mm	730030ADP	15cm
32 mm	730032ADP	15cm
34 mm	730034ADP	15cm

RANGE: GELWEAVE VALSALVA (VS115)		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
24 mm	730024ADPS	15cm
26 mm	730026ADPS	15cm
28 mm	730028ADPS	15cm
30 mm	730030ADPS	15cm
32 mm	730032ADPS	15cm
34 mm	730034ADPS	15cm

Declaration of Conformity Record

RANGE: GELWEAVE VALSALVA WITH SHORTER SKIRT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
24mm	730024ADP20	15cm
26mm	730026ADP22	15cm
28mm	30028ADP23	15cm
30mm	730030ADP24	15cm
32mm	730032ADP26	15cm
34mm	730034ADP28	15cm

RANGE: GELWEAVE VALSALVA WITH LONGER SKIRT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
24mm	730024ADP50	15cm
26mm	730026ADP50	15cm
28mm	730028ADP50	15cm
30mm	730030ADP50	15cm
32mm	730032ADP50	15cm
34mm	730034ADP50	15cm

RANGE: GELWEAVE VALSALVA WITH COLLAR		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
16mm	730016ADPS5	15cm
18mm	730018ADPS5	15cm
20mm	730020ADPS5	15cm
22mm	730022ADPS5	15cm
24 mm	730024ADPS5	15 cm
26 mm	730026ADPS5	15 cm
28 mm	730028ADPS5	15 cm
30 mm	730030ADPS5	15 cm
32 mm	730032ADPS5	15 cm
34 mm	730034ADPS5	15 cm

RANGE: GELWEAVE VALSALVA ANTE-FLO		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
24 x 8 mm	734024/8ADP	40cm + 1 @ 15cm
26 x 8 mm	734026/8ADP	40cm + 1 @ 15cm
28 x 8 mm	734028/8ADP	40cm + 1 @ 15cm
30 x 8 mm	734030/8ADP	40cm + 1 @ 15cm
32 x 8 mm	734032/8ADP	40cm + 1 @ 15cm
34 x 8 mm	734034/8ADP	40cm + 1 @ 15cm
24 x 10 mm	734024/10ADP	40cm + 1 @ 15cm
26 x 10 mm	734026/10ADP	40cm + 1 @ 15cm
28 x 10 mm	734028/10ADP	40cm + 1 @ 15cm
30 x 10 mm	734030/10ADP	40cm + 1 @ 15cm
32 x 10 mm	734032/10ADP	40cm + 1 @ 15cm
34 x 10 mm	734034/10ADP	40cm + 1 @ 15cm

RANGE: GELWEAVE VALSALVA ANTE-FLO WITH SHORTER SKIRT		
SIZE	CATALOGUE NUMBER	USABLE LENGTH
24 x 8 mm	734024/8ADP20	40cm + 1 @ 15cm
26 x 8 mm	734026/8ADP22	40cm + 1 @ 15cm
28 x 8 mm	734028/8ADP23	40cm + 1 @ 15cm
30 x 8 mm	734030/8ADP24	40cm + 1 @ 15cm
32 x 8 mm	734032/8ADP26	40cm + 1 @ 15cm
34 x 8 mm	734034/8ADP28	40cm + 1 @ 15cm
24 x 10 mm	734024/10ADP20	40cm + 1 @ 15cm
26 x 10 mm	734026/10ADP22	40cm + 1 @ 15cm
28 x 10 mm	734028/10ADP23	40cm + 1 @ 15cm
30 x 10 mm	734030/10ADP24	40cm + 1 @ 15cm
32 x 10 mm	734032/10ADP26	40cm + 1 @ 15cm
34 x 10 mm	734034/10ADP28	40cm + 1 @ 15cm

