
26183TZ

Electrode shape

Type of usage

Connector pin for unipolar coagulation



Ball end

Reusable

Bipolar

HOPKINS® Telescopes

Diameter 3.3 mm



26007BA



HOPKINS® Forward-Oblique Telescope 30°, enlarged view, diameter 3.3 mm, length 25 cm, autoclavable,

fiber optic light transmission incorporated, color code: red

Diameter 10 mm



26003BA




HOPKINS® Forward-Oblique Telescope 30°, enlarged view, diameter 10 mm, length 31 cm, autoclavable,

fiber optic light transmission incorporated, color code: red



26003 EC





26003EC ENDOCAMELEON® HOPKINS® Telescope,

diameter 10 mm Ø, length 31 cm, **autoclavable**, variable direction of view from 0° – 90° , adjustment knob for selecting the desired direction of view, fiber optic light transmission incorporated, color code: gold



533TVB

Adaptor, autoclavable, permits telescope changing under sterile conditions

It is recommended to check the suitability of the product for the intended procedure prior to use.

Fiber Optic Light Cables

for Cold Light Fountains



	Fiber Optic 495 NE	Light Cable with Straight Connector Fiber Optic Light Cable, diameter 4.8 mm, length 300 cm	
	495 NCS	Fiber Optic Light Cable, extremely heat-resistant, diameter 4.8 mm, length 250 cm	
	NEW 495 NCSC	Fiber Optic Light Cable, extremely heat-resistant, with safety locking device, diameter 4.8 mm, length 250 cm	
	495 NB	Fiber Optic Light Cable, diameter 4.8 mm, length 180 cm	
	495 ND	Fiber Optic Light Cable, diameter 3.5 mm, length 300 cm	
	495 NA	Fiber Optic Light Cable, diameter 3.5 mm, length 230 cm	
	NEW 495 NAC	Fiber Optic Light Cable, with safety locking device, diameter 3.5 mm, length 230 cm	
	495 NL	Fiber Optic Light Cable, diameter 3.5 mm, length 180 cm	
	495 NTA	Fiber Optic Light Cable, diameter 2.5 mm, length 230 cm	
	495 NT	Fiber Optic Light Cable, diameter 2.5 mm, length 180 cm	
	NEW 495 TIP	Fiber Optic Light Cable, highly heat resistant, diameter 4.8 mm, length 300 cm, optimized for 3D TIPCAM®1 and TIPCAM®	
	Fiber Optic 495 NWL	Light Cable with 90° Deflection to the Light Source Fiber Optic Light Cable, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 300 cm	
	495 NWM	Fiber Optic Light Cable, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 230 cm	
	495 NW	Fiber Optic Light Cable, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 180 cm	
1	495 NTX	Fiber Optic Light Cable, with 90° deflection to the cold light fountain on the fountain side, diameter 2.5 mm, length 230 cm	
	495 NTW	Fiber Optic Light Cable, with 90° deflection to the cold light fountain on the fountain side, diameter 2.5 mm, length 180 cm	
	Fiber Optic Light Cable with 90° Deflection to the Instrument		
	NEW 495 NVC	Fiber Optic Light Cable, with 90° deflection to the instrument, very narrow radius of curvature, diameter 4.8 mm, length 300 cm	
	495 NVB	Fiber Optic Light Cable, with 90° deflection to the instrument, diameter 4.8 mm, length 300 cm	
	495 NVL	Fiber Optic Light Cable, with 90° deflection to the instrument, diameter 3.5 mm, length 300 cm	
495 NVC 495 NV/NVL/NVB	495 NV	Fiber Optic Light Cable, with 90° deflection to the instrument, diameter 3.5 mm, length 230 cm	
4	Light Adapt	tor, angled 90°	
	490 EVV	Light Adaptor, angled 90°, diameter 4.8 mm, free rotatable, to connect with standard telescopes	

Adaptor for connecting KARL STORZ fiber optic light cables with endoscopes and light sources from other manufacturers, see page TP 132

3-012

ENDOMAT® SELECT VET





	031524-10*	Tubing Set for Irrigation, FC, sterile, for single use, package of 10, for use with HAMOU® ENDOMAT® SCB 26 3311 20-1 and ENDOMAT® SELECT	STERILE
	UP 007	Tubing Set, irrigation, FC, reusable, sterilizable	
For use with program	ns ART SA and	ART LA	
	031523-10*	Tubing Set for Irrigation, PC, sterile, for single use, package of 10, for use with HAMOU® ENDOMAT® SCB 26 331120-1 and ENDOMAT® SELECT	STERILE
6	UP 008	Tubing Set, irrigation, PC, reusable, sterilizable	
For use with program	SURG Direct	Suction	
99	030647-10*	Tubing Set for Suction, DS, sterile, for single use, package of 10, for use with ENDOMAT® SELECT	STERILE (
00	UP 009	Tubing Set, irrigation, DS, reusable, sterilizable	
For use with program	SURG Bottle	Suction	
	031647-10*	Tubing Set for Suction, BS, sterile, for single use, package of 10, for use with ENDOMAT® SELECT	STERILE (
00)	UP 010	Tubing Set, irrigation, BS, reusable, sterilizable	

Suction and Irrigation Systems

HAMOU® ENDOMAT® SCB



Components / Spare Parts

Catalog page

26331101-1 HAMOU® ENDOMAT® SCB

26331120-1 HAMOU® ENDOMAT® SCB

20090170 SCB Connecting Cable, length 100 cm
031523-10* Cassette Tubing Set PC, sterile, for single use

031524-10* Cassette Tubing Set FC, sterile, for single use





U 26

26331120

Spare Part



2027590 Mains Fuse, T 2.0 AL (SB), package of 10

Components / Spare Parts

Catalog page

26331142 Silicone Tubing Set, for suction, sterilizable



Spare Parts



27500 LUER-Lock Tube Connector, male, tube diameter 9 mm

59351111018 **LUER-Lock Connector,** male



20300180 Tubing Connector Set

2-15

