














	
Schiller nellcor module		ADP-10023A	Schiller Argus TM-7		Din 7P-DB9F
Schiller Masimo		ADP-10023B	Argus Pro, Argus LCM		Din 7P-DB9F
		ADP-10023C			
Schiller (Kontrolon)		ADP-10031A	DEFIGARD 5000; use Masimo LNOP spo2 sensor; or RSP-10031; DSP-10031		6/12 pins to Masimo LNOP
		ADP-10031B	DEFIGARD 5000; use Masimo LNCS spo2 sensors; or RSP-10032; DSP-10032		6/12 pins to Masimo LNCS
Schiller PB1000 masimo cable	2.310261	ADP-10037A	Schiller PB1000, metal connector 7 pins metal connector (60 degree) to masimo LNOP spo2 connector; use with masimo LNOP spo2 sensor		
Schiller Innomed SpO2 cable		ADP-10037B	redel 7 pins 40degree dual keyings Innomed innocare T; use with BCI spo2 sensor		
Schiller PB1000 masimo cable		ADP-10037C	Schiller PB1000, metal connector 7 pins metal connector to masimo LNCS spo2 connector; use with masimo LNCS spo2 sensor		
Siemens		ADP-10024A	Siemens Medical : Sirecust 700, Sirecust 630, Sirem Spo2 Module		Siemens 10P-DB9F
	3368433; 3375834	ADP-10024B	Siemens Medical: SC 6002XL, SC7000, SC8000, SC9000XL, use with Nellcor spo2 sensor or RSP-10041A etc		Siemens 7P-DB9F

Drager		ADP-10024C	Narkomed – Oxysat 2		12pins connector to DB9F
	M35370	ADP-10024D	Vitalert 1000, Cicero, Vitara, use with Nellcor spo2 sensor or RSP-10041A etc		Lemo 8P-DB9F
	MS18683 MS17330	ADP-10024E	compatible Drager MS18683 MS17330 spo2 cable; Nellcor Oximax module		Siemens7P-DB9F
	MS24303 MS17522	ADP-10024F	drager to masimo LNCS spo2 extension cable		Siemens7P-DB9F
	MS18680 MS17041	ADP-10024G	drager to masimo LNOP spo2 extension cable		Siemens7P-masimo 6P connector
Spacelabs	700-0020-0	ADP-10025A	Spacelab : 90351-D/Gm, 90465, 90466,90467, 90489 use with Nellcor spo2 sensor or RSP-10041A etc		HYP12P-DB9F
	700-0030-00	ADP-10025B	Spacelab : New modules 90651-08, 90496, Ultraview, use with Nellcor spo2 sensor or RSP-10041A etc		Rectangle 10P-DB9F
	700-0792-00	ADP-10025C	Nellcor Oximax cable use with Nellcor Oximax spo2 probes;		Rectangle 10P-DB9F
		ADP-10025D	masimo cable use with masim spo2 probes;		Rectangle 10P-DB9F
Omnimed		ADP-10026	Omni 200		AMP14P-DB9F