

# ordering information

One of the most complete coil offers on the market

## Prestige Plus Complex Xtra Pack

Model No.	Implant size
PRES0202CXPPPL	2mm x 2cm
PRES0204CXPPPL	2mm x 4cm
PRES0206CXPPPL	2mm x 6cm
PRES0307CXPPPL	3mm x 7cm
PRES0310CXPPPL	3mm x 10cm

## Prestige Plus Complex Pack

Model No.	Implant size
PRES0204HXPPL	2mm x 4cm
PRES0206HXPPL	2mm x 6cm
PRES0208HXPPL	2mm x 8cm
PRES0310HXPPL	3mm x 10cm
PRES0315HXPPL	3mm x 15cm

## Prestige Plus Helical Xtra Pack

Model No.	Implant size
PRES1240COM18	12mm x 40cm
PRES1447COM18	14mm x 47cm
PRES1655COM18	16mm x 55cm
PRES1860COM18	18mm x 60cm
PRES2065COM18	20mm x 65cm
PRES2465COM18	24mm x 65cm

## Prestige Plus Helical Pack

Model No.	Implant size
PRES0407CPKPL	4mm x 7cm
PRES0410CPKPL	4mm x 10cm
PRES0413CPKPL	4mm x 13cm
PRES0510CPKPL	5mm x 10cm
PRES0517CPKPL	5mm x 17cm
PRES0620CPKPL	6mm x 20cm
PRES0720CPKPL	7mm x 20cm
PRES0830CPKPL	8mm x 30cm
PRES0930CPKPL	9mm x 30cm
PRES1030CPKPL	10mm x 30cm
PRES1040CPKPL	10mm x 40cm

Model No.	Implant size
PRES0410HPKPL	4mm x 10cm
PRES0420HPKPL	4mm x 20cm
PRES0515HPKPL	5mm x 15cm
PRES0520HPKPL	5mm x 20cm
PRES0620HPKPL	6mm x 20cm
PRES0720HPKPL	7mm x 20cm
PRES0830HPKPL	8mm x 30cm
PRES0930HPKPL	9mm x 30cm
PRES1040HPKPL	10mm x 40cm

## ULTRA Detachment Controller

Model No.	Box content
ULTRA	1 detachment controller

\*Data on file

The Prestige Coil System is intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/EEC amended by 2007/47/EC). Manufactured by BALT USA LLC, 29 Parker, Irvine, CA 92618. Carefully read the instructions for use before use. First CE marking: 2021. For a list of indications, warnings, precautions and contraindications, refer to IFU-060.

The content of this document, in particular data, information, trademarks and logos is BALT S.A.S and affiliates' sole property. Consequently, all representation and/ or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT S.A.S and affiliates' copyrights and other intellectual proprietary rights ©2021 BALT S.A.S and affiliates all rights reserved. This document with associated pictures is non-contractual and is solely dedicated to healthcare professionals and BALT S.A.S and affiliates' distributors. The products commercialized by BALT S.A.S and affiliates shall exclusively be used in accordance with the package inserts included in the boxes. DC057GB (05/2021)

Balt  
10, rue de la Croix Vigneron, 95160 Montmorency, France  
Tél. : +33 (0)1 39 89 46 41  
Fax : +33 (0)1 34 17 03 46  
www.balt-corp.com




# prestige plus

volume and versatility



# prestige plus

Intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels.



## Design

**Premium coil softness** through variations of filar and primary wind for deliverability and conformability

Body coil/hypotube pusher design provides **increased ease of delivery**, even in tortuous anatomy

**Smooth, precise embolization** thanks to implant softness and enhanced pusher technology

## Optimal detachment

**Instantaneous detachment** in under 1 second\* allows for rapid treatment & efficient case times

**Reliable detachment\*** system provides confidence during coiling cases

## Complete offer

With a full range of Complex and Helical coils to accommodate a variety of clinical application

**Unique large volume coils\*** (Complex 18 and Complex Pack) for dense packing



entire size range  
deliverable through

.0165"-.021"/.022" catheters


peripheral treatment

## Implant design

- ✓ Platinum, non fibered coil → fully visible
- ✓ Complete set of configurations: Helical & Complex
- ✓ Deliverable through microcatheters with ID = 0.0165" - 0.027"
- ✓ Second radiopaque mark: possibility to use it with double or single mark microcatheters.
- ✓ New 185 cm total length.



helical




complex

The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. ©2018 BALT and affiliates, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The distribution, the sale and the use of BALT's products is under distributor's liability.

EXTERNAL USE -CONFIDENTIAL



# Pusher design




The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. ©2018 BALT and affiliates, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The distribution, the sale and the use of BALT's products is under distributor's liability.

EXTERNAL USE -CONFIDENTIAL



## Components



**Prestige coil system**

To **use with** (not included in the package):

- ✓ ULTRA Detachment Controller
- ✓ Compatible guidewire
- ✓ Delivery microcatheter I.D. between .0165" and .027"
- ✓ Compatible sheath / guiding catheter

The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. ©2018 BALT and affiliates, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The distribution, the sale and the use of BALT's products is under distributor's liability.

EXTERNAL USE –CONFIDENTIAL

