





Home ► Shop All Products ► PCR Equipment And Supplies ► Thermal Cycler Accessories ►

Region Of Interest (ROI) And Bac...







Q Certificates



Applied Biosystems™

Region of Interest (ROI) and Background Plates, 96-well

These Region-of-Interest (ROI) and Background Plates are used to maintain the ViiA™ 7, QuantStudio™ 3, QuantStudio™ 5, QuantStudio™ 6, Read more

Have Questions? Contact Us

Catalog Number	Quantity
4432364	2 x 96-well plates

Catalog number 4432364 ☆ e (MDL)



Q

These Region-of-Interest (ROI) and Background Plates are used to maintain the ViiA™ 7, QuantStudio™ 3, QuantStudio™

Product Overview

FAQ

Documents

For Research Use Only. Not for use in diagnostic procedures.

Specifications	
Calibration Type	Spectral Calibration
For Use With (Application)	Real Time PCR (qPCR)
For Use With (Equipment)	ViiA™ 7 System
Quantity	2 x 96-well plates
Shipping Condition	Dry Ice
Product Type	Background Plate
Unit Size	Each

Contents & Storage

Store at -15°C to -25°C

Frequently asked questions (FAQs)

I want to purchase the ROI/Uniformity, background calibration and RNase P verification plates for the QuantStudio 5 Real Time PCR System, 96-well, 0.2 ml (Cat. No. A28139). What are the Catalog numbers for these plates?

^

Catalog numbers are as follows:



Q

visit the following webpage.

I cannot open the QuantStudio 3/5 Design & Analysis 2 (DA2) software on my Mac. What are the software Find additional tips, troubleshooting help, and resources within our Real-Time PCR and Digital PCR Instruments requirements for Mac?

Support Center.

^

To open the software download on your Mac, you will need to adjust your security settings using the following instructions from Apple support:

Safely open apps on your mac

Please note, the software should run on Mac OS 10.10 or later.

Find additional tips, troubleshooting help, and resources within our QuantStudio 3 and QuantStudio 5 Real Time What is the expiration date for the Region of Interest (ROI) and Background Plates, 96-well (Cat. No. 4432364)? PCR Systems Support Center.

The expiration date for the Region of Interest (ROI) and Background Plates, 96-well (Cat. No. 4432364) is 12 months from the date of manufacturer, which can be found on the Certificate of Analysis (COA).

You can find the COA for a specific lot number by using the Documents and Certificatessearch page or by filling out a Certificate of Analysis (COA) or Certificate of Origin (COO) Request Form.

Find additional disammulty of the ingular and resources within our landstudio 6 and 74 few find plants one performed within our landstudio 6 and 74 few find plants one systems of the contractions.

^

You can get a summary of the calibrations from the QuantStudio Software by going to Instrument > Instrument Console > Manage Instrument > Information > Print Calibration Log.

Eind additional tips, troubleshooting help, and resources within our Real-Time PCR and Digital PCR Instruments Do I need to back up my calibration files before I upgrade the software for the ViiA7 Real-Time PCR System? Support Center.

^

No, you do not need to back up the calibration files as all of the calibration files are stored on the ViiA7 Real-Time PCR instrument. Once the computer software is properly installed and the computer is connected to the instrument, the calibration information will automatically update in the QuantStudio Software.

Find additional tips, troubleshooting help, and resources within our Real-Time PCR and Digital PCR Instruments Support Center.


View all Region of Interest (ROI) and Background Plates, 96-well FAQs >

Documents & Downloads

Certificates

Search by lot number or partial lot number

Search



5 results displayed, search above for a specific certificate

Request a Certificate

Safety Data Sheets



SDS

Product Information

Manuals

Product Information Sheet: Region of Interest (ROI) and Background Calibration Plates