

BLOOD BAG WITH LEUCOCYTE FILTER CPDA-1 SOLUTION		
STERILE • NON PYROGENIC • NON TOXIC For collection of 450ml of Blood. 63ml of Anticoagulant Citrate Phosphate Dextrose Adenine Solution (USP) Each 100ml of CPDA-1 Contains: Citric Acid monohydrate 0,327g Sodium Citrate dihydrate 2,63g Sodium Biphosphate monohydrate 0,222g Dextrose monohydrate 3,19g Adenine 0,0275g Water for injection q,s	Donor Number (Name)	
	<table border="1"> <tr> <td>Blood Group</td> <td>Rh - Type</td> </tr> </table>	Blood Group
Blood Group	Rh - Type	
CAUTION 1. Do not use unless the anticoagulant is clear. 2. Store blood between + 2°C and + 6°C. 3. Crossmatch before transfusion. 4. Mix blood thoroughly immediately before use. 5. Do not add medication to this blood. 6. Transfusion set must have a filter.	Collection Date	
	Do not infuse after	
DO NOT VENT	Serology: non-reactive to	
Lot: Exp: Mfg:	Test	
Shandong Weigao Group Medical Polymer Co.,Ltd. Add:No.312 Shichang Road, Weihai, Shandong, China.		

95

95

Collection Date	Donation Number	EXPIRES
Place Product Label here after bag is filled.		Store Platelets 3 days at 20°C - 24°C Store Red Blood Cells 35 days at 2°C - 6°C
DO NOT VENT		
STERILE, NON-PYROGENIC FLUID PATH This product may transmit infectious agents. PROPERLY IDENTIFY INTENDED RECIPIENT		LOT NO.:

95

75