

non-woven IV cannula dressing elastopor IV

features

- centrally situated absorbent pad protects the non-woven material from sticking to the injection site
- non-woven lining should be placed under cannula’ s wings, which prevents it from causing damage to the skin

indication

- dressing application around the inserted cannula
- stabilization, protection and fixation of intravenous cannulas



Ref No.	dimensions	Packaging unit paper-paper	Bulk packaging box
811013	6 cm x 8 cm	1 pc	100 x 1 pc
811012	5.1 cm x 7.6 cm	1 pc	100 x 1 pc
811023	6 cm x 8 cm	1 pc	50 x 1 pc
811022	5.1 cm x 7.6 cm	1 pc	50 x 1 pc