



EVOLUTIONARY PROCESS, REVOLUTIONARY SOLUTION

Cantel is dedicated to developing and providing innovative solutions in infection prevention and control worldwide.

The ISA™ Endoscope Reprocessor represents the unique new generation response. It is a rapid, efficient and flexible reprocessing solution designed for endoscopy units demanding the most advanced technology and the highest quality standards.





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

REPROCESS High-level disinfection is the cornerstone of infection prevention. Reprocessing technologies from Cantel are designed to ensure patient safety by adhering to the strictest of standards for high level disinfection, optimize workflow efficiency, and deliver versatility by supporting a wide range of scopes.

ISA™ Endoscope Reprocessor The Fast And Flexible Reprocessing Solution

FAST & SIMPLE OPERATION FOR IMPROVED EFFICIENCY

- Full cycle in 20 minutes with channel monitoring
- Simple, user-friendly handling and operation
- Insertion of endoscopes and operation information through RFID system – helps to avoid errors and saves time




COMPACT & ERGONOMIC DESIGN FOR MAXIMUM USER COMFORT

- Large and easy to load basin accommodates all types of flexible endoscopes
- Intuitive interface assures easy and fast access to all the cycle information
- Compact size contributes to easy installation, operation and maintenance



ENSURING SAFETY FOR OPERATORS, PATIENTS & ENDOSCOPES

- Hands free operation and RFID-based technology for complete infection control
- Full traceability through printed and electronic cycle documentation
- Validated process fully compliant to the standards UNI EN ISO 15883-1/4 and UNI CEN ISO/TS 15883-5 (including Austrian, German and French test soil)



THE RIGHT FORMULA FOR THE BEST PERFORMANCE & RESULT

- Dedicated room temperature chemicals compatible with all types of flexible endoscopes
- Improved water and chemical consumption
- Designed specifically for use with ISA Endoscope Reprocessor



PRODUCT SPECIFICATIONS ISA™ ENDOSCOPE REPROCESSOR





Dimensions

700 mm (I) x 1400 mm (h) x 650 mm (d) - monitor included

Weight

75 kg

Electrical Requirements

100 V-240 V, 50-60Hz

Operating Room Temperature

5° - 40° C

Humidity

20 - 80%

Reprocessing Time

20 minutes (full cycle), 12 minutes only (HLD cycle)

Process Capacity

1 endoscope per cycle

High Level Disinfectant/Sterilant

ISASPOR™ Single Shot Disinfectant/Sterilant, two parts, Peracetic acid and ISAZONE™ Activator. RAPICIDE™ PA Single Shot Disinfectant/Sterilant, two part, Peracetic acid.

Detergent

ISACLEAN™ Detergent, multienzymatic detergent INTERCEPT™ PLUS Detergent, Non-Enzymatic detergent

ISA™ ENDOSCOPE REPROCESSOR ORDERING INFORMATION		
PART NUMBER	DESCRIPTION	QUANTITY
500187	ISASPOR™ Single Shot High-Level Disinfectant/Sterilant	1 x 10L Sol. A+B
500188	ISASPOR Single Shot High-Level Disinfectant/Sterilant	1 x 5L Sol. A+B
ML02-0158	RAPICIDE™ PA Single Shot High-Level Disinfectant/Sterilant	1 x 5L A+B
ML02-0115	RAPICIDE PA Single Shot High-Level Disinfectant/Sterilant	2 x 5L AA+BB
500189	ISACLEAN™ Detergent	2 x 5L
ML02-0145	INTERCEPT™ PLUS Detergent	2 x 5L
93250-222	ISA™ CPC Filters	0,45µm + 0,1µm
ISA4400	ISA™ Air Filter	1
ISAROTS2	ISA™ Printer Roll	1

ISATM, ISASPORTM, ISAZONETM and ISACLEANTM are trademarks of Cantel Medical (Italy) S.r.l. RAPICIDETM and INTERCEPTTM are trademarks of Medivators Inc.

MEDIVATORS IS NOW CANTEL

www.cantelmedical.eu

Medivators BV | p: +31 45 5 471 471 (customer service)





REGISTRUL DE STAT AL DISPOZITIVELOR MEDICALE

