

Safe AQ Angel

Blood Glucose Meter



- FAD-GDH System
- No Coding
- 5s testing time
- 0.6µl sample volume
- Virtually pain free

Safe AQ Angel

Blood Glucose Meter

SPECIFICATION

FAD-GDH System

Blood sample : fresh capillary whole blood;
venous whole blood

Calibration sample : venous plasma

Test range : 1.1-33.3mmol/L(20-600mg/dL)

Test time : 5s

Blood volume : 0.6 μ l


HCT : 30-60%

Memory : 200 test results

Expiry date : 24 months

Test enviroment : 10°C~35°C

Battery lifetime : Up to 1000 tests



1 Insert Safe AQ Angel blood glucose test strip into the meter until it stops



2 Apply the blood sample



3 Display the test result after a 5-second count down

Sinocare Inc.

Add: No.265 Guyuan Road Hi-tech Zone, Changsha, Hunan, China, 410205
Tel: 86-731-89935929 Fax: 86-731-89825189
Email: info@sinocare.com
www.sinocare.com

Sinocare

• Caring for love