

**Product Name: One-piece Ostomy bag, with the opening for drainage (emptying)**  
**Item No.: A1O2270**



**Purpose of use**

- It is used for stoma cleaning, nursing, excreta collection and skin care around the stoma.
- It is a product used to stick on the stoma to collect excreta.
- The product is in good contact with skin and intestinal cavity. (Non-sterile ).

**Specifications:**

Bag size: 29\*14,5 cm;

Cutting size: 15-75 mm;

Color: Skin color;

Filter type: Built-in round type activated carbon filter, for gas exhaust, to eliminate completely the odors;

Flange (sucker): Flexible, is an integral part of a bag, with special hypoallergenic hydrocolloid adhesive;

Drainage closure mechanism: With adhesive tape;

**Key benefits:**

- Made of layers of high-strength raw materials
- Provides safety and protection against odors
- Hypoallergenic hydrocolloid adhesive support
- Secure adhesion
- Maximum skin protection

Address:No. 599 Haijing Rd. ,Longgang City, Wenzhou, Zhejiang, China.

Email:strong@zjailebao.com

Mobile Phone:+0086 18868272121

## **Classification and Conformity Routes**

### **Classification**

According to the Annex VIII of MDR, Rules 1, the product is subject to: Class I

### **Conformity Routes**

Conformity as per terms of Regulation (EU) 2017/745 Annex IV (Declaration of Conformity)

### **Certification**

- ISO 13485:2016
- CE
- FDA, CFDA

Address: No. 599 Haijing Rd., Longgang City, Wenzhou, Zhejiang, China.

Email: [strong@zjailebao.com](mailto:strong@zjailebao.com)

Mobile Phone: +0086 18868272121

## LABELING:

**QTY: 10 pieces/bag**



### **Ostomy Collection Bags UDI-DI:**

Model No./Specification : A1O2270

Product : One piece system open


### **Precautions:**

1. Patients with purulent stoma should be used under the guidance of doctors or professional stoma specialist.
2. Valid for 3 years, please use within the validity period.
3. The packaging is damaged and it is forbidden to use it.
4. Single-use, destroy after use.



Zhejiang Ailebao Medical Technology Co.,Ltd.

2nd Floor, No.192 Xingxian Road, Longgang City, Wenzhou City, Zhejiang



Name: MedPath GmbH.

Mies-van-der-Rohe-Steasse 8,80807 Munich, Germany

Tel:+49(0)89 189174474

Address:No. 599 Haijing Rd. , Longgang City, Wenzhou, Zhejiang, China.

Email:strong@zjailebao.com

Mobile Phone:+0086 18868272121

**Product Name: One-piece Urostomy bag**

**Item No.: A1U5145**



**Purpose of use**

- It is used for stoma cleaning, nursing, excreta collection and skin care around the stoma.
- It is a product used to stick on the stoma to collect excreta.
- The product is in good contact with skin and intestinal cavity. (Non-sterile ).

**Specifications:**

Bag size: 160\*242cm; (Excluding urine mouth)

Cutting size: 13-45mm;

Color: White color;

From a single component (one-piece). The adhesive part made of the hydrocolloid.

Collector bag with built-in anti-reflux valve. Possibility of emptying. With non-woven material on the back.

Non-allergic.;

**Key benefits:**

- Made of layers of high-strength raw materials
- Hypoallergenic hydrocolloid adhesive support
- Secure adhesion
- Maximum skin protection

Address: No. 599 Haijing Rd., Longgang City, Wenzhou, Zhejiang, China.

Email: strong@zjailebao.com

Mobile Phone: +0086 18868272121

# AILEBAO Zhejiang Ailebao Medical Technology CO., LTD

## **Classification and Conformity Routes**

### **Classification**

According to the Annex VIII of MDR, Rules 1, the product is subject to: Class I

### **Conformity Routes**

Conformity as per terms of Regulation (EU) 2017/745 Annex IV (Declaration of Conformity)

### **Certification**

- ISO 13485:2016
- CE
- FDA, CFDA

Address: No. 599 Haijing Rd., Longgang City, Wenzhou, Zhejiang, China.

Email: strong@zjailebao.com

Mobile Phone: +0086 18868272121

**LABELING:**

**QTY: 10 pieces/bag**



**Urostomy Collection Bags UDI-DI:**

Model No./Specification : A1U5145

Product : One piece


**Precautions:**

1. Patients with purulent stoma should be used under the guidance of doctors or professional stoma specialist.
2. Valid for 3 years, please use within the validity period.
3. The packaging is damaged and it is forbidden to use it.
4. Single-use, destroy after use.



Zhejiang Ailebao Medical Technology Co.,Ltd.

2nd Floor, No.192 Xingxian Road, Longgang City, Wenzhou City, Zhejiang



Name: MedPath GmbH

Mies-van-der-Rohe-Steasse 8,80807 Munich, Germany

Tel:+49(0)89 189174474

Address: No. 599 Haijing Rd., Longgang City, Wenzhou, Zhejiang, China.

Email: strong@zjailebao.com

Mobile Phone: +0086 18868272121