


Sterling™

PTA Balloon Dilatation Catheter



Demand

from your 0.018" / 0.46 mm Balloon Dilatation Catheter

More

The #1 0.018" / 0.46 mm** PTA Balloon Dilatation Catheter

with the most comprehensive
size offering available

2-10 mm diameters 10-220 mm lengths*

The most comprehensive
size offering available

10 seconds*

Best-in-class
deflation time

0.020" / 0.51 mm*

Lowest in-class
lesion entry profile

The Sterling™ family† is designed to offer best-in-class lesion entry
profile and deflation times with the most clinically comprehensive
size offering and indications available.



Demand Exceptional Crossability

- Exceptionally low
0.020" / 0.51 mm lesion entry profile
on a 0.018" / 0.46 mm wire
- 7 x 220 through 5 F / 0.65 mm




Demand the Broadest Size Matrix

- The most comprehensive size
offering available in both
Over-The-Wire and Monorail™
rapid-exchange
- Indicated for wider acute use
including carotid‡
and post stent dilatation

* Average measurements taken by Boston Scientific (n=3, 5 x 150 mm balloons). Data on file.
Bench testing may not be representative of clinical performance.
** US and EU market share data per internal and external analysis.
† The Sterling family includes Sterling and Sterling SL.
‡ Monorail only


Sterling™ Monorail™

PTA Balloon Dilatation Catheter



Sterling™ Over-the-Wire

PTA Balloon Dilatation Catheter



Diameter (mm)	Length (mm)												
	Balloon	10	15	20	30	40	60	80	100	120	150	200	220
2.0 mm													
2.5 mm													
3.0 mm													
3.5 mm													
4.0 mm													
4.5 mm													
5.0 mm													
5.5 mm													
6.0 mm													
6.5 mm													
7.0 mm													
8.0 mm													
9.0 mm													
10.0 mm													

Platform	Shaft Length (cm)			
	40	80/90	135	150
OTW				
Rapid Exchange				

Boston Scientific

Advancing science for life™

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome.

PI-195601-AB OCT2014 Printed in Germany by medicalvision.

www.bostonscientific-international.com

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved. DINPER4642EB