



Test report No. 10619hd

EVALUATION OF FUNGICIDAL OR YEASTICIDAL ACTIVITY IN THE MEDICAL AREA (EN 13624)

Name of the product: CHEMISEPT GEL

Batch number: 198030619

Date of test report: 24.07.2019

Client, representative:

Chemi-Pharm Ltd.

Tänassilma tee 11

Tänassilma küla

Saku vald, 76406

ESTONIA

Test report No. 10619hd

EVALUATION OF FUNGICIDAL OR YEASTICIDAL ACTIVITY IN THE MEDICAL AREA
(EN 13624)

Name of the product: CHEMISEPT GEL
Batch number: 198030619
Order number: 19023
Manufacturer: Chemi-Pharm Ltd.
Client, representative: Chemi-Pharm Ltd., Tännassilma tee 11, Tännassilma küla, Saku vald,
76406, ESTONIA
Maris Millner, +372-51-77-090
Date of delivery: 03.07.2019
Test material conditions: No specific features, sample in the manufacturers tare
Storage conditions: In room temperature, dark
Active substance – conc.: Ethyl alcohol 72.5 (w/w) %; isopropyl alcohol (w/w) 7.5%
Appearance of the product: Transparent gel
Test concentration: Ready to use
Contact time: 60s
Interfering substance: 3.0 g/l bovine albumin = clean conditions
Rinsing liquid: Tryptone 1 g/l + NaCl 9 g/l
Neutralizer: -
Test organisms: *Aspergillus brasiliensis* ATCC 16404
Testing method: EVS-EN 13624:2013
Quantitative suspension test for the evaluation of fungicidal or
yeastcidal activity in the medical area.
Testing date: 16.07.2019 – 18.07.2019
Results: look appendix 1-2



Allar Laaneleht
Chief specialist

Date of test report: 24.07.2019

Appendix 1

TEST RESULTS (fungicidal suspension test)

EVS-EN 13624:2013; Phase 2, step 1;
Membrane filtration method;
Rinsing liquid: Tryptone 1 g/l + NaCl 9 g/l;
Test organism: *Aspergillus brasiliensis* ATCC 16404;
Test temperature: +20° C; Incubation temperature: +30° C
Interfering substance: 3.0 g/l bovine albumin = clean conditions;
Nordic Tersus Laboratory LLC.; Date of test: 16.07.2019
Responsible person: Allar Laaneleht

Validation and controls

Clean conditions

Validation suspension N_{vo}			Experimental conditions (A)			Filtration control (B)			Method validation (C)		
V_{C1}	V_{C2}	\bar{x}	V_{C1}	V_{C2}	\bar{x}	V_{C1}	V_{C2}	\bar{x}	V_{C1}	V_{C2}	\bar{x}
61	74	67.5	33	39	36	42	30	36	41	49	45
$30 \leq \bar{x} N_{vo} \leq 160$? yes X; no <input type="checkbox"/>			$\bar{x} A \geq 0.5 \bar{x} N_{vo}$? yes X; no <input type="checkbox"/>			$\bar{x} B \geq 0.5 \bar{x} N_{vo}$? yes X; no <input type="checkbox"/>			$\bar{x} C \geq 0.5 \bar{x} N_{vo}$? yes X; no <input type="checkbox"/>		

Test suspension and test

Testsuspension: N and N_0	N	V_{C1}	V_{C2}	$\bar{x}_{wm} = 1.79 \times 10^8$; $\log N = 8.25$ $N_0 = N/100$; $\log N_0 = 6.25$ $6.17 \leq \log N_0 \leq 6.70$; yes X; no <input type="checkbox"/>
	10^{-6}	194	164	
	10^{-7}	21	15	

Experimental results

Concentration of the product	Dilution step	V_{C1}	V_{C2}	Na ($=\bar{x} \cdot 10$)	$\log Na$	$\log R$	Contact time	Conditions
Ready to use	-	17	15	160	2.20	4.05	60 s	Clean

Explanations:

V_C = count per ml (one plate or more)
 \bar{x} = average of V_{C1} and V_{C2} (1. + 2. duplicate)
 N = cfu/ml microbes in testsuspension
 N_0 = cfu/ml at the start of the contact time ($t=0$)
 N_{vo} = cfu/ml in the validation suspension ($t=0$)
 Na = surviving microbes after the test
 R = reduction factor ($R = N_0 / Na$; $\log R = \log N_0 - \log Na$)

Appendix 2

Interpretation:

The ready to use product CHEMISEPT GEL (batch no. 198030619) was tested according to the test method EVS-EN 13624:2013. The test was performed at $20\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$, under clean conditions during contact time of 60s. The membrane filtration method was used for testing the product's effectiveness against the reference strain: *Aspergillus brasiliensis* ATCC 16404. Under clean conditions the tested product was effective against the reference strain within 60 s.

Conclusion:

The surviving count of the reference strain showed at least 4 lg reduction meaning that the ready to use product CHEMISEPT GEL has a fungicidal effect under clean conditions within 60 s.



Allar Laaneleht
Chief specialist
24.07.2019