

# **TECHNICAL DATA SHEET**

## STERILIZATION FLAT/GUSSET REEL, 60gsm

PRODUCT CODE	:	4AFR (Flat Reel) - 4AGR (Gusset Reel)
RELATED STANDARDS	:	4A Medical flat and gusset reels are in compliance with the provisions of the Medical Device Regulation 2017/745 for medical devices.
INTENDED USE	:	4A Medical reels are designed to provide excellent barrier for sterile medical device packaging purposes. Reels meet the Medical Industry's requirement for high quality, hygiene and safety levels. 4A Medical reels are used for the packing of the Medical Devices before sterilization and it maintains the sterility of the products after sterilization. 4A Medical reels are ease to package and indicator used for the sterilization provides information with the sterilization status of the medical devices. Sterilized material in sterilization package is storable for up to 5 years after the sterilization in a clean room conditions. Sealing temparature is 180-200C.
FILLING VOLUME	:	According to DIN 58953-7 (German national standard): Seleable paper bags and seleable pouches and reals shall be filled up to 3/4 max of their volume. This limits avoids too high stress / burdon to the seals. Belowe the buttom seal there must be space of minimum 30 mm between load and sealing area. There is no known limit bu standarts for weight and it is based on experience. However, it is suggested to follow the 3/4 rule. In terms of solid, heavy load (metal instruments), it is the weight itself, in case of any sharp objects which could demage the material, and the volume of condensate (directly related to mass=weight) which is created during the steam sterilization.
OPERATION CONDITIONS	:	4A Medical reels are used for steam, ethylene oxide, formaldehyde sterilization methods. The sterilization conditions should be determined by the end user regarding to material to be sterilized.

# TEHNICALDATA SHEET

	PROPERTIES	UNITS	STANDARD	TYPICAL VALUE
SUBSTANCE	3	g/m²	ISO 536	60,0
THICKNESS		μm	ISO 534	76
BENDTSEN	POROSITY	ml/min	ISO 5636-3	650
AIR PERMEA	ANCE	μm (Pa*s)	ISO 5636-3	7,4
BENDTSEN	ROUGHNESS FS	ml/min	ISO 8791-2	250
BENDTSEN	ROUGHNESS WS	ml/min	ISO 8791-2	250
PORE SIZE		μm	EN 868-2 (app.D)	17
TENSILE ST TENSILE ST WET TENSII WET TENSII	RENGTH, MD	kN/m	EN ISO 1924-2	5,9
TENSILE ST	RENGTH, CD	kN/m	EN ISO 1924-2	2,9
WET TENSII	LE STRENGTH, MD	kN/m	ISO 3781	1,6
WET TENSII	LE STRENGTH, CD	kN/m	ISO 3781	0,9
BURST STR	ENGHT	kPa	ISO 2758	300
WET BURST	<b>1</b>	kPa	ISO 3689	100
TEARING ST	TRENGHT MD	mN	ISO 1974	600
TEARING ST	RENGHT CD	mN	ISO 1974	650
COBB TEST	(60s)	g/m²	ISO 535	15,0
WATER REP	ELLENCY	S	EN 868-2 (app.E)	30,0
FLOUORESC	CENCE	pts/dm <sup>2</sup>	EN 868-2 (app.B)	0,0



## **TECHNICAL DATA SHEET**

### STERILIZATION FLAT/GUSSET REEL, 60gsm

	PROPERTIES		UNIT	TYPICAL VALUE		
LAMINATED FILM				LAM.FILM	$52\pm5$	
	THICKNESS		Micron	PE	$40\pm5$	
				PET	$12\pm 2$	
	BREAKING FACTOR		N/15mm	35		
LAN	TEAR STRENGTH MD/CD		mN	300		
'WW	BOND STRENGHT		N/m	265		
F	ELONGATION		%	70		
	TEMPERATURE DURABILIT	Υ	min/degrees	30 min/140C		
	FREE FROM LEAD AND HEAVY METALS AND TOXIC MATERIALS					
PRODUCT SPECIFICATION :						
	PROPERTIES	UNIT		TYPICAL VALUE		
	SEAL STRENGTH		N/15mm	2,5		
POUCH	BUBBLE TEST			None		
	PINHOLE DETERMINATION			None		
	DIMENSION CONTROL		cm or mm	Desired dimensions		
	LEAKAGE TEST			None		

E

INDICATOR	STERILIZATION METHOD	BEFORE STERILIZATION	AFTER STERILIZATION
	STEAM	BLUE	GREEN
	ЕО	PINK	ORANGE
	FORMALDEHYDE	PURPLE	GREEN
		4A Medical reels are packaged in carton boxes as below: Flat Reels are 200m long and Gusset Reels are 100m long.	

Product dimensions	Pieces in carton
50mm X 200m / 50mm X 100m	10
75mm X 200m / 75mm X 100m	8
100mm X 200m / 100mm X 100m	8
150mm X 200m / 150mm X 100m	4
200mm X 200m / 200mm X 100m	4
250mm X 200m / 250mm X 100m	2
300mm X 200m / 300mm X 100m	2
350mm X 200m / 350mm X 100m	2
400mm X 200m / 400mm X 100m	2
450mm X 200m / 450mm X 100m	2
500mm X 200m /500mm X 100m	2

None Returns the

specified color

MANIFACTURER

:

PEEL DIRECTION

INDICATOR CONTROL

**DORT-A TIP MALZEMELERI SAN. ITH. IHR. TIC. LTD. STI.** Balıkhisar Mahallesi, Köyiçi Serpmeleri, No:795 Akyurt, Ankara, Turkey Tel: +90 312 363 50 52 / 53 Fax: +90 312 363 52 10 Web: http://www.4amedical.com E-mail: export@4amedical.com