



Steelco detergent and disinfection range

Steelco Alkal

Powerful detergent for difficult to clean soiling

Steelco Alkal is a powerful alkaline cleaner suitable for heavily soiled instruments that are challenging to clean:



- Suitable for stainless steel and glass
- Effective in removing electrosurgery deposits, dried blood and fatty tissue
- Used in combination with Steelco Acid

Steelco Acid

Neutralizing agent and descaler

Steelco Acid is a mild neutralizing agent designed to neutralize alkaline residues following the main wash phase of a washer disinfectant:



- Dosed after using an alkaline detergent
- Mild acid compatible with many materials
- Combats the formation of limescale
- Suitable as an acid detergent in laboratory applications

Rinse aid- lubricant

Steelco RL

Rinse aid and lubricant in one

Steelco RL is a rinse aid added to the final rinse to speed up drying times by reducing the surface tension of water for a stain free shiny finish.



- Reduces water stains
- Compatible with most materials
- Speeds up drying times
- Has lubricant properties
- Environmentally safe and bio compatible

Steelco Zyme

Neutral enzyme based detergent

Steelco Zyme is a near neutral PH complex multi enzyme based detergent formulated to break down and remove protein



- Suitable for use in manual cleaning, ultrasonic baths, washers and cart washers
- Excellent material compatibility
- Environmentally friendly
- PH neutral