

Specificații Tehnice

Numărul procedurii de achiziție:						ocds-b3wdp1-MD-1696579691013 din 12.10.2023		
Obiectul achiziției:						Ace pentru biopsie osoasă vertebrală		
Nr. Lot	Denumire Lot	Denumirea bunurilor și/sau a serviciilor	Modelul articolului	Țara de origine	Producătorul	Specificarea tehnică deplină solicitată de către autoritatea contractantă	Specificarea tehnică deplină propusă de către ofertant	Standarde de referință
1	2	3	4	5	6	7	8	9
1	Ace pentru biopsie osoasă vertebrală	Ace pentru biopsie osoasă vertebrală	F18A BONE BIOPSY F18A BONE BIOPSY KIT	SUA	Medtronic	Sistem pentru prelevarea intraoperatorie a biopsiei vertebrale - Diametru 8G. Lungimea de la 150mm. Ac cu gradare. Stilet solid. Canula pentru biopsie. Set steril. - Catalogul producătorului pentru produsele oferite - catalog și/sau pliant confirmat cu semnătura electronică a ofertantului. - Instrucțiunea de utilizare tradusă în limba de stat sau altă limbă de circulație internațională (Engleză/Rusă) confirmată cu semnătura electronică a ofertantului - la livrare - Instruirea specialiștilor în centre de referință specializate a producătorului pentru produsele care nu au fost folosite anterior, în termen de o lună de la solicitare.	Sistem pentru prelevarea intraoperatorie a biopsiei vertebrale - Diametru 8G. Lungimea de la 150mm. Ac cu gradare. Stilet solid. Canula pentru biopsie. Set steril. - Catalogul producătorului pentru produsele oferite - catalog și/sau pliant confirmat cu semnătura electronică a ofertantului. - Instrucțiunea de utilizare tradusă în limba de stat sau altă limbă de circulație internațională (Engleză/Rusă) confirmată cu semnătura electronică a ofertantului - la livrare - Instruirea specialiștilor în centre de referință specializate a producătorului pentru produsele care nu au fost folosite anterior, în termen de o lună de la solicitare.	ISO, CE

Semnat: _____ Numele, Prenumele: Dmitrii Kojevnikov În calitate de: Administrator

Ofertantul: S.C. OXIVIT-MED S.R.L Adresa: mun. Chişinău, str. Decebal 82-90

IDNO: 1007600044280;
adresa: mun. Chişinău, str. Decebal 82-90
tel. 022 808 002, 022 808 003
e-mail: oxivit.medical@gmail.com



Anexa nr. 23
la Documentația standard
conform Ordin MF Nr. 115
din 15.09.2021

Specificații de Preț

Numărul procedurii de achiziție:					ocds-b3wdp1-MD-1696579691013 din 12.10.2023							
Obiectul achiziției:					Ace pentru biopsie osoasă vertebrală							
Nr. Lot	Denumire Lot	Denumirea bunurilor și/sau a serviciilor	Unitatea de măsură	Cantitatea	Preț unitar (fără TVA)	Preț unitar (cu TVA)	Suma fără TVA	Suma cu TVA	Termenul de livrare/ prestare	Cod CPV	Clasificație bugetară (IBAN)	Discount %
1	2	3	4	5	6	7	8	9	10	11	12	13
1	Ace pentru biopsie osoasă vertebrală	Ace pentru biopsie osoasă vertebrală	Bucăți	30	2 622.00	3 146.40	78 660.00	94 392.00	Livrarea bunurilor se va efectua în termen de maxim 10 zile după semnarea contractului, la sediul beneficiarului.	33100000-1	-	-
TOTAL Oferta							78 660.00	94 392.00				

Semnat: _____ Numele, Prenumele: Dmitrii Kojevnikov În calitate de: Administrator

Ofertantul: S.C. OXIVIT-MED S.R.L Adresa: mun. Chişinău, str. Decebal 82-90

DECLARAȚIE DE ELIGIBILITATE

Către IMSP Institutul Oncologic

(denumirea autorității/entității contractante și adresa completă)

Stimați domni,

Subsemnatul, reprezentant împuternicit al S.C. OXIVIT-MED S.R.L, în calitate de ofertant, declar pe propria răspundere, sub sancțiunea excluderii din procedură și sub sancțiunile aplicate faptei de fals în acte publice, că nu mă aflu în una dintre situațiile prevăzute la art. 19 din Legea nr. 131/2015 privind achizițiile publice.

Mă oblig, la solicitarea autorității/entității contractante, în scopul verificării și confirmării declarației, să prezint orice document doveditor de care dispun.

Data completării 11.10.2023

Ofertant/candidat
S.C. OXIVIT-MED S.R.L
Dmitrii Kojevnikov

(semnătura autorizată)

OXIVIT MED

IDNO: 1007600044280;
adresa: mun. Chișinău, str. Decebal 82-90
tel. 022 808 002, 022 808 003
e-mail: oxivit.medical@gmail.com

**Către Grupul de lucru pentru evaluarea
Procedurii de achiziție Nr:
ocds-b3wdp1-MD-1696579691013 din 12.10.2023
din cadrul IMSP Institutul Oncologic**

DECLARAȚIE

Prin prezenta, SRL „Oxivit-med”, declara ca:

- Va oferi instrucțiunea de utilizare tradusă în limba de stat sau altă limbă de circulație internațională (Engleză/Rusă)
- Va oferi instruirea specialiștilor în centre de referință specializate a producătorului pentru produsele care nu au fost folosite anterior, în termen de o lună de la solicitare.

Data completării: 11.10.2023


Cu stimă,
Dmitrii Kojevnikov
Administrator

S.C. OXIVIT-MED S.R.L. _____
(semnătura)



Kyphon® Interventional Therapies

PRODUCT CATALOG



2016
edition

Featuring Kyphon® Balloon Kyphoplasty and Premium Vertebroplasty

**MICHELSON
TECHNOLOGY
AT WORK**

Table of Contents

2015

Fracture Reduction

Kyphon® Inflatable Bone Tamps (IBTs)

Kyphon XPANDER® II IBT	3
Kyphon XPANDER® IBT	5
Kyphon Express™ II IBT	7
Kyphon Express™ IBT	9
Kyphon Exact® IBT	11

Syringes

Kyphon® Inflation Syringe	13
Kyphon XPANDER® Inflation Syringe	14

Curettes

Kyphon® Latitude™ II	15
Kyphon Express™ Curette	16

KyphoPak® Trays

Kyphon XPANDER® II First Fracture Osteo Introducer® System	17
Kyphon XPANDER® First Fracture Osteo Introducer® System	18
Kyphon XPANDER® II Additional Fracture	19
Kyphon XPANDER® Additional Fracture	20
Kyphon Express™ II First Fracture	21
Kyphon Express™ First Fracture with CDS	22
Kyphon Express™ II Additional Fracture	23

Bone Cement

Kyphon® Xpede™	24
Kyphon® Xpede™ Bone Cement and Kyphon® Mixer Pack with Kyphon Multi-tap™ Adapter	25
Kyphon® ActivOs™ 10	27
Kyphon® HV-R®	28


Kyphon® Cement Delivery System	30
Kyphon® CDS Cement Cartridges	32
Kyphon® Bone Filler Device (BFD) & Kyphon® Directional BFD	33
Kyphon Express™ BFD & Kyphon Express™ Directional BFD	34
Kyphon® Mixer	35
Kyphon Multi-tap™ Adapter	36

Bone Access

Kyphon® Osteo Introducer® System	38
Kyphon Advanced™ Osteo Introducer® System	39
Kyphon One-Step™ Osteo Introducer® System - Trocar & Diamond	40
Kyphon One-Step™ Osteo Introducer® Device - Bevel	41
Kyphon Advanced™ Osteo Introducer® Device	42
Kyphon® 11 Gauge Bone Access Needles	43
Kyphon Express™ Osteo Introducer® System - Diamond & Bevel	44
Kyphon Express™ Osteo Introducer® Device - Trocar	45
Kyphon Express™ Advanced Osteo Introducer® Device	46
Kyphon Express™ Osteo Introducer® Device - Diamond	47
Kyphon Express™ Osteo Introducer® Device - Bevel	48
Kyphon Express™ Osteo Introducer® Device - Drill	49
Kyphon® Bone Biopsy Device	50
Kyphon Express™ Bone Biopsy Device	51
Kyphon® Bone Biopsy Kit	52
Kyphon V™ Premium Vertebroplasty	53
Kyphon V™ Osteo Introducer 11Ga Kits	54
Kyphon V™ Osteo Introducer 13Ga Kits	55
Kyphon V™ Bone Biopsy Kits	56
Vertebroplasty Kits and A La Carte Packages	57
SKY Bundles	59
KyphoPak® Trays	61
Kyphon® Instrument Size Chart	67

Fracture Reduction | Balloons | Size 3

Kyphon XPANDER® II Inflatable Bone Tamp




Fracture Reduction | Balloons | Size 3**Kyphon XPANDER® II Inflatable Bone Tamp - Balloon Dimensions**

	IBT Size	Initial Length (mm)	Max Volume (cc)	Pressure (psi)
Kyphon XPANDER® II Inflatable Bone Tamp	10/3	10	4	700
	15/3	15	5	700
	20/3	20	6	700

Fracture Reduction | Balloons | Size 3

Kyphon XPANDER® Inflatable Bone Tamp



Fracture Reduction | Balloons | Size 3

Kyphon XPANDER® Inflatable Bone Tamp - Balloon Dimensions

	IBT Size	Initial Length (mm)	Max Volume (cc)	Pressure (psi)
Kyphon XPANDER® Inflatable Bone Tamp	10/3	10	4	400
	15/3	15	4	400
	20/3	20	6	400

Fracture Reduction | Balloons | Size 2

Kyphon Express™ II Inflatable Bone Tamp



10 mm (KE102)



15 mm (KE152)




20 mm (KE202)

Fracture Reduction | Balloons | Size 2**Kyphon Express™ II Inflatable Bone Tamp - Balloon Dimensions**


	IBT Size	Initial Length (mm)	Max Volume (cc)	Pressure (psi)
Kyphon Express™ II Inflatable Bone Tamp	10/2	10	3	700
	15/2	15	4	700
	20/2	20	5	700

Fracture Reduction | Balloons | Size 2

Kyphon Express™ Inflatable Bone Tamp



10 mm (K15A)



15 mm (K15B)

Fracture Reduction | Balloons | Size 2**Kyphon Express™ Inflatable Bone Tamp - Balloon Dimensions**

	IBT Size	Initial Length (mm)	Max Volume (cc)	Pressure (psi)
Kyphon Express™ Inflatable Bone Tamp	10/2	10	4	400
	15/2	15	5	400

Fracture Reduction | Balloons | Size 3

Kyphon Exact® Inflatable Bone Tamp



10 mm - Directional Guide Included (K17A)




15 mm - Directional Guide Included (K17B)

Fracture Reduction | Balloons | Size 3**Kyphon Exact® Inflatable Bone Tamp - Balloon Dimensions**

	IBT Size	Initial Length (mm)	Max Volume (cc)	Pressure (psi)
Kyphon Exact® Inflatable Bone Tamp	10/3	10	3	300
	15/3	15	4	300

Fracture Reduction | Syringes


Kyphon® Inflation Syringe with locking syringe (A08E)



Compatible with all Kyphon® Inflatable Bone Tamps

Fracture Reduction | Syringes


Kyphon XPANDER® Inflation Syringe with locking syringe (A08A)



Not compatible with Kyphon XPANDER® II or Kyphon® Express™ II Inflatable Bone Tamps

Fracture Reduction | Curettes | Size 3

Kyphon® Latitude™ II Curette



Fracture Reduction | Curettes | Size 2

Kyphon Express™ Curette




7.0 mm T-Tip (A13A)

KyphoPak® Trays

Kyphon XPANDER® II First Fracture with Osteo Introducer® System (KPX153RB)




Kyphon XPANDER® II Inflatable Bone Tamps, 15/3 (KX153)



Kyphon® 11g Bone Access Needle



Kyphon® Bone Filler Devices (F04B)



Kyphon® Osteo Introducer® System (T15D)




Kyphon® Inflation Syringes (A08E)

KyphoPak® Trays


Kyphon XPANDER® First Fracture with Osteo Introducer® System (KPT1502)




Kyphon XPANDER® Inflatable Bone Tamps, 15/3 (K09A)




Kyphon® 11g Bone Access Needle



Kyphon® Bone Filler Devices (F04B)



Kyphon® Osteo Introducer® System (T15D)




Kyphon XPANDER® Inflation Syringes (A08A)

KyphoPak® Trays

Kyphon XPANDER® II Additional Fracture (KPX153NB)



Kyphon XPANDER® II Inflatable Bone Tamps, 15/3 (KX153) and Cannulas




Kyphon® Bone Filler Devices (F04B)

KyphoPak® Trays

Kyphon XPANDER® Additional Fracture (KPT1504)



Kyphon XPANDER® Inflatable Bone Tamps, 15/3 (K09A) and Cannulas



Kyphon® Bone Filler Devices (F04B)

KyphoPak® Express™ II Trays

Kyphon Express™ II First Fracture with Osteo Introducer® System (KEX152EB)



Kyphon Express™ II Inflatable Bone Tamps, 15/2 (KE152)



Kyphon Express™ Osteo Introducer® System - Diamond and Bevel (T34A)



Kyphon Express™ Bone Filler Devices (F06A)



Kyphon® Inflation Syringes (A08E)

KyphoPak® Trays

KyphoPak® Trays

Kyphon Express™ First Fracture Osteo Introducer® with Cement Delivery System (KPE1003-CDS)




Kyphon Express™ Inflatable Bone Tamps, 15/2 (K15B)




Kyphon® CDS Cement Cartridges (CC02A)



Kyphon® Cement Delivery Gun and Bone Filler Device
Size 2 (CDS2A)



Kyphon Express™ Osteo Introducer® System - Diamond
and Bevel (T34A)



Kyphon XPANDER® Inflation Syringes (A08A)

KyphoPak® Trays

Kyphon Express™ II Additional Fracture (KEX152NB)



Kyphon Express™ II Inflatable Bone Tamps, 15/2 (KE152) and Cannulas



Kyphon Express™ Bone Filler Devices (F06A)

Bone Cement

Kyphon® Xpede™ Bone Cement and Kyphon® Mixer Pack (CX01B)



Kyphon® Xpede™ Bone Cement and Kyphon® Mixer Pack with Kyphon Multi-tap™ Adapter (CX01B-MT2)



Kyphon® Xpede™ Bone Cement and Kyphon® Mixer Pack with Kyphon Multi-tap™ Adapter (CX01B-MT3)



Kyphon Multi-tap™ Adapter (A23)

Bone Cement

Kyphon® ActivOs™ 10 Bone Cement (C10A)

*Mixer sold separately (A07A)



Bone Cement

Kyphon® HV-R® Bone Cement (C01A)



Cement & Mixer Pack


Kyphon® HV-R® Bone Cement and Kyphon® Mixer Pack (C01B)



Delivery Systems | Size 3

Kyphon® Cement* Delivery Gun and Bone Filler Device (CDS3A)


*Cement cartridges sold separately (CC02A)



Delivery Systems | Size 2

Kyphon® Cement* Delivery Gun and Bone Filler Device (CDS2A)

*Cement cartridges sold separately (CC02A)



Delivery Systems

Kyphon® CDS Cement Cartridges (CC02A)




Delivery Systems | Size 3

Kyphon® Bone Filler Device (F04B)



Kyphon® Directional Bone Filler Device (F04C)




Delivery Systems | Size 2

Kyphon Express™ Bone Filler Device (F06A)



Kyphon Express™ Directional Bone Filler Device (F06B)




Mixer

Kyphon® Mixer (A07A)



Delivery Systems | Size 2

Kyphon Multi-tap™ Adapter (A22)



Delivery Systems | Size 3


Kyphon Multi-tap™ Adapter (A23)



Bone Access | Size 3 | **4.2 mm** outer diameter/**8** gauge

Kyphon® Osteo Introducer® System (T15D)


*Replaces Kyphon® Osteo Introducer® System (T05D)



Bone Access | Size 3 | **4.2 mm** outer diameter/**8** gauge

Kyphon Advanced™ Osteo Introducer® System (T15E)


*Replaces Kyphon® Advanced Osteo Introducer® System (T05E)



Bone Access | Size 3 | **4.2 mm** outer diameter/**8** gauge

Kyphon One-Step™ Osteo Introducer® System Trocar and Diamond (T15J)


*Replaces Kyphon® One-Step Osteo Introducer® System (T05J)



Bone Access | Size 3 | *4.2 mm outer diameter/8 gauge*

Kyphon One-Step™ Osteo Introducer® Device Bevel (T15K)


*Replaces Kyphon One-Step™ Osteo Introducer® Device (T05K)



Bone Access | Size 3 | **4.2 mm** outer diameter/**8** gauge




Kyphon Advanced™ Osteo Introducer® Device (A20A)

*Replaces Kyphon Advanced™ Osteo Introducer® Device (A10A)



Bone Access | Size 3



Kyphon® 11 Gauge Bone Access Needles



Bone Access | Size 2 | 3.4 mm outer diameter/10 gauge

Kyphon Express™ Osteo Introducer® System Diamond and Bevel (T34A)

Replaces Kyphon Express™ Osteo Introducer® System (T24A)



Bone Access | Size 2 | **3.4 mm** outer diameter/**10** gauge

Kyphon Express™ Osteo Introducer® Device Trocar (T34B)

Replaces Kyphon Express™ Osteo Introducer® Device (T24B)



Bone Access | Size 2 | **3.4 mm** outer diameter/**10** gauge


Kyphon Express™ Advanced Osteo Introducer® Device (T34C)

Replaces Kyphon Express™ Advanced Osteo Introducer® Device (T24C)




Bone Access | Size 2 | **3.4 mm** outer diameter/**10** gauge

Kyphon Express™ Osteo Introducer® Device Diamond (T34D)



Bone Access | Size 2 | **3.4 mm** outer diameter/**10** gauge

Kyphon Express™ Osteo Introducer® Device Bevel (T34E)



Bone Access | Size 2 | **3.4 mm** outer diameter/**10** gauge

Kyphon Express™ Osteo Introducer® Device Drill (T34F)



Drill

Bone Biopsy | Size 3 | *4.2 mm outer diameter/8 gauge/22.5 cm length*

Kyphon® Bone Biopsy Device (F05A)

Nozzle and Plunger Separated



Nozzle and Plunger Together



Bone Biopsy | Size 2 | *3.4 mm outer diameter/10 gauge*

Kyphon Express™ Bone Biopsy Device (F07A)




Bone Biopsy | Size 3 | *4.2 mm outer diameter/8 gauge*

Kyphon® Bone Biopsy Kit (F18A)

*Replaces Kyphon® Bone Biopsy Kit (F08A)

Bone Access, Biopsy Devices



Vertebroplasty

Kyphon V™ Osteo Introducer® 10Ga Trocar Kit (VPK1002)



VPT01B



Kyphon Kurve™ Kit (VPK1003)




VPT02A

10g cannula, curved BFD, and tamp




Vertebroplasty

Kyphon V™ Osteo Introducer® 11Ga Trocar/Bevel Kit (VPK011D)



Kyphon V™ Osteo Introducer® 11Ga Bevel Kit (VPK011B)



Vertebroplasty


Kyphon V™ Osteo Introducer® 13Ga Trocar/Bevel Kit (VPK013D)



VPT13D



Kyphon V™ Osteo Introducer® 13Ga Bevel Kit (VPK013B)



VPT13B



Vertebroplasty

Kyphon V™ Bone Biopsy Device 13Ga Kit (VPB11D) 13Ga Biopsy Device only (VPB11)



Kyphon V™ Bone Biopsy Device 15Ga Kit (VPB13D) 15Ga Biopsy Device only (VPB13)



Vertebroplasty Kits and A La Carte Packages

Catalog No.	Description	Kurve BFD and Tamp	Trocar	Trocar Tapered		Bevel Stylet		Matched Bevel		Bone Biopsy and Tamp		CDS Gun	CDS Cartridges
		13ga	10ga	11ga	13ga	11ga	13ga	11ga	13ga	13ga	15ga		
Kits													
VPK1002	Kyphon V™ Osteo Introducer® 10Ga Trocar Kit		1									1	2
VPK011D	Kyphon V™ Osteo Introducer® 11Ga Trocar/Bevel Kit			1		1						1	2
VPK011B	Kyphon V™ Osteo Introducer® 11Ga Bevel Kit							1				1	2
VPK013D	Kyphon V™ Osteo Introducer® 13Ga Trocar/Bevel Kit				1		1					1	2
VPK013B	Kyphon V™ Osteo Introducer® 13Ga Bevel Kit								1			1	2
VPK1003	Kyphon Kurve Kit	1	1									1	2
A La Carte													
Curved Bone Filler													
VPT02A	Kyphon Kurve Curved Bone Filler Device	1											
10 Gauge													
VPT01B	Kyphon V™ Osteo Introducer® Trocar 10Ga		1										

Vertebroplasty Kits and A La Carte Packages

Catalog No.	Description	Kurve BFD and Tamp	Trocar	Trocar Tapered		Bevel Stylet		Matched Bevel		Bone Biopsy and Tamp		CDS Gun	CDS Cartridges
		13ga	10ga	11ga	13ga	11ga	13ga	11ga	13ga	13ga	15ga		
11 Gauge													
VPT11D	Kyphon V™ Osteo Introducer® 11Ga Trocar/Bevel			1		1							
VPT11B	Kyphon V™ Osteo Introducer® 11Ga Bevel							1					
13 Gauge													
VPT13D	Kyphon V™ Osteo Introducer® 13Ga Trocar/Bevel				1		1						
VPT13B	Kyphon V™ Osteo Introducer® 13Ga Bevel								1				
Bone Biopsy													
VPB11	Kyphon V™ Bone Biopsy Device 13Ga									1			
VPB11D	Kyphon V™ Bone Biopsy Device 13Ga Kit			1		1				1			
VPB13	Kyphon V™ Bone Biopsy Device 15Ga										1		
VPB13D	Kyphon V™ Bone Biopsy Device 15Ga Kit				1		1				1		

SKY Bundles

with Kyphon® Balloon Kyphoplasty SKY Bundles with BFD Configurations

Single Pack BFD Configurations						
Catalog #	Balloon	Introducer	Drill	Syringe	Biopsy	BFDs
BUNASC10ADV	K15A (10/2)	T34C (Advanced)	-	A08A	F07A (Express Biopsy)	F06A (BFD) – qty 4
BUNASC15ADV	K15B (15/2)					
BUNASC10BEV	K15A (10/2)	T34E (Bevel)	T34F (Size 2 Drill)			
BUNASC15BEV	K15B (15/2)					
BUNASC10DIA	K15A (10/2)	T34D (Diamond)				
BUNASC15DIA	K15B (15/2)					
BUNASC10TRO	K15A (10/2)	T34B (Trocars)				
BUNASC15TRO	K15B (15/2)					
3-Pack BFD Configurations						
Catalog #	Balloon	Introducer	Drill	Syringe	Biopsy	BFDs
BUNASC10ADV-03	K15A (10/2)- qty 3	T24C (Advanced) – qty 3	-	A08A – qty 3	F07A (Express Biopsy) – qty 3	F06A (BFD) – qty 12
BUNASC15ADV-03	K15B (15/2)- qty 3					
BUNASC10BEV-03	K15A (10/2)- qty 3	T34E (Bevel) – qty 3	T34F (Size 2 Drill) – qty 3			
BUNASC15BEV-03	K15B (15/2)- qty 3					
BUNASC10DIA-03	K15A (10/2)- qty 3	T34D (Diamond) – qty 3				
BUNASC15DIA-03	K15B (15/2)- qty 3					
BUNASC10TRO-03	K15A (10/2)- qty 3	T34B (Trocars) – qty 3				
BUNASC15TRO-03	K15B (15/2)- qty 3					

SKY Bundles

with Kyphon® Balloon Kyphoplasty SKY Bundles with CDS Configurations

Single Pack CDS Configurations														
Bundle (Catalog#)	Balloon	Qty	Syringe	Qty	Introducer	Qty	Drill	Qty	CDS	Qty	Cartridges	Qty	Biopsy	Qty
BUNASC10ADVC	K15A	1	A08A	1	T34C	1	-	-	CDS2A	1	CC02A	1	F07A	1
BUNASC10TROC	K15A	1	A08A	1	T34B	1	T34F	1	CDS2A	1	CC02A	1	F07A	1
BUNASC10DIAC	K15A	1	A08A	1	T34D	1	T34F	1	CDS2A	1	CC02A	1	F07A	1
BUNASC10BEVC	K15A	1	A08A	1	T34E	1	T34F	1	CDS2A	1	CC02A	1	F07A	1
BUNASC15ADVC	K15B	1	A08A	1	T34C	1	-	-	CDS2A	1	CC02A	1	F07A	1
BUNASC15TROC	K15B	1	A08A	1	T34B	1	T34F	1	CDS2A	1	CC02A	1	F07A	1
BUNASC15DIAC	K15B	1	A08A	1	T34D	1	T34F	1	CDS2A	1	CC02A	1	F07A	1
BUNASC15BEVC	K15B	1	A08A	1	T34E	1	T34F	1	CDS2A	1	CC02A	1	F07A	1
3-Pack CDS Configurations														
Bundle (Catalog#)	Balloon	Qty	Syringe	Qty	Introducer	Qty	Drill	Qty	CDS	Qty	Cartridges	Qty	Biopsy	Qty
BUNASC10ADVC-03	K15A	3	A08A	3	T34C	3	-	-	CDS2A	3	CC02A	3	F07A	3
BUNASC10TROC-03	K15A	3	A08A	3	T34B	3	T34F	3	CDS2A	3	CC02A	3	F07A	3
BUNASC10DIAC-03	K15A	3	A08A	3	T34D	3	T34F	3	CDS2A	3	CC02A	3	F07A	3
BUNASC10BEVC-03	K15A	3	A08A	3	T34E	3	T34F	3	CDS2A	3	CC02A	3	F07A	3
BUNASC15ADVC-03	K15B	3	A08A	3	T34C	3	-	-	CDS2A	3	CC02A	3	F07A	3
BUNASC15TROC-03	K15B	3	A08A	3	T34B	3	T34F	3	CDS2A	3	CC02A	3	F07A	3
BUNASC15DIAC-03	K15B	3	A08A	3	T34D	3	T34F	3	CDS2A	3	CC02A	3	F07A	3
BUNASC15BEVC-03	K15B	3	A08A	3	T34E	3	T34F	3	CDS2A	3	CC02A	3	F07A	3

KyphoPak® Trays | Components & Options

with Kyphon XPANDER® II Inflatable Bone Tamps

Catalog No.	Description	First Fracture	Additional Fracture	Kyphon® 11 ga Bone Access Needle	T15D	T15E	T15J	KX103	KX153	KX203	F04B	A08E	CC02A	CDS3A
KPX103RB	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® Osteo Introducer® System	✓		2	1			2			6	2		
KPX103RB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® Osteo Introducer® with CDS System	✓		2	1			2				2	1	1
KPX103AB	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® Advanced Osteo Introducer® System	✓		2		1		2			6	2		
KPX103AB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® Advanced Osteo Introducer® with CDS System	✓		2		1		2				2	1	1
KPX103NB	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® Additional Fracture		✓					2			6			
KPX103PB	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® One-Step Osteo Introducer® System	✓					1	2			6	2		
KPX103PB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 10/3 IBT, Kyphon® One-Step Osteo Introducer® with CDS System	✓					1	2				2	1	1
KPX153RB	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® Osteo Introducer® System	✓		2	1				2		6	2		
KPX153RB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® Osteo Introducer® with CDS System	✓		2	1				2			2	1	1
KPX153AB	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® Advanced Osteo Introducer® System	✓		2		1			2		6	2		
KPX153AB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® Advanced Osteo Introducer® with CDS System	✓		2		1			2			2	1	1

KyphoPak® Trays | Components & Options

with Kyphon XPANDER® II Inflatable Bone Tamps

Catalog No.	Description	First Fracture	Additional Fracture	Kyphon® 11 ga Bone Access Needle	T15D	T15E	T15J	KX103	KX153	KX203	F04B	A08E	CC02A	CDS3A
KPX153NB	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® Additional Fracture	✓							2	6				
KPX153PB	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® One-Step Osteo Introducer® System	✓					1		2	6	2			
KPX153PB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 15/3 IBT, Kyphon® One-Step Osteo Introducer® with CDS System	✓					1		2		2	1	1	
KPX203RB	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® Osteo Introducer® System	✓	2	1						2	6	2		
KPX203RB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® Osteo Introducer® with CDS System	✓	2	1						2		2	1	1
KPX203AB	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® Advanced Osteo Introducer® System	✓	2		1					2	6	2		
KPX203AB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® Advanced Osteo Introducer® with CDS System	✓	2		1					2		2	1	1
KPX203NB	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® Additional Fracture	✓								2	6			
KPX203PB	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® One-Step Osteo Introducer® System	✓					1			2	6	2		
KPX203PB-CDS	KyphoPak® Tray, Kyphon XPANDER® II 20/3 IBT, Kyphon® One-Step Osteo Introducer® with CDS System	✓					1			2		2	1	1

KyphoPak® Trays | Components & Options

with Kyphon XPANDER® Inflatable Bone Tamps

Catalog No.	Description	First Fracture	Additional Fracture	Kyphon™ T1 ga Bone Access Needle	T15D	T15E	T15J	K13A	K09A	K08A	F04B	A08A	CC02A	CDS3A
KPT1002	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon® Osteo Introducer® System	✓		2	1			2			6	2		
KPT1002-CDS	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon® Osteo Introducer® with CDS System	✓		2	1			2				2	1	1
KPT1003	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon Advanced™ Osteo Introducer® System	✓		2		1		2			6	2		
KPT1003-CDS	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon Advanced™ Osteo Introducer® with CDS System	✓		2		1		2				2	1	1
KPT1004	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon® Additional Fracture		✓					2			6			
KPT1005	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon One-Step™ Osteo Introducer® System	✓					1	2			6	2		
KPT1005-CDS	KyphoPak® Tray, Kyphon XPANDER® 10/3 IBT, Kyphon One-Step™ Osteo Introducer® with CDS System	✓					1	2				2	1	1
KPT1502	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon® Osteo Introducer® System	✓		2	1				2		6	2		
KPT1502-CDS	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon® Osteo Introducer® with CDS System	✓		2	1				2			2	1	1
KPT1503	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon Advanced™ Osteo Introducer® System	✓		2		1			2		6	2		
KPT1503-CDS	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon Advanced™ Osteo Introducer® with CDS System	✓		2		1			2			2	1	1

KyphoPak® Trays | Components & Options

with Kyphon XPANDER® Inflatable Bone Tamps

Catalog No.	Description	First Fracture	Additional Fracture	Kyphon® T1 ga Bone Access Needle	T15D	T15E	T15J	K13A	K09A	K08A	F04B	A08A	CC02A	CD3A
KPT1504	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon® Additional Fracture	✓							2		6			
KPT1505	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon One-Step™ Osteo Introducer® System	✓					1		2		6	2		
KPT1505-CDS	KyphoPak® Tray, Kyphon XPANDER® 15/3 IBT, Kyphon One-Step™ Osteo Introducer® with CDS System	✓					1		2			2	1	1
KPT2002	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon® Osteo Introducer® System	✓		2	1					2	6	2		
KPT2002-CDS	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon® Osteo Introducer® with CDS System	✓		2	1					2		2	1	1
KPT2003	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon Advanced™ Osteo Introducer® System	✓		2		1				2	6	2		
KPT2003-CDS	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon Advanced™ Osteo Introducer® with CDS System	✓		2		1				2		2	1	1
KPT2004	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon® Additional Fracture		✓							2	6			
KPT2005	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon One-Step™ Osteo Introducer® System	✓					1			2	6	2		
KPT2005-CDS	KyphoPak® Tray, Kyphon XPANDER® 20/3 IBT, Kyphon One-Step™ Osteo Introducer® with CDS System	✓					1			2		2	1	1

KyphoPak® Trays | Components & Options

with Kyphon Express™ II Inflatable Bone Tamps

Catalog No.	Description	First Fracture	Additional Fracture	T34A	KE102	KE152	KE202	F06A	A08E	CC02A	CDS2A
KEX102EB	KyphoPak® Tray, 10/2 IBT, Kyphon Express™ II Osteo Introducer® System	✓		1	2			8	2		
KEX102EB-CDS	KyphoPak® Tray, 10/2 IBT, Kyphon Express™ II Osteo Introducer® System with CDS System Bundle	✓		1	2				2	1	1
KEX102NB	KyphoPak® Tray, 10/2 IBT, Kyphon® Additional Fracture		✓		2			8			
KEX152EB	KyphoPak® Tray, 15/2 IBT, Kyphon Express™ II Osteo Introducer® System	✓		1		2		8	2		
KEX152EB-CDS	KyphoPak® Tray, 15/2 IBT, Kyphon Express™ II Osteo Introducer® System with CDS System Bundle	✓		1		2			2	1	1
KEX152NB	KyphoPak® Tray, 15/2 IBT, Kyphon® Additional Fracture		✓			2		8			
KEX202EB	KyphoPak® Tray, 20/2 IBT, Kyphon Express™ II Osteo Introducer® System	✓		1			2	8	2		
KEX202EB-CDS	KyphoPak® Tray, 20/2 IBT, Kyphon Express™ II Osteo Introducer® System with CDS System Bundle	✓		1			2		2	1	1
KEX202NB	KyphoPak® Tray, 20/2 IBT, Kyphon® Additional Fracture		✓				2	8			

KyphoPak® Trays | Components & Options

with Kyphon Express™ Inflatable Bone Tamps

Catalog No.	Description	First Fracture	Additional Fracture	T34A	K15A	K15B	F06A	A08A	CC02A	CDS2A
KPE1001	KyphoPak® Tray, 10/2 IBT, Kyphon Express™ Osteo Introducer® System	✓		1	2		8	2		
KPE1001-CDS	KyphoPak® Tray, 10/2 IBT, Kyphon Express™ Osteo Introducer® with CDS System	✓		1	2			2	1	1
KPE1002	KyphoPak® Tray, 10/2 IBT, Kyphon® Additional Fracture		✓		2		8			
KPE1003	KyphoPak® Tray, 15/2 IBT, Kyphon Express™ Osteo Introducer® System	✓		1		2	8	2		
KPE1003-CDS	KyphoPak® Tray, 15/2 IBT, Kyphon Express™ Osteo Introducer® with CDS System	✓		1		2		2	1	1
KPE1004	KyphoPak® Tray, 15/2 IBT, Kyphon® Additional Fracture		✓			2	8			

Kyphon® Instrument Size Chart

Access	Gauge	O.D.	I.D.
		Millimeters	
Working Cannula, standard size 3	7 - 8	4.293	3.733
Kyphon® Advanced™, Osteo Introducer®, and Kyphon One-Step™ Stylet	9 - 10	3.505	-
Kyphon Express™ Working Cannula (Proximal End)	7	4.572	4.064
Kyphon Express™ Working Cannula (Distal End)	9 - 10	3.429	2.921
Kyphon Express™ Osteo Introducer® Stylet (Proximal End)	9 - 10	3.505	-
Kyphon Express™ Osteo Introducer® Stylet (Distal End)	12 - 13	2.692	-
Fracture Reduction			
Kyphon XPANDER™ II and XPANDER™ Inflatable Bone Tamp Shaft	12 - 13	2.591	-
Kyphon® Exact® Inflatable Bone Tamp Shaft	10	3.404	-
Kyphon® Latitude® II Curette Shaft	10	3.404	-
Kyphon Express™ II and Express™ Inflatable Bone Tamp Shaft	14 - 15	2.032	-
Kyphon Express™ Curette Shaft (Proximal/Distal End)	10/12 - 13	3.404/2.743	-
Fill			
Kyphon® Bone Filler Device/Bone Biopsy Device Nozzle	10	3.404	2.896
Kyphon® Bone Filler Device/Bone Biopsy Device Plunger	12	2.769	-
Kyphon Express™ Bone Filler Device/Bone Biopsy Device Nozzle (Proximal End)	9	3.759	3.251
Kyphon Express™ Bone Filler Device/Bone Biopsy Device Nozzle (Distal End)	13	2.413	1.905
Kyphon Express™ Proximal Plunger	11	3.048	-
Kyphon Express™ Distal Plunger	15	1.829	-

www.medtronic.com

Medtronic

Spinal and Biologics Business
Worldwide Headquarters

2600 Sofamor Danek Drive
Memphis, TN 38132

1800 Pyramid Place
Memphis, TN 38132

Telephone: 408.548.6500
Customer Service: 866.959.7466

The complication rate with Kyphon® Balloon Kyphoplasty has been demonstrated to be low. There are risks associated with the procedure (e.g., cement extravasation), including serious complications, and though rare, some of which may be fatal. For complete information regarding indications for use, contraindications, warnings, precautions, adverse events, and methods of use, please reference the devices' Instructions for Use included with the product.

UC201608691 EN ©2016 Medtronic Sofamor Danek USA, Inc. All Rights Reserved.
PMD011792-3.0



ACCESSORIES

BALLOON KYPHOPLASTY



PRODUCT DETAILS

CURETTES

Curettes are used to scrape and score bone in the vertebral body during balloon kyphoplasty procedures.

KYPHON LATITUDE II CURETTE

T-tip or wedge-tip, compatible only with 8ga access tools



KYPHON EXPRESS CURETTE

T-tip, compatible with both 8 and 10ga access tools



Curette with T-tip (A13A)

CURVED NEEDLE

The curved needle is designed to optimize flexibility and stiffness, with an open distal tip for precisely targeted cement flow across the vertebral body.

KYPHON KURVE CURVED BONE FILLER DEVICE

13ga, use with Kyphon™ Cement Delivery System



Curved Bone Filler Device

BONE BIOPSY

These tools let you easily collect a biopsy sample after accessing the vertebral body.

KYPHON EXPRESS II BIOPSY DEVICE


Size 2, 10ga



Biopsy Device (F07A)

KYPHON XPANDER II BIOPSY DEVICE AND KIT

Size 3, 8ga



Biopsy Device (F05A)

Biopsy Kit (F18A)

IMPORTANT SAFETY INFORMATION

The complication rate with Kyphon™ Balloon Kyphoplasty has been demonstrated to be low. There are risks associated with the procedure (e.g., cement extravasation), including serious complications, and through rare, some of which may be fatal.

Risks of acrylic bone cements include cement leakage, which may cause tissue damage, nerve or circulatory problems, and other serious adverse events, such as:

- Cardiac arrest
- Cerebrovascular accident
- Myocardial infarction
- Pulmonary embolism
- Cardiac embolism

For complete information regarding indications for use, contraindications, warnings, precautions, adverse events, and methods of use, please reference the devices' Instructions for Use included with the product.



Our Company
Careers
Contact & Support
Locations
ID Cards

Order Products
Diabetes Supplies
MRI Tools
Recycling
Product Manuals

Product Security
Electromagnetic Compatibility for Cardiac Devices
California Transparency in Supply Chains Act
Information About Proposition 65 for California Customers
Digital Millennium Copyright Act (DMCA) Notification