

| Procedure Name |                  | TOMOS Head,<br>TOMOS Hand,<br>TOMOS Knee,<br>TOMOS Abdomen | TOMOS C-Spine<br>TOMOS L-Spine,<br>TOMOS Hip Joint,<br>TOMOS Chest | TOMO      |           |
|----------------|------------------|--|--|-----------|-----------|
| General        | Image Processing | Table Action Mode  | Normal   | Normal    | Normal    |
|                |                  | TOUCH Disp. Mode   | RF   | RF        | RF        |
|                |                  | Auto Shutter   | MED  | MED       | MED       |
|                |                  | VReverse   | OFF  | OFF       | OFF       |
|                |                  | HReverse   | OFF  | OFF       | OFF       |
|                |                  | Field of View  | 12 inch  | 17 inch   | HOLD      |
|                |                  | Marker   | OFF  | OFF       | OFF       |
|                |                  | Position Memory  | no select  | no select | no select |
|                | Correction Step  | RkV(kV)  | 2  | 2         | 2         |
|                |                  | Acq.Time (step)  | 3  | 3         | 3         |
|                | Destination      | Send Dest 1  | no select  | no select | no select |
|                |                  | Send Dest 2  | no select  | no select | no select |
|                |                  | Auto Send Timing   | Next ACQ   | Next ACQ  | Next ACQ  |
|                |                  | Print Preset   | no select  | no select | no select |