

4- **Guedel Air Way** **Code No.: 2005**



❖ **Description:**

- The Guedel Airway is used to maintain an unobstructed Oropharyngeal Airway and it prevents biting of tongue and airway occlusion.
- Suitable for maintaining an unobstructed oropharyngeal airway during general anesthesia and in unconscious patients.
- Bite block to prevent biting of tongue and airway occlusion.
- Rounded traumatic edges.
- Smooth airway path for easy cleaning.
- Color coded for easy identification of sizes.
- Sizes: 00, 0, 1, 2, 3, 4, 5.

Size	00	0	1	2	3	4	5
Code	2005-00	2005-0	2005-1	2005-2	2005-3	2005-4	2005-5
Colour	Blue	Black	White	Green	Orange	Red	Yellow

❖ **Material:**

- H.D.P.E.
- P.P.
- Master Batch (Colored).