



Home » Products » Medical Imaging Products » Angiographic Syringes



# Angiographic Syringes

· Angiographic syringe is intended to deliver contrast media, which is compatible with high-pressure injection equipment during CT, DSA and MRI.

INQUIRY

## FEATURES

- Used for CT, MR, DSA-- It could be used for computed tomographic scanning, angiography, Nuclear magnetic resonance and imaging. The contrast media carried is the gold standard for identifying cancerous tissue, presenting high-definition images for lesions;
- Great variety of Syringes; it could match so many injectors of nearly all manufacturers on the market;



special configurations;

- Good materials to provide more resistant to high temperature and high pressure, good permeability;
- Clear body provides excellent clarification and smooth as glass feeling; it could see air bubble and contrast medium.

## USE METHOD

This video tutorial shows you how to use our angiographic syringe step by step -



**D A P A M E T E D**



Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for Medrad	Single	Dual
42.16.10000	42.16.10001	200	Used for CT: MCT & plus; OP-100; VCT610; ECT710	CT-200-A1	CT-200/200-A1
42.16.10002	42.16.10003	200	Used for CT: STELANT	CT-200-A2	CT-200/200-A2
42.16.10004		150	Used for DSA: Mark V&Mark V Provis	DSA-150-A1	
42.16.10005		200	Used for DSA: Mark V& Mark V Provis	DSA-200-A1	
42.16.10007		130	Used for DSA: Mark III & Mark IV	DSA-130-A1	
	42.16.10008	65	Used for MRI: Spectris		MRI-65/65-A1
	42.16.10009	115	Used for MRI: Spectris		MRI-115/115-A1
	42.16.10010	65/115	Used for MRI: Spectris		MRI-65/115-A1
42.16.10011	42.16.10012	190	Used for CT: SALENT	CT-190-A1	CT-190/190-A1
42.16.10016		150	Used for DSA: Mark 7	DSA-150-A2	



Single		Dual		Used for Libel-Flarsheim	Single	Dual
42.16.20000	42.16.20001	200		Used for CT: CT 9000& CT 9000 ADV	CT-200-B1	CT-200/200-B
42.16.20002		150		Used for DSA: Angiomat 6000	DSA-150-B1	
42.16.20003		150		Used for CT & DSA: Angiomat Illumena	CT/DSA-150-B1	
42.16.20005	42.16.20006	60		Used for MRI: Optistar	MRI-60-B1	MRI-60/60-B:

Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for EZEM	Single	Dual
42.16.50000	42.16.50001	200	Used for CT: Empower	CT-200-EM	CT-200/200-EI



Single		Dual		Used for Nemoto Kyorindo	
42,16,30000	42.16.30001	100		Used for CT: A-25	CT-100-NE CT-100/100-N
42.16.30002	42,16.30003	200		Used for CT: A-25 A-60	CT-200-NE CT-200/200-N
42.16.30005		120		Used for DSA: Nemoto 120S	DSA-120-NE
42.16.30006	42.16.30007	60		Nemoto Dual Shot. Sonic Shot 50	MRI-60-NE MRI-60/60-N
	42,16,30008	100/200		Used for CT.A-25	CT-100/200-N
	42.16.30012	60/200		Used for CT: A-25.A-60	CT-60/200-N

Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for Metron	Single	Dual



	42.16.40004	65/200	Used for MRI: ACCUTRON		MRI-65/200-M
--	-------------	--------	------------------------	--	--------------



## CE CERTIFICATION

This is a CE certified product. Here is the picture for your reference:



**Manufacturer:**

**SCW MEDICATH LTD.**

NO.4, Baolong 6th Road, Baolong Industrial  
Town, Longgang District, Shenzhen,  
518116, Guangdong, P.R. China

We, the manufacturer, herewith declare that the products  
**Angiographic Syringes**

Models: DSA-60-A1, DSA-150-A1, CT-100-X1, CT-200-A1, MRI-115-A1,

DSA-100-B1, DSA-150-B1, CT-200-B1, MRI-65-A1, CT-200-A2,  
DSA-200-A1, DSA-130-A1, CT/DSA-150-B1, MRI-60-B1, CT-200-NE,  
DSA-90-NE, DSA-120-NE, CT-100-NE, CT-190-A1, CT-190/190-A1,  
CT-200/200-A2, CT-200/200-B1, CT-100/100-NE, CT-100/200-NE,  
CT-200/200-NE, CT-200-EM, CT-200/200-EM, CT-200-MED,  
CT-200/200-MED, MRI-65/115-A1, MRI-65/65-A1, MRI-115/115-A1,  
MRI-60/60-B1, MRI-60/60-NE, MRI-65-MED, MRI-65/200-MED

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class IIa according to Annex IX Rule 2 of the  
Directive 93/42/EEC. It bears the mark

**CE 0197**

The product concerned has been manufactured under a quality management system  
according to Annex II of Directive 93/42/EEC.

Compliance of the designated product with the Directive 93/42/EEC has been assured via  
assessment of the quality management system by the Notified Body

**TÜV Rheinland LGA Products GmbH**  
Tillystraße 2  
90431 Nürnberg  
Deutschland

Certificate No.: HD 60144232 0001

Issue date: 26.05.2020

Expiry date: 26.05.2024

following the procedure relating to the EC Declaration of Conformity set out in Annex II of  
Directive 93/42/EEC.

This Declaration of conformity is valid in connection with the release document for the  
respective batch of produced devices.

The above mentioned declaration of conformity is exclusively under the responsibility of

**SCW MEDICATH LTD.**  
NO.4, Baolong 6<sup>th</sup> Road, Baolong Industrial Town,  
Longgang District, Shenzhen, 518116, Guangdong, P.R. China

Shenzhen, 2020/07/16

Place, date

Miriam Xie, RA

Legally binding signature, Function





# PRODUCT RECOMMENDATION



Angiographic Syringes Accessories