COMPATIBLE COMPONENTS REFERENCE CARD

Model 53401 Single Chamber Temporary Pacer

•



SURGICAL CABLES

Model		Туре	Channel	Length
	5487	Disposable, short	AorV	6
	5487L		A or V	12
		Disposable, long		
	5832S	Reusable, small clip	A or V	12
	5833S	Disposable, small clip,	A or V	6
	5833SL	short	A or V	12
		Disposable, small clip, long		

Medtronic

PATIENT CABLES

5433A Reusable, A 6 short 5433AL Reusable, A 12 5433V Reusable, V 6
long
EAZZV Bousship V/
short
5433VL Reusable, V 12 long
5846A Disposable, A 6 short
5846AL Disposable, A 12 long
5846V Disposable, V 6 short
5846VL Disposable, V 12 long

L = Long (12'), S = Small Clip, A = Atrial, V = Ventricle

OTHER

Model		Туре
	5409	Disposable Pouch
	53408	Replacement Covers (Main Cover & Rapid Atrial Pacing Cover)
	53407	IV Pole Replacement Hook & Replacement Knobs
	53409	Carrying Case

TEMPORARY PACING LEADS StreamLine[™] Pacing Leads (Heart Wires)

Model		Туре	Channel	Length
	6491F	Unipolar, Pediatric	A or V	525 mm
	6492F	Unipolar, Atrial	A	525 mm
ð	6494F	Unipolar (2 per pack)	A or V	610 mm
	6495F	Bipolar, Coaxial	A or V	475 mm
	6500F	Unipolar, Premium	A or V	525 mm

Brief Statement

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

medtronic.ca

L6Y OR3 Canada

 99 Hereford Street
 Tel: (905) 460-3800

 Brampton, Ontario
 Fax: (905) 460-3998

Medtronic and the Medtronic logo are trademarks of Medtronic. *Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. UC201703700 EC © 2017 Medtronic. All Rights Reserved.

Medtronic