


Multi-Drug Test

Tests for choices	Specimen	Catalog No.	Format	Package
AMP/BAR/BUP/	Saliva	GBDSA-715A to GBDSA-7105A	Multi-Drug Cassette	25T/Box
BZO/ COC/COT/	Saliva	GBDSA-715B to GBDSA-7105B	Multi-Drug Cassette	25T/Box
EDDP/MTD/MET/	Saliva	GBDSA-715C to GBDSA-7105C	Multi-Drug Tube	25T/Box
MDMA/6-MAM/OPI/	Saliva	GBDSA-715D to GBDSA-7105D	Multi-Drug Cube	25T/Box
OXY/PCP/PPX/THC	Saliva	GBDSA-715E to GBDSA-7105E	Multi-Drug Cassette	25T/Box

Rapid Test Reader



Report

- Can print report in small piece of paper
- Can be applied to A4 format report (PC)
- Provide small database, can batch processing and also save real-time image, automatic traceability data (PC)
- Exported in several formats : PDF and Excel Format (PC)

- Storage ▶ More than 1000 results
- Screen ▶ 5-inch resistive touch screen
- Power ▶ 5V, 5A, 25W
- Size/Wt ▶ 22 × 21 × 14 cm, 1.5kg



Methylphenidate Test Dip Card




Cotinine Test Dip Card



Methaqualone Test Dip Card



Saliva Alcohol Test Strip



Healgen Scientific Limited Liability Company
Add: 3818 Fuqua street, Houston, USA. TX 77047
Tel: +1 713-733-8088
Toll free: 866 982 3818
Fax: +1 713-733-8848
E-mail: healgen_sales@healgen.us. (For South America and North America)
Web:<http://www.healgen.com>

ZheJiang Orient Gene Biotech Co., Ltd
Add:679#, East Yangguang Avenue, Dipu Town,Anji
P.C.: 313300
Tel: +86-572-5303755/5303756
Fax: +86-572-5226222
E-mail: sales@orientgene.com (For Rest of world)
Web:<http://www.orientgene.com>



Rev.07/2016 L00114-01




DOA Rapid Tests
Say No to Drugs

DOA Rapid Tests-Urine


Say No to Drugs

Convenient Procedure


Strip




Cassette




Multi-Drug Cassette/Panel



Generic V Cup



Split Cup (A cup)



5 Minutes



Drugs of Abuse

Product Name	Specimen	Catalog No.	Format	Label	Sensitivity	Package
Amphetamine (AMP) Test	Urine	GBAMP-101 ^{V+}	Strip	Gold	1000 ⁺ /2000ng/ml	100T/Box
		GBAMP-U101 ^{V+}	Strip	Gold	500 ⁺ ng/ml	100T/Box
		GBAMP-A101 ^V	Strip	Gold	300ng/ml	100T/Box
		GBAMP-102 ^{V+}	Device	Gold	1000 ⁺ /2000ng/ml	25T/Box
		GBAMP-A102 ^V	Device	Gold	300ng/ml	25T/Box
Secobarbital(SEB) Test	Urine	GBSEB-101 ^{V+}	Strip	Gold	500 ⁺ ng/ml	25T/Box
		GBSEB-A101 ^V	Strip	Gold	300 ⁺ /600ng/ml	100T/Box
		GBSEB-102 ^{V+}	Device	Gold	200ng/ml	100T/Box
		GBSEB-A102 ^V	Device	Gold	300 ⁺ /600ng/ml	25T/Box
Buprenorphine (BUP) Test	Urine	GBBUP-101 ^{V+}	Strip	Gold	10 ⁺ ng/ml	100T/Box
		GBBUP-102 ^{V+}	Device	Gold	10 ⁺ ng/ml	25T/Box
Oxazepam(OZP) Test	Urine	GBOZP-101 ^{V+}	Strip	Gold	300 ⁺ /600ng/ml	100T/Box
		GBOZP-A101 ^V	Strip	Gold	200ng/ml	100T/Box
		GBOZP-102 ^{V+}	Device	Gold	300 ⁺ /600ng/ml	25T/Box
		GBOZP-A102 ^V	Device	Gold	200ng/ml	25T/Box
		GBCOC-101 ^{V+}	Strip	Gold	300 ⁺ /600ng/ml	100T/Box
Cocaine (COC) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	150 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	300 ⁺ /600ng/ml	25T/Box
		GBMET-102 ^{V+}	Device	Gold	150 ⁺ ng/ml	25T/Box
		GBMET-A102 ^V	Device	Gold	300 ⁺ /600ng/ml	25T/Box
Cotinine (COT) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
Methadone Metabolite (EDDP) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
Fentanyl (FEN) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
Synthetic Marijuana (K2) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
Ketamine (KET) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
Ecstasy (MDMA) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-102 ^{V+}	Device	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A102 ^V	Device	Gold	100 ⁺ ng/ml	100T/Box
Morphine (MOP) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-102 ^{V+}	Device	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A102 ^V	Device	Gold	100 ⁺ ng/ml	100T/Box
Methaqualone (MQL) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Methadone (MTD) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Oxycodone (OXY) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Phencyclidine (PCP) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Propoxyphene (PPX) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Nortriptyline (NTL) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Cannabinoids(CAB) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-102 ^{V+}	Device	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A102 ^V	Device	Gold	100 ⁺ ng/ml	100T/Box
Tramadol (TRA) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Ethyl Glucuronide (EtG) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
D-lysergic acid diethylamide (LSD) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box
Methylphenidate(MPD) Test	Urine	GBMET-101 ^{V+}	Strip	Gold	100 ⁺ ng/ml	100T/Box
		GBMET-A101 ^V	Strip	Gold	100 ⁺ ng/ml	100T/Box


Multi-Drug Test

Tests for choices	Specimen	Catalog No.	Format	Combination	Format
AMP/SEB/BUP/OZP/COC/COT/EDDP/FEN/	Urine	GBDOA124 to GBDOA1124 ^V	Multi-Drug Panel	2-12Drugs	20T/Box
K2/KET/MDMA/MET/MOP/MQL/MTD/OXY/	Urine	GBDOA125 to GBDOA1125 ^V	Multi-Drug Cassette	2-12Drugs	20T/Box
CAB/NTL/TRA/PCP/PPX/LSD/MPD	Urine	GBDOA126 to GBDOA1126 ^V	Multi-Drug Cup	2-12Drugs	20T/Box


✓ CE Market + Cleared for US 510(k)

DOA Rapid Tests-Saliva


Oral Fluid Drug Test



3 Minutes



One Step Oral Fluid Drug Test



Oral Fluid Test

Product Name	Specimen	Catalog No.	Label	Sensitivity
Amphetamine (AMP) Test	Saliva	GBAMP-702a	Gold	50ng/ml
Barbiturates (BAR) Test	Saliva	GBBAR-702a	Gold	50ng/ml
Buprenorphine (BUP) Test	Saliva	GBBUP-702a	Gold	300ng/ml
Benzodiazepines (BZO) Test	Saliva	GBBZO-702a	Gold	5ng/ml
Cocaine (COC) Test	Saliva	GBBUP-702b	Gold	10ng/ml
Cotinine (COT) Test	Saliva	GBBZO-702a	Gold	10ng/ml
Methadone Metabolite (EDDP) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Methadone (MTD) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Methamphetamine (MET) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Ecstasy (MDMA) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Opiates (OPI) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Oxycodone (OXY) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Phencyclidine (PCP) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Propoxyphene (PPX) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Marijuana (THC) Test	Saliva	GBBZO-702a	Gold	50ng/ml
6-Monoacetylmorphine(6-MAM) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Ketamine(KET) Test	Saliva	GBBZO-702a	Gold	50ng/ml
Saliva Alcohol Test	Saliva	GBBZO-702a	Gold	50ng/ml



E



D



C



Device B



Device A

