



## **Bone Cements**



The **Cement Range** includes bone cements for primary arthroplasty surgeries (**CEMFIX® 1 & CEMFIX® 3**) and revision surgeries of infected prostheses (**GENTAFIX® 1**, **GENTAFIX® 3 & GENTAFIX® 3** MV).

- Self-hardening PMMA cements.
- Mechanical properties far above the ISO 5833 standard.
- Low monomer release.
- Rapid self-hardening solution to reduce prosthetis micromovements.
- Stable implant fixation.
- Barium sulfate, an excellent radiopaque agent.
- Easily prepared within a bowl & spatula or with vacuum mixing system.

	LOCALIZATION	PRODUCT	VISCOSITY	GENTAMICIN (local antibiotherapy for revision surgery)	MIXING/INJECTION SYSTEM		ORDERING			
					Bowl & Spatula	Vacuukit	Weight	References	REGISTRATION	
	KNEE HIP SHOULDER	CEMFIX® 1	STANDARD		<b>*</b>		40g	T040140	CE (class IIb)	
		GENTAFIX® 1	STANDARD	G	<b>=</b>		40g	T040140G	CE (class III)	
		GENTAFIX® 3 MV	MEDIUM	G	<b>=</b>	<b>7</b>	40g	T040341G	CE (class III)	
	HIP (Stem) SHOULDER (Stem)	CEMFIX® 3	LOW		<b>=</b>	7	40g	T040340	CE (class IIb)	
		GENTAFIX® 3	LOW	G	<b>=</b>	<del>'</del>   <b>7</b>	40g	T040340G	CE (class III)	

**GENTAMICIN**, issued from aminoglycosides family, has a broad spectrum and is effective against most bacteria encountered. The quantity of gentamicin contained in one dose has been optimised to guarantee a long term release and then prevent infections.



## Sales and engineering departments:

Z.I. de Montredon - 12, rue Apollo 31240 L'UNION (France) Tel (33) 5 34 25 10 60 Fax (33) 5 34 25 27 39

## Administration and supply-chain offices:

Z.I de la Herray 65500 VIC en BIGORRE (France) Tel (33) 5 62 96 88 38 Fax (33) 5 62 96 28 72

