



S	M	L
34-36	36-38	38-40

ARN 104

The product has been designed to provide higher stabilization than the sponge and the Nelson collars following a neck trauma or cervical disc hernia. Made of plastizote and suitable for MR monitoring as it does not contain any metal parts.



S	M	L
34-36	36-38	38-40

ARN 104B

Functions like the classical Philadelphia collar. It has a tracheal opening in the front. Made of plastizote and suitable for MR monitoring as it does not contain any metal parts.

ARMOR CLASSICS



S	M	L	XL	XXL
70-80	80-90	90-100	100-110	110-120



trunk

Available colors;



ARC 320

This product is actively used for facet joint problems of the lumbago, unbalanced postural muscles of the lumbar spine, ligament problems, osteokondrosis, spondiliosis, spondiliartosis and chronic back pains. Air permeable elastic filamented woven textile. 32 cm back height.



S	M	L	XL	XXL
70-80	80-90	90-100	100-110	110-120




Available colors;



ARC 320K

This product is actively used for facet joint problems of the lumbago, unbalanced postural muscles of the lumbar spine, ligament problems, osteokondrosis, spondiliosis, spondiliartosis and chronic back pains. Air permeable elastic filamented woven textile. 32 cm back height. Comes with built-in additional belt.

ARMOR CLASSICS



S	M	L	XL	XXL
70-80	80-90	90-100	100-110	110-120



trunk

Available colors;



ARC 330K

This product is applicable to correct the dorsal vertebra and reduce kyphosis, in traumas of lumbar and thoracic vertebra and postural abnormalities. Contains a pair of steel inserts.



S	M	L	XL	XXL
70-80	80-90	90-100	100-110	110-120




Available colors;



ARC 340

This product has been designed to correct dorsal vertebra, reduce kyphosis and for thoracic vertebral traumas and postural deficiencies. Contains a «Y» shaped light weight durable aluminum insert.

ANTI-EMBOLISM STOCKINGS



MEDICAL STOCKINGS

	cB	cD	cF	cG
1 - Small	19-21	29-32	39-47	42-55
2 - Medium	21-23	30-35	41-50	44-58
3 - Large	23-25	32-38	45-54	48-62
4 - X Large	25-27	33-42	46-57	50-67
5 - XX Large	27-29	35-45	47-60	53-73

Ciorapii compresivi antitromboza, pana la coapsa ARS08 previn aparitia formarii unor mici cheaguri de sange care sunt extrem de periculoase pentru circulatia sanguina a picioarelor. Aceste afectiuni apar de obicei la persoanele care sustin o activitate fizica redusa la nivelul musculaturii membrelor inferioare.

Ciorapii antitromboza se recomanda:

- ✓ In preventia afectiunilor ce tin de circulatia sangelui, cum ar fi venele varicoase si trombozele;
- ✓ Pacientilor care sunt imobilizati la pat din cauza unor probleme de sanatate la membrele inferioare;
- ✓ Pe perioada internarii in spital, inainte sau dupa interventii chirurgicale, care tin bolnavul la pat;
- ✓ Persoanelor care se stiu deja cu afectiuni ale circulatiei venoase la nivelul membrelor inferioare.

Nu se recomanda folosirea in caz de: edeme ale membrelor, picioare diforme, neuropatie periferica, boala arteriala periferica ocluziva, sensibilitati ale pielii cum ar fi dermatite si alergii la material.

Marimea potrivita se va selecta in urma masuratorilor indicate in imagini conform cu dimensiunile din tabel.