

LANCETS AND LANCET DEVICE

23919 AUTO LANCET DEVICE

Reusable auto lancet plastic device, adjustable tip, allows precise configuration of lancet penetration, with 10 adjustable depths. Lancet ejector for safe disposal of the lancet. Instructions: GB, FR, IT, ES.



23919

- 23917 GIMA LANCETS 28G - sterile - box of 25 pcs
 - 23916 GIMA LANCETS 28G - sterile - box of 100 pcs
 - 23918 GIMA LANCETS 30G - sterile - box of 100 pcs
- Fine needles usable with code 23919 and most of lancet devices in the market to draw a blood sample for glucose testing.



Multilingual box: GB, FR, IT, DE, ES, PT, GR, PL, CZ, Arabic

ABBOTT GLUCOSE STRIPS

ABBOTT FREESTYLE OPTIUM BLOOD GLUCOSE TEST STRIPS		
GIMA code	Compatibility	Minimum order
23942	FreeStyle Optium Neo, FreeStyle Libre, FreeStyle Optium, Optium Xceed, Optium, Optium Easy, Boots	box of 25
23943		box of 50
ABBOTT FREESTYLE LITE BLOOD GLUCOSE TEST STRIPS		
23944	FreeStyle Lite, FreeStyle Freedom Lite, FreeStyle InsuLinx, FreeStyle Navigator II	box of 25
23945		box of 50



Abbott


23944-5

23942-3

GIMALANCET SQUARE AND ROUND SAFETY LANCETS - STERILE

Easy to use, comfortable, safety lancets to meet various capillary blood requirements.

- **single-use:** the build-in-self-destructing design can prevent reuse of the device
- **no accidental injury risk:** needle is hidden before and after use to prevent accidental injuries
- **sterilization guaranteed:** the lancet is Gamma sterilized and the needle tip is protected with a sterility tab until its use allowing a long shelf life
- **no blood contamination**



Dual-Spring Push-Pull design to eliminate needle vibrations. Ultra-sharp tri-bevel needle

- 1 Remove the protective cap.
- 2 Push gently against the testing site.
- 3 Dispose of safety lancets in sharps container.

Easy to use pre-loaded and pressure activated device

High quality tri-bevel needle high-speed penetration with single spring

GIMA ROUND SAFETY LANCETS - NEEDLES									
GIMA code	Colour	Depth of penetration	Gauge	Blood Flow	Finger Stick	Glucose	Hematocrit	Bleeding	Minimum order
24451	Yellow	2.2 mm	G21	40-60 µl	•	•	•	High	box of 100
24453	Green	2.2 mm	G23	40-60 µl	•	•	•	High	box of 100
24454	Light blue	1.8 mm	G26	5-10 µl	•	•	•	Medium	box of 100
24455	Purple	1.8 mm	G28	5-10 µl	•	•	•	Low	box of 100
GIMA SQUARE SAFETY LANCETS - NEEDLES									
24461	Orange	2.2 mm	G21	40-50 µl	•	•	•	High	box of 100
24463	White	1.8 mm	G23	20-40 µl	•	•	•	Medium	box of 100
24465	Purple	1.8 mm	G28	5-10 µl	•	•	•	Low	box of 100



24461 24465 24451 24453 24454 24455

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, CZ, SE, RO, SK, NL, DK, FI, EE, HR, HU, LT, LV, NO, SI, BG, RU, GR, Arabic.

SURGILANCE SAFETY AUTOMATIC LANCETS - STERILE

SURGILANCE SAFETY AUTOMATIC LANCETS - NEEDLES

Multi purpose capillary blood sampling devices for all testing requirements.

The one-step safety lancet is easier to use, no arming is required, and safer, once lancet is used it is rendered inoperable.

High-speed delivery and penetration method minimizes patient pain and operator error.

Latex-free, hypo-allergenic, disposable, sterile.

Multilingual manual on www.medipurpose.com



CE 1008

Multilingual box: GB, FR, IT, ES, DE, PT, NL, SE, DK, RU, GR.

24515

SURGILANCE™ SAFETY LANCETS - NEEDLES									
GIMA code	Colour	Depth of penetration	Entire needle Gauge	Needle tip Gauge	Blood Flow	Finger Stick	Glucose	Hematocrit	Minimum order
24513	Yellow	1.0 mm	G21	G26	5-10 µl	•	•	•	box of 100
24514	Grey	1.8 mm	G21	G23	10-20 µl	•	•	•	box of 100
24515	Orange	2.2 mm	G21	G22	20-40 µl	•	•	•	box of 100
24517	Pink	2.8 mm	G21	G21	40-60 µl	•	•	•	box of 100

