



Viale G.di Vittorio, 30  
20090 Pantigliate (MI) Italy  
Phone 02.55230061  
[www.histoline.com](http://www.histoline.com)  
histoline@histoline.com

## Data Sheet Water Bath WB 2800



Version 1.2

Histo-Line Laboratories

18/10/2019



## INDEX

1.Features.....	3
1.2 Information and specified use .....	3
2.Technical specifications .....	4
2.1 Electrical data.....	4
2.2 Capacity.....	4
2.3 Dimensions and weight.....	5
2.4 Instrument equipment.....	5

## 1.Features



### *1.2 Information and specified use*

The Water Bath WB 2800 is a small water bath and stretch slices is designed to lay the sections of histological samples.

The Teflon coating allows easy cleaning.

The instrument is provided with a control panel to set and monitor the desired water temperature, and is equipped with overheating protection.

The water temperature is kept constant thanks to the internal microprocessor.

After switching off the instrument, or if there is a power surge, the re-power displays the last set temperature (memory function).

The board allows to support the slides around to be dried.

At the front is a practical place "invitation" for the emptying.



## 2. Technical specifications

### 2.1 Electrical data

Certification: CE, IVD

Protection class: I

Fuses type Fst: 10 A

Pollution class: 2

Over-voltage Category II Installation

Power supply: 220V ~ 50Hz and 110V ~ 60Hz

Maximum power draw: 300 V/A

### 2.2 Capacity

Water Bath: ambient °C ~ 75°C increments 1°C

Internal dimensions of the bath.

Width: 230 mm

Depth: 180 mm

Height: 60 mm

### *2.3 Dimensions and weight*

Width: 328 mm

Depth: 300 mm

Height: 96 mm

Weight: 2,4 kg

### *2.4 Instrument equipment*

Power cable

Manual instructions