



Angiographic Syringes

- Clear body Provides excellent clarification and smooth as glass feeling ,it could see air bubble and contrast medium
- It could be used for computed tomographic scanning, angiography, Nuclear magnetic resonance and imaging
- Great Variety of Syringes, it could match so many injectors



| Product No. | | Volume(ml) | Description | Models | |
|-------------|-------------|------------|--|------------|----------------|
| Single | Dual | | Used for Medrad | Single | Dual |
| 42.16.10000 | 42.16.10001 | 200 | Used for CT: MCT+ OP-100;VCT610; ECT710 | CT-200-A1 | CT-200/200-A1 |
| 42.16.10002 | 42.16.10003 | 200 | Used for CT: STELLANT | CT-200-A2 | CT-200/200-A2 |
| 42.16.10004 | | 150 | Used for DSA: Mark V & Mark V Provis | DSA-150-A1 | |
| 42.16.10005 | | 200 | Used for DSA: Mark V & Mark V Provis | DSA-200-A1 | |
| 42.16.10007 | | 130 | Used for DSA: Mark III & Mark IV | DSA-130-A1 | |
| | 42.16.10008 | 65 | Used for MRI: Spectris | | MRI-65/65-A1 |
| | 42.16.10009 | 115 | Used for MRI: Spectris | | MRI-115/115-A1 |
| | 42.16.10010 | 65/115 | Used for MRI: Spectris | | MRI-65/115-A1 |
| 42.16.10011 | 42.16.10012 | 190 | Used for CT: SALIENT | CT-190-A1 | CT-190/190-A1 |
| 42.16.10016 | | 150 | Used for DSA:Mark 7 | DSA-150-A2 | |
| | 42.16.10019 | 200 | Used for GT STELLANT(12HOURS) | | CT-200/200-A2 |