

Anexa tehnica, Lot 10

Model EBA 200

Parametri solicitati	Parametri oferiti
Centrifugă, de laborator (8-12 tuburi) Cod 150910 Descriere Centrifugele de laborator sunt destinate pentru a centrifuga sângele. Parametrul Specificația Viteza de rotație Minimală, rpm reglabilă Maximală, rpm ≥ 6000 Setări setarea vitezei Capacitatea Tipul tuburilor tuburi de 10 - 15 ml Numărul de tuburi 8-12 tuburi Timer Gama de timp 0 - ≥ 60 min Incrementarea ≤ 1 min Securitatea Blocarea capacului în timpul lucrului da Indicatori Indicarori vizual și acustic da Debalansare da Pornire/oprire da Capac deschis da Display Digital da Nivelul de zgomot ≤ 70 dB da Fereastra, orificiu sau alt acces, necesar pentru efectuarea procedurii de verificare periodica, conform normelor si standartelor in vigoare da	Centrifugă, de laborator (8 tuburi) Cod 150910 Descriere Centrifugele de laborator sunt destinate pentru a centrifuga sângele. Parametrul Specificația Viteza de rotație Minimală, rpm reglabilă Maximală, rpm 6000 Setări setarea vitezei Capacitatea Tipul tuburilor tuburi de 10 - 15 ml Numărul de tuburi 8 tuburi Timer Gama de timp 0 - 99 min Increment 1 min Securitatea Blocarea capacului în timpul lucrului da Indicatori Indicarori vizual și acustic da Debalansare da Pornire/oprire da Capac deschis da Display Digital da Nivelul de zgomot 50 dB da Fereastra, orificiu sau alt acces, necesar pentru efectuarea procedurii de verificare periodica, conform normelor si standartelor in vigoare da

Certificate

mdc medical device certification GmbH
certifies that



Andreas Hettich GmbH & Co. KG
Föhrenstrasse 12
78532 Tuttlingen
Germany

for the scope

**design, development, manufacturing, distribution and servicing of
laboratory centrifuges, centrifuges for separation of blood components for
transfusion purposes and microbiological incubators.**

has introduced and applies a

Quality Management System

The mdc audit has proven that this quality management system
meets all requirements of the following standard

EN ISO 13485

Medical devices – Quality management systems –
Requirements for regulatory purposes

EN ISO 13485:2016 + AC:2016 - ISO 13485:2016

Valid from	2020-04-01
Valid until	2023-03-31
Registration no.	D1459300001
Report no.	P19-01336-156285
Stuttgart	2020-03-26

A blue ink signature, likely of the Head of Certification Body, written in a cursive style.

Head of Certification Body



CUSTOMER INFORMATION

+49 7461 705 1318
christian.vondergruen@hettichlab.com
www.hettichlab.com

TRANSITION PHASE - MDR 2017/745 AND IVDR 2017/746

Tuttlingen, 09.04.2020

Dear valued Customer,

Andreas Hettich GmbH & Co. KG as legal manufacturer is required to comply with EU Medical Device Regulation 2017/745 and EU In-Vitro Diagnostic Medical Device Regulation 2017/746 for the following medical devices / in vitro diagnostic medical devices:


Medical Devices	ROTO SILENTA 630 RS, ROTIXA 500 RS, ROTANTA 460; ROTANTA 460R, ROTANTA 460RC; ROTANTA 460RF
IVD Medical Devices	EBA 200; EBA 200S, EBA 270, EBA 280; EBA 280S, MIKRO 220; MIKRO 220R, UNIVERSAL 320; UNIVERSAL 320R, MIKRO 185, MIKRO 200; MIKRO 200R, ROTOFIX 32A, ROTINA 380; ROTINA 380R, ROTINA 420; ROTINA 420R, HAEMATOKRIT 200; HETTCUBE ff.

We started already our transition phase in the beginning of 2019 and are currently under implementation of the new requirements.

Our current EN ISO 13485:2016 Certificate, MDSAP Certificate and CE-Certificate MDD 93/42/EEC Annex II remains valid until the end of the transition phase for our in vitro diagnostic devices under IVDR 2017/746 until **26.05.2022** and for our medical devices under MDR 2017/745 until **26.05.2024**.

We will publish the successful transition of our product portfolio on our web page as soon as available.

Mit freundlichen Grüßen / Best regards,
Andreas Hettich GmbH & Co. KG

i.A. 

Christian von der Grün
Head of Regulatory & Quality Affairs

SMALL CENTRIFUGES

— EBA 200 | 200 S

TOP PERFORMANCE FOR SMALL LABORATORIES

The EBA 200 and EBA 200 S are practical, compact centrifuges for small sample sizes. The high speed S model can deliver reliable PPP results in 3 minutes or less. Included an 8-place fixed angle rotor to hold standard blood and urine tubes up to 15 ml in capacity.



— HIGHLIGHTS

- max. RPM: 200 - 6,000 min⁻¹ | 8,000 min⁻¹
Adjustable in increments of 10
- max. RCF: 3,461 | 6,153
- max. capacity: 8 x 15 ml
- small centrifuge including 8-place angle rotor
- compact design
- IVD-conform according to directive 98/79/EC
- maximum noise level of ≤ 50 dB(A) | ≤ 55 dB(A)
- impulse key for short cycle mode
- 2 individual acceleration and deceleration stages
- easy operation with keypad

— FEATURES

- plastic housing
- metal lid
- viewing port in the lid
- one-hand lid lock
- lid dropping protection
- emergency lid lock release
- light alloy chamber
- brushless drive
- error display
- imbalance monitoring and switch-off
- backlit panel with actual values of all parameters
- 3 adjustable messages after completion of the centrifugation run (audio + visual)

TECHNICAL DATA

	EBA 200	EBA 200 S
voltage *)	200 – 240 V 1 ~	200 – 240 V 1 ~
frequency	50 – 60 Hz	50 – 60 Hz
consumption	100 VA	160 VA
emission, immunity	EN/IEC 61326-1, class B	EN/IEC 61326-1, class B
max. capacity	8 x 15 ml	8 x 15 ml
max. RPM	6,000 min ⁻¹	8,000 min ⁻¹
max. RCF	3,461	6,153
running time	1 – 99 min, ∞ continuous run, short cycle mode (impulse button)	1 – 99 min, ∞ continuous run, short cycle mode (impulse button)
dimensions (WxDxH)	261 x 353 x 228 mm	261 x 353 x 228 mm
weight	approx. 9 kg	approx. 11 kg
max. noise level	≤ 50 dB (A)	≤ 55 dB (A)
Cat. No.	1800	1802
100 – 127 V 1 ~ / 50–60 Hz	1800-01	1802-01
emission, immunity	FCC class B	FCC class B

*) Other voltages on request.

AVAILABLE ROTORS

ANGLE ROTORS	angle	max. RPM	max. capacity	Cat. No	page
 angle rotor, 8-place	33°	6,000 min ⁻¹ 8,000 min ⁻¹	8 x 15 ml	INCLUSIVE	3

— ANGLE ROTOR, 8-PLACE



Rotor

max. RPM	EBA 200 200 S	6,000 min ⁻¹ 8,000 min ⁻¹
max. RCF		3,461 6,153
max. capacity		8 x 15 ml
run up / run down, braked in sec		17 / 12 37 / 17
angle		33°
Cat. No.		INCLUSIVE



Vessels

capacity in ml		0.5	1.5	2	4	5	6	15	1.1–1.4	2.6–3.4	2.7–3	4.5–5	4.9	7.5–8.5	9–10	5
Ø x L in mm		10.7x36	11 x 38	11 x 38	10x88	12x75	12x82	17x100	8x66	13x65	11x66	11x92	13x90	15x92	16x92	17x59
max. RCF ²⁾	EBA 200	2,214	2,173	2,173	2,817	2,697	2,697	3,461	2,697	2,697	2,697	3,461	3,461	3,461	3,461	2,576
max. RCF ²⁾	EBA 200 S	3,935	3,864	3,864	5,009	4,794	4,794	6,153	4,794	4,794	4,794	6,153	6,153	6,153	6,153	4,579
radius in mm		55	54	54	70	67	67	86	67	67	67	86	86	86	86	64
Cat. No.		Pediatric	microliter tubes		tubes ²⁾			blood collection / urine tubes							-	



Adapter

boring Ø x L in mm	11 x 35	11 x 35	11 x 35	11.5 x 67.5	13.5 x 60	13.5 x 60	17.7 x 88	13.5 x 60	13.5 x 60	13.5 x 60	17.7 x 88	17.7 x 88	17.7 x 88	17.7 x 88	17 x 25
vessels per rotor	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8
Cat. No.	1063-8 (8 pcs.)			6305	1054-A	1054-A	-	1054-A	1054-A	1054-A	-	-	-	-	1064

Vessels

capacity in ml	15	10	1.6–5	4–7	8	8.5–10	12
Ø x L in mm	17 x 120	15 x 102	13 x 75	13 x 100	16 x 125	16 x 100	17 x 102
max. RCF ²⁾	EBA 200	3,461	3,461	2,697	3,461	3,461	3,461
max. RCF ²⁾	EBA 200 S	6,153	6,153	4,794	6,153	6,153	6,153
radius in mm	86	86	86	86	86	86	86
Cat. No.	-	blood collection / urine tubes					



Adapter

boring Ø x L in mm	17.7 x 88	17.7 x 88	13.5 x 60	13.5 x 79	17.7 x 88	17.7 x 88	17.7 x 88
vessels per rotor	4	8	8	8	4	8	4
Cat. No.	-	-	1054-A	1058	-	-	-

2) Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote 2) is 4,000.

12) For the EBA 200 only! Model EBA 200 S is supplied as standard with an adapter for these tubes.

CERTIFICATIONS / REGISTRATIONS

Product certification:



Product registration:



Company certifications:



DOWNLOADS

[↓ Operating Manual – EBA 200 | 200 S](#)
[↓ General Catalog](#)

DIMENSIONS – EBA 200 | 200 S

