


NON-INVASIVE MASK



ELBOW CONNECTOR

- **QUICK RELEASE**
Quick release design allows patients to disconnect and connect to the tubing with minimum effort, increasing patient mobility.
- **MONITORING**
Standard 4.3mm port is located on the connector to fulfill the demands of patient monitoring.
- **SWIVEL CONNECTION**
Swivel connector allows clinician to adjust the angle according to the setting of the ventilation circuit.

FRAME & MASK

- **OPTIMAL FIT**
The reinforcing frame incorporates with the silicone mask to prevent deformation. The soft-touch cushion helps seal the mask without irritating the face.
- **EASY CLEAN**
The connection between elbow connector and frame can be disassembled easily to allow effective cleansing.
- **STERILIZATION**
High quality medical grade materials are autoclavable to 121 °C.

HEADGEAR

- **SIMPLE CLIP**
The clip-in design allows patients to lock the headgear into the frame and remove it efficiently.
- **GOOD FIT**
The adjustable strap made with flexible fabric is adaptable to various patients. The permeable mesh on the backside of headgear increases patient comfort.

ORDER INFORMATION

REF	DESCRIPTION	SIZE	PACKAGE
1701	NIV Face Mask	L	1set/box, 20sets/carton
1702	NIV Face Mask	M	1set/box, 20sets/carton
1703	NIV Face Mask	S	1set/box, 20sets/carton