

To whom it may concern

Manufacturer's Declaration

Date: 30/05/2025

We hereby declare that PRO-300 Tissue Processor tissue processor is equipped with the following features:

- S8** • The processing chamber is resistant to routine processing reagents such as alcohol, xylene, toluene, etc., but not to highly acidic reagents (which are never used in histology).
- S14** • The dehydration process is automatic. Then the tissues must be embedded in embedding stations (to create paraffin blocks).
- S10** • Various types of tissues, up to 10 mm thick, can be processed according to the processing procedure (time and reagents)



Histo-Line Laboratories S.r.l.

Name: Mr. Fabrizio Illuminati
Position: Director
Histo Line Laboratories Srl