





HONJIN



Chimb	≤3%	
Open Force	≤35%	
Sterilization Force	Sterilization about 30min under 121±0.2°C, not change color, and will not crack	
Fastness of aluminum coat	Brush the aluminum Surface by 80% Ethanol about 30 seconds, then brush the surface by 70% isopropanol about 30 seconds, the aluminum coat not appear abrasion	
Material	Aluminum	
Size	Inner Diameter	34±0.1 mm
	Thickness	0.30±0.02 mm
	Inner Height	10.30±0.10 mm
Produce basis: YBB00082005-2015		

Shanxi Hongjin Medical Glass Co.,Ltd				Title: Aluminum seals with three-bridge tabs Ø34mm			
DRAWING BY:		UNIT:	GEN. TOLERANCE:		SIZE: A4	REV:0	
Hou Mingzhen	APPROVAL	DEPARTMENT	NAME - SURNAME	SIGN	DATE	DWG NO: HJ-AS-34	
		FACTORY MANAGER	Yang Wen			FILE:	
		QA / RA	Pei Zhuo			SCALE:1/1	PAGE:1
Date: 2022/09/20							