



PRODUCT DESCRIPTION

1. Product Name and Product Code:

AL-080189.V001 Bacterial Viral Filter with Soft Cap (Sterile EO)





2. Manufacturer:

Meditera Tibbi Malzeme Sanayi ve Ticaret A.Ş.

3. Sterilization Method:

Ethylene oxide sterilization.

4. Packaging:

Tyvek packaging.

5. Intended Use and Functional Description:

6. AL-080189.V001 manufactured by Meditera, is a disposable, single-use barrier type, bi-directional fitter fabricated from a plastic resin that is supplied to the customer packaged and non-sterile. Fabricated with a filtering medium that is highly effective in reducing the numbers of bath bacterial and viral contaminates that may be present in a patient's exhaled gas, the device's design also minimizes the resistance to air flow. The product is intended to protect both the patient and machine from potential transmission of pathogens by droplets and aerosols.

The filter is transparent and has SBC raw material.

It has been manufactured for single patient use and maximum recommended hours of use is 24 hours.



PRODUCT DESCRIPTION

7. Shelf life:

5 years from the date of production .

8. Product Description:

8.1. Classification: Class IIa Rule 2

8.2. This product does not contain any metalic parts

8.3. Spesification of the product :

Filter Code	Viral efficiency	Bacterial efficiency	Resistance to flow			Weight (gr)	Dead Space (ml)	Moisture output at 500 ml Vt	Tidal Volume (ml)
			30 lt/min	60 lt/min	90 lt/min				
AL-080189	>99.999%	>99.9999%	4,3mm H ₂ O	12 mm H ₂ O	25 mm H ₂ O	19	33	NA	150- 1500



BACTERIAL	-VIRAL	FII	TER

	SBC raw material			
Code	AL-08018			
Sterile Code	AL-080189			
Qnty/Box	50			
Tidal volume (ml)	150-1500			
Dead Space (ml)	33			
Bacterial Efficiency	>99,9999%			
Viral Efficiency	>99,999%			
Resistance to Flow	30 L/min 4,3 mm H ₂ O			
	60 L/min 12 mm H ₂ O			
	90 L/mîn 25 mm H ₂ O			
Weight (gr)	19			
Fittings	22mmM - 15/22mmMF			



PRODUCT DESCRIPTION

9. Standarts:

EN ISO 13485:2003, 93/42/EEC Medical Device Directive, EN ISO 9360-1:2009, EN ISO 9360-2:2009, EN ISO 11607-1:2009, EN ISO 11607-2:2009, EN ISO 1135-1:2007, EN ISO 23328-1, EN ISO 23328-2

10. Waste Method:

After use product will be "contaminated medical waste" and package will be "packaging waste" so they should be handled according to relevant national and international standards and regulations.

MEDITER TIC. A.S.

Ibni Meke OSB M. TO SI TOL 4 Sk. No.29 35900 Tire-IZMIR

Tire V.D: 960 029 8703 Tic. Sicil No: 2346

Mersis No: 0060 0298 7030 0016

Tel: 0 232 5/3 51 10 Foxi 0 232 913 91 14