

BALLOON INFLATION DEVICES – GUN TYPE

Features

- Easy to inflate even at high temperatures
- Pressure is released quickly
- The modified humanized design fits comfortably when operation



Product No.	Volume (ml)	Contents
42.04.40053	20	Balloon inflation device, Stopcock("Right OFF")+Hemostasis Valve Set(Push pull with 25cm Extension line & Stopcocks)

CONNECTING TUBINGS (1200PSI)

www.scwmed.com

Features

- Protects against microbiological infect
- ganoid; smooth-faced
- succinct design
- High pressure resistant
- Easy operation

Product No.	Description
-------------	-------------

42.02.20005 180 cm Female Luer Lock to Male Rotator



HEMOSTASIS VALVE SETS

Features

- Various types, to meet the clinician's habits
- Pressure resistance and good sealing valve to reduce the blood loss.

Contents

- Hemostasis Valve
- Guide Wire Insertion Tool
- Torque Device



Product No.	Description
42.03.10000	Twist Y Connector

HEMOSTASIS VALVE SETS

Features

- Various types to meet the clinician's habits
- Pressure resistance and good sealing valve to reduce the blood loss.

Contents

- Hemostasis Valve
- Guide Wire Insertion Tool
- Torque Device



Product No.	Description
42.03.10005	Twist Y Connector with 25 cm extension line and stopcock

MANIFOLDS-THREE PORTS

Features

- Large internal lumen (2.3mm diameter)
- Easy to rotation
- Latex free

Product No.	Description	Pressure
42.05.40003	"OFF" Right-Handed w/Rotating Adapter, 3 Port	250 PSI



BALLOON INFLATION DEVICES

Features

- Accurate reading, reliable inflation and deflation.
- Modified humanized design handle facilitates doctor's operation.
- Locking design enables one handed operation in deflation.



Product No.	Volume (ml)	Contents
42.04.40007	30	Balloon inflation device, Stopcock (RIGHT OFF)+Hemostasis valve set(Push Pull with 25cm extension line&stopcock)

Angiographic Syringes

- Clear body Provides excellent clarification and smooth as glass feeling, it could see air bubble and contrast medium
- It could be used for computed tomographic scanning, angiography, Nuclear magnetic resonance and imaging
- Great Variety of Syringes, it could match so many injectors.



Product No.	Volume(ml)	Description	Models
Single		Used for Medrad	Single
42.16.10000	200	Used for CT: MCT+ OP-100;VCT610; ECT710	CT-200-A1

Features

- Selective compression of the radial artery to allow blood return and preserve patency.
- Transparent structure for visual control of the puncture side
- Two band lengths available for a precise fit and optimal patient comfort

Product No.	Size (cm)
42.17.30001	24cm
42.17.30003	29cm



HEMOSTASIS VALVE SETS

Features

- Various types to meet the clinician's habits
- Pressure resistance and good sealing valve to reduce the blood loss.

Contents

- Hemostasis Valve
- Guide Wire Insertion Tool
- Torque Device



Product No.	Description
42.03.10003	Push pull Y Connector with 45 cm extension line and stopcock

Pressure Line

1200 PSI

- ▶ For pressures up to 83 bar (1200 PSI)
- ▶ Ø 1.7 x 4.9 mm tubing
- ▶ Male luer lock and female luer lock connectors



Rotating Male Connector



Ref				Length
		Rotating		
M/F	M/M	M/F	M/M	
135 3025 1	135 4025 1	135 1025 1	135 2025 1	25 cm
135 3050 1	135 4050 1	135 1050 1	135 2050 1	50 cm
135 3075 1	135 4075 1	135 1075 1	135 2075 1	75 cm
135 3100 1	135 4100 1	135 1100 1	135 2100 1	100 cm
135 3125 1	135 4125 1	135 1125 1	135 2125 1	125 cm

Angiographic Syringes

- Clear body Provides excellent clarification and smooth as glass feeling ,it could see air bubble and contrast medium
- It could be used for computed tomographic scanning, angiography, Nuclear magnetic resonance and imaging
- Great Variety of Syringes, it could match so many injectors



Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for Medrad	Single	Dual
42.16.10005		200	Used for DSA: Mark V & Mark V Provis	DSA-200-A1	

Lot 213

DOSE-CONTROL SYRINGES



Features

- Clear body provides excellent clarification and smooth as glass feeling
- Solid piston barrel Maintains stable and persistent under pressure
- Rotating and Fixed Luers , in luer connection. It allows for flexibility and dependence
- Safety of space ,Intended to reduce the potential of air bubbles getting into the catheter
- Safety of reservoir,available with 0.5 mL taking precautions against bottoming out during injection

Product No.	Description	Tip Style
42.07.10002	12ml Thumb Ring&Finger Rings	Rotator Tip

