





### GUANGDONG BAIHE MEDICAL TECHNOLOGY CO.,LTD.

Address: No. 89, Taoyuan East Road, Nanhai, Foshan, Guangdong, China

Tel: +86-757-81207300 ext. 881; or ext 276 Fax: +86-757-81207311 Email: raymond@baihemedical.com

info@baihemedical.eu http://www.baihemedical.eu







Provide comprehensive medical solutions.



# Company Profile

Established in 1999, Baihe Medical is one of the leading medical device manufacturers in the South of China, specialized in medical disposables for use in anesthesia, blood purification, infusion and wound care.

With total 2000 employees and land occupation over 40,000 square meters, Baihe Medical's headquarter locates in Southern Medical Industrial park, Foshan City, Guangdong Province, 200 kilometers away from Hong Kong. Baihe Medical currently operates 6 manufacturing plants, meanwhile she has 4 wholly-owned subsidiary companies, including Biosun Medical, United Medical, Jiakang Medical and Dihua Precision Moulding.

The Quality Assurance System is certified and correspondent with ISO-13485, ISO 9001. More than 20 products are certified with CE mark. Central Venous Catheter and Hemodialysis Catheter can be marketed in USA. Baihe is the biggest manufacturer of two above mentioned products in terms of production capacity in China.

The company upholds "Focusing on Quality, Respect for Life", dedicating to becoming one of the leading suppliers of medical device.

























百合医疗2019-8-1.indd 3-4







Movable clamp allows anchorage at puncture site regardless of depth of catheter, which minimizes trauma and irritation to the puncturing site. Depth marking assists in accurate placement of central venous catheter from the right or left subclavian or jugular vein. Soft tip reduces trauma to vessel, minimizing vessel erosion, hemothorax and cardiac tamponade. Single, double, triple and quad lumen is available for choice.

Radiopacity facilitates confirmation of catheter placement. The multi-lumen version's tips are more radiopaque, which makes it easy to confirm the placement of fluoroscopic tip. Vessel dilator ensures "super soft" catheters to be easily placed percutaneously.

### Optional

Triangle Hub

Puncture Needle.

Alternative



Straight or Y-shaped Type Nitinol or Stainless Steel



Guide-wire with Flexible J-tip



Sharp Scalpel





Compound-packaging Tray

Various Accessories for Choice

#### Standard Kits include:

- I. Central Venous Catheter
- 2. Guide-wire
- 3. Vessel Dilator
- 4. Clamp
- 5. Fastener: Catheter Clamp
- 6. Introducer Needle
- 7. Introducer Syringe
- 8. Injection Needle
- 9. Injection Cap

### Optional Compound Kits include:

- 1. Central Venous Catheter Standard Kit Accessories
- 2. 5ml Syringe
- 3. Surgical Gloves
- 4. Surgical Pledget
- 5. Surgery Sheet
- 6. ECG cable
- 7. Surgery Towel
- 8. Sterile Brush
- 9. Gauze Pad
- 10. Suture of Needle
- 11. Wound Dressing
- 12. Scalpel

**Anesthetic Products** 

百合医疗2019-8-1.indd 7-8 2020/1/18 14:16:36

### Order Information

Order	Intol	rmati	on			
Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV-1121	Single	24Ga	5	4F*5cm	0.018"×43.5cm	Straight/Y
FV-1122	Single	24Ga	8	4F*5cm	0.018"×43.5cm	Needle
FV-11299	Single	24Ga	9	4F*5cm	0.018"×43.5cm	
FV-1123	Single	24Ga	10	4F*5cm	0.018"×43.5cm	20Ga
FV-112912	Single	24Ga	12	4F*5cm	0.018"×43.5cm	
FV-1124	Single	24Ga	13	4F*5cm	0.018"×43.5cm	24Ga
FV-112915	Single	24Ga	15	4F*5cm	0.018"×43.5cm	
FV-1126	Single	24Ga	20	4F*5cm	0.018"×53.5cm	
FV-1221	Single	22Ga	5	4F×5cm	0.018"×50cm	Straight/Y
FV-1222	Single	22Ga	8	4F×5cm	0.018"×50cm	-
FV-1223	Single	22Ga	10	4F×5cm	0.018"×50cm	Needle
FV-1224	Single	22Ga	13	4F×5cm	0.018"×50cm	20Ga
FV-122915	Single	22Ga	15	4F×5cm	0.018"×50cm	
FV-1225	Single	22Ga	16	4F×5cm	0.018"×50cm	22Ga
FV-1226	Single	22Ga	20	4F×5cm	0.018"×60cm	
FV-1228	Single	22Ga	30	4F×5cm	0.018"×70cm	
FV-1321	Single	20Ga	5	4F×5cm	0.021"×50cm	Ctraight/V
FV-1322	Single	20Ga	8	4F×5cm	0.021"×50cm	Straight/Y
FV-1323	Single	20Ga	10	4F×5cm	0.021"×50cm	Needle
FV-1324	Single	20Ga	13	4F×5cm	0.021"×50cm	20Ga
FV-132915	Single	20Ga	15	4F×5cm	0.021"×50cm	
FV-1325	Single	20Ga	16	4F×5cm	0.021"×50cm	20Ga
FV-1326	Single	20Ga	20	4F×5cm	0.021"×60cm	
FV-1328	Single	20Ga	30	4F×5cm	0.021"×70cm	
FV-1421	Single	18Ga	5	5.5F×6cm	S 0.025"×50cm	
1 4-1421	Sirigic	1000	3	5.51 AGCITI	Y 0.021"×50cm	
FV-1422	Single	18Ga	8	5.5F×6cm	\$ 0.025"×50cm	
1 4-1422	311 IGIC	1000	O	5.51 ×6CIII	Y 0.021"×50cm	
FV-1423	Single	18Ga	10	5.5F×6cm	\$ 0.025"×50cm	Straight/Y
1 4-1425	311 IGIC	1000	10	3.31 ×6CIII	Y 0.021"×50cm	Needle
FV-1424	Single	18Ga	13	5.5F×6cm	\$ 0.025"×50cm	19Ga/20Ga
1 4-1424	311 IGIC	1000	10	3.31 ×6CIII	Y 0.021"×50cm	1960/2060
FV-142915	Single	18Ga	15	5.5F×6cm	\$ 0.025"×50cm	
1 4-142713	311 IGIC	1000	10	5.51 ×6CIII	Y 0.021"×50cm	(100
FV-1425	Single	18Ga	16	5.5F×6cm	\$ 0.025"×50cm	18Ga
1 7-1425	311 IGIC	1000	10	5.51 ×6CIII	Y 0.021"×50cm	
FV-1426	Single	18Ga	20	5.5F×6cm	\$ 0.025"×60cm	
1 4-1420	Jii igie	1000	20	0.01 ~00111	Y 0.021"×60cm	
FV-1428	Single	18Ga	30	5.5F×6cm	\$ 0.025"×70cm	
1 1 1420	Jirigic	1000	50	3.51 · OCITI	Y 0.021"×70cm	

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introduce Needle
FV-1521	Single	16Ga	5	6F×10cm	0.032"×50cm	Ctrain a lat ()
FV-1522	Single	16Ga	8	6F×10cm	0.032"×50cm	Straight/Y
FV-1523	Single	16Ga	10	6F×10cm	0.032"×50cm	Needle
FV-1524	Single	16Ga	13	6F×10cm	0.032"×50cm	18Ga
FV-152915	Single	16Ga	15	6F×10cm	0.032"×50cm	
FV-1525	Single	16Ga	16	6F×10cm	0.032"×50cm	16Ga
FV-1526	Single	16Ga	20	6F×10cm	0.032"×60cm	
FV-1528	Single	16Ga	30	6F×10cm	0.032"×70cm	
FV-1624	Single	14Ga	13	8F×10cm	0.032"×50cm	Straight/Y
FV-162915	Single	14Ga	15	8F×10cm	0.032"×50cm	Needle 18Ga
FV-1625	Single	14Ga	16	8F×10cm	0.032"×50cm	
FV-1626	Single	14Ga	20	8F×10cm	0.032"×60cm	14Ga
FV-1628	Single	14Ga	30	8F×10cm	0.032"×70cm	
FV/FC-2421	Double	4Fr	5	5.5F×6cm	0.018"×50cm	
FV/FC-2422	Double	4Fr	8	5.5F×6cm	0.018"×50cm	Straight/Y
FV/FC-2423	Double	4Fr	10	5.5F×6cm	0.018"×50cm	Needle
FV/FC-2424	Double	4Fr	13	5.5F×6cm	0.018"×50cm	20Ga
FV/FC-242915	Double	4Fr	15	5.5F×6cm	0.018"×50cm	
FV/FC-2425	Double	4Fr	16	5.5F×6cm	0.018"×50cm	22
FV/FC-2426	Double	4Fr	20	5.5F×6cm	0.018"×60cm	22
FV/FC-2428	Double	4Fr	30	5.5F×6cm	0.018"×70cm	
FV/FC-2521	Double	5Fr	5	6F×7cm	0.021"×50cm	
FV/FC-2522	Double	5Fr	8	6F×7cm	0.021"×50cm	Straight/Y
FV/FC-2523	Double	5Fr	10	6F×7cm	0.021"×50cm	Needle
FV/FC-2524	Double	5Fr	13	6F×7cm	0.021"×50cm	20Ga
FV/FC-252915	Double	5Fr	15	6F×7cm	0.021"×50cm	
FV/FC-2525	Double	5Fr	16	6F×7cm	0.021"×50cm	18
FV/FC-2526	Double	5Fr	20	6F×7cm	0.021"×60cm	20
FV/FC-2528	Double	5Fr	30	6F×7cm	0.021"×70cm	
FV/FC-2724	Double	7Fr	13	8F×10cm	0.035"×50cm	Straight/Y
FV/FC-272915	Double	7Fr	15	8F×10cm	0.035"×50cm	Needle 18Ga
FV/FC-2725	Double	7Fr	16	8F×10cm	0.035"×50cm	
FV/FC-2726	Double	7Fr	20	8F×10cm	0.035"×60cm	14
FV/FC-2728	Double	7Fr	30	8F×10cm	0.035"×70cm	18

## Order Information

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV-2824	Double	8Fr	13	9F×10cm	0.035"×50cm	Straight/Y
FV-282915	Double	8Fr	15	9F×10cm	0.035"×50cm	Needle 18Ga
FV-2825	Double	8Fr	16	9F×10cm	0.035"×50cm	14
FV-2826	Double	8Fr	20	9F×10cm	0.035"×60cm	14
FV-2828	Double	8Fr	30	9F×10cm	0.035"×70cm	14
FV-2921	Double	8.5Fr	5	9F×10cm	0.035"×50cm	
FV-2922	Double	8.5Fr	8	9F×10cm	0.035"×50cm	Straight/Y
FV-2923	Double	8.5Fr	10	9F×10cm	0.035"×50cm	Needle
FV-2924	Double	8.5Fr	13	9F×10cm	0.035"×50cm	18Ga
FV-2925	Double	8.5Fr	16	9F×10cm	0.035"×50cm	14
FV-2926	Double	8.5Fr	20	9F×10cm	0.035"×60cm	14
FV-2927	Double	8.5Fr	25	9F×10cm	0.035"×60cm	14
FV-2928	Double	8.5Fr	30	9F×10cm	0.035"×70cm	
FV/FC-3421	Triple	4.5Fr	5	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3422	Triple	4.5Fr	8	6F×5cm/6F×7cm	0.021"×50cm	Straight/Y
FV/FC-3423	Triple	4.5Fr	10	6F×5cm/6F×7cm	0.021"×50cm	Needle
FV/FC-3424	Triple	4.5Fr	13	6F×5cm/6F×7cm	0.021"×50cm	20Ga
FV/FC-342915	Triple	4.5Fr	15	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3425	Triple	4.5Fr	16	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3426	Triple	4.5Fr	20	6F×5cm/6F×7cm	0.021"×60cm	23 23
FV/FC-3428	Triple	4.5Fr	30	6F×5cm/6F×7cm	0.021"×70cm	
FV/FC-3521	Triple	5.5Fr	5	6.5F×10cm	0.021"×50cm	
FV/FC-3522	Triple	5.5Fr	8	6.5F×10cm	0.021"×50cm	Straight/Y
FV/FC-3523	Triple	5.5Fr	10	6.5F×10cm	0.021"×50cm	Needle
FV/FC-3524	Triple	5.5Fr	13	6.5F×10cm	0.021"×50cm	20Ga
FV/FC-352915	Triple	5.5Fr	15	6.5F×10cm	0.021"×50cm	
FV/FC-3525	Triple	5.5Fr	16	6.5F×10cm	0.021"×50cm	20
FV/FC-3526	Triple	5.5Fr	20	6.5F×10cm	0.021"×60cm	22 22
FV/FC-3528	Triple	5.5Fr	30	6.5F×10cm	0.021"×70cm	

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV/FC-3724	Triple	7Fr	13	8F×10cm	0.035"×50cm	Straight/Y
FV/FC-372915	Triple	7Fr	15	8F×10cm	0.035"×50cm	Needle 18Ga
FV/FC-3725	Triple	7Fr	16	8F×10cm	0.035"×50cm	16
FV/FC-3726	Triple	7Fr	20	8F×10cm	0.035"×60cm	
FV/FC-3728	Triple	7Fr	30	8F×10cm	0.035"×70cm	18 18
FV/FC-3924	Triple	8.5Fr	13	9F×10cm	0.035"×50cm	Straight/Y
FV/FC-392915	Triple	8.5Fr	15	9F×10cm	0.035"×50cm	Needle 18Ga
FV/FC-3925	Triple	8.5Fr	16	9F×10cm	0.035"×50cm	
FV/FC-3926	Triple	8.5Fr	20	9F×10cm	0.035"×60cm	16 16
FV/FC-3928	Triple	8.5Fr	30	9F×10cm	0.035"×70cm	14
FV-4921	Quad	8.5Fr	5	9F×10cm	0.035"×50cm	
FV-4922	Quad	8.5Fr	8	9F×10cm	0.035"×50cm	Straight/Y Needle 18Ga
FV-4923	Quad	8.5Fr	10	9F×10cm	0.035"×50cm	
FV-4924	Quad	8.5Fr	13	9F×10cm	0.035"×50cm	
FV-492915	Quad	8.5Fr	15	9F×10cm	0.035"×50cm	1/
FV-4925	Quad	8.5Fr	16	9F×10cm	0.035"×50cm	(B) (B)
FV-4926	Quad	8.5Fr	20	9F×10cm	0.035"×60cm	14
FV-4928	Quad	8.5Fr	30	9F×10cm	0.035"×70cm	

#### Remarks:

- 1. The product codes in this chart are for CVC kits with straight introducer needleconfiguration.
- 2. CVC kits with Y-shaped introducer needle configuration are also available in adulttypes. The product code of CVC kit with Y-shaped introducer needle configurationis to add letter "Y" after the product code of CVC kit with straight introducer needleconfiguration, for example, "FC-2726Y". Before you order, please clarify which configuration you request, either straight or Y-shape introducer needle.
- 3. Pediatric types of CVC include single lumen 22Ga,20Ga,18Ga;double lumen 4Fr,5Fr;triple lumen 4. 5Fr and triple lumen 5.5Fr
- 4. Adult types of CVC include single lumen 16Ga,14Ga;double lumen 7Fr,8Fr,8.5Fr;triple lumen 7Fr,8. 5Fr and quad lumen 8.5Fr.
- 5. "FV"in the product code for single lumen CVC means with cylinder hub.
- 6. "FC"in the product code for double lumen, triple lumen, quad lumen CVC means with triangular hub." FV"in the product code for double lumen, triple lumen, quadlumen CVC means with round hub.
- 7. CVC kits with ECG cable are also available, the product code is to add letter " ECG" after the codes.

8. The Guidewire and introducer needle can be customized.

百合医疗2019-8-1.indd 9-10 2020/1/18 14:16:37