

PP Box

Use medical-grade transparent polypropylene material. ·temperature range from- 90 °C to 121 °C, can be used for low temperature storage or autoclave. ·Perfect match low temperature freezer, avoid appearing sequence disorder, greatly improve space utilization . ·Transparent cover designed for easy observation of the sample in the box. ·Cover with grid and code for easy inventory management and error free access. ·Available in five colors: transparent, pink, blue, green and yellow. ·There are gaps for the surrounding of the box, which facilitates rapid cooling of the sample and liquid discharge, and is not easy to frost. ·Angle special design to avoid mis-operation of box cover. ·Standard size design, suitable for all standard storage racks.

Model ZKPP-100-2 (100*2 inch freezer tubes)

