


*Take a walk with us...*





## table of content

### hip stem systems

Aida® short stem	4
EcoFit®	6
DiaLoc®	16
ic-straight stem Type I and Type II	18
Actinia®	21
BethaLoc®	23
AJS®	25
Cortina	27
Bicana®	30
MUTARS® RS revision system	32

### hip cup systems

EcoFit®	40
EcoFit® SC	44
cup inserts	48
EcoFit® 2M	50
EcoFit® 2M SP	51
MUTARS® RS cup / PRS	54
ic-reconstructions implants	56
MUTARS® LUMiC®	57
PE-hip cup Mueller II	60
ic-bipolar head	64
spongiosa screw flat head	66
ic-head	68
EPORE® acetabular spacer	71
trochanter plate	72

### accessories


intramedullary plug	73
ic-cement restrictor resorbable	73
ic-cerclage / ic-titanium cerclage wire	73
implabond®, ic Lavage, cempadic®	75

### single instruments

extractor adapter	77
broach handles	78
stem extractor	80
alignment guide	81
stem impactor	82

---





**hip stems**



## Aida® hip stem cementless

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® Duo, cpTi- and HA-coating  
acc. to ISO 13779-2

### REF


standard:

8007-0000	size 0
8007-0001	size 1
8007-0002	size 2
8007-0003	size 3
8007-0004	size 4
8007-0005	size 5
8007-0006	size 6
8007-0007	size 7
8007-0008	size 8


lateralized:

8007-0010	size 0
8007-0011	size 1
8007-0012	size 2
8007-0013	size 3
8007-0014	size 4
8007-0015	size 5
8007-0016	size 6
8007-0017	size 7
8007-0018	size 8


# instruments



**Aida® easy lock container**  
8007-7999



**Aida® easy lock GIS® container**  
8007-7997



**Aida® container II**  
8007-7993

## container without image:

**Aida® container III**  
8007-7994

**Aida® broach handle GIS® left/right container**  
8007-7995

**Aida® broach handle double offset left/right container**  
8007-7992



## EcoFit® short stem cementless cpTi

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® cpTi-coating

### REF

standard:	
3032-0625	size 6,25
3032-0075	size 7,5
3032-0875	size 8,75
3032-0100	size 10
3032-1125	size 11,25
3032-0125	size 12,5
3032-1375	size 13,75
3032-1500	size 15
3032-0175	size 17,5
3032-0200	size 20
lateralized:	
3042-0625	size 6,25
3042-0075	size 7,5
3042-0875	size 8,75
3042-0100	size 10
3042-1125	size 11,25
3042-0125	size 12,5
3042-1375	size 13,75
3042-1500	size 15
3042-0175	size 17,5
3042-0200	size 20

## EcoFit® short stem cementless cpTi and HA

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® Duo, cpTi-and HA coating  
acc. to ISO 13779-2

### REF

standard:	
3043-0625	size 6,25
3043-0075	size 7,5
3043-0875	size 8,75
3043-0100	size 10
3043-1125	size 11,25
3043-0125	size 12,5
3043-1375	size 13,75
3043-1500	size 15
3043-0175	size 17,5
3043-0200	size 20
lateralized:	
3044-0625	size 6,25
3044-0075	size 7,5
3044-0875	size 8,75
3044-0100	size 10
3044-1125	size 11,25
3044-0125	size 12,5
3044-1375	size 13,75
3044-1500	size 15
3044-0175	size 17,5
3044-0200	size 20





## EcoFit® hip stem cementless cpTi

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® cpTi-coating



### REF

standard:

3039-3062	size 6,25
3039-3075	size 7,5
3039-3087	size 8,75
3039-3100	size 10
3039-3112	size 11,25
3039-3125	size 12,5
3039-3137	size 13,75
3039-3150	size 15
3039-3175	size 17,5
3039-3200	size 20*

lateralized:

3039-4062	size 6,25
3039-4075	size 7,5
3039-4087	size 8,75
3039-4100	size 10
3039-4112	size 11,25
3039-4125	size 12,5
3039-4137	size 13,75
3039-4150	size 15
3039-4175	size 17,5
3039-4200	size 20*



## EcoFit® hip stem 133° cementless cpTi

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® cpTi-coating

### REF

standard:

3037-3062	size 6,25
3037-3075	size 7,5
3037-3087	size 8,75
3037-3100	size 10
3037-3112	size 11,25
3037-3125	size 12,5
3037-3137	size 13,75
3037-3150	size 15
3037-3162	size 16,25
3037-3175	size 17,5
3037-3187	size 18,75
3037-3200	size 20*

lateralized:

3037-4062	size 6,25
3037-4075	size 7,5
3037-4087	size 8,75
3037-4100	size 10
3037-4112	size 11,25
3037-4125	size 12,5
3037-4137	size 13,75
3037-4150	size 15
3037-4162	size 16,25
3037-4175	size 17,5
3037-4187	size 18,75
3037-4200	size 20*



## EcoFit® hip stem 133° cementless HA

Mat.: implatan®,  $TiAl_6V_4$  acc. to ISO 5832-3 with implaFix® Duo, cpTi- and HA coating acc. to ISO 13779-2

### REF

standard:

3036-3062	size 6,25
3036-3075	size 7,5
3036-3087	size 8,75
3036-3100	size 10
3036-3112	size 11,25
3036-3125	size 12,5
3036-3137	size 13,75
3036-3150	size 15
3036-3162	size 16,25
3036-3175	size 17,5
3036-3187	size 18,75
3036-3200	size 20*

lateralized:

3036-4062	size 6,25
3036-4075	size 7,5
3036-4087	size 8,75
3036-4100	size 10
3036-4112	size 11,25
3036-4125	size 12,5
3036-4137	size 13,75
3036-4150	size 15
3036-4162	size 16,25
3036-4175	size 17,5
3036-4187	size 18,75
3036-4200	size 20*



## EcoFit® hip stem coxa vara cementless cpTi

Mat.: implatan®,  $TiAl_6V_4$  acc. to ISO 5832-3 with implaFix® cpTi coating

### REF

standard:

3038-8062	size 6,25
3038-8075	size 7,5
3038-8087	size 8,75
3038-8100	size 10
3038-8112	size 11,25
3038-8125	size 12,5
3038-8137	size 13,75
3038-8150	size 15
3038-8175	size 17,5
3038-8200	size 20



## EcoFit® hip stem coxa vara cementless HA

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3 with implaFix® Duo, cpTi - and HA-coating acc. to ISO 13779-2

### REF

standard:	
3038-2062	size 6,25
3038-2075	size 7,5
3038-2087	size 8,75
3038-2100	size 10
3038-2112	size 11,25
3038-2125	size 12,5
3038-2137	size 13,75
3038-2150	size 15
3038-2175	size 17,5
3038-2200	size 20*



## EcoFit® hip stem cementless HA

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3 with implaFix® Duo, cpTi- and HA-coating acc. to ISO 13779-2

### REF

standard:	
3038-3062	size 6,25
3038-3075	size 7,5
3038-3087	size 8,75
3038-3100	size 10
3038-3112	size 11,25
3038-3125	size 12,5
3038-3137	size 13,75
3038-3150	size 15
3038-3175	size 17,5
3038-3200	size 20*
lateralized:	
3038-4062	size 6,25
3038-4075	size 7,5
3038-4087	size 8,75
3038-4100	size 10
3038-4112	size 11,25
3038-4125	size 12,5
3038-4137	size 13,75
3038-4150	size 15
3038-4175	size 17,5
3038-4200	size 20*



## EcoFit® hip stem cementless T-HA

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3  
with *implaFix*® Duo, *cpTi* - and Total HA-coating  
acc. to ISO 13779-2

### REF

standard:

3038-6062	size 6,25
3038-6075	size 7,5
3038-6087	size 8,75
3038-6100	size 10
3038-6112	size 11,25
3038-6125	size 12,5
3038-6137	size 13,75
3038-6150	size 15
3038-6175	size 17,5
3038-6200	size 20*

lateralized:

3038-7062	size 6,25
3038-7075	size 7,5
3038-7087	size 8,75
3038-7100	size 10
3038-7112	size 11,25
3038-7125	size 12,5
3038-7137	size 13,75
3038-7150	size 15
3038-7175	size 17,5
3038-7200	size 20*



## EcoFit® hip stem cemented\*\*

Mat.: *implavit*®, *CoCrMo* acc. to ISO 5832-4

### REF

standard:

3039-6062	size 6,25
3039-6075	size 7,5
3039-6100	size 10
3039-6125	size 12,5
3039-6150	size 15
3039-6175	size 17,5*

lateralized:

3039-8062	size 6,25
3039-8075	size 7,5
3039-8100	size 10
3039-8125	size 12,5
3039-8150	size 15
3039-8175	size 17,5*



## EcoFit® hip stem cemented highly polished

Mat.: implavit®, CoCrMo acc. to ISO 5832-4 highly polished

### REF

standard:	
3038-0062	size 6,25
3038-0075	size 7,5
3038-0100	size 10
3038-0125	size 12,5
3038-0150	size 15
3038-0175	size 17,5*
lateralized:	
3038-1062	size 6,25
3038-1075	size 7,5
3038-1100	size 10
3038-1125	size 12,5
3038-1150	size 15
3038-1175	size 17,5*



## EcoFit® hip stem cemented 133° highly polished

Mat.: implavit®, CoCrMo acc. to ISO 5832-4 highly polished

### REF

standard:	
3033-0062	size 6,25
3033-0075	size 7,5
3033-0100	size 10
3033-0125	size 12,5
3033-0150	size 15
3033-0175	size 17,5*
lateralized:	
3033-1062	size 6,25
3033-1075	size 7,5
3033-1100	size 10
3033-1125	size 12,5
3033-1150	size 15
3033-1175	size 17,5*

\* article only available on request

\*\*article with TiN coating available (reference code + N)



## EcoFit® hip stem coxa vara cemented

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF

standard:

3039-9062	size 6,25
3039-9075	size 7,5
3039-9100	size 10
3039-9125	size 12,5
3039-9150	size 15
3039-9175	size 17,5*



## EcoFit® hip stem coxa vara cemented highly polished

Mat.: implavit®, CoCrMo acc. to ISO 5832-4 highly polished

### REF

standard:

3038-9062	size 6,25
3038-9075	size 7,5
3038-9100	size 10
3038-9125	size 12,5
3038-9150	size 15
3038-9175	size 17,5*



## EcoFit® hip stem cemented 133°

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF

standard:

3034-6062	size 6,25
3034-6075	size 7,5
3034-6100	size 10
3034-6125	size 12,5
3034-6150	size 15
3034-6175	size 17,5*

lateralized:

3034-8062	size 6,25
3034-8075	size 7,5
3034-8100	size 10
3034-8125	size 12,5
3034-8150	size 15
3034-8175	size 17,5*



## Subsider

Mat.: PMMA (polymethylmethacrylate)  
for the use in conjunction with highly polished  
cemented EcoFit® stem / ic-straight stem Type I


0299-0750	size 2
0299-1000	size 3
0299-1250	size 4
0299-1500	size 5
0299-1650	size 6

\* article only available on request


\*\*article with TiN coating available (reference code + N)




# instruments




**EcoFit® easy lock  
container**  
7999-7043




**EcoFit® easy lock GIS®  
container**  
7999-7049




**EcoFit® container II**  
7999-7074




**EcoFit® basic container**  
7999-6050




**EcoFit® container III**  
7999-7075



**EcoFit® broach handle GIS®  
left/right container**  
7999-7076




**EcoFit® 133° GIS® container II**  
7999-7077




**EcoFit® 133° container II**  
7999-7078


# instruments




**EcoFit® short stem easy lock container**  
7999-7057



**EcoFit® short stem easy lock GIS® container**  
7999-7058



**EcoFit® short stem container with double offset handle**  
7999-7052



**EcoFit® short stem GIS® container II**  
7999-7079

## container without images:

**EcoFit® easy lock GIS® container 133°**

7999-7011

**EcoFit® coxa vara GIS® container II**

7999-7081

**EcoFit® container coxa vara easy lock GIS®**

7999-7070

**EcoFit® easy lock container 133°**

7999-7010

**EcoFit® coxa vara container II**

7999-7082

**EcoFit® coxa vara double offset container II**

7999-7083

**EcoFit® hip stem container 133° with double offset handle**

7999-7056

**EcoFit® 133° double offset container II**

7999-7084

**EcoFit® short stem double offset container II**

7999-7085



## DiaLoc® hip stem cementless

Mat.:  $TiAl_6Nb_7$ , acc. to ISO 5832-11

### REF

standard:

8500-00-3	size -3
8500-00-2	size -2
8500-00-1	size -1
8500-0000	size 0
8500-0001	size 1
8500-0002	size 2
8500-0003	size 3
8500-0004	size 4
8500-0005	size 5
8500-0006	size 6
8500-0007	size 7
8500-0008	size 8
8500-0009	size 9
8500-0010	size 10

lateralized:

8500-0101	size 1
8500-0102	size 2
8500-0103	size 3
8500-0104	size 4
8500-0105	size 5
8500-0106	size 6
8500-0107	size 7
8500-0108	size 8
8500-0109	size 9



## DiaLoc® RS revision stem cementless


Mat.:  $TiAl_6Nb_7$ , acc. to ISO 5832-11

### REF


standard:

8500-0203	size 3
8500-0204	size 4
8500-0205	size 5
8500-0206	size 6
8500-0207	size 7
8500-0208	size 8
8500-0209	size 9
8500-0210	size 10

# instruments



**DiaLoc® container**  
8501-0035



**DiaLoc® RS container**  
8501-0039



**DiaLoc® container -2/-3**  
8501-0050

# ic-straight stem



## ic-straight stem Type I cemented

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF

standard:

3513-0750	size 7,5**
3513-0875	size 8,75
3513-1000	size 10**
3513-1125	size 11,25**
3513-1250	size 12,5**
3513-1375	size 13,75**
3513-1500	size 15**
3513-1625	size 16,25**
3513-1750	size 17,5**
3513-2000	size 20**

lateralized:

3514-0750	size 7,5**
3514-0875	size 8,75
3514-1000	size 10**
3514-1125	size 11,25**
3514-1250	size 12,5**
3514-1375	size 13,75
3514-1500	size 15**
3514-1625	size 16,25
3514-1750	size 17,5
3514-2000	size 20

## ic-straight stem Type II cemented

Mat.: implavit®, CoCrMo acc. to ISO 5832-4


### REF

standard:

3516-0750	size 7,5
3516-0875	size 8,75
3516-1000	size 10
3516-1125	size 11,25
3516-1250	size 12,5
3516-1375	size 13,75
3516-1500	size 15
3516-1625	size 16,25
3516-1750	size 17,5
3516-2000	size 20

lateralized:

3518-0750	size 7,5
3518-0875	size 8,75
3518-1000	size 10
3518-1125	size 11,25
3518-1250	size 12,5
3518-1375	size 13,75
3518-1500	size 15
3518-1625	size 16,25
3518-1750	size 17,5
3518-2000	size 20



# ic-straight stem



## ic-straight stem Type I cemented highly polished incl. Subsider

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF


standard:

3513-0750HP	size 7,5
3513-0875HP	size 8,75
3513-1000HP	size 10
3513-1125HP	size 11,25
3513-1250HP	size 12,5
3513-1375HP	size 13,75
3513-1500HP	size 15
3513-1625HP	size 16,25
3513-1750HP	size 17,5
3513-2000HP	size 20


lateralized:

3514-0750HP	size 7,5
3514-1125HP	size 11,25
3514-1250HP	size 12,5
3514-1375HP	size 13,75


# instruments




**ic-straight stem  
container**  
7516-9907



**ic-straight stem  
container III**  
7516-9910



**ic-straight stem  
container II**  
7516-9909



**ic-straight stem  
broach handle GIS® right/left  
container**  
7516-9911



**ic-straight stem  
broach handle double offset right/left  
container**  
7516-9912

## Actinia® hip stem cementless

Mat.: *implatan®*,  $TiAl_6V_4$  acc. to ISO 5832-3  
*implaFix®* HA-coating acc. to ISO 13779-2



### REF

standard:	
8004-0008	size 8
8004-0009	size 9
8004-0010	size 10
8004-0011	size 11
8004-0012	size 12
8004-0013	size 13
8004-0014	size 14
8004-0015	size 15
8004-0016	size 16
8004-0018	size 18
8004-0020	size 20
lateralized:	
8004-0109	size 9
8004-0110	size 10
8004-0111	size 11
8004-0112	size 12
8004-0113	size 13
8004-0114	size 14
8004-0115	size 15
8004-0116	size 16
8004-0118	size 18
8004-0120	size 20

## Actinia® hip stem cemented

Mat.: *implavit®* CoCrMo acc. to ISO 5832-4




### REF


standard:	
8004-0208	size 8
8004-0209	size 9
8004-0210	size 10
8004-0211	size 11
8004-0212	size 12
8004-0213	size 13
8004-0214	size 14
8004-0215	size 15
8004-0216	size 16
8004-0218	size 18
8004-0220	size 20
lateralized:	
8004-0309	size 9
8004-0310	size 10
8004-0311	size 11
8004-0312	size 12
8004-0313	size 13
8004-0314	size 14
8004-0315	size 15
8004-0316	size 16
8004-0318	size 18
8004-0320	size 20




# instruments




**Actinia®  
container**  
8004-9000




**Actinia® GIS®  
container**  
8004-9001




**Actinia® broach  
container**  
8004-9002




**Actinia® GIS®  
container II**  
8004-9007




**Actinia® broach/broach handle straight  
container**  
8004-9074



**Actinia® broach/broach handle GIS®  
left/right container**  
8004-9073



**Actinia®  
container III**  
8004-9072



**Actinia® broach/broach handle double  
offset left/right container**  
8004-9075



## BethaLoc<sup>®</sup> hip stem cementless HA

Mat.: *implatan<sup>®</sup>*,  $TiAl_6V_4$  acc. to ISO 5832-3  
*implaFix<sup>®</sup>* HA-coating acc. to ISO 13779-2

### REF


standard:

8700-1001	sizee 1
8700-1002	size 2
8700-1003	size 3
8700-1004	size 4
8700-1005	size 5
8700-1006	size 6
8700-1007	size 7
8700-1008	size 8
8700-1009	size 9
8700-1010	size 10


lateralized:

8700-1024	size 4
8700-1025	size 5
8700-1026	size 6
8700-1027	size 7
8700-1028	size 8
8700-1029	size 9
8700-1030	size 10

# instruments



**BethaLoc®  
container**  
8700-1088



**BethaLoc®  
container I**  
8700-1089

## container without images:

**BethaLoc® container II**  
8700-1091



## AJS® hip stem cementless

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® HA-coating  
acc. to ISO 13779-2

### REF

8006-1105	size 1	left
8006-1205	size 2	left
8006-1305	size 3	left
8006-1405	size 4	left
8006-1505	size 5	left
8006-1605	size 6	left
8006-1705	size 7	left
8006-1805	size 8	left
8006-1905	size 9	left
8006-1110	size 1	right
8006-1210	size 2	right
8006-1310	size 3	right
8006-1410	size 4	right
8006-1510	size 5	right
8006-1610	size 6	right
8006-1710	size 7	right
8006-1810	size 8	right
8006-1910	size 9	right




## AJS® hip stem cemented\*\*

Mat.: implavit®, CoCrMo acc. to ISO 5832-4


### REF

8006-2105	size 1	left
8006-2305	size 3	left
8006-2405	size 4	left
8006-2505	size 5	left
8006-2605	size 6	left
8006-2705	size 7	left
8006-2805	size 8	left
8006-2110	size 1	right
8006-2310	size 3	right
8006-2410	size 4	right
8006-2510	size 5	right
8006-2610	size 6	right
8006-2710	size 7	right
8006-2810	size 8	right


# instruments



**AJS® snap lock GIS®  
container**  
8006-7996



**AJS® snap lock  
container**  
8006-7999



**AJS® broach  
container left**  
8006-7997

**AJS® broach  
container right**  
8006-7998

## Cortina hip stem cemented\*\*

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF

standard:

8205-1130	size 1/130	left
8205-2135	size 2/135	left
8205-3140	size 3/140	left
8205-4145	size 4/145	left
8205-5150	size 5/150	left
8205-1170	size 1/170	left
8205-2170	size 2/170	left
8205-3170	size 3/170	left
8205-4170	size 4/170	left
8205-5170	size 5/170	left
8200-1130	size 1/130	right
8200-2135	size 2/135	right
8200-3140	size 3/140	right
8200-4145	size 4/145	right
8200-5150	size 5/150	right
8200-1170	size 1/170	right
8200-2170	size 2/170	right
8200-3170	size 3/170	right
8200-4170	size 4/170	right
8200-5170	size 5/170	right

lateralized:

8206-1130	size 1/130	left
8206-2135	size 2/135	left
8206-3140	size 3/140	left
8206-4145	size 4/145	left
8206-5150	size 5/150	left
8206-1170	size 1/170	left
8206-2170	size 2/170	left
8206-3170	size 3/170	left
8206-4170	size 4/170	left
8206-5170	size 5/170	left
8201-1130	size 1/130	right
8201-2135	size 2/135	right
8201-3140	size 3/140	right
8201-4145	size 4/145	right
8201-5150	size 5/150	right
8201-1170	size 1/170	right
8201-2170	size 2/170	right
8201-3170	size 3/170	right
8201-4170	size 4/170	right
8201-5170	size 5/170	right



\* article available on request

\*\*article with TiN coating available (reference number + N)

## Cortina hip stem (revision) cemented\*\*

Mat.: *implavit*<sup>®</sup>, CoCrMo acc. to ISO 5832-4

### REF

standard:


8205-2200	size 2/200	left
8205-3200	size 3/200	left
8205-4200	size 4/200	left
8205-5200	size 5/200	left
8205-3250	size 3/250	left
8205-4250	size 4/250	left
8205-5250	size 5/250	left
8205-4300	size 4/300	left
8205-5300	size 5/300	left
8200-2200	size 2/200	right
8200-3200	size 3/200	right
8200-4200	size 4/200	right
8200-5200	size 5/200	right
8200-3250	size 3/250	right
8200-4250	size 4/250	right
8200-5250	size 5/250	right
8200-4300	size 4/300	right
8200-5300	size 5/300	right

lateralized:

8206-2200	size 2/200	left
8206-3200	size 3/200	left
8206-4200	size 4/200	left
8206-5200	size 5/200	left
8206-3250	size 3/250	left
8206-4250	size 4/250	left
8206-5250	size 5/250	left
8206-4300	size 4/300	left
8206-5300	size 5/300	left
8201-2200	size 2/200	right
8201-3200	size 3/200	right
8201-4200	size 4/200	right
8201-5200	size 5/200	right
8201-3250	size 3/250	right
8201-4250	size 4/250	right
8201-5250	size 5/250	right
8201-4300	size 4/300	right
8201-5300	size 5/300	right




# instruments




**Cortina basic  
container**  
7200-1100




**Cortina broach  
container left**  
7200-1105



**Cortina broach  
container right**  
7200-1110



**Cortina broach  
container 170 left**  
7200-1705



**Cortina broach  
container 170 right**  
7200-1710

## container without images:

**Cortina broach handle GIS® left/right container**  
7200-1115

**Cortina broach handle double offset left/right container**  
7200-1116

**Cortina basic container size 3**  
7200-1101

**Cortina basic container size 3 and 4**  
7200-1102

**Cortina container size 2/3/4**  
7200-1103

**Cortina broach container 200-300mm right**  
7200-1106

**Cortina broach container 200-300mm left**  
7200-1107

**Cortina broach container 170mm right**  
7200-1108

**Cortina broach container 170mm left**  
7200-1109





## Bicana® hip stem cemented

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF

#### 117°:

640101	size 1/130 left
640103	size 2/130 left
640105	size 3/130 left
640150	size 1/150 left
640107	size 2/150 left
640109	size 3/150 left
640111	size 4/150 left
640113	size 5/150 left
640102	size 1/130 right
640104	size 2/130 right
640106	size 3/130 right
640151	size 1/150 right
640108	size 2/150 right
640110	size 3/150 right
640112	size 4/150 right
640114	size 5/150 right

#### 126°:

640115	size 1/130 left
640117	size 2/130 left
640119	size 3/130 left
640152	size 1/150 left
640121	size 2/150 left
640123	size 3/150 left
640125	size 4/150 left
640127	size 5/150 left
640143	size 6/150 left
640116	size 1/130 right
640118	size 2/130 right
640120	size 3/130 right
640153	size 1/150 right
640122	size 2/150 right
640124	size 3/150 right
640126	size 4/150 right
640128	size 5/150 right
640144	size 6/150 right



## Bicana® hip stem cemented

Mat.: implavit®, CoCrMo acc. to ISO 5832-4

### REF

#### 126° + 10,5:

640129	size 1/130 left
640131	size 2/130 left
640133	size 3/130 left
640154	size 1/150 left
640135	size 2/150 left
640137	size 3/150 left
640139	size 4/150 left
640141	size 5/150 left
640130	size 1/130 right
640132	size 2/130 right
640134	size 3/130 right
640155	size 1/150 right
640136	size 2/150 right
640138	size 3/150 right
640140	size 4/150 right
640142	size 5/150 right

#### 135°:

640301	size 1/130 left
640303	size 2/130 left
640305	size 3/130 left
640320	size 1/150 left
640307	size 2/150 left
640309	size 3/150 left
640311	size 4/150 left
640313	size 5/150 left
640301	size 1/130 right
640304	size 2/130 right
640306	size 3/130 right
640321	size 1/150 right
640308	size 2/150 right
640310	size 3/150 right
640312	size 4/150 right
640314	size 5/150 right



## Bicana<sup>®</sup> hip stem revision cemented

Mat.: implavit<sup>®</sup>, CoCrMo acc. to ISO 5832-4

### REF


#### 126°:

640614	size 3/200 left
640616	size 4/200 left
640618	size 5/200 left
640602	size 3/250 left
640604	size 4/250 left
640606	size 5/250 left
640608	size 3/300 left
640610	size 4/300 left
640612	size 5/300 left
640615	size 3/200 right
640617	size 4/200 right
640619	size 5/200 right
640603	size 3/250 right
640605	size 4/250 right
640607	size 5/250 right
640609	size 3/300 right
640611	size 4/300 right
640613	size 5/300 right


#### 135°:

640626	size 3/250 left
640628	size 4/250 left
640630	size 5/250 left
640632	size 3/300 left
640634	size 4/300 left
640636	size 5/300 left
640627	size 3/250 right
640629	size 4/250 right
640631	size 5/250 right
640633	size 3/300 right
640635	size 4/300 right
640637	size 5/300 right

# instruments



**Bicana® container (top)**  
640-000



**Bicana® container (bottom)**  
640-000

## container without images:

**Bicana® container I**  
6400-0001

**Bicana® container II**  
6400-0002

# MUTARS® RS



## RS proximal component \*S incl. counter screw

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3

### REF

6710-1527	127°/32 mm
6710-1535	135°/32 mm
6710-1627	127°/42 mm
6710-1635	135°/42 mm



## RS metaphyseal component HA

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3  
with *implaFix*® HA-coating acc. to ISO  
13779-2

### REF

6730-4121	40 mm
6730-5121	50 mm
6730-4221	40 mm small
6730-5221	50 mm small



## RS extension piece HA

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3  
with *implaFix*® HA-coating acc. to ISO  
13779-2

### REF

6730-0125	25 mm
-----------	-------



## RS screw

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3

### REF

6720-4008	M8x40 mm
6720-5008	M8x50 mm
6720-6508	M8x65 mm
6720-7508	M8x75 mm
6720-9008	M8x90 mm
6720-1008	M8x100 mm
6720-1158	M8x115 mm
6720-1258	M8x125 mm

\*S: available with silver coating (reference number + S)



## RS stem cementless HA

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3  
with *implaFix*® HA-coating acc. to ISO 13779-2

### REF

6762-1512	12/150 mm
6762-1513	13/150 mm
6762-1514	14/150 mm
6762-1515	15/150 mm
6762-1516	16/150 mm
6762-1517	17/150 mm
6762-1518	18/150 mm
6762-1519	19/150 mm
6762-1520	20/150 mm
6762-1521	21/150 mm
6762-1522	22/150 mm
6762-2012	12/200 mm
6762-2013	13/200 mm
6762-2014	14/200 mm
6762-2015***	15/200 mm
6762-2016***	16/200 mm
6762-2017***	17/200 mm
6762-2018***	18/200 mm
6762-2019***	19/200 mm
6762-2020***	20/200 mm
6762-2021***	21/200 mm
6762-2022***	22/200 mm
6762-2514	14/250 mm
6762-2515	15/250 mm
6762-2516	16/250 mm
6762-2517***	17/250 mm
6762-2518***	18/250 mm
6762-2519***	19/250 mm
6762-2520***	20/250 mm
6762-2521***	21/250 mm
6762-2522***	22/250 mm

\*\*\* marked sizes of length have distal interlocking screw holes for  $\varnothing 4,5$  mm screws

## RS stem extra small cementless HA

Mat.: *implatan*®,  $TiAl_6V_4$  acc. to ISO 5832-3  
with *implaFix*® HA-coating acc. to ISO 13779-2

### REF

6764-1514HA	14/190 mm
6764-2014HA	14/240 mm

# MUTARS® RS



## RS stem cemented \*\*

Mat.: implavit®, CoCrMo acc. to ISO 5832-4


### REF

6760-1212	12/120 mm*
6760-1412	14/120 mm*
6760-1612	16/120 mm*
6760-1812	18/120 mm*
6760-1215	12/150 mm
6760-1415	14/150 mm
6760-1615	16/150 mm
6760-1815	18/150 mm
6761-1220	12/200 mm
6761-1420	14/200 mm
6761-1620	16/200 mm
6761-1820	18/200 mm


\* article available on request

\*\*article with TiN coating available (reference number + N)


# instruments




**MUTARS® RS broach container (bottom)**  
7999-6721




**MUTARS® RS broach container (top)**  
7999-6721




**MUTARS® RS ES container 2**  
7999-6715



**MUTARS® RS ES container 3**  
7999-6716



**MUTARS® RS trial stem container**  
7999-6724



**MUTARS® RS container 5**  
7999-6705




**flexible reamer nitinol\* container**  
7999-7000


**\*flexible reamer nitinol container**  
always to be ordered with guide rod nitinol  
7702-0001 (single use)




# instruments




**MUTARS® RS ES extracting container**  
7999-6717




**MUTARS® RS broach container 150mm (intermediate sizes)**  
7999-6722




**MUTARS® RS broach container 200mm (intermediate sizes)**  
7999-6723



**MUTARS® trial stem container 150mm + 200mm**  
7999-6727



**MUTARS® RS trial stem cemented container**  
7999-6725



**MUTARS® RS trial stem 120mm cemented container**  
7999-6726

## container without images:

**MUTARS® trial stem container 250mm (for intermediate sizes)**  
7999-6731



**hip cups**



## EcoFit® hip cup EPORE® cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub>

### REF

0220-0742	Ø 42 mm
0220-0744	Ø 44 mm
0220-0746	Ø 46 mm
0220-0748	Ø 48 mm
0220-0750	Ø 50 mm
0220-0752	Ø 52 mm
0220-0754	Ø 54 mm
0220-0756	Ø 56 mm
0220-0758	Ø 58 mm
0220-0760	Ø 60 mm
0220-0762	Ø 62 mm
0220-0764	Ø 64 mm
0220-0766	Ø 66 mm
0220-0768	Ø 68 mm
0220-0770	Ø 70 mm
0220-0772	Ø 72 mm



## EcoFit® hip cup EPORE® NH cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub>

### REF

0220-0942	Ø 42 mm
0220-0944	Ø 44 mm
0220-0946	Ø 46 mm
0220-0948	Ø 48 mm
0220-0950	Ø 50 mm
0220-0952	Ø 52 mm
0220-0954	Ø 54 mm
0220-0956	Ø 56 mm
0220-0958	Ø 58 mm
0220-0960	Ø 60 mm
0220-0962	Ø 62 mm
0220-0964	Ø 64 mm
0220-0966	Ø 66 mm
0220-0968	Ø 68 mm
0220-0970	Ø 70 mm
0220-0972	Ø 72 mm

## EcoFit® hip cup EPORE® Multi-Hole cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub>



### REF

0220-0642	Ø 42 mm
0220-0644	Ø 44 mm
0220-0646	Ø 46 mm
0220-0648	Ø 48 mm
0220-0650	Ø 50 mm
0220-0652	Ø 52 mm
0220-0654	Ø 54 mm
0220-0656	Ø 56 mm
0220-0658	Ø 58 mm
0220-0660	Ø 60 mm
0220-0662	Ø 62 mm
0220-0664	Ø 64 mm
0220-0666	Ø 66 mm
0220-0668	Ø 68 mm
0220-0670	Ø 70 mm
0220-0672	Ø 72 mm

## EcoFit® hip cup EPORE® / TCP Multi- Hole cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> with TCP coating



### REF

0220-6042	Ø 42 mm
0220-6044	Ø 44 mm
0220-6046	Ø 46 mm
0220-6048	Ø 48 mm
0220-6050	Ø 50 mm
0220-6052	Ø 52 mm
0220-6054	Ø 54 mm
0220-6056	Ø 56 mm
0220-6058	Ø 58 mm
0220-6060	Ø 60 mm
0220-6062	Ø 62 mm
0220-6064	Ø 64 mm
0220-6066	Ø 66 mm
0220-6068	Ø 68 mm
0220-6070	Ø 70 mm
0220-6072	Ø 72 mm

## EcoFit® hip cup EPORE® / TCP cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> with TCP coating



### REF

0220-2042	Ø 42 mm
0220-2044	Ø 44 mm
0220-2046	Ø 46 mm
0220-2048	Ø 48 mm
0220-2050	Ø 50 mm
0220-2052	Ø 52 mm
0220-2054	Ø 54 mm
0220-2056	Ø 56 mm
0220-2058	Ø 58 mm
0220-2060	Ø 60 mm
0220-2062	Ø 62 mm
0220-2064	Ø 64 mm
0220-2066	Ø 66 mm
0220-2068	Ø 68 mm
0220-2070	Ø 70 mm
0220-2072	Ø 72 mm

## EcoFit® hip cup cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® cpTi-coating



### REF

0220-0046	Ø 46 mm
0220-0048	Ø 48 mm
0220-0050	Ø 50 mm
0220-0052	Ø 52 mm
0220-0054	Ø 54 mm
0220-0056	Ø 56 mm
0220-0058	Ø 58 mm
0220-0060	Ø 60 mm
0220-0062	Ø 62 mm
0220-0064	Ø 64 mm
0220-0066	Ø 66 mm
0220-0068	Ø 68 mm



## EcoFit® hip cup HA cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® Duo, cpTi- and  
HA-coating acc. to ISO 13779-2

### REF

0220-0346	Ø 46 mm
0220-0348	Ø 48 mm
0220-0350	Ø 50 mm
0220-0352	Ø 52 mm
0220-0354	Ø 54 mm
0220-0356	Ø 56 mm
0220-0358	Ø 58 mm
0220-0360	Ø 60 mm
0220-0362	Ø 62 mm
0220-0364	Ø 64 mm
0220-0366	Ø 66 mm
0220-0368	Ø 68 mm



## EcoFit® NH hip cup cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® cpTi-coating

### REF

0220-0146	Ø 46 mm
0220-0148	Ø 48 mm
0220-0150	Ø 50 mm
0220-0152	Ø 52 mm
0220-0154	Ø 54 mm
0220-0156	Ø 56 mm
0220-0158	Ø 58 mm
0220-0160	Ø 60 mm
0220-0162	Ø 62 mm
0220-0164	Ø 64 mm
0220-0166	Ø 66 mm
0220-0168	Ø 68 mm



## EcoFit® NH hip cup HA cementless incl. central hole cover

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with implaFix® Duo, cpTi- and HA-coating acc. to ISO 13779-2

### REF

0220-0446	Ø 46 mm
0220-0448	Ø 48 mm
0220-0450	Ø 50 mm
0220-0452	Ø 52 mm
0220-0454	Ø 54 mm
0220-0456	Ø 56 mm
0220-0458	Ø 58 mm
0220-0460	Ø 60 mm
0220-0462	Ø 62 mm
0220-0464	Ø 64 mm
0220-0466	Ø 66 mm
0220-0468	Ø 68 mm



## EcoFit® SC cementless screw cup cementless

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3


### REF

0220-0246	Ø 46 mm
0220-0248	Ø 48 mm
0220-0250	Ø 50 mm
0220-0252	Ø 52 mm
0220-0254	Ø 54 mm
0220-0256	Ø 56 mm
0220-0258	Ø 58 mm
0220-0260	Ø 60 mm
0220-0262	Ø 62 mm
0220-0264	Ø 64 mm
0220-0266	Ø 66 mm
0220-0268	Ø 68 mm


# instruments




**acetabular reamer container**  
42-64mm in 1mm steps  
2950-1048




**GIS® retractor container**  
7999-7514




**EcoFit® cup container 36mm**  
(top)  
0220-1036




**EcoFit® cup container 36mm**  
(bottom)  
0220-1036




**EcoFit® cup GIS® container 36mm**  
(top)  
0220-1136



**EcoFit® cup GIS® container 36mm**  
(bottom)  
0220-1136




**trial container for ceramic insert**  
40mm  
0220-0040




**GIS® acetabular reamer container**  
with curved reamer handle  
7999-0800




# instruments




**EcoFit® cup GIS® container  
(top)**  
0220-2136




**EcoFit® cup GIS® container  
(bottom)**  
0220-2136



**EcoFit® NH cup container 36mm  
(top)**  
0220-1236




**EcoFit® NH cup container 36mm  
(bottom)**  
0220-1236




**EcoFit® cup EPORE® container 36mm  
(top)**  
0220-0336



**EcoFit® cup EPORE® container 36mm  
(bottom)**  
0220-0336




**EcoFit® cup EPORE® container  
(top)**  
0220-0036




**EcoFit® cup EPORE® container  
(bottom)**  
0220-0036


# instruments




**EcoFit® cup EPORE® container GIS®**  
**(top)**  
0220-0236




**EcoFit® cup EPORE® container GIS®**  
**(bottom)**  
0220-0236




**EcoFit® cup EPORE® container GIS®**  
**(top)**  
0220-0436




**EcoFit® cup EPORE® container GIS®**  
**(bottom)**  
0220-0436




**EcoFit® cup EPORE® container**  
**42-44mm**  
0220-0136



**EcoFit® cup EPORE® container**  
**70-72mm**  
0220-0138




**extension container GIS®**  
0220-1030




**container for EcoFit® cup insert**  
0220-2020


# instruments




**EcoFit® SC container GIS®**  
0220-1039



**EcoFit® SC container 36mm**  
0220-1037



**EcoFit® cup EPORE® NH container  
(top)**  
0220-0037



**EcoFit® cup EPORE® NH container  
(bottom)**  
0220-0037

# PE cup insert




## PE cup insert 0°

Mat.: UHMWPE acc. to ISO 5834-2

### REF

0280-2039	Ø 28/39 mm
0280-2044	Ø 28/44 mm
0280-2048	Ø 28/48 mm
0280-2052	Ø 28/52 mm
0280-3444	Ø 32/44 mm
0280-3448	Ø 32/48 mm
0280-4448	Ø 36/48 mm
0280-3452	Ø 32/52 mm
0280-4452	Ø 36/52 mm



## PE cup insert 10°

Mat.: UHMWPE acc. to ISO 5834-2

### REF

0280-2139	Ø 28/39 mm
0280-2144	Ø 28/44 mm
0280-2148	Ø 28/48 mm
0280-2152	Ø 28/52 mm
0280-3144	Ø 32/44 mm
0280-3148	Ø 32/48 mm
0280-4548	Ø 36/48 mm
0280-3152	Ø 32/52 mm
0280-4552	Ø 36/52 mm



## implacross®


## PE cup insert 0°

Mat.: crosslinked UHMWPE

### REF

0223-2839	Ø 28/39 mm
0223-2844	Ø 28/44 mm
0223-2848	Ø 28/48 mm
0223-2852	Ø 28/52 mm
0223-3239	Ø 32/39 mm
0223-3244	Ø 32/44 mm
0223-3644	Ø 36/44 mm
0223-3248	Ø 32/48 mm
0223-3648	Ø 36/48 mm
0223-3252	Ø 32/52 mm
0223-3652	Ø 36/52 mm

# PE cup insert



## implacross® PE cup insert 10°

*Mat.: crosslinked UHMWPE*

### REF

0224-2835	Ø 28/35 mm
0224-2839	Ø 28/39 mm
0224-2844	Ø 28/44 mm
0224-2848	Ø 28/48 mm
0224-2852	Ø 28/52 mm
0224-3239	Ø 32/39 mm
0224-3244	Ø 32/44 mm
0224-3644	Ø 36/44 mm
0224-3248	Ø 32/48 mm
0224-3648	Ø 36/48 mm
0224-3252	Ø 32/52 mm
0224-3652	Ø 36/52 mm



## implacross® E PE cup insert 10°

*Mat.: implacross® E, crosslinked UHMWPE with Vitamin E*

### REF

0226-2835	Ø 28/35 mm
0234-3239	Ø 32/39 mm
0234-3244	Ø 32/44 mm
0234-3644	Ø 36/44 mm
0234-3248	Ø 32/48 mm
0234-3648	Ø 36/48 mm
0234-3252	Ø 32/52 mm
0234-3652	Ø 36/52 mm



## implacross® E PE cup insert 0°

*Mat.: implacross® E, crosslinked UHMWPE with Vitamin E*

### REF

0233-3239	Ø 32/39 mm
0233-3244	Ø 32/44 mm
0233-3644	Ø 36/44 mm
0233-3248	Ø 32/48 mm
0233-3648	Ø 36/48 mm
0233-3252	Ø 32/52 mm
0233-3652	Ø 36/52 mm

# cup insert



## BIOLOX® delta cup-insert

Mat.:  $Al_2O_3$  and  $ZrO_3$  acc. to ISO 6474-2

### REF

0220-2835	Ø 28/35 mm
0220-3239	Ø 32/39 mm
0220-3244	Ø 32/44 mm
0220-3644	Ø 36/44 mm
0220-3248	Ø 32/48 mm
0220-3648	Ø 36/48 mm
0220-3252	Ø 32/52 mm
0220-3652	Ø 36/52 mm
0220-4048	Ø 40/48 mm
0220-4052	Ø 40/52 mm

## EcoFit® insert

Mat.: *implavit*® CoCrMo acc. to ISO 5832-12

### REF

0223-3844	Ø 38/44 mm
0223-4248	Ø 42/48 mm
0223-4652	Ø 46/52 mm



## EcoFit® insert TiN\*

Mat.: *implavit*® CoCrMo acc. to ISO 5832-12 with TiN coating

### REF

0220-3844	Ø 38/44 mm
0220-4248	Ø 42/48 mm
0220-4652	Ø 46/52 mm

\* available on request

## EcoFit® 2M hip cup, cementless

Mat.: implavit®, CoCrMo acc. to ISO 5832-4  
with implaFix® Duo, cpTi- and  
HA-coating acc. to ISO 13779-2



### REF

0220-1042	38/42 mm
0220-1044	40/44 mm
0220-1046	42/46 mm
0220-1048	44/48 mm
0220-1050	46/50 mm
0220-1052	48/52 mm
0220-1054	50/54 mm
0220-1056	52/56 mm
0220-1058	54/58 mm
0220-1060	56/60 mm
0220-1062	58/62 mm
0220-1064	58/64 mm

## EcoFit® 2M hip cup, cemented

Mat.: implavit®, CoCrMo acc. to ISO 5832-4



### REF

0220-1144	38/44 mm
0220-1146	40/46 mm
0220-1148	42/48 mm
0220-1150	44/50 mm
0220-1152	46/52 mm
0220-1154	48/54 mm
0220-1156	50/56 mm
0220-1158	52/58 mm
0220-1160	54/60 mm
0220-1162	56/62 mm
0220-1164	58/64 mm

## EcoFit® 2M hip cup, cemented TiN

Mat.: implavit®, CoCrMo acc. to ISO 5832-4 with  
TiN-coating



### REF

0220-1144N	38/44 mm
0220-1146N	40/46 mm
0220-1148N	42/48 mm
0220-1150N	44/50 mm
0220-1152N	46/52 mm
0220-1154N	48/54 mm
0220-1156N	50/56 mm
0220-1158N	52/58 mm
0220-1160N	54/60 mm
0220-1162N	56/62 mm
0220-1164N	58/64 mm

## EcoFit® 2M hip cup, cementless SP

Mat.: implavit®, CoCrMo acc. to ISO 5832-4  
with implaFix® Duo, cpTi- and  
HA-coating acc. to ISO 13779-2



### REF

0220-1244	40/44 mm
0220-1246	42/46 mm
0220-1248	44/48 mm
0220-1250	46/50 mm
0220-1252	48/52 mm
0220-1254	50/54 mm
0220-1256	52/56 mm
0220-1258	54/58 mm
0220-1260	56/60 mm
0220-1262	58/62 mm





## EcoFit® 2M head

Mat.: UHMWPE acc. to ISO 5834-2

### REF

2906-2238	22/38 mm
2906-2240	22/40 mm
2906-2842	28/42 mm
2906-2844	28/44 mm
2906-2846	28/46 mm
2906-2848	28/48 mm
2906-2850	28/50 mm
2906-2852	28/52 mm
2906-2854	28/54 mm
2906-2856	28/56 mm
2906-2858	28/58 mm
2906-3248	32/48 mm
2906-3250	32/50 mm
2906-3252	32/52 mm
2906-3254	32/54 mm
2906-3256	32/56 mm
2906-3258	32/58 mm




## 2M implacross® E head

Mat.: crosslinked UHMWPE with Vitamin E


### REF

2905-2238	22/38 mm
2905-2240	22/40 mm
2905-2842	28/42 mm
2905-2844	28/44 mm
2905-2846	28/46 mm
2905-2848	28/48 mm
2905-2850	28/50 mm
2905-2852	28/52 mm
2905-2854	28/54 mm
2905-2856	28/56 mm
2905-2858	28/58 mm
2905-3248	32/48 mm
2905-3250	32/50 mm
2905-3252	32/52 mm
2905-3254	32/54 mm
2905-3256	32/56 mm
2905-3258	32/58 mm


# instruments




**EcoFit® 2M container (top)**  
0220-1081




**EcoFit® 2M container (bottom)**  
0220-1081




**EcoFit® 2M container II**  
0220-1082



**EcoFit® 2M cup SP basic container (top)**  
0220-1086



**EcoFit® 2M cup SP basic container (bottom)**  
0220-1086



**EcoFit® 2M cup SP container**  
0220-1085

# MUTARS® RS cup / PRS




## MUTARS® RS cup

Mat.: *implatan*®,  $TiAl_6V_4$

### REF

5712-0546	Ø 46 mm left
5712-0550	Ø 50 mm left
5712-0554	Ø 54 mm left
5712-0558	Ø 58 mm left
5712-0562	Ø 62 mm left

5712-0046	Ø 46 mm right
5712-0050	Ø 50 mm right
5712-0054	Ø 54 mm right
5712-0058	Ø 58 mm right
5712-0062	Ø 62 mm right



## MUTARS® PRS incl. central hole cover

Mat.: *implatan*®,  $TiAl_6V_4$

### REF

0296-0052	52 mm
0296-0056	56 mm
0296-0060	60 mm
0296-0064	64 mm
0296-0068	68 mm
0296-0072	72 mm

## central hole cover threaded M10x1


Mat.: *implatan*®,  $TiAl_6V_4$

### REF

0280-3100	M10x1
-----------	-------



# MUTARS® RS cup



## implacross® PE-insert 15° neutral 0mm

Mat.: crosslinked UHMWPE

### REF

0227-3239	Ø 32/39 mm
0227-3644	Ø 36/44 mm
0227-3648	Ø 36/48 mm
0227-3652	Ø 36/52 mm



## implacross® PE-insert 15° Offset 4mm

Mat.: crosslinked UHMWPE

### REF

0228-3239	Ø 32/39 mm
0228-3644	Ø 36/44 mm
0228-3648	Ø 36/48 mm
0228-3652	Ø 36/52 mm



## central plug M16x1

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3

### REF

0220-0216	M16 x 1
-----------	---------




## 2M Inlay 15° for MUTARS® RS cup and LUMiC® TiN

Mat.: implavit®, CoCrMo acc. to ISO 5832-4  
with TiN coating / UHMWPE acc. to ISO 5834-2

### REF

0242-3839	Ø 38/39 mm
0242-4244	Ø 42/44 mm
0242-4448	Ø 44/48 mm
0242-4652	Ø 46/52 mm



## fastening bolt

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3

### REF

5712-0000
-----------

# ic-reinforcement cage/acetabular ring



## ic-acetabular ring

Mat.: pure titanium acc. to ISO 5832-2 Grade 1

### REF

7000-1242	42 mm
7000-1244	44 mm
7000-1246	46 mm
7100-1248	48 mm
7000-1250	50 mm
7000-1252	52 mm
7000-1254	54 mm
7000-1256	56 mm
7000-1258	58 mm




## ic-reinforcement cage

Mat.: pure titanium acc. to ISO 5832-2 Grade 1

### REF

7000-0244	44 mm right
7000-0250	50 mm right
7000-0256	56 mm right
7000-0262	62 mm right
7001-0244	44 mm left
7001-0250	50 mm left
7001-0256	56 mm left
7001-0262	62 mm left

# instruments



ic-reinforcement container  
7999-7002



## MUTARS® LUMiC® cup incl. counter screw

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3

### REF

5711-0050	50mm
5711-0054	54mm
5711-0160	60mm
5711-0060S	60mm silver
5711-0250	50mm HA
5711-0254	54mm HA
5711-0260	60mm HA



## MUTARS® LUMiC® stem cementless HA

Mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
implaFix® HA-coating acc. to ISO 13779-2

### REF

5711-1865	8x65 mm
5711-1875	8x75 mm
5711-1885	8x85 mm
5711-1065	10x65 mm
5711-1075	10x75 mm
5711-1085	10x85 mm



## MUTARS® LUMiC® stem cemented

Mat.: implavit®; CoCrMo acc. to ISO 5832-4

### REF

5711-2865	8x65 mm
5711-2875	8x75 mm
5711-2885	8x85 mm




## MUTARS® LUMiC® screw

Mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3


### REF

5711-1002	M6x28
-----------	-------


# instruments



**MUTARS® LUMiC® container 1 (top)**  
7999-5701




**MUTARS® LUMiC® container 1 (bottom)**  
7999-5701




**MUTARS® LUMiC® container 2**  
7999-5702



**trial insert 15° container MUTARS® RS cup / LUMiC®**  
2950-1067




**MUTARS® RS cup container (top)**  
7999-7712



**MUTARS® RS cup container (bottom)**  
7999-7712



**MUTARS® RS cup trial cup container**  
7999-7713 left  
7999-7713 right




**MUTARS® RS cup PE trial insert container**  
7999-7715


# instruments




**MUTARS® LUMiC® container 2M trials**  
7999-5704




**MUTARS® RS cup container 2M trials**  
7999-7717




**container MUTARS® PRS (top)**  
0296-0001




**container MUTARS® PRS (bottom)**  
0296-0001




**acetabular reamer container  
MUTARS® PRS (top)**  
0296-0002



**acetabular reamer container  
MUTARS® PRS (bottom)**  
0296-0002



**container extracting instrument for  
2M insert 15° (top)**  
2951-3002



**container for extracting instrument  
for 2M insert 15° (bottom)**  
2951-3002



# PE-cup Mueller II



## PE-cup Mueller II cemented

*Mat.: UHMWPE acc. to ISO 5834-2; Stainless  
Steel X-ray wire acc. to ISO 5832-1*

### REF

1021-2240	22/40 mm
1021-2242	22/42 mm
1021-2244	22/44 mm
1021-2246	22/46 mm
1021-2248	22/48 mm
1021-2844	28/44 mm
1021-2846	28/46 mm
1021-2848	28/48 mm
1021-2850	28/50 mm
1021-2852	28/52 mm
1021-2854	28/54 mm
1021-2856	28/56 mm
1021-2858	28/58 mm
1021-2860	28/60 mm
1021-2862	28/62 mm
1021-2864	28/64 mm
1021-3244	32/44 mm
1021-3246	32/46 mm
1021-3248	32/48 mm
1021-3250	32/50 mm
1021-3252	32/52 mm
1021-3254	32/54 mm
1021-3256	32/56 mm
1021-3258	32/58 mm
1021-3260	32/60 mm
1021-3262	32/62 mm
1021-3264	32/64 mm



## PE-cup Mueller II snap cemented

*Mat.: UHMWPE acc. to ISO 5834-2; Stainless  
Steel X-ray wire acc. to ISO 5832-1*

### REF

1041-3244	32/44 mm
1041-3246	32/46 mm
1041-3248	32/48 mm
1041-3250	32/50 mm
1041-3252	32/52 mm
1041-3254	32/54 mm
1041-3256	32/56 mm
1041-3258	32/58 mm
1041-3260	32/60 mm
1041-3262	32/62 mm
1041-3264	32/64 mm

# PE-cup Mueller II



## PE-cup Mueller II dysplasia 10° cemented

Mat.: UHMWPE acc. to ISO 5834-2; Stainless Steel X-ray wire acc. to ISO 5832-1

### REF

1071-2844	28/44 mm
1071-2846	28/46 mm
1071-2848	28/48 mm
1071-2850	28/50 mm
1071-2852	28/52 mm
1071-2854	28/54 mm
1071-2856	28/56 mm
1071-2858	28/58 mm
1071-2860	28/60 mm
1071-2862	28/62 mm
1071-3244	32/44 mm
1071-3246	32/46 mm
1071-3248	32/48 mm
1071-3250	32/50 mm
1071-3252	32/52 mm
1071-3254	32/54 mm
1071-3256	32/56 mm
1071-3258	32/58 mm
1071-3260	32/60 mm
1071-3262	32/62 mm
1071-3264	32/64 mm

## PE-cup Mueller II implacross® E cemented


Mat.: crosslinked UHMWPE acc. to ISO 5834-2; Stainless Steel X-ray wire acc. to ISO 5832-1

### REF


1010-3244	32/44 mm
1010-3246	32/46 mm
1010-3248	32/48 mm
1010-3250	32/50 mm
1010-3252	32/52 mm
1010-3254	32/54 mm
1010-3256	32/56 mm
1010-3258	32/58 mm
1010-3260	32/60 mm
1010-3262	32/62 mm
1010-3264	32/64 mm
1010-3648	36/48 mm
1010-3650	36/50 mm
1010-3652	36/52 mm
1010-3654	36/54 mm
1010-3656	36/56 mm
1010-3658	36/58 mm
1010-3660	36/60 mm
1010-3662	36/62 mm
1010-3664	36/64 mm




# instruments




**acetabulum reamer container**  
0282-0001




**Mueller PE cup container 28/32mm  
incl. acetabulum reamer**  
7960-0273




**Mueller PE cup container**  
7960-0270



**Mueller PE cup container 32/36mm**  
7960-0272




**Mueller PE cup container I**  
7960-0274




**Mueller PE cup container 32/36mm II**  
7960-0276

# instruments



**Mueller PE cup container 28/32mm  
incl. acetabulum reamer III (top)**  
7960-0277



**Mueller PE cup container 28/32mm  
incl. acetabulum reamer III (bottom)**  
7960-0277

# ic-bipolar head

## ic-Duokopf CoCrMo

Mat.: outer shell: implavit®, CoCrMo acc. to ISO 5832-4, PE-insert and PE-locking ring: UHMWPE acc. to ISO 5834-2




### REF


2151-2238	22/38 mm
2151-2239	22/39 mm*
2151-2240	22/40 mm
2151-2241	22/41 mm*
2151-0042	28/42 mm
2151-0043	28/43 mm*
2151-0044	28/44 mm
2151-0045	28/45 mm*
2151-0046	28/46 mm
2151-0047	28/47 mm*
2151-0048	28/48 mm
2151-0049	28/49 mm*
2151-0050	28/50 mm
2151-0051	28/51 mm*
2151-0052	28/52 mm
2151-0053	28/53 mm*
2151-0054	28/54 mm
2151-0055	28/55 mm*
2151-0056	28/56 mm
2151-0057	28/57 mm*
2151-0058	28/58 mm
2151-0059	28/59 mm*
2151-0060	28/60 mm
2151-0062	28/62 mm

\*implants are only available on request.

# instruments



**ic-bipolar head  
container**  
7960-9999



**ic-bipolar head container 44-62mm  
in 1mm steps**  
7960-9997



**adapter for bipolar head sizing shell  
container 28mm S-XL**  
7966-0001

## instruments/container without image:

**ic-bipolar head container 38-42mm**  
7960-9998

**ic-bipolar head container 38-62mm in 1mm steps**  
7960-9991

**femoral head template 42-62mm in 2mm steps**  
7960-7000

**femoral head template 38-43mm in 1mm steps**  
7960-7001

**femoral head template 44-49mm in 1mm steps**  
7960-7002

**femoral head template 50-55mm in 1mm steps**  
7960-7003

**femoral head template 56-62mm in 1mm steps**  
7960-7004

# screws



## spongiosa screw flat head Ø 6.5mm

Mat.: *implatan*<sup>®</sup>, *TiAl<sub>6</sub>V<sub>4</sub>* acc. to ISO 5832-3

### REF

0280-1015	length 15 mm
0280-1020	length 20 mm
0280-1025	length 25 mm
0280-1030	length 30 mm
0280-1035	length 35 mm
0280-1040	length 40 mm
0280-1045	length 45 mm
0280-1050	length 50 mm
0280-1055	length 55 mm
0280-1060	length 60 mm
0280-1065	length 65 mm
0280-1070	length 70 mm
0280-1075	length 75 mm
0280-1080	length 80 mm




## cancellous screw Ø 8mm

Mat.: *implatan*<sup>®</sup>, *TiAl<sub>6</sub>V<sub>4</sub>* acc. to ISO 5832-3

### REF

0287-0070	length 70 mm
0287-0080	length 80 mm
0287-0090	length 90 mm
0287-0100	length 100 mm



## shim for EPORE<sup>®</sup> acetabular spacer 8mm

Mat.: *implatan*<sup>®</sup>, *TiAl<sub>6</sub>V<sub>4</sub>* acc. to ISO 5832-3

### REF

0296-0006



## shim for spongiosa screw flat head 6.5mm

Mat.: *implatan*<sup>®</sup>, *TiAl<sub>6</sub>V<sub>4</sub>* acc. to ISO 5832-3

### REF

0296-0000

# screws




## cortical screw Ø 4.5mm

Mat.: implatan®, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3


### REF

5792-4525	length 25 mm
5792-4530	length 30 mm
5792-4535	length 35 mm
5792-4540	length 40 mm
5792-4545	length 45 mm
5792-4550	length 50 mm
5792-4555	length 55 mm
5792-4560	length 60 mm

# instruments



container screw instruments  
7999-7001



container screw instruments  
modular for MUTARS® PRS  
0296-0003



# ic-heads

## ic-head CoCrMo taper 12/14mm

Mat.: *implavit*<sup>®</sup>, CoCrMo acc. to ISO 5832-12  
taper 12/14



### REF

2312-2200	22 mm, S
2312-2205	22 mm, M
2312-2210	22 mm, L
2387-2800	28 mm, S
2387-2805	28 mm, M
2387-2810	28 mm, L
2387-2815	28 mm, XL
2387-2820	28 mm, XXL
2387-2825	28 mm, XXXL
2387-3200	32 mm, S
2387-3205	32 mm, M
2387-3210	32 mm, L
2387-3215	32 mm, XL
2387-3220	32 mm, XXL
2387-3225	32 mm, XXXL
2387-3600	36 mm, S
2387-3605	36 mm, M
2387-3610	36 mm, L
2387-3615	36 mm, XL
2387-3620	36 mm, XXL
2387-3625	36 mm, XXXL

## ic-head titanium taper 12/14mm

Mat.: *implatan*<sup>®</sup>, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3  
with TiN-coating; taper 12/14



### REF

2787-2800	28 mm, S
2787-2805	28 mm, M
2787-2810	28 mm, L
2787-2815	28 mm, XL
2787-2820	28 mm, XXL
2787-2825	28 mm, XXXL
2787-3200	32 mm, S
2787-3205	32 mm, M
2787-3210	32 mm, L
2787-3215	32 mm, XL
2787-3220	32 mm, XXL
2787-3225	32 mm, XXXL
2787-3600	36 mm, S
2787-3605	36 mm, M
2787-3610	36 mm, L
2787-3615	36 mm, XL
2787-3620	36 mm, XXL
2787-3625	36 mm, XXXL



## ic-head CoCrMo TiN taper 12/14mm

Mat.: *implavit*<sup>®</sup>, CoCrMo acc. to ISO 5832-12 with TiN-coating taper 12/14

### REF

2322-2200	22 mm, S
2322-2205	22 mm, M
2322-2210	22 mm, L

## ic-head BIOLOX<sup>®</sup> delta taper 12/14mm

Mat.:  $Al_2O_3$  and  $ZrO_3$  acc. to ISO 6474-2 taper 12/14

### REF

2586-2800	28 mm, S
2586-2805	28 mm, M
2586-2810	28 mm, L
2586-3200	32 mm, S
2586-3205	32 mm, M
2586-3210	32 mm, L
2586-3215	32 mm, XL
2586-3600	36 mm, S
2586-3605	36 mm, M
2586-3610	36 mm, L
2586-3615	36 mm, XL
2586-4000	40 mm, S
2586-4005	40 mm, M
2586-4010	40 mm, L
2586-4015	40 mm, XL



## ic-head BIOLOX<sup>®</sup> delta taper 12/14mm\*

Mat.:  $Al_2O_3$  and  $ZrO_3$  acc. to ISO 6474-2

### REF

2589-2800	28 mm
2589-3200	32 mm
2589-3600	36 mm
2589-4000	40 mm



## taper adapter for ic-head revision BIOLOX<sup>®</sup> delta 12/14


Mat.:  $TiAl_6V_4$  acc. to ISO 5832-3 and ASTM F138

### REF


2588-0000	S
2588-0005	M
2588-0010	L
2588-0015	XL




# instruments




**trial head 28mm, 32mm S-XL  
container**  
7999-7960




**trial head container 28mm+32mm  
S-XXL**  
7999-7961




**trial head 28mm+32mm+36mm S-XL  
container**  
7999-7963



**trial head container 22/28/32mm, S-XL**  
7999-7967



**trial head container 28mm, 32mm,  
36mm S-XXXL**  
7999-7964



**trial head container 36mm S-XL**  
7999-7969

\* Only in combination with taper adapter for ic-head revision BIOLOX® delta.

# EPORE<sup>®</sup> acetabular spacer



## EPORE<sup>®</sup> acetabular spacer

Mat.: implatan<sup>®</sup>, TiAl<sub>6</sub>V<sub>4</sub>

### REF

02955410	54/10 mm
02955415	54/15 mm
02955420	54/20 mm
02955810	58/10 mm
02955815	58/15 mm
02955820	58/20 mm
02956210	62/10 mm
02956215	62/15 mm
02956220	62/20 mm
02956610	66/10 mm
02956615	66/15 mm
02956620	66/20 mm

## EPORE<sup>®</sup> acetabular spacer 8mm

Mat.: implatan<sup>®</sup>, TiAl<sub>6</sub>V<sub>4</sub>

### REF

0296-4610	46/10 mm
0296-4615	46/15 mm
0296-4620	46/20 mm
0296-5010	50/10 mm
0296-5015	50/15 mm
0296-5020	50/20 mm
0296-5410	54/10 mm
0296-5415	54/15 mm
0296-5420	54/20 mm
0296-5430	54/30 mm*
0296-5810	58/10 mm
0296-5815	58/15 mm
0296-5820	58/20 mm
0296-5830	58/30 mm*
0296-6210	62/10 mm
0296-6215	62/15 mm
0296-6220	62/20 mm
0296-6230	62/30 mm*
0296-6610	66/10 mm
0296-6615	66/15 mm
0296-6620	66/20 mm
0296-6630	66/30 mm*
0296-7010	70/10 mm
0296-7015	70/15 mm
0296-7020	70/20 mm
0296-7030	70/30 mm*



\* article only available on request

# EPORE<sup>®</sup> Buttress/Shims



## EPORE<sup>®</sup> Buttress Augment

Mat.: implatan<sup>®</sup>, TiAl<sub>6</sub>V<sub>4</sub>

### REF

0296-0100	left posterior/right anterior
0296-0101	right posterior/left anterior
0296-0154	straight 54mm
0296-0158	straight 58mm



## EPORE<sup>®</sup> Shim Augment

Mat.: implatan<sup>®</sup>, TiAl<sub>6</sub>V<sub>4</sub>

### REF

0296-1005	5°
0296-1010	10°
0296-1015	15°

# trochanter plate




## trochanter plate

Mat.: *implatan*<sup>®</sup>,  $TiAl_6V_4$


### REF

84000120	120mm
84000160	160mm
84000200	200mm
84000240	240mm

# instruments



**container trochanter plate**  
7999-7004




**EPORE<sup>®</sup> acetabulum spacer container**  
7295-2000

## container without image

**EPORE<sup>®</sup> acetabular spacer container 8mm**  
7295-2060

**EPORE<sup>®</sup> Buttress and Shims Augment container**  
0297-0002



## intramedullary plug

Mat.: UHMWPE acc. to ISO 5834-2; X-ray wire  
TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3

### REF

0299-4000 small medullary cavity ≥ Ø: 9 mm  
0299-4010 large medullary cavity ≥ Ø: 14 mm



## CEMSTOP / ic-cement restrictor resorbable

Mat.: Gelatine, Glycerol, water und  
Methylhydroxybenzoate

### REF

T770008 / 0299-3008 Ø 8mm  
T770010 / 0299-3010 Ø 10mm  
T770012 / 0299-3012 Ø 12mm  
T770014 / 0299-3014 Ø 14mm  
T770016 / 0299-3016 Ø 16mm  
T770018 / 0299-3018 Ø 18mm



## ic-cerclage

Mat.: pure titanium acc. to ISO 5832-2 Grade 1

### REF

0060-1008 8mm



## ic-cerclage band titanium 2R

Mat.: pure titanium acc. to ISO 5832-2 Grade 1

### REF

0060-1018 8mm




## ic-titanium cerclage wire

Mat.: pure titanium acc. to ISO 5832-2 Grade 1


### REF

0060-2014 1.4mm x 5m  
0060-2013 1.4mm x 50cm

# instruments



**ic-cerclage  
container**  
7999-0060



**ic-cement restrictor resorbable  
container**  
0299-3000



**handle for intramedullary plug**  
7512-4001



# accessories



## implabond bone cement

### REF

D041140	standard viscosity 40g
D041140G	standard viscosity 40g with Gentamicine
D041340	low viscosity 40g
D041340G	low viscosity 40g with Gentamicine



## Cempadic® spacer mold for EcoFit® stem cemented

### REF

0298-0841	for EcoFit® stem 6,25 mm
0298-0842	for EcoFit® stem 7,5 mm
0298-0843	for EcoFit® stem 10 mm
0298-0844	for EcoFit® stem 12,5 mm
0298-0845	for EcoFit® stem 15 mm
0298-0846	for EcoFit® stem 17,5 mm

## Cempadic® acetabulum cement spacer mold

### REF

0298-0944	for 44 mm
0298-0946	for 46 mm
0298-0948	for 48 mm
0298-0850	for 50 mm
0298-0852	for 52 mm
0298-0854	for 54 mm
0298-0856	for 56 mm
0298-0858	for 58 mm
0298-0960	for 60 mm
0298-0962	for 62 mm


## Cempadic® spacer rod EcoFit®

Mat.: implan®<sub>6</sub>, TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-2

### REF

0298-0861	for EcoFit® stem 6,25-10 mm
0298-0862	for EcoFit® stem 12,5-17,5 mm

# accessories




## ic-Lavage

### REF

- 0067-8000 ic-Lavage set with nozzles and splash shield  
1bo = 10 pieces
- 0067-8002 ic-Lavage set femoral brush with suction  
1bo = 5 pieces

# extractor adapter




**ic-extractor adapter  
container**  
6500-0020

## container without images

**ic-extractor adapter container incl. slide hammer with snap mechanism**  
6500-0010

## content:



# broach handles



description	REF	description
Aida® broach handle easy lock double offset left	7512-8014 / 7512-8015	Aida® broach handle easy lock double offset right
Cortina broach handle easy lock double offset left	7200-2011 / 7200-2012	Cortina broach handle easy lock double offset right
Actinia® broach handle easy lock double offset left	8004-9047 / 8004-9046	Actinia® broach handle easy lock double offset right
BethaLoc® broach handle easy lock double offset left	8700-1171 / 8700-1172	BethaLoc® broach handle easy lock double offset right
ic-straight stem broach handle easy lock double offset left	7512-0261 / 7512-0260	ic-straight stem broach handle easy lock double offset right
EcoFit® broach handle with double offset left	7512-0073 / 7512-0074	EcoFit® broach handle with double offset right



description	REF	description
Aida® broach handle GIS® easy lock left	7512-8004 / 7512-8005	Aida® broach handle GIS® easy lock right
Cortina broach handle GIS® easy lock left	7200-2013 / 7200-2014	Cortina broach handle GIS® easy lock right
Actinia® broach handle GIS® easy lock left	8004-9135 / 8004-9134	Actinia® broach handle GIS® easy lock right
BethaLoc® broach handle GIS® easy lock left	8700-1175 / 8700-1180	BethaLoc® broach handle GIS® easy lock right
ic-straight stem broach handle GIS® easy lock left	7512-0263 / 7512-0262	ic-straight stem broach handle GIS® easy lock right
EcoFit® broach handle GIS® easy lock left	7512-0067 / 7512-0068	EcoFit® broach handle GIS® easy lock right

# broach handles



REF	description
7512-0112	EcoFit® broach handle easy lock
7512-0002	Aida® broach handle easy lock
7200-2010	Cortina broach handle easy lock
8004-9044	Actinia® broach handle easy lock
8700-1077	BethaLoc® broach handle easy lock
7512-0259	ic-straight stem broach handle easy lock



REF	description
7512-0085	DAA broach handle
7512-0084	DAA broach handle (bevelled striker plate)

# stem extractor



REF	description
7512-2033	stem extractor revision
7512-203302	stem extractor revision part 2



REF	description
7512-2034	stem extractor




REF	description
6500-0012	slide hammer with snap mechanism
6500-0008	MUTARS® RS broach handle (without image)

# alignment guide



REF	description
7512-2203	alignment guide for curved handle



REF	description
7512-2204	alignment guide for straight handle

# stem impactor



REF	description
7512-2022	EcoFit® stem impactor unguided curved
8004-9037	Actinia® stem impactor unguided curved (image similar)

REF	description
3039-0109	EcoFit® stem impactor (unguided)
8004-9031	stem impactor (image similar)



# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---





# implantcast

implantcast GmbH  
Lueneburger Schanze 26  
D-21614 Buxtehude  
Germany

phone: +49 4161 744-0  
fax: +49 4161 744-200  
e-mail: [info@implantcast.de](mailto:info@implantcast.de)  
internet: [www.implantcast.de](http://www.implantcast.de)



your local distributor:

HIPOVERE-291221

